



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6044 Cornerstone Ct. West, Ste. E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet

HSP90 α (C-Terminal) Inhibitor Screening Assay Kit

Catalog # 50317

Size: 384 reactions

DESCRIPTION: The *HSP90 α (C-Terminal) Inhibitor Screening Assay Kit* is designed to measure the inhibition of HSP90 α binding to its protein target PPID (also known as Cyclophilin D). The assay kit comes in a convenient AlphaLISA[®] format, with enough HSP90 α (535-732), assay buffer, detection buffer, and purified GST-tagged PPID to perform a total of 384 enzyme reactions. With this kit, only three simple steps on a microtiter plate are required. First, a sample containing HSP90 α , PPID, and an inhibitor of choice is incubated for thirty minutes. Next, acceptor beads are added, then donor beads, followed by reading the Alpha-counts.

COMPONENTS:

Catalog #	Component	Amount	Storage	
50316	HSP90 α (C-terminal), Biotin-labeled	50 μ g	-80 $^{\circ}$ C	(Avoid freeze/thaw cycles!)
71095	PPID, GST-tag	100 μ g	-80 $^{\circ}$ C	
50324	3x HSP90 Assay Buffer 2	4 ml	-20 $^{\circ}$ C	
	3x HSP90 Detection Buffer	3 ml	-20 $^{\circ}$ C	

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

Glutathione AlphaLISA[®] Acceptor Beads, 5 mg/ml (PerkinElmer #AL109C)
AlphaScreen[®] Streptavidin-conjugated donor beads, 5 mg/ml (PerkinElmer #6760002S)
Optiplate-384 (PerkinElmer #6007290)
AlphaScreen[®] microplate reader
Adjustable micropipettor and sterile tips

APPLICATIONS: Useful for screening inhibitors of HSP90 α and for HSP90 α binding assays

CONTRAINDICATIONS: Only limited amounts of DMSO (< 1%) can be included as it has been shown to disrupt HSP90 α :PPID interaction. Green and blue dyes that absorb light in the AlphaScreen[®] signal emission range (520-620 nm), such as Trypan Blue. Avoid the use of the potent singlet oxygen quenchers such as sodium azide (NaN₃) or metal ions (Fe²⁺, Fe³⁺, Cu²⁺, Zn²⁺ and Ni²⁺). The presence of >1% RPMI 1640 culture medium leads to a signal reduction due to the presence of excess biotin and iron in this medium. MEM, which lacks these components, does not affect AlphaScreen[®] assays.

STABILITY: At least one year from date of receipt when stored as directed.

REFERENCE(S): Allan, R.K. *et al. J. Biol.Chem* 2006 **281(11)**: 7161-71.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082**, Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com

ASSAY PROTOCOL:

All samples and controls should be tested in duplicate. Use slow shaking for all incubations.

Step 1:

- 1) Thaw **HSP90 α** and **PPID** on ice. Upon first thaw, briefly spin tube containing protein to recover full content of the tube. Aliquot protein into single use aliquots. Store remaining undiluted protein in aliquots at -80°C immediately. *Note: the HSP90 α and PPID proteins are very sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted protein.*
- 2) Dilute **3x HSP90 Assay Buffer 2** with water to **1x HSP90 Assay Buffer 2**.
- 3) Dilute **HSP90 α** in **1x HSP90 Assay Buffer 2** at 1.5 ng/ μ l. Add 4 μ l of diluted HSP α protein to each well designated for the “Positive Control”, “Test Inhibitor”, and “Blank”. To the wells labeled “Substrate Control”, add 4 μ l of **1x HSP90 Assay Buffer 2**. Discard any unused diluted protein after use.
- 4) Add 2 μ l of inhibitor solution to each well designated “Test Inhibitor”. For the “Positive Control”, “Substrate Control” and “Blank”, add 2 μ l of the same solution without inhibitor (inhibitor buffer). *Note: Keep DMSO concentration below 1%.*
- 5) Dilute **PPID** in **1x HSP90 Assay Buffer 2** at 10 ng/ μ l. Keep diluted protein on ice until use.
- 6) Add 4 μ l of **1x HSP90 Assay Buffer 2** to the well designated “Blank”.
- 7) Initiate reaction by adding 4 μ l of diluted **PPID** to each well designated for the “Substrate control”, “Positive Control”, and “Test Inhibitor”. Incubate at room temperature for 30 minutes.

	Blank	Substrate Control	Positive Control	Test Inhibitor
Diluted HSP90 α (1.5 ng/ μ l)	4 μ l	–	4 μ l	4 μ l
Test Inhibitor/Activator	–	–	–	2 μ l
Inhibitor buffer (no inhibitor)	2 μ l	2 μ l	2 μ l	–
1x HSP90 Assay Buffer 2	4 μ l	4 μ l	–	–
Diluted PPID (10 ng/ μ l)	–	4 μ l	4 μ l	4 μ l
Total	10 μl	10 μl	10 μl	10 μl

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082**, Fax **1.858.481.8694**
 Or you can Email us at: info@bpsbioscience.com
 Please visit our website at: www.bpsbioscience.com

Step 2:

Note: Protect your samples from direct exposure to light!


- 1) Dilute **3x Detection Buffer** with water to **1x Detection Buffer**. Dilute **Glutathione AlphaLISA® Acceptor Beads** (PerkinElmer #AL109C) 250-fold with **1x Detection Buffer**. Add 10 µl per well. Shake plate briefly. Incubate at room temperature for 30 minutes.

Step 3:

- 1) Dilute **Streptavidin-conjugated donor beads** (PE #6760002S) 125-fold with **1x Detection Buffer**. Add 10 µl per well. Incubate at room temperature for 1 hour.
- 2) Read Alpha-counts.

Due to lot to lot variability in AlphaScreen® bead performance, it may be necessary to optimize assay conditions. For example, slight adjustments to HSP90α or PPID concentrations may improve signal-to-noise ratio.

EXAMPLE OF ASSAY RESULTS:



HSP90α (C-terminal):PPID binding activity measured using *HSP90α (C-Terminal) Inhibitor Screening Assay Kit*, BPS Bioscience, Catalog # 50317. *Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.com.*

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082**, Fax **1.858.481.8694**
 Or you can Email us at: info@bpsbioscience.com
 Please visit our website at: www.bpsbioscience.com



6044 Cornerstone Ct. West, Ste. E
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

RELATED PRODUCTS:

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
HSP90 α (C-terminal), Biotin-labeled	50316	100 μ g
HSP90 β (C-terminal), Biotin-labeled	50313	100 μ g
HSP90 α recombinant enzyme	50290	200 μ g
HSP90 β recombinant enzyme	50292	200 μ g
Aha1 recombinant enzyme	50291	200 μ g
Geldanamycin inhibitor	27008	5 mg
MS-275 (Entinostat) inhibitor	27011	25 mg
Novobiocin inhibitor	27501	250 μ l
HSP90 α Assay Kit (96 well)	50293	96 rxns
HSP90 β Assay Kit (96 well)	50294	96 rxns
HSP90 α Assay Kit (384 well)	50298	384 rxns
HSP90 β Assay Kit (384 well)	50299	384 rxns

AlphaScreen[®] and AlphaLISA[®] are registered trademarks of PerkinElmer, Inc.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082**, Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com