



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
Tel: 1.858.202.1401  
Fax: 1.858.481.8694  
Email: info@bpsbioscience.com

## Data Sheet

### **TET1 Chemiluminescent Assay Kit**

Catalog #50651

**DESCRIPTION:** The *TET1 Chemiluminescent Assay Kit* is designed to measure TET1 activity for screening and profiling applications. TET1 belongs to the Ten Eleven Translocation (TET) family proteins that catalyze 5-methylcytosine oxidation and generate 5-methylcytosine derivatives, including 5-hydroxymethylcytosine. TET1 may repress gene transcription in a catalytic activity-independent manner. The *TET1 Chemiluminescent Assay Kit* comes in a convenient format, with a 96-well strip plate precoated with methylated DNA substrate, primary antibody, a secondary HRP-labeled antibody, TET1 assay buffer, and purified TET1 for 100 enzyme reactions. The key to the *TET1 Chemiluminescent Assay Kit* is a highly specific antibody that recognizes hydroxymethylated substrate. With this kit, only three simple steps on a microtiter plate are required for detection of TET1 activity. First, the TET1 enzyme is incubated with the methylated substrate for two hours. Then the primary antibody is added. Next, the plate is treated with an HRP-labeled secondary antibody. Finally, the HRP substrate is added to produce chemiluminescence that can be measured using a chemiluminescence reader.

#### **COMPONENTS:**

Catalog #	Component	Amount	Storage	
50161	TET1, FLAG-tag	20 µg	-80°C	<b>Avoid freeze/ thaw cycles!</b>
52140Z2	Primary antibody 27	100 µl	-80°C	
52131H	Secondary HRP-labeled antibody 2	10 µl	-80°C	
52132	4x TET1 assay buffer 1	3 ml	-80°C	
79556	Blocking buffer 1	50 ml	+4°C	
79670	ELISA ECL substrate A (transparent bottle)	6 ml	+4°C	
79670	ELISA ECL substrate B (brown bottle)	6 ml	+4°C	
	8-well strip plate module precoated with DNA substrate	1	+4°C	

#### **MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:**

TBST buffer (1 x TBS, pH 8.0, containing 0.05% Tween-20)  
Luminometer or fluorescent microplate reader capable of reading chemiluminescence  
Adjustable micropipettor and sterile tips  
Rotating or rocker platform

**APPLICATIONS:** Great for studying enzyme kinetics and HTS applications.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.202.1401  
**Fax:** 1.858.481.8694  
**Email:** info@bpsbioscience.com

**CONTRAINDICATIONS:** DMSO >1%, strong acids or bases, ionic detergents, high salt  
**STABILITY:** One year from date of receipt when stored as directed.

**REFERENCE:** Li Tan, Yujiang Geno Shi. *Development* 2012; **139(11):** 1895-1902.

#### ASSAY PROTOCOL:

All samples and controls should be tested in duplicate.

#### Step 1:

- 1) Rehydrate the microwells by adding 200  $\mu$ l of TBST buffer (1 x TBS, pH 8.0, containing 0.05% Tween-20) to every well. Incubate 15 minutes at room temperature. Tap the strip plate onto clean paper towels to remove liquid.
- 2) Add 12.5  $\mu$ l **H<sub>2</sub>O** and 7.5  $\mu$ l **4x TET1 Assay Buffer 1** to each well. Dilute **4x TET1 Assay Buffer 1** with water to make **1x TET1 Assay Buffer 1**.
- 3) Thaw **TET1** on ice. Upon first thaw, briefly spin tube containing enzyme to recover full content of the tube. Aliquot **TET1** enzyme into single use aliquots. Store remaining undiluted enzyme in aliquots at -80°C. *Note: TET1 is very sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted enzyme.*
- 4) Dilute **TET1** in **1x TET1 Assay Buffer 1** at 5-10 ng/ $\mu$ l (100-200 ng/reaction). Keep diluted enzyme on ice until use. Discard any unused diluted enzyme after use.
- 5) Add 10  $\mu$ l of inhibitor solution of each well designated "Test Inhibitor." For the "Positive Control" and "Blank," add 10  $\mu$ l of 5% DMSO in Water (inhibitor buffer). *Note: Keep final DMSO concentration  $\leq$ 1%.*

	Positive Control	Test Inhibitor	Blank
H <sub>2</sub> O	12.5 $\mu$ l	12.5 $\mu$ l	12.5 $\mu$ l
4x TET1 Assay Buffer 1	7.5 $\mu$ l	7.5 $\mu$ l	7.5 $\mu$ l
1x TET1 Assay Buffer 1	-	-	20 $\mu$ l
Test Inhibitor	-	10 $\mu$ l	-
5% DMSO in Water (inhibitor buffer)	10 $\mu$ l	-	10 $\mu$ l
TET1 (5-10 ng/ $\mu$ l)	20 $\mu$ l	20 $\mu$ l	-
<b>Total</b>	<b>50 <math>\mu</math>l</b>	<b>50 <math>\mu</math>l</b>	<b>50 <math>\mu</math>l</b>

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.202.1401  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

- 6) Add 20  $\mu$ l of **1x TET1 Assay Buffer 1** to wells designated as "Blank." Initiate reaction by adding 20  $\mu$ l of diluted **TET1** (prepared as described above) to wells designated "Positive Control" and "Test Inhibitor." Incubate at room temperature for two hours.
- 7) Wash the strip plate three times with TBST buffer. Blot dry onto clean paper towels.
- 8) Add 100  $\mu$ l of Blocking buffer 1 to every well. Shake on a rotating platform for 10 min. Remove supernatant as above.

### **Step 2:**

- 1) Dilute "**Primary antibody 27**" 100-fold with **Blocking buffer 1**.
- 2) Add 100  $\mu$ l per well. Incubate 1 hour at room temperature with slow shaking.
- 3) Wash strip plate three times with TBST buffer and incubate in Blocking buffer 1 as in steps 1-7 and 1-8.

### **Step 3:**

- 1) Dilute "**Secondary HRP-labeled antibody 2**" 1,000-fold with **Blocking buffer 1**.
- 2) Add 100  $\mu$ l per well. Incubate for 30 min. at room temperature with slow shaking.
- 3) Wash strip plate three times with TBST buffer and incubate in Blocking buffer 1 as in steps 1-7 and 1-8.
- 4) Just before use, mix on ice 50  $\mu$ l **ELISA ECL substrate A**, 50  $\mu$ l **ELISA ECL substrate B**, and add 100  $\mu$ l per well. Discard any unused chemiluminescent reagent after use.
- 5) Immediately read sample in a luminometer or microtiter-plate capable of reading chemiluminescence.

### **Reading Chemiluminescence:**

Chemiluminescence is the emission of light (luminescence) which results from a chemical reaction. The detection of chemiluminescence requires no wavelength selection because the method used is emission photometry and is not emission spectrophotometry.

To properly read chemiluminescence, make sure the plate reader is set for LUMINESCENCE mode. Typical integration time is 1 second, delay after plate movement is 100 msec. Do not use a filter when measuring light emission. Typical settings for the Synergy 2 BioTek plate reader are: use the "hole" position on the filter wheel; Optics position: Top; Read type: endpoint. Sensitivity may be adjusted based on the luminescence of a control assay without enzyme (typically we set this value as 100).


OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)

### Examples of Assay Results:



TET1 enzyme inhibition by 2,4-pyridinecarboxylic acid (2,4-PC), measured using the *TET1 Chemiluminescent Assay Kit*, BPS Bioscience #50651. Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

### RELATED PRODUCTS

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
TET1 recombinant protein	50161	50 $\mu\text{g}$
TET2 recombinant protein	50162	50 $\mu\text{g}$
TET2 Chemiluminescent Assay Kit	50652	96 reactions
Anti-5-hmC polyclonal antibody	25205	100 $\mu\text{l}$
Anti-5-mC monoclonal antibody 33D3	25207	100 $\mu\text{g}$
Anti-5-mC polyclonal antibody	25200	50 $\mu\text{g}$
MECP2 recombinant protein	50250	50 $\mu\text{g}$
DNMT1 recombinant protein	51101	10 $\mu\text{g}$
DNMT2 recombinant protein	51102	10 $\mu\text{g}$
DNMT3A/DNMT3L protein complex	51106	10 $\mu\text{g}$
DNMT3B/DNMT3L protein complex	51104	10 $\mu\text{g}$
DNMT1 Assay Kit	52050L	96 reactions
DNMT3A Assay Kit	52033	96 reactions
DNMT3B Assay Kit	52034	96 reactions

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court W, Ste B  
 San Diego, CA 92121  
**Tel:** 1.858.202.1401  
**Fax:** 1.858.481.8694  
**Email:** info@bpsbioscience.com

## TROUBLESHOOTING GUIDE

Problem	Possible Cause	Solution
Luminescence signal of positive control reaction is same as "blank" value.	TET1 has lost activity	Enzyme loses activity upon repeated freeze/thaw cycles. Use fresh TET1, BPS Bioscience #50161. Store enzyme in single-use aliquots. Increase time of enzyme incubation. Increase enzyme concentration.
	Antibody reaction is insufficient	Increase time for primary antibody incubation. Avoid freeze/thaw cycles of antibodies.
	Incorrect settings on instruments	Refer to instrument instructions for settings to increase sensitivity of light detection. See section on "Reading Chemiluminescence" above.
	Chemiluminescent reagents mixed too soon	Chemiluminescent solution should be used within 15 minutes of mixing. Ensure both reagents are properly mixed.
Luminescent signal is erratic or varies widely among wells	Inaccurate pipetting/technique	Run duplicates of all reactions. Use a multichannel pipettor. Use master mixes to minimize errors.
	Bubbles in wells	Pipette slowly to avoid bubble formation. Tap plate lightly to disperse bubbles; be careful not to splash between wells.
Background (signal to noise ratio) is high	Insufficient washes	Be sure to include blocking steps after wash steps. Increase number of washes. Increase wash volume. Increase Tween-20 concentration to 0.1% in TBST.
	Sample solvent is inhibiting the enzyme	Run negative control assay including solvent. Maintain DMSO level at <1% Increase time of enzyme incubation.
	Results are outside the linear range of the assay	Use different concentrations of TET1, BPS Bioscience #50161, to create a standard curve.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.  
 To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**  
 Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)  
 Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)