



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet

SRE Reporter Kit (MAPK/ERK Signaling Pathway) Catalog #: 60511

Background

The MAPK/ERK signaling pathway is a major participant in the regulation of cell growth and differentiation. It can be activated by various extracellular stimuli including mitogens, growth factors, and cytokines. Upon stimulation, MEK1/2 phosphorylate and activate ERK1/2. The activated ERK translocates to the nucleus where it phosphorylates and activates transcription factors. The TCFs (Ternary Complex Factors), including Elk1, are among the best-characterized transcription factor substrates of ERK. When phosphorylated by ERK, Elk1 forms a complex with Serum Response Factor (SRF) and binds to Serum Response Element (SRE), resulting in the expression of numerous mitogen-inducible genes.

Description

The *SRE Reporter Kit* is designed for monitoring the activity of the ERK signaling pathway and the transcriptional activity of SRF in cultured cells. The kit contains a transfection-ready SRE luciferase reporter vector, which is an ERK pathway-responsive reporter. This reporter contains the firefly luciferase gene under the control of multimerized SRE responsive elements located upstream of a minimal promoter. The SRE reporter is premixed with a constitutively-expressing *Renilla* luciferase vector that serves as an internal control for transfection efficiency.

The kit also includes a non-inducible firefly luciferase vector premixed with constitutively-expressing *Renilla* luciferase vector as a negative control. The non-inducible luciferase vector contains the firefly luciferase gene under the control of a minimal promoter, without any additional response elements. The negative control is critical for determining pathway-specific effects and the background luciferase activity.

Applications

- Monitor MAPK/ERK signaling pathway activity and SRF-mediated activity.
- Screen for activators or inhibitors of the MAPK/ERK signaling pathway.
- Study effects of RNAi or gene overexpression on the activity of the MAPK/ERK pathway.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Components

Component	Specification	Amount	Storage
Reporter (Component A)	SRE luciferase reporter vector + constitutively expressing <i>Renilla</i> luciferase vector	500 µl (55 ng DNA/ µl)	-20°C
Negative Control Reporter (Component B)	Non-inducible luciferase vector + constitutively expressing <i>Renilla</i> luciferase vector	500 µl (55 ng DNA/ µl)	-20°C

Note: These vectors are designed for transient transfection. They are NOT SUITABLE for transformation and amplification in bacteria.

Materials Required but Not Supplied

- Mammalian cell line and appropriate cell culture medium
- 96-well tissue culture plate or 96-well tissue culture-treated white clear-bottom assay plate
- Transfection reagent for mammalian cell line [We use Lipofectamine™ 2000 (Invitrogen # 11668027). However, other transfection reagents work equally well.]
- Opti-MEM I Reduced Serum Medium (Invitrogen #31985-062)
- Dual Luciferase Assay System:
Dual Luciferase (Firefly-Renilla) Assay System (BPS Cat. #60683): This system assays cells directly in growth medium. It can be used with any luminometer. Automated injectors are not required.
- Luminometer

Generalized Transfection and Assay Protocols

The following procedure is designed to transfect the reporter into HEK293 cells using Lipofectamine 2000 in a 96-well format. To transfect cells in different tissue culture formats, adjust the amounts of reagents and cell number in proportion to the relative surface area. If using a transfection reagent other than Lipofectamine 2000, follow the manufacturer's transfection protocol. Transfection conditions should be optimized according to the cell type and study requirements.

All amounts and volumes in the following setup are given on a per-well basis.

1. One day before transfection, seed cells at a density of ~ 30,000 cells per well in 100 µl of growth medium so that cells will be 90% confluent at the time of transfection.
2. The next day, for each well, prepare complexes as follows:
 - a. Dilute DNA mixtures in 15 µl of Opti-MEM I medium (antibiotic-free). Mix gently.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Depending upon the experimental design, the DNA mixtures may be any of following combinations:

- **1 μ l of Reporter** (component A); in this experiment, the control transfection is **1 μ l of Negative Control Reporter** (component B).
- **1 μ l of Reporter** (component A) + experimental vector expressing gene of interest; in this experiment, the control transfections are: **1 μ l of Reporter** (component A) + negative control expression vector, **1 μ l of Negative Control Reporter** (component B) + experimental vector expressing gene of interest, and **1 μ l of Negative Control Reporter** (component B) + negative control expression vector.
- **1 μ l of Reporter** (component A) + specific siRNA; in this experiment, the control transfections are: **1 μ l of Reporter** (component A) + negative control siRNA, **1 μ l of Negative Control Reporter** (component B) + specific siRNA, and **1 μ l of Negative Control Reporter** (component B) + negative control siRNA.

Note: we recommend setting up each condition in at least triplicate, and preparing transfection cocktail for multiple wells to minimize pipetting errors.

b. Mix Lipofectamine 2000 gently before use, then dilute 0.35 μ l of Lipofectamine 2000 in 15 μ l of Opti-MEM I medium (antibiotic-free). Incubate for 5 minutes at room temperature.

Note: Prepare this dilution cocktail in volumes sufficient for the whole experiment.

c. After the 5 minute incubation, combine the diluted DNA with diluted Lipofectamine 2000. Mix gently and incubate for 25 minutes at room temperature.

3. Add the 30 μ l of the complexes to each well containing cells and medium. Mix gently by tapping the plate.
4. Incubate cells at 37°C in a CO₂ incubator. After ~5 to 6 hours of transfection, change medium to fresh medium with 0.5% serum. Incubate cells at 37°C in a CO₂ incubator overnight.
5. The next day, induce the SRE reporter with medium containing activators of the ERK pathway such as high percentage of serum or growth factors. Incubate cells at 37°C in a CO₂ incubator for ~ 6 hours. After 6-hour treatment, perform the Dual Luciferase Assay System following the protocol on the BPS data sheet (BPS Cat. #60683).

To study the effect of inhibitors on the ERK pathway, after ~5-6 hours of transfection, treat cells with inhibitors in medium containing 0.5% serum. The next day, treat cells with activators for 6 hours, then perform the luciferase assay.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Sample protocol to determine the effect of serum or EGF on SRE reporter activity in HEK293 cells

1. One day before transfection, seed HEK293 cells at a density of 30,000 cells per well into white clear-bottom 96-well plate in 100 μ l of growth medium (MEM/EBSS (Hyclone #SH30024.01), 10% FBS, 1% non-essential amino acids, 1 mM Na-pyruvate, 1% Pen/Strep). Incubate cells overnight at 37°C in a CO₂ incubator.
2. The next day, transfect 1 μ l of SRE reporter (component A) into cells following the procedure in **Generalized Transfection and Assay Protocols**.
3. After ~ 6 hours of transfection, change medium to 50 μ l of medium containing 0.5% FBS (MEM, 0.5% FBS, with non-essential amino acids, Na-pyruvate, and 1% Pen/Strep). Incubate cells at 37°C in a CO₂ incubator for ~ 16 to 18 hours.
4. The next day after transfection, treat cells with 50 μ l of medium containing a high percentage of FBS, with or without EGF, or medium containing 0.5% FBS with EGF. For unstimulated control wells, use cells in medium with 0.5% FBS. Add 50 μ l of growth medium to cell-free control wells to determine the background luminescence). Set up each treatment in at least triplicate.
5. Incubate cells at 37°C in a CO₂ incubator for ~ 6 hours.
6. After ~48 hours of transfection, perform dual luciferase assay using BPS Dual Luciferase (Firefly-Renilla) Assay System (BPS Cat. #60683): Dilute 100x Firefly Luciferase Reagent Substrate (Component B) into Firefly Luciferase Reagent Buffer (Component A). Add 50 μ l of Firefly Luciferase reagent per well and rock at room temperature for ~15 minutes, then measure firefly luminescence using a luminometer. Dilute 100x Renilla Luciferase Reagent Substrate (Component D) into Renilla Luciferase Reagent Buffer (Component C). Add 50 μ l of Renilla Luciferase reagent per well, rock at room temperature for ~1 minute and measure Renilla luminescence.
7. To obtain the normalized luciferase activity for the SRE reporter, subtract the background luminescence, then calculate the ratio of firefly luminescence from SRE reporter to *Renilla* luminescence from the control *Renilla* luciferase vector.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

Figure 1. Serum induced the expression of SRE reporter. The results are shown as fold induction of normalized SRE reporter activity. Fold induction is determined by comparing values against the mean value for control cells with 0.5% FBS treatment.





Figure 2. EGF induced the expression of SRE reporter. The results are shown as fold induction of normalized SRE reporter activity. Fold induction are determined by comparing values against the mean value for control cells with 0.5% FBS treatment only.




OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
 Or you can Email us at: info@bpsbioscience.com
 Please visit our website at: www.bpsbioscience.com

Figure 3. Dose response of SRE reporter activity to EGF in the presence of 0.5% FBS. The results are shown as fold induction of normalized SRE reporter activity. Fold induction is determined by comparing values against the mean value for control cells without EGF treatment.

The EC₅₀ of EGF is ~0.97 ng/ml



Sample protocol to determine the effect of inhibitors of the ERK pathway on SRE reporter activity in HEK293 cells

1. One day before transfection, seed HEK293 cells at a density of 30,000 cells per well into white clear-bottom 96-well plate in 100 μ l of growth medium. Incubate cells overnight at 37°C in a CO₂ incubator.
2. The next day, transfect 1 μ l of SRE reporter (component A) into cells following the procedure in **Generalized Transfection and Assay Protocols**.
3. After ~6 hours of transfection, treat transfected cells with three-fold serial dilution of U0126 (MEK inhibitor) in 50 μ l of medium containing 0.5% FBS. For wells without U0126, treat cells with medium containing 0.5% FBS only. Incubate cells at 37°C in a CO₂ incubator for ~ 16 to 18 hours.
4. The next day after transfection, treat the cells with recombinant EGF (final concentration 10 ng/ml) in 50 μ l of medium containing 0.5% FBS with U0126. For unstimulated control wells, determine the basal activity using cells in medium with 0.5% FBS. To determine background luminescence, ; add 50 μ l of medium to cell-free control wells. Set up each treatment in at least triplicate.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**


Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

5. Incubate cells at 37°C in a CO₂ incubator for ~ 6 hours.
6. After ~48 hours of transfection, perform dual luciferase assay using BPS Dual Luciferase (Firefly-Renilla) Assay System (BPS Cat. #60683): Dilute 100x Firefly Luciferase Reagent Substrate (Component B) into Firefly Luciferase Reagent Buffer (Component A). Add 50 µl of Firefly Luciferase reagent per well and rock at room temperature for ~15 minutes, then measure firefly luminescence using a luminometer. Dilute 100x Renilla Luciferase Reagent Substrate (Component D) into Renilla Luciferase Reagent Buffer (Component C). Add 50 µl of Renilla Luciferase reagent per well, rock at room temperature for ~1 minute and measure Renilla luminescence.
7. To obtain the normalized luciferase activity for the SRE reporter, subtract the background luminescence, then calculate the ratio of firefly luminescence from the SRE reporter to *Renilla* luminescence from the control *Renilla* luciferase vector.

Figure 4. Inhibition of EGF-induced SRE reporter activity by ERK pathway inhibitor, U0126. The results are shown as percentage of SRE reporter activity. The normalized luciferase activity for cells stimulated with EGF in the absence of U0126 is set at 100%.

The IC₅₀ of U0126 is ~ 0.7 µM



References

- Wong, K.K. (2009) Recent developments in anti-cancer agents targeting the Ras/Raf/ MEK/ERK pathway. *Recent Pat Anticancer Drug Discov.* **4(1)**:28-35.
Treisman, R. (1992) The serum response element. *Trends Biochem Sci.* **17(10)**: 423-426.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Related Products

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
Dual Luciferase (Firefly-Renilla) Assay System	60683-1	10 mL
Dual Luciferase (Firefly-Renilla) Assay System	60683-2	100 mL
Dual Luciferase (Firefly-Renilla) Assay System	60683-3	1 L
EGF, human	90201-1	100 µg
EGF, human	90201-2	500 µg
EGF, mouse	90200-1	100 µg
EGF, mouse	90200-2	500 µg
ERK1	40055	10 µg
ERK2	40299	10 µg
MAP3K14 (NIK)	40090	10 µg
MAPKAPK2 (MK2)	40088	100 µg
MAPK10 (JNK3)	40092	10 µg
MEK1 (K97R)	40075	100 µg
MEK1, mouse	40121	10 µg
MEK1, human	40123	10 µg
MEK1, GST-tag	40527	50 µg
MEK2	40125	10 µg
MEKK2	40122	10 µg
MEKK3	40124	10 µg
U0126	27012	5 mg

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com