



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet

Fluorogenic MMP12 Assay Kit

Catalog #78017
Size: 96 reactions

BACKGROUND: MMP12 (matrix metalloproteinase 12) is a member of the matrix metalloproteinase (MMP) family involved in the degradation of the extracellular matrix. MMP12 is also associated with the regulation of cytokines and chemokines, suggesting a role for MMP12 in inflammation. MMP12 is also essential for arterial stiffening in mice and is a highly prognostic marker of arterial stiffness, making MMP12 a therapeutic target for cardiovascular disease.

DESCRIPTION: The *Fluorogenic MMP12 Assay Kit* is designed to measure MMP12 activity for screening and profiling applications, in a homogeneous assay with no time-consuming washing steps. The kit comes in a convenient 96-well format, with purified MMP12 enzyme, fluorogenic substrate, and MMP assay buffer for 100 enzyme reactions.

COMPONENTS:

Catalog #	Component	Amount	Storage	
	MMP12	3 µg	-80°C	Avoid freeze/ thaw cycles!
79919	MMP Substrate (1 mM)	10 µl	-80°C	
79917	1X MMP Assay Buffer 1	25 ml	-20°C	
79685	Black, low binding black microtiter plate	1	Room Temperature	

APPLICATIONS: Great for studying enzyme kinetics and HTS applications.

STABILITY: At least one year from date of receipt when stored as directed.

REFERENCE(S):

1. Morris, D.G., *et al.* 2003. "Loss of integrin $\alpha\beta6$ -mediated TGF- β activation causes MMP12-dependent emphysema." *Nature* **422(6928)**: 169-173.
2. Hunninghake, G.M., *et al.* 2009. "MMP12, lung function, and COPD in high-risk populations." *New England Journal of Medicine* **361(27)**: 2599-2608.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

ASSAY PROTOCOL:

All samples and controls should be tested in duplicate.

- 1) Dilute 1 mM MMP substrate 1:100 in 1X assay buffer, to make a 10 μ M solution. Dilute only enough as is required for the assay. Store remaining 1 mM substrate in aliquots at -80°C.
- 2) Prepare the substrate solution: N wells \times (20 μ l 1X assay buffer + 5 μ l diluted (10 μ M) MMP Substrate).
- 3) Add 25 μ l of the substrate solution to each well (Final concentration of the MMP substrate in a 50 μ l reaction is 1 μ M).

Component	Positive Control	Test Sample	Blank
Substrate solution	25 μ l	25 μ l	25 μ l
Test Inhibitor	–	5 μ l	–
Inhibitor buffer (usually 10% DMSO in assay buffer)	5 μ l	–	5 μ l
MMP12 (1.25 ng/ μ l)	20 μ l	20 μ l	–
1X Assay Buffer	–	–	20 μ l
Total	50 μl	50 μl	50 μl

- 4) Prepare the inhibitor solution.

The final concentration of DMSO in the assay should not exceed 1%. If the inhibitor compound is dissolved in DMSO, make a 100-fold higher concentration of the compound than the highest concentration you want to test in DMSO. Then make a 10-fold dilution in 1X assay buffer (at this step the compound concentration is 10-fold higher than the final concentration in 10% DMSO). To determine an IC₅₀ or to test lower concentrations of the compound, prepare a series of further dilutions in 1X assay buffer containing 10% DMSO (the final concentration of the DMSO will be 1% in all samples).

If the inhibitor compound is dissolved in water, make a solution of the compound 10-fold higher than the final concentration in 1X assay buffer.

- 5) Add 5 μ l inhibitor solution to each well designated “Test Sample.” Add 5 μ l of the same solution without the inhibitor (inhibitor buffer, usually 10% DMSO in assay buffer) to “Blank” and “Positive Control” wells.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

- 6) Thaw MMP12 on ice. Upon first thaw, briefly spin tube containing enzyme to recover the full content of the tube. Aliquot MMP12 into single use aliquots. Store remaining undiluted enzyme in aliquots at -80°C. Note: MMP12 enzyme is sensitive to freeze/thaw cycles. Do not re-use diluted enzyme.
- 7) Dilute MMP12 in 1x assay buffer at 1.25 ng/μl (25 ng per reaction).
- 8) Add 20 μl diluted MMP12 enzyme solution to wells designated as “Positive Control” and “Test Sample.” Add 20 μl 1X assay buffer to the “Blank” wells.
- 9) Incubate at room temperature for 2 hours. Measure the fluorescence intensity in a microtiter plate-reading fluorimeter capable of excitation at wavelength 328 nm ±10 nm and detection of emission at wavelength 393 nm ±10 nm. The fluorescence intensity can also be measured kinetically. “Blank” value is subtracted from all other values.

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:


Fluorescent microplate reader capable of reading exc/em=328 nm/393 nm

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

EXAMPLE OF ASSAY RESULTS:


Inhibition of MMP12 enzyme activity by MMP408 (Sigma #444291), measured using the *Fluorogenic MMP12 Assay Kit (BPS Bioscience #78017)*. Fluorescence intensity was measured using a Tecan fluorescent microplate reader. *Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.com*

RELATED PRODUCTS

<u>Product</u>	<u>Cat. #</u>	<u>Size</u>
MMP1, His-Tag (Human)	80214	20 μ g
MMP2, His-Tag (Human)	80213	20 μ g
MMP3(K45E), His-Tag (Human)	11346	100 μ g
MMP8, His-Tag (Human)	100552	100 μ g
MMP9(Q279R), His-Tag (Human)	80215	20 μ g
Fluorogenic MMP3 (K45E) Assay Kit	79907	384 rxns.
Fluorogenic MMP2 Assay Kit	79918	96 rxns.
Fluorogenic MMP9 (Q279R) Assay Kit	79915	96 rxns.
Fluorogenic MMP10 Assay Kit	79986	96 rxns.
Fluorogenic MMP13 Assay Kit	79991	96 rxns.
Fluorogenic MMP14 Assay Kit	79993	96 rxns.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com