



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court West, Ste. B  
San Diego Ca 92121  
Tel: 1.858.202.1401  
Fax: 1.858.481.8694  
Email: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## **Data Sheet**

### ***CRISPRa (SAM) Jurkat Cell Line***

**Catalog #: 78080**

#### **Descriptions**

This cell line has been engineered for use with the CRISPR Synergistic Activation Mediator (SAM) system to induce transcriptional activation and expression of any gene of interest. Cells stably express a mutated dCas9 (*Streptococcus pyogenes* CRISPR associated protein 9), lacking any endonuclease activity, fused to VP64, a transcriptional activator. Stable dCas9-VP64 expression is maintained with Blasticidin resistance. Cells also stably express P65 (Transcription Factor p65, or Nuclear Factor NF- $\kappa$ -B p65) and HSF1 (Heat Shock Factor 1) fused with an MS2 tag, which is maintained with Hygromycin resistance. When these cells are transfected with an MS2-tagged sgRNA targeting the promoter region of the gene of interest, dCas9-VP64 and MS2-P65-HSF1 are recruited to the site in the genomic DNA and begin transcription, inducing expression of the desired gene.

#### **Application**

1. Quickly generating CRISPR-activated cell pools or cell lines in Jurkat cells.
2. Implementing sgRNA CRISPRa screens in Jurkat cells.

#### **Format**

Each vial contains  $\sim 2 \times 10^6$  cells in 1 ml of FBS with 10% DMSO.

#### **Storage**

Store in liquid nitrogen immediately upon receipt.

#### **Culture conditions**

**Thaw Medium 2 (BPS Bioscience, #60184):** RPMI 1640 medium (ThermoFisher, #A1049101) supplemented with 10% FBS (ThermoFisher, #26140079), 1% Penicillin/Streptomycin (Hyclone, #SV30010.01).

**Growth Medium 2L (BPS Bioscience, #78094):** Thaw Medium 2 (BPS Bioscience, #60184) plus 200  $\mu$ g/ml of Hygromycin B (Invivogen, #ant-hg) and 5  $\mu$ g/ml Blasticidin (Invivogen, #ant-bl) to ensure recombinant expression.

Cells should be grown at 37°C with 5% CO<sub>2</sub> using Growth Medium 2L to ensure recombinant expression is maintained.

It is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, and then transfer the entire contents of the vial to a tube containing 10 ml of Thaw Medium 2 (**no Hygromycin or Blasticidin**). Then spin the cells down, remove the supernatant, and resuspend the cells in 5 ml of pre-warmed Thaw Medium 2 (**no Hygromycin or Blasticidin**). Transfer the resuspended cells to a T25 flask and incubate at 37°C in a 5% CO<sub>2</sub> incubator. After 24 hours of culture, add an additional 3-4 ml of Thaw Medium 2 (**no Hygromycin or Blasticidin**).

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)

**Blasticidin**). At first passage, switch to Growth Medium 2L (contains Hygromycin and Blasticidin). Cells should be split before they reach  $2 \times 10^6$  cells/ml.


**Cryopreservation**: When cells reach 90% confluency, spin cells, and remove medium from the pellet. Resuspend the cells in freezing medium (10% DMSO in FBS). Freeze cells using a reduced rate freezing box ( $-0.5^{\circ}\text{C}$  to  $-1^{\circ}\text{C}$  per minute) down to  $-80^{\circ}\text{C}$ , then move cells to liquid nitrogen for long term storage. Cells have been demonstrated to be stable for at least 15 passages; BPS recommends preparing frozen stocks at a low passage number so cells are not used beyond passage 20.

### **Mycoplasma Testing**

This cell line has been screened using the MycoAlert™ Mycoplasma Detection Kit (Lonza, #LT07-118) to confirm the absence of Mycoplasma contamination. MycoAlert Assay Control Set (Lonza, #LT07-518) was used as a positive control.

### **Validation**

Expression of dCas9-VP64 and MS2-P65-HSF1 was functionally validated by transfecting the cells with sgRNA-MS2 targeting PD-1.



### **Figure 1. Induction of PD-1 in CRISPRa (SAM) Jurkat cells.**

CRISPRa (SAM) Jurkat cells were electroporated with sgRNA-MS2 targeting PD-1 (Programmed Cell Death protein 1, or CD279, BPS Bioscience #78091) to induce PD-1 expression. Cells were stained with PE-labeled anti-PD-1 antibody (BioLegend, #637309) and analyzed by FACS. Parental CRISPRa (SAM) Jurkat cells are shown in green, and the transfected CRISPRa (SAM) Jurkat cells are shown in blue.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court West, Ste. B  
San Diego Ca 92121  
Tel: 1.858.202.1401  
Fax: 1.858.481.8694  
Email: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

### Sequence

dCas9 (nuclease deficient *Streptococcus pyogenes* Cas9, in blue) fused with a linker (black) and VP64 (red):

DKKYSIGLAIGTNSVGVAVITDEYKVPSSKFKVLGNTDRHSIKKNLIGALLFDSGETAEATRLKRT  
ARRRYTRRKNRICYLQEIFSNEMAKVDDSSFFHRLEESFLVEEDKKHERHPIFGNIVDEVAYHEK  
YPTIYHLRKKLV DSTDKADLRILIYLAHAMIKFRGHFLIEGDLNPDNSDVKLFIQLVQTYNQLFE  
ENPINASGVDAKAILSARLSKSRRLLENLIAQLPGEKKNLFGNLIASLGLTPNFKSNFDLAEDAK  
LQLSKD TYDDLDNLLAQIGDQYADLFLAAKNLSDAILLSDILRVNTEITKAPLSASMIKRYDEHH  
QDLTLLKALVRQQLPEKYKEIFFDQSKNGYAGYIDGGASQEEFYKFIKPILEKMDGTEELLVKNL  
REDLLRKQRTFDNGSIPHQIHLGELHAILRRQEDFYFPLKDNREKIEKILTRIPYYVGPLARGNS  
RFAWMTRKSEETITPWNFEEVVDK GASAQSFIERMTNFDKNLPNEKVLPHSLLEYFTVYNEL  
TKVKYVTEGMRKPAFLSGEQKKAIVDLLFKTNRKVTVKQLKEDYFKKIECFDSVEISGVEDRFNA  
SLGTYHDLLKIIKDKDFLDNEENEDILEDIVLTLTFEDREMIEERLKYAHLFDDKVMKQLKRRR  
YTGWGRLSRKLINGIRDKQSGKTILDFLKSDGFANRNFMLIHDDSLTFKEDIQKAQVSGQGDS  
LHEHIANLAGSPAIKKGIQTVKVDELVKVMGRHKPENIVIAMARENQTTQKGQKNSRERMKRI  
EEGKELGSQILKEHPVENTQLQNEKLYLYLQNGRDMYVDQELDINRLSDYDVDHIVPQSFLK  
DSDIDNKVLRSDKARGKSDNVPSEEVVKKMKNYWRQLLNAKLITQRKFDNLTKAERGGLSEL  
DKAGFIKRQLVETRQITKHVAQILDSRMNTKYDENDKLIREVKVITLKSCLVSDFRKDFQFYKVR  
EINNYHHAHDAYLNAVVG TALIKKYPKLESEFVYGDYKVVYDVRKMIKSEQEIGKATAKYFFYSN  
IMNFFKTEITLANGEIRKRPLIETNGETGEIVWDKGRDFATVRKVL SMPQVNIKKTEVQTGGFS  
KESILPKRNSDKLIARKKDWDPKKYGGFDSPTVAYSVLVAKVEK GKSKKLKSVKELLGITIMER  
SSFENPIDFLEAKGYKEVKKDLIKLPKYSLFELENGRKRMLASAGELQKGNELALPSKYVNFL  
YLASHYEKLGKSPEDNEQKQLFVEQHKHYLDEIIEQISEFSKRVLADANLDKVL SAYNKHRDKPI  
REQAENIIHLFTLNLGAPAAFKYFDTTIDRKRYTSTKEVLDATLIHQSI TGLYETRIDLSQLGGDS  
AGGGGSGGGGSGGGGSGPKKKR KVAAAGSGRADALDDFDL DMLGSDALDDFDL DMLGSDA  
LDDFDL DMLGSDALDDFDL DML

MS2 (in blue) fused with a linker (black), P65 (red), and HSF1 (green):

ASNFTQFVLVDNNGTGDVTVAPSNFANGVAEWISSNSRSQAYKVTCSVRQSSAQKRKYTIKVE  
VPKVATQTVGGVELPVAAWRSYLN MELTIPIFATNSDCELIVKAMQGLLDGNPIPSAIAANSI  
YSAGGGGSGGGGSGGGGSGPKKKR KVAAAGSPSGQISNQALALAPSSAPVLAQTMVPSSAM  
VPLAQPPAPAPVLTGPPQSL SAPVPKSTQAGEGTLSEALLHLQFDADEDLGALLGNSTDPGV  
FTDLASVDNSEFQQLLNQGVSM SHSTAEPMLMEYPEAITRLVTGSQRPPDPAPTPLGTSGLPN  
GLSGDEDFSSIADMDFSALLSQISSGQGGGGSGF SVDT SALLDLFSPSVTPDMSLPDL DSSL  
ASIQELLSPQEPPRPEAENSSPDSGKQLVHYTAQPLFLLDPGSVDTGSNDLPVLFELGEGSYF  
SEGDGFAEDPTISLLTGSEPPKAKDPTVS

### Reference

Konermann, S., *et al.* (2015). Genome-scale transcriptional activation by an engineered CRISPR-Cas9 complex. *Nature*. **517(7536)**: 583-588.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.  
To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**  
Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)  
Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court West, Ste. B  
San Diego Ca 92121  
Tel: 1.858.202.1401  
Fax: 1.858.481.8694  
Email: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## Related Products

<u>Product</u>	<u>Cat. #</u>	<u>Size</u>
Cas9-Expressing Raji cells	78071	2 vials
Cas9-Expressing MDA-MB-231 cells	78069	2 vials
Cas9-Expressing A549 cells	78072	2 vials
Cas9-Expressing HCT116 cells	78073	2 vials
Cas9 Lentivirus (puromycin selection)	78066	500 µl x 2
Cas9, His-tag (S. pyogenes)	100206-1	50 µg

## Notes

*The CRISPR/CAS9 technology is covered under numerous patents, including U.S. Patent Nos. 8,697,359 and 8,771,945, as well as corresponding foreign patents applications, and patent rights.*

## License Disclosure

Purchase of this cell line grants you with a 10-year license to use this cell line in your immediate laboratory, for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make the cell line available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit use of this cell line in humans or for therapeutic or drug use. The license does not permit modification of the cell line in any way. Inappropriate use or distribution of this cell line will result in revocation of the license and result in an immediate cease of sales and distribution of BPS products to your laboratory. BPS does not warrant the suitability of the cell line for any particular use, and does not accept any liability in connection with the handling or use of the cell line. Modifications of this cell line, transfer to another facility, or commercial use of the cells may require a separate license and additional fees; contact [sales@bpsbioscience.com](mailto:sales@bpsbioscience.com) for details. Publications using this cell line should reference BPS Bioscience, Inc., San Diego.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.  
To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)