



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Description

The Mouse TIGIT: Mouse CD155 (mTIGIT:mCD155) Homogeneous Assay Kit is designed to measure the inhibition of mTIGIT binding to mCD155. The mTIGIT:mCD155 Homogeneous Assay Kit comes in a convenient AlphaLISA® format with purified biotinylated Mouse TIGIT, His-tagged Mouse CD155, and assay buffer to perform a total of 384 reactions. With this kit, only three simple steps on a microtiter plate are required. First, a sample containing mTIGIT and an inhibitor of choice is incubated with the mCD155 for 60 minutes. Next, acceptor beads are added, then donor beads, followed by reading the Alpha-counts.

Background

T-cell immunoreceptor with Ig and ITIM domains (TIGIT) is a receptor that is expressed on the surface of T cells and NK cells that binds to CD155 and CD112 on the surface of dendritic cells. Binding of TIGIT with CD155 or CD112 results in inhibition of T cell and NK cell activation. Antibodies and other agents that inhibit this signaling pathway have been shown to increase the immune response, especially in the case of certain cancers.

Applications

Useful for screening for inhibitors of TIGIT binding to CD155.

Supplied Materials

Catalog #	Name	Amount	Storage
79269	Mouse TIGIT-Fc-biotin	15 µg	-80°C
71167	Mouse CD155-His	30 µg	-80°C
79311	3x Immuno Buffer 1	4 ml	-20°C

Materials Required but Not Supplied

Name	Catalog #
AlphaLISA Ni Chelate Acceptor beads, 5 mg/ml	PerkinElmer #AL108C
AlphaScreen Streptavidin-conjugated Donor beads, 5 mg/ml	PerkinElmer #6760002S
Optiplate-384	PerkinElmer #6007290
AlphaScreen microplate reader	
Adjustable micropipettor and sterile tips	

Storage Conditions

This assay kit will perform optimally for up to one year from date of receipt when the materials are stored as directed.

Safety

This product is for research purposes only and not for human or therapeutic use. This product should be considered hazardous and is harmful by inhalation, in contact with skin, eyes, clothing, and if swallowed. If contact occurs, wash thoroughly.

Contraindications

Only limited amounts of DMSO can be included, as it has been shown to disrupt mTIGIT:mCD155 interaction. Avoid green and blue dyes that absorb light in the AlphaScreen signal emission range (520-620 nm), such as Trypan Blue. Avoid the use of the potent singlet oxygen quenchers such as sodium azide (NaN₃) or metal ions (Fe²⁺, Fe³⁺, Cu²⁺, Zn²⁺ and Ni²⁺). The presence of >1% RPMI 1640 culture medium leads to a signal reduction due to the presence of excess biotin and iron in this medium. MEM, which lacks these components, does not affect AlphaScreen assays.

Assay Protocol

All samples and controls should be tested in duplicate. Use slow shaking for all incubations.

Preparing Your Reagents

1. Dilute one part 3x Immuno Buffer 1 with 2 parts of distilled water (3-fold dilution) to make 1x Immuno Buffer 1. Make only a sufficient quantity needed for the assay; store remaining stock solution in aliquots at -20°C.
2. Thaw mTIGIT-biotin on ice. Upon first thaw, briefly spin tube containing protein to recover full contents of the tube. Aliquot the protein into single use aliquots. Store remaining undiluted protein in aliquots at -80°C immediately. Note: mTIGIT-biotin is very sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted protein.
3. Dilute mTIGIT-biotin in 1x Immuno Buffer 1 to 16 ng/μl. Keep diluted proteins on ice until use. Discard any remaining unused diluted protein after use.
4. Prepare the master mixture: N wells × (2 μl 3x Immuno Buffer 1 + 2 μl diluted mTIGIT-biotin + 2 μl distilled water). Add 6 μl of master mixture to every well.

Component	Blank	Positive Control	Test Inhibitor
3x Immuno Buffer 1	2 μl	2 μl	2 μl
mCD155-His (32 ng/μl)	2 μl	2 μl	2 μl
Distilled water	2 μl	2 μl	2 μl
Test Inhibitor	-	-	2 μl
Inhibitor buffer (no inhibitor)	2 μl	2 μl	-
1x Immuno Buffer 1	2 μl	-	-
mTIGIT-biotin (16 ng/μl)	-	2 μl	2 μl
Total	10 μl	10 μl	10 μl

5. Add 2 μl of inhibitor solution to each well designated "Test Inhibitor". For the "Positive Control" and "Blank", add 2 μl of the same solution without inhibitor (inhibitor buffer). Note: If possible, keep final DMSO concentration below 0.5%.
6. Add 2 μl of 1x Immuno Buffer 1 to the well designated "Blank".
7. Incubate at room temperature for 30 minutes with slow shaking.
8. Thaw mCD155-His on ice. Upon first thaw, briefly spin tube containing protein to recover full contents of the tube. Aliquot the protein into single use aliquots. Store remaining undiluted protein in aliquots at -80°C immediately. Note: mCD155-His is very sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted protein.

- Dilute mCD155-His in 1x Immuno Buffer 1 to 32 ng/μl. Keep diluted protein on ice until ready to use. Discard any remaining unused diluted protein after use.

Initiating the reaction

Initiate reaction by adding 2 μl of diluted mCD155-His prepared as described above to each well designated “Positive Control” and “Test Inhibitor”. Incubate at room temperature for 60 minutes.



Protect your samples from direct exposure to light. Photobleaching will occur.


- Dilute Ni Chelate Acceptor beads (PerkinElmer #AL108C) 250-fold with 1x Immuno Buffer 1. Add 10 μl per well. Shake plate briefly. Incubate at room temperature for 30 minutes.
- Dilute Streptavidin-conjugated donor beads (PE #6760002S) 125-fold with 1x Immuno Buffer 1. Add 10 μl per well. Incubate at room temperature for 30 minutes.
- Read Alpha-counts.



Due to lot to lot variability in AlphaScreen® bead performance, it may be necessary to optimize assay conditions. For example, slight adjustments to mTIGIT-biotin or mCD155-His concentrations may improve signal-to-noise ratio.

Example Results

mTIGIT-mCD155 Binding Activity



Mouse TIGIT:Mouse CD155 inhibition, measured using the mTIGIT:mCD155 Inhibitor Screening Assay Kit (BPS Bioscience Bioscience #78141) and an anti-TIGIT antibody (BPS Bioscience #71218). Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at support@bpsbioscience.com

General considerations

“Blank” Control: The “Blank” control is important to determine the background absorbance in the assay.

Visit bpsbioscience.com/assay-kits-faq for detailed troubleshooting instructions. For all further questions, please email support@bpsbioscience.com

References

1. Yu, X., *et al.*, *Nat. Immunol.* 2009; **10(1)**: 48-57.
2. Stanietsky, N., *et al.*, *Proc. Natl. Acad. Sci.* 2009; **106(42)**: 17858-17863.

Related Products

<i>Products</i>	<i>Catalog #</i>	<i>Size</i>
Mouse TIGIT, Fc-fusion (IgG1), Avi-tag, Biotin-labeled	#79269	25 µg; 50 µg
Mouse CD155 (PVR), His-tag	#71167	100 µg
Mouse CD155 (PVR), His-tag, Biotin-labeled	#71168	50 µg
Human TIGIT	#71218	100 µg
Human TIGIT, Fc fusion, Biotin-labeled	#71251	50 µg
Human TIGIT, Fc fusion	#71186	100 µg
Human CD112, His-tag	#71197	100 µg
Human CD112, His-tag, Biotin-labeled	#71234	50 µg
Human CD155 (PVR), His-tag	#71181	100 µg
Human CD112, His-tag, Biotin-labeled	#71234	50 µg
Human CD155 (PVR), His-tag	#71181	100 µg
Human TIGIT:CD155 Homogenous Assay Kit	#72029	384 rxns
Human TIGIT:CD112 Homogenous Assay Kit	#72030	384 rxns

AlphaScreen® and AlphaLISA® are registered trademarks of PerkinElmer, Inc.