



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet

NFAT Reporter – Hek293 Cell line *PKC/ Ca²⁺ Pathway* Catalog #: 79298

Background

The protein kinase c (PKC)/ Ca²⁺ response pathway leads to activation of the transcription factor nuclear factor of activator T cells (NFAT). NFAT is regulated by Ca²⁺ and the Ca²⁺/calmodulin-dependent serine phosphatase calcineurin. NFAT proteins are phosphorylated and reside in the cytoplasm in resting cells; upon stimulation, they are dephosphorylated by calcineurin, translocate to the nucleus, and induce gene expression.

Product Description

The NFAT Reporter – Hek293 cell line contains a firefly luciferase gene under the control of NFAT response element stably integrated into Hek293 cells. This cell line is validated for the response to the stimulation of phorbol 12-myristate 13-acetate (PMA) with ionomycin.

Application

- Monitor intracellular calcium levels.
- Screen for activators or inhibitors of the PKC/ Ca²⁺ pathway.

Format

Each vial contains ~2 X 10⁶ cells in 1 ml of 10% DMSO.

Storage

Immediately upon receipt, store in liquid nitrogen.

Functional Validation and Assay Performance

The following assays are designed for 96-well format. To perform the assay in different tissue culture formats, the cell number and reagent volumes should be scaled appropriately.

Mycoplasma testing

The cell line has been screened using the PCR-based VenorGeM[®] Mycoplasma Detection kit (Sigma-Aldrich, #MP0025) to confirm the absence of Mycoplasma species.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Culture Conditions:

Thaw Medium 1 (BPS Bioscience, #60187): Optimized cell culture medium for thawing and plating HEK293 cells. Includes 10% FBS, non-essential amino acids, sodium pyruvate and 1% Penicillin/Streptomycin.

Growth Medium 1B (BPS Cat. #79531): Thaw Medium 1 (BPS Bioscience, #60187) with 400 µg/ml of Geneticin (Life Technologies #11811031).

Cells should be grown at 37°C with 5% CO₂ using Growth Medium 1B. It may be necessary to adjust the percentage of CO₂ in the incubator depending on the NaHCO₃ level in the basal medium.

It is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, transfer to a tube containing 10 ml of Thaw Medium 1, spin down cells, remove supernatant and resuspend cells in pre-warmed Thaw Medium 1. Transfer the resuspended cells to a T25 flask and incubate at 37°C in a 5% CO₂ incubator. At first passage, switch to Growth Medium 1B. Cells should be split before they reach 100% confluence. To passage the cells, rinse cells with phosphate buffered saline (PBS), detach cells from culture vessel with Trypsin/EDTA, add Growth Medium 1B and transfer to a tube, spin down cells, resuspend cells and seed appropriate aliquots of cell suspension into new culture vessels.

Note: Just after thawing and at low density, the cells may grow at a slower rate. It is recommended to split the cells with ~ 1:4 ratio at the beginning of culturing. After several passages, the cell growth rate increases and the cells can be split with 1:8 -1:20 ratio weekly.

Functional Validation and Assay Performance

The following assays are designed for 96-well format. To perform the assay in different tissue culture formats, the cell number and reagent volumes should be scaled appropriately.

Materials Required but Not Supplied

- PMA (LC Laboratories, #P1680) Prepare stock solution in DMSO.
- Ionomycin (Sigma, #I3909) Prepare stock solution in DMSO.
- Assay medium: Thaw Medium 1 (no Geneticin)
- Growth Medium 1B (BPS Cat. #79531)
- 96-well tissue culture plate or 96-well tissue culture-treated white clear-bottom assay plate
- RO-31-8220 (EMD Millipore, #557520) or other PKC inhibitor
- ONE-Step Luciferase Detection Reagents (BPS Bioscience, #60690) for measuring firefly luciferase activity.
- Luminometer

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**


Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

A. Response of NFAT Reporter – Hek293 cells to TCR crosslinkers.

1. Harvest NFAT Reporter – Hek293 cells from culture in Growth Medium 1B and seed cells at a density of ~ 30,000 cells per well into white clear-bottom 96-well microplate in 45 µl of assay medium.
2. Incubate the plate at 37°C in a CO₂ incubator overnight (~18-24 hours).
3. Dilute TCR crosslinker (PMA with ionomycin) into assay medium at 10x desired final concentration and add 5 µl of dilution to each well. (We recommend a starting concentration around 30 nM PMA and 1 µM Ionomycin.) The final DMSO concentration can be up to 0.5%.
 Add 5 µl of assay medium with same concentration of DMSO but without the crosslinker to the unstimulated control wells.
 Add 50 µl of assay medium with DMSO to cell-free control wells (for determining background luminescence).
 Set up each treatment in at least triplicate.
4. Incubate cells at 37°C in a CO₂ incubator overnight (~18 hours).
5. The next day, perform luciferase assay using the ONE-Step luciferase assay system: Add 50 µl of ONE-Step Luciferase reagent per well and rock at room temperature for ~15 minutes. Measure luminescence using a luminometer.
6. Data Analysis: Subtract the average background luminescence (cell-free control wells) from the luminescence reading of all wells.
 The fold induction of NFAT luciferase reporter expression = background-subtracted luminescence of stimulated well / average background-subtracted luminescence of unstimulated control wells.

Figure 1. NFAT Reporter (Luc) – Hek293 cell response to TCR Crosslinker Ionomycin with PMA.



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
 To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
 Or you can Email us at: info@bpsbioscience.com
 Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

B. Response of NFAT Reporter – Hek293 cells to Inhibitors of TCR crosslinkers.

1. Harvest NFAT Reporter – Hek293 cells from culture in Growth Medium 1B and seed cells at a density of ~ 30,000 cells per well into white clear-bottom 96-well microplate in 45 μ l of assay medium.
2. Incubate the plate at 37°C in a CO₂ incubator overnight (~18-24 hours).
3. Dilute inhibitor (RO-31-8220) into assay medium to 10x desired final concentration and add 5 μ l of dilution to each well. The final DMSO concentration can be up to 0.5%. Add 5 μ l of assay medium with same concentration of DMSO but without the inhibitor to the un-inhibited control wells.
4. Dilute TCR crosslinker (PMA with ionomycin) into assay medium to 10x desired final concentration (We recommend a starting concentration around 30 nM PMA and 1 μ M Ionomycin) and add 5 μ l of dilution to each well. The final DMSO concentration can be up to 0.5%.
Add 5 μ l of assay medium with same concentration of DMSO but without the crosslinker to the unstimulated control wells.
Add 55 μ l of assay medium with DMSO to cell-free control wells (for determining background luminescence).
Set up each treatment in at least triplicate.
5. Incubate cells at 37°C in a CO₂ incubator overnight (~18 hours).

The next day, perform luciferase assay using the ONE-Step luciferase assay system: Add 55 μ l of One-Step Luciferase reagent per well and rock at room temperature for ~15 minutes. Measure luminescence using a luminometer.

6. Data Analysis: Subtract the average background luminescence (cell-free control wells) from the luminescence reading of all wells.

The percent luminescence of NFAT luciferase reporter expression is calculated by dividing each condition by the background-subtracted luminescence of the stimulated well. The background-subtracted luminescence of cells stimulated with 1 μ M Ionomycin and 30 nM PMA in the absence of RO-31-8220 was set at 100%.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.


To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

Figure 2. Inhibition of PKC by RO-31-8220 in PMA and Ionomycin induced NFAT reporter (Luc)-HEK293 cells.

The results are shown as percentage of luminescence. The background-subtracted luminescence of cells stimulated with 1 μ M Ionomycin and 30 nM PMA in the absence of RO-31-8220 was set at 100%.



Related Products

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
NFAT Reporter (Luc) – Jurkat Recombinant Cell Line	60621	2 vials
TIGIT / NFAT Reporter - Jurkat Cell Line	60538	2 vials
PD-1 / NFAT Reporter - Jurkat Recombinant Cell Line	60535	2 vials
LAG3 / NFAT Reporter - Jurkat Recombinant Cell Line	71278	2 vials
PKC α (PKCalpha), GST-tag	40157	10 μ g
PKC β I (PKCbeta1), GST-tag	40158	10 μ g
PKC β II (PKCbeta1), GST-tag	40159	10 μ g
PKC γ (PKCgamma), GST-tag	40160	10 μ g
ERK Signaling Pathway SRE Reporter – HEK293 Cell Line	60406	2 vials
Hedgehog Pathway Gli Reporter – NIH3T3 Cell Line	60409	2 vials
JAK/STAT Signaling Pathway ISRE Reporter – HEK293 Cell Line	60510	2 vials
JNK Signaling Pathway AP1 Reporter – HEK293 Cell Line	60405	2 vials
NK-kB Reporter (Luc) – HEK293 Cell Line	60650	2 vials
RARalpha Reporter (Luc) – HEK293 Cell Line	60503	2 vials
Wnt Signaling Pathway TCF/LEF Reporter – HEK293 Cell Line	60501	2 vials

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com