



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet

Transfection Collection™ - NF-κB Reporter Cellular Assay Pack (HEK293) Catalog #: 79327

Product Description

The NF-κB Reporter Cellular Assay Pack provides all the key reagents required to monitor the activity of the nuclear factor Kappa B (NF-κB) signal transduction pathways. The pack contains the NF-κB Reporter (Luc)-HEK293 Recombinant Cell Line, a luciferase reporter cell line that contains a firefly luciferase gene under the control of four copies of the NF-κB response element located upstream of the minimal TATA promoter. After activation by pro-inflammatory cytokines or stimulants of lymphokine receptors, endogenous NF-κB transcription factors bind to the DNA response elements, inducing transcription of the luciferase reporter gene. This cell line is validated for the response to TNFalpha and to treatment with NF-κB inhibitor, evodiamine.

Additionally, the pack includes cell culture medium (Thaw Medium 1) that has been optimized for use with HEK293 cells. Thaw Medium 1 includes 10% fetal bovine serum, non-essential amino acids, sodium pyruvate, and 1% Pen/Strep. Finally, the pack provides the ONE-Step™ Luciferase Detection System. This reagent provides highly sensitive, stable detection of firefly (*Photinus pyralis*) luciferase activity. The ONE-Step luciferase reagent can be used directly in cells in growth medium and can be detected with any luminometer; automated injectors are not required.

Application

The NF-κB reporter cell line is designed for screening inhibitors of NF-κB and for monitoring NF-κB signaling pathway activity.

Components

Cat. #	Component	Amount	Storage
60650	NF-κB Reporter (Luc) - HEK293 Cell Line	2 vials*	liquid nitrogen
60690-1	ONE-Step Luciferase Buffer (Component A)	10 ml	-20°C
	ONE-Step Luciferase Reagent Substrate, 100x (Component B)	100 µl	-20°C <i>Protect from light</i>
60187	Thaw Medium 1	100 ml	+4°C

*Each vial contains ~2 X 10⁶ cells in 1 ml of 10% DMSO in FBS.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

General Culture conditions

Thaw Medium 1 (BPS Bioscience, #60187): Medium optimized for culturing HEK293 cells.

Growth Medium 1C (BPS Bioscience, #79532): Thaw Medium 1 plus 50 µg/ml of Hygromycin B (Hyclone, #SV30070.01)

Cells should be grown at 37°C with 5% CO₂ using Growth Medium 1C (Thaw Medium 1 plus Hygromycin). NF-κB reporter (Luc)-HEK293 cells should exhibit a typical cell division time of 24 hours.

To thaw the cells, it is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, transfer to a tube containing 10 ml of Thaw Medium 1 (**no Hygromycin**), spin down cells, resuspend cells in pre-warmed Thaw Medium 1 (**no Hygromycin**), transfer resuspended cells to a T25 flask and culture in a 37°C CO₂ incubator. At first passage switch to Growth Medium 1C (**contains Hygromycin**). Cells should be split before they reach complete confluence.

To passage the cells, rinse cells with phosphate buffered saline (PBS), and detach cells from culture vessel with 0.05% Trypsin/EDTA. Add Growth Medium 1C and transfer to a tube, spin down the cells, then resuspend cells and seed appropriate aliquots of cell suspension into new culture vessels. Subcultivation ration: 1:10 to 1:20 weekly.

To freeze down the cells, rinse cells with phosphate buffered saline (PBS), and detach cells from culture vessel with Trypsin/EDTA. Add Growth Medium 1C and transfer to a tube, spin down cells, and resuspend in freezing medium (10% DMSO + 90% FBS). Place at -80°C overnight and place in liquid nitrogen the next day. Alternatively, vials may be placed directly in liquid nitrogen.

Mycoplasma testing

The cell line has been screened using the PCR-based VenorGeM Mycoplasma Detection kit (Sigma-Aldrich, #MP0025) to confirm the absence of *Mycoplasma* species.

Assay performance

The following assays are designed for 96-well format. To perform the assay in different tissue culture formats, cell number and reagent volume should be scaled appropriately.

Materials Required but Not Supplied

- TNFα (BPS Bioscience, #90245)
- IL-1β (BPS Bioscience, #90168)
- Growth Medium 1C (BPS Bioscience, #79532)
- Evodiamine (Abcam, #142427): inhibitor of NF-κB activation

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

- 96-well tissue culture treated white clear-bottom assay plate (Corning, #3610)
- Luminometer

A. TNF α dose response

1. Harvest NF- κ B reporter (Luc)-HEK293 cells and seed cells at a density of 30,000 cells per well into white opaque 96-well microplate in 45 μ l of Thaw Medium 1 (without Hygromycin).
2. Incubate cells at 37°C with 5% CO₂ overnight.
3. The next day, prepare threefold serial dilution of TNF α in Thaw Medium 1 and add 5 μ l to TNF α -stimulated wells.
Add 5 μ l of Thaw Medium 1 to the unstimulated control wells (for measuring uninduced level of NF- κ B reporter activity).
Add 50 μ l of Thaw Medium 1 to cell-free control wells (for determining background luminescence).
4. Incubate at 37°C with 5% CO₂ for 5-6 hours.
5. Perform the luciferase detection assay using the ONE-Step™ Luciferase Assay System according to the protocol below:

Luciferase Detection Procedure

6. Thaw Luciferase Reagent Buffer (**Component A**) by placing the reagent in a room temperature water bath. Equilibrate the buffer to room temperature and mix well before use. Note: Luciferase Reagent Buffer must be at ~ room temperature (20- 25°C) before use.
7. Calculate the amount of Luciferase Reagent needed for the experiment (**Component A + Component B**). Immediately prior to performing the experiment, prepare the luciferase assay working solution by diluting Luciferase Reagent Substrate (**Component B**) into Luciferase Reagent Buffer (**Component A**) at a 1:100 ratio and mix well. Avoid exposing to excessive light. *Only use enough of each component for the experiment, remaining **Component A** and **Component B** should be stored separately at -20°C.*
8. Remove multi-well plate containing mammalian cells from incubator. *Note: plates must be compatible with luminescence measurement with luminometer being used.*
9. Add 50 μ l of luciferase assay working solution (**Component A + Component B**) to the culture medium in each well.
10. Gently rock the plates for \geq 15 minutes at room temperature. Measure firefly luminescence using a luminometer.

The signal under these conditions will be stable for more than 2 hours at room temperature. For maximal light intensity, measure samples within 1 hour of reagent addition.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com




6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Analysis: Background luminescence is a characteristic of luminometer performance, therefore, background luminescence must be subtracted from all readings for accuracy. Obtain the background-subtracted luminescence by subtracting the average background luminescence (cell-free control wells) from the luminescence reading of all wells.

The fold induction of NF- κ B luciferase reporter expression = background-subtracted luminescence of TNF α -stimulated well / average background-subtracted luminescence of unstimulated control wells

Figure 1. TNF α dose response in NF- κ B reporter (Luc)-HEK293 cells. The results were shown as fold induction of luciferase reporter expression. Fold induction was determined by comparing values against the mean value for control cells without TNF α treatment.

The EC50 of TNF α in this cell line is ~2ng/ml.



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

B. IL-1 β dose response

1. Harvest NF- κ B reporter (Luc)-HEK293 cells from culture in Growth Medium 1C and seed cells at a density of 30,000 cells per well into white opaque 96-well microplate in 45 μ l of Thaw Medium 1 (no hygromycin).
2. Incubate cells at 37°C with 5% CO₂ overnight.
3. The next day, prepare threefold dilutions of IL-1 β in Thaw Medium 1 and add 5 μ l to IL-1 β -stimulated wells.

Add 5 μ l of Thaw Medium 1 to the unstimulated control wells (for measuring uninduced level of NF- κ B reporter activity).

Add 50 μ l of Thaw Medium 1 to cell-free control wells (for determining background luminescence).

5. Incubate at 37°C with 5% CO₂ for 5-6 hours.

Perform the luciferase detection assay using the ONE-Step™ Luciferase Assay System according to the protocol below:

Luciferase Detection Procedure

6. Thaw Luciferase Reagent Buffer (**Component A**) by placing the reagent in a room temperature water bath. Equilibrate the buffer to room temperature and mix well before use. Note: Luciferase Reagent Buffer must be at ~ room temperature (20- 25°C) before use.
7. Calculate the amount of Luciferase Reagent needed for the experiment (**Component A + Component B**). Immediately prior to performing the experiment, prepare the luciferase assay working solution by diluting Luciferase Reagent Substrate (**Component B**) into Luciferase Reagent Buffer (**Component A**) at a 1:100 ratio and mix well. Avoid exposing to excessive light. *Only use enough of each component for the experiment, remaining **Component A** and **Component B** should be stored separately at -20°C.*
8. Remove multi-well plate containing mammalian cells from incubator. *Note: plates must be compatible with luminescence measurement with luminometer being used.*
9. Add 50 μ l of luciferase assay working solution (**Component A + Component B**) to the culture medium in each well.
10. Gently rock the plates for \geq 15 minutes at room temperature. Measure firefly luminescence using a luminometer.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com


The signal under these conditions will be stable for more than 2 hours at room temperature. For maximal light intensity, measure samples within 1 hour of reagent addition.

Data Analysis: Background luminescence is a characteristic of luminometer performance, therefore, background luminescence must be subtracted from all readings for accuracy. Obtain the background-subtracted luminescence by NF- κ B subtracting the average background luminescence (cell-free control wells) from the luminescence reading of all wells.

The fold induction of luciferase reporter expression = background-subtracted luminescence of IL-1 β -stimulated well / average background-subtracted luminescence of unstimulated control wells

Figure 2. IL-1 β dose response in NF- κ B reporter (Luc)-HEK293 cells. The results were shown as fold induction of luciferase reporter expression. Fold induction was determined by comparing values against the mean value for control cells without IL-1 β treatment.

The EC50 of IL-1 β in this cell line is ~0.077 ng/ml.



OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

C. Inhibition of TNF α -induced NF- κ B activity

1. Harvest NF- κ B reporter (Luc)-HEK293 cells from culture in Growth Medium 1C and seed cells at a density of 30,000 cells per well into white opaque 96-well microplate in 45 μ l of Thaw Medium 1 (no hygromycin).
2. Incubate cells at 37°C with 5% CO₂ overnight.
3. Add 5 μ l of Thaw Medium 1 with or without NF- κ B inhibitor to wells. Incubate cells overnight at 37°C with 5% CO₂. [Alternately, inhibitor may be added to cells and incubated at 37°C with 5% CO₂ for 2-4 hours before addition of TNF α .]
4. The next day, prepare threefold serial dilution of TNF α in Thaw Medium 1 and add 5 μ l to TNF α -stimulated wells.
Add 5 μ l of Thaw Medium 1 to the unstimulated control wells.
Add 55 μ l of Thaw Medium 1 to cell-free control wells.

Incubate at 37°C with 5% CO₂ for 5-6 hours.

5. Perform the luciferase detection assay using the ONE-Step™ Luciferase Assay System according to the protocol below:

Luciferase Detection Procedure

6. Thaw Luciferase Reagent Buffer (**Component A**) by placing the reagent in a room temperature water bath. Equilibrate the buffer to room temperature and mix well before use. Note: Luciferase Reagent Buffer must be at ~ room temperature (20- 25°C) before use.
7. Calculate the amount of Luciferase Reagent needed for the experiment (**Component A + Component B**). Immediately prior to performing the experiment, prepare the luciferase assay working solution by diluting Luciferase Reagent Substrate (**Component B**) into Luciferase Reagent Buffer (**Component A**) at a 1:100 ratio and mix well. Avoid exposing to excessive light. *Only use enough of each component for the experiment, remaining **Component A** and **Component B** should be stored separately at -20°C.*
8. Remove multi-well plate containing mammalian cells from incubator. *Note: plates must be compatible with luminescence measurement with luminometer being used.*
9. Add 55 μ l of luciferase assay working solution (**Component A + Component B**) to the culture medium in each well.
10. Gently rock the plates for \geq 15 minutes at room temperature. Measure firefly luminescence using a luminometer.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com


Please visit our website at: www.bpsbioscience.com

The signal under these conditions will be stable for more than 2 hours at room temperature. For maximal light intensity, measure samples within 1 hour of reagent addition.

Data Analysis: Background luminescence is a characteristic of luminometer performance, therefore, background luminescence must be subtracted from all readings for accuracy. Obtain the background-subtracted luminescence by subtracting the average background luminescence (cell-free control wells) from the luminescence reading of all wells.

The fold induction of NF- κ B luciferase reporter expression = background-subtracted luminescence of TNF α -stimulated well / average background-subtracted luminescence of unstimulated control wells

Figure 3. Inhibition of TNF α -induced NF- κ B activity by NF- κ B inhibitor, evodiamine, in NF- κ B reporter (Luc)-HEK293 cells



References

1. Pessara U, Koch N (1990) Tumor necrosis factor alpha regulates expression of the major histocompatibility complex class II-associated invariant chain by binding of an NF- κ B-like factor to a promoter element. *Mol. Cell. Biol.* **10(8)**:4146-4154.
2. Baeuerle PA (1998) Pro-inflammatory signaling: last pieces in the NF- κ B puzzle? *Curr. Biol.* **8(1)**:R19-R22.
3. Takada Y, Kobayashi Y, Aggarwal BB (2005) Evodiamine abolishes constitutive and inducible NF- κ B activation by inhibiting I κ B α kinase activation, thereby suppressing NF- κ B-regulated antiapoptotic and metastatic gene expression, up-regulating apoptosis, and inhibiting invasion. *J. Biol. Chem.* **280(17)**:17203-17212

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.
 To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court West, Suite B
San Diego, CA 92121
Tel: 1.858.829.3082
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

License Disclosure

Purchase of this cell line grants you with a 10-year license to use this cell line in your immediate laboratory, for research use only. This license does not permit you to share, distribute, sell, sublicense, or otherwise make the cell line available for use to other laboratories, departments, research institutions, hospitals, universities, or biotech companies. The license does not permit use of this cell line in humans or for therapeutic or drug use. The license does not permit modification of the cell line in any way. Inappropriate use or distribution of this cell line will result in revocation of the license and result in an immediate cease of sales and distribution of BPS products to your laboratory. BPS does not warrant the suitability of the cell line for any particular use, and does not accept any liability in connection with the handling or use of the cell line. Modifications of this cell line, transfer to another facility, or commercial use of the cells may require a separate license and additional fees; contact sales@bpsbioscience.com for details. Publications using this cell line should reference BPS Bioscience, Inc., San Diego.

Refills

<u>Product</u>	<u>Cat. #</u>	<u>Size</u>
NF- κ B Reporter (Luc)-HEK293 Recombinant Cell Line	60650	2 vials
ONE-Step Luciferase Assay Detection System	60690-1	10 ml
ONE-Step Luciferase Assay Detection System	60690-2	100 ml
ONE-Step Luciferase Assay Detection System	60690-3	1 L
Thaw Medium 1	60187	100 ml

Related Products

<u>Product</u>	<u>Cat. #</u>	<u>Size</u>
NF- κ B Reporter (Luc) - Jurkat Recombinant Cell Line	60651	2 vials
NF- κ B Reporter (Luc) - CHO-K1 Recombinant Cell Line	60622	2 vials
NF- κ B Reporter (Luc) - A549 Recombinant Cell Line	60625	2 vials
NF- κ B Reporter (Luc) - HCT116 Recombinant Cell Line	60623	2 vials
Transfection Collection™ : NF- κ B Transient Pack	79268	500 rxns
NF- κ B Reporter Kit	60614	500 rxns
TLR9/ NF- κ B Reporter – HEK293 Recombinant Cell Line	60485	2 vials
OX40/ NF- κ B Reporter – HEK293 Recombinant Cell Line	60482	2 vials
GITR/ NF- κ B Reporter – HEK293 Recombinant Cell Line	60546	2 vials
CD40/ NF- κ B Reporter – HEK293 Recombinant Cell Line	60626	2 vials
Interleukin-1 beta (IL-1 β), human	90168-B	10 μ g
TNF α , human	90244-A	10 μ g
TNF α , mouse	90246-B	20 μ g

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com