



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com


www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet **Sortase Sampling Kit** Catalog #79709

DESCRIPTION:

Sortases catalyze a transpeptidation reaction, where the enzyme cleaves the amide bond between the threonine and glycine of the sortase recognition sequence (LPXTG for *S. aureus* Sortase A), generating a thioacyl intermediate. Subsequently, this intermediate is resolved by the N-terminus of an oligoglycine nucleophile, creating a new peptide bond that links the substrate to the incoming nucleophile. Sortase variants are available that may exhibit enhanced activity with the substrate of interest. This kit provides wildtype Sortase and four different mutants to evaluate which sortase is best for the chosen application.



COMPONENTS:

Catalog #	Component	Amount	Storage
71086	Sortase wildtype enzyme	10 µg	-80°C Avoid freeze/thaw cycles
79393	Sortase 5Δ enzyme	10 µg	
71047	Sortase 6Δ enzyme	10 µg	
71048	Sortase 7Δ enzyme	10 µg	
72518	Sortase 8Δ enzyme	10 µg	
79394	Reaction Buffer	100 µl	4°C

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

Target protein (1 mg/mL, 10 µg per reaction)
 Poly-glycine label (5 mM, 1 µl per reaction, <10 kDa)
 37°C incubator
 Microcentrifuge
 Electrophoresis sample buffer (for SDS-PAGE) or 1 M EDTA solution (for mass spectroscopy)

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**
 Or you can Email us at: info@bpsbioscience.com
 Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

APPLICATIONS: Great for evaluating substrate specificity or optimizing the sortase ligation reaction.

STABILITY: At least 6 months from date of receipt, when stored as directed.

REFERENCE: Popp M. *Methods Mol Biol.* 2015; **1266**:185-98.

LABELING PROTOCOL:

- 1) Thaw a single vial of Sortase on ice.
- 2) While on ice, set up the labeling reaction as follows (add all reagents to the vial containing Sortase):

Reagent	Volume to add
Target Protein (1 mg/mL)	10 μ l
Reaction Buffer	2 μ l
Poly-Glycine Label (5 mM)	1 μ l
Sortase	10 μ l
Total	23 μl

- 3) Mix reaction well by gently pipetting up and down.
- 4) Incubate reaction at 37°C for 30 minutes and mix every 5-10 minutes (or place on shaker if available).
- 5) Stop reaction by adding SDS-PAGE sample buffer, or 2.3 μ l of 1 M EDTA (10 mM final).
- 6) Evaluate labeling by SDS-PAGE or mass spectroscopy.

For questions or concerns, please contact BPS Bioscience, Inc. at info@bpsbioscience.com

BPS is a leader in the custom expression of recombinant proteins. Contact us about custom projects for the expression of target proteins containing the sortase recognition sequence.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

RELATED PRODUCTS

<u>Product</u>	<u>Cat. #</u>	<u>Size</u>
Sortase wildtype enzyme	71086	10 µg
Sortase 5Δ enzyme	79393	10 µg
Sortase 6Δ enzyme	71047	10 µg
Sortase 7Δ enzyme	71048	10 µg
Sortase 8Δ enzyme	72518	10 µg
SiMPLe Protein Labeling Kit (Sortase Mediated Protein Ligation)	79392	3 units

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: [**info@bpsbioscience.com**](mailto:info@bpsbioscience.com)

Please visit our website at: [**www.bpsbioscience.com**](http://www.bpsbioscience.com)