



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Ct. West, Ste. B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet
PRMT7 Homogeneous Assay Kit (Histone H2b)
Catalog #79720
Size: 384 reactions

DESCRIPTION: The *PRMT7 Homogeneous Assay Kit (Histone H2b)* is designed to measure PRMT7 activity for screening and profiling applications. PRMT7 methylates arginine residues on various protein substrates and is involved in DNA transcription, RNA splicing, DNA repair, cell differentiation, and metastasis. The *PRMT7 Homogeneous Assay Kit* comes in a convenient AlphaLISA® format, with His-tagged PRMT7 substrate, primary antibody, methylation assay buffer, and purified PRMT7 for 384 enzyme reactions. The key to the *PRMT7 Homogeneous Assay Kit* is a highly specific antibody that recognizes methylated substrate. With this kit, only three simple steps on a microtiter plate are required for methyltransferase detection. First, a sample containing PRMT7 enzyme is incubated with the substrate. Next, acceptor beads and primary antibody are added, then donor beads, followed by reading the Alpha-counts.

COMPONENTS:

Catalog #	Component	Amount	Storage	
51054	PRMT7	40 µg	-80°C	<i>Avoid freeze/ thaw cycles!</i>
52120-A	20 µM S-adenosylmethionine	400 µl	-80°C	
52140Z3	Primary antibody 28	20 µl	-80°C	
	PRMT7 Histone substrate	200 µl	-80°C	
	4x PRMT7 assay buffer (add DTT before experiment)	3x1 ml	-20°C	
	4x PRMT7 Detection buffer	2 ml	-20°C	

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

AlphaLISA® anti-rIgG acceptor beads, 5 mg/ml (PerkinElmer #AL104C)
AlphaScreen® Nickel donor beads, 5 mg/ml (PerkinElmer #AS101D)
Optiplate-384 (PerkinElmer #6007290)

APPLICATIONS: Great for studying enzyme kinetics and HTS applications.

CONTRAINDICATIONS: Green and blue dyes that absorb light in the AlphaScreen® signal emission range (520-620 nm), such as Trypan Blue. Avoid the use of the potent singlet oxygen quenchers such as sodium azide (NaN₃) or metal ions (Fe²⁺, Fe³⁺, Cu²⁺, Zn²⁺ and Ni²⁺). The presence of >1% RPMI 1640 culture medium leads to a signal reduction due to the presence of excess biotin and iron in this medium. MEM, which lacks these components, does not affect AlphaScreen® assays.

STABILITY: At least one year from date of receipt when stored as directed.

REFERENCE(S): You Feng, *et al.* 2014. *J Biol Chem* **289**:32604.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Ct. West, Ste. B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

ASSAY PROTOCOL:

All samples and controls should be tested in duplicate.

Step 1:

- 1) Add 10 μ l 0.5M DTT (not provided) to 1-ml tube with **4x PRMT7 Assay Buffer**. Prepare **1x PRMT7 Assay Buffer** by adding 1 part of **4x PRMT7 Assay Buffer** to 3 parts water (v/v).
- 2) Prepare the master mixture: N wells x (2.0 μ l **4x PRMT7 Assay buffer** + 0.5 μ l **S-adenosylmethionine** (20 μ M) + 0.5 μ l **PRMT7 substrate** + 2 μ l water). Add 5 μ l to wells designated "Positive Control", "Test Sample", and "Blank". To wells labeled "Substrate Control", add 2 μ l **4x PRMT7 Assay buffer** + 0.5 μ l **S-adenosylmethionine** (20 μ M) + 2.5 μ l water.
- 3) Add 3 μ l of Inhibitor solution of each well labeled as "Test Inhibitor". We recommend to use **1x PRMT7 Assay buffer** for intermediate dilution of inhibitors. For the "Positive Control", "Substrate Control" and "Blank", add 3 μ l of the same solution without inhibitor (Inhibitor buffer).
- 4) Prepare **1x PRMT7 buffer** by adding 1 part of **4x PRMT7 Assay buffer** to 3 parts water (v/v).
- 5) Thaw **PRMT7** on ice. Upon first thaw, briefly spin tube containing enzyme to recover full content of the tube. Aliquot **PRMT7** enzyme into single use aliquots. Store remaining undiluted enzyme in aliquots at -80°C. *Note: PRMT7 is very sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted enzyme. Depending on the nature of inhibitor, pre-incubation with the enzyme may be necessary.*
- 6) Dilute **PRMT7** in **1X PRMT7 Assay buffer** at 50 ng/ μ l (100 ng/2 μ l). Keep diluted enzyme on ice until use. Discard any unused diluted enzyme after use.

	Positive Control	Test Sample	Substrate Control	Blank
4x PRMT7 Assay buffer	2 μ l	2 μ l	2 μ l	2 μ l
20 μ M S-adenosylmethionine	0.5 μ l	0.5 μ l	0.5 μ l	0.5 μ l
PRMT7 Histone substrate	0.5 μ l	0.5 μ l	-	0.5 μ l
H ₂ O	2 μ l	2 μ l	2.5 μ l	2 μ l
Test Inhibitor/Activator	-	3 μ l	-	-
Inhibitor Buffer (1xPRMT7 Assay buffer)	3 μ l	-	3 μ l	3 μ l
PRMT7 (50 ng/ μ l)	2 μ l	2 μ l	2 μ l	-
1x PRMT7 Assay buffer	-	-	-	2 μ l
Total	10 μl	10 μl	10 μl	10 μl

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Ct. West, Ste. B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

- 7) To the wells designated as "Blank", add 2 μ l of **1x PRMT7 Assay buffer**.
- 8) Initiate reaction by adding 2 μ l of diluted **PRMT7** enzyme to the wells designated "Positive Control", "Substrate Control", and "Test Sample". Incubate two hours at room temperature with slow shaking.

Protect your samples from direct exposure to light for steps 2 and 3!

(Alternatively, preincubate diluted enzyme with inhibitor first for 30 minutes, then initiate 2 hours reaction by addition of S-adenosylmethionine and Histone substrate.)

Step 2:

- 1) Dilute anti-Rabbit Acceptor beads (PerkinElmer #AL104C) 1:250-fold with **1x Detection buffer** (made by diluting **4x Detection buffer** 1:4 in distilled water). Add 5 μ l per well. Shake plate briefly.
- 2) Dilute **Primary antibody 28** 100-fold with **1x Detection buffer**. Add 5 μ l per well. Shake plate. Incubate 30 min at room temperature.
(Alternatively, dilute anti-Rabbit Acceptor beads (1:500) and Primary antibody 28 (1:200) with 1x Detection buffer in one step. Add 10 μ L of acceptor beads/antibody mixture per well.)


Step 3:

- 1) Dilute Nickel donor beads (PerkinElmer #AS101D) 125-fold with **1x Detection buffer**. Add 10 μ l per well. Incubate for 30-60 min. at room temperature.
- 2) Read Alpha-counts.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com

Examples of Assay Results:



PRMT7 enzyme activity, measured using the *PRMT7 Homogeneous Assay Kit (Histone H2b)*, BPS Bioscience Cat. #79720. Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.com

RELATED PRODUCTS:

PRMT1 Homogeneous Assay Kit	#52054	384 reactions
PRMT3 Homogeneous Assay Kit	#52055	384 reactions
PRMT6 Homogeneous Assay Kit	#52056	384 reactions
PRMT8 Homogeneous Assay Kit	#52058	384 reactions
PRMT1 Chemiluminescent Assay Kit	#52004L	96 reactions
PRMT3 Chemiluminescent Assay Kit	#52005L	96 reactions
PRMT4 Chemiluminescent Assay Kit	#52041L	96 reactions
PRMT5 Chemiluminescent Assay Kit	#52002L	96 reactions
PRMT6 Chemiluminescent Assay Kit	#52046	96 reactions
PRMT5 recombinant protein (HEK293)	#51045	20 µg
PRMT5/MEP50 recombinant protein (Sf9)	#51048	20 µg
PRMT1 recombinant protein (<i>E. coli</i>)	#51040	50 µg
PRMT1 recombinant protein (Sf9)	#51041	20 µg
PRMT3 recombinant protein	#51043	50 µg
PRMT4 (CARM 1) recombinant protein	#51047	20 µg
PRMT6 recombinant protein	#51049	20 µg
PRMT7 recombinant protein	#51054	20 µg
PRMT8 recombinant protein	#51052	20 µg
PRMT9 recombinant protein	#51053	20 µg

AlphaScreen® and AlphaLISA® are registered trademarks of PerkinElmer, Inc.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com