



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Data Sheet

Fluorogenic PTP1B (Catalytic Region) Assay Kit

Catalog #79764
Size: 96 reactions

DESCRIPTION: Protein phosphorylation is one of the most important post-translational modification processes. Phosphorylation is reversibly regulated by Protein Kinases (PKs) and Protein Phosphatases (PTPs). PTP1B (PTPN1) is known to catalyze dephosphorylation of insulin receptor kinases and plays a critical role in insulin signaling. The *Fluorogenic PTP1B (Catalytic Region) Assay Kit* is designed to inhibitors of the catalytic region of PTP1B in a homogeneous assay with no time-consuming washing steps. The PTP1B assay kit comes in a convenient 96-well format, with purified PTP1B enzyme, fluorogenic substrate, and PTP assay buffer for 100 enzyme reactions. Note: To identify inhibitors of the regulatory region of PTP1B, please use our *Fluorogenic PTP1B (Full Length) Assay Kit, #79766*

COMPONENTS:

| Catalog # | Component | Amount | Storage | |
|-----------|---|--------|------------------|---|
| 30010 | Recombinant Human PTP1B (1-321) | ≥1 µg | -80°C | Avoid freeze/ thaw cycles! |
| | 0.5 mM PTP Substrate | 50 µl | -80°C | |
| 79716 | 5X PTP Assay Buffer | 20 ml | -20°C | |
| | 0.5 M DTT | 20 µl | -20°C | |
| 79685 | Black, low binding black microtiter plate | 1 | Room Temperature | |

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

Fluorescent microplate reader capable of reading exc/em=360nm/460nm

APPLICATIONS: Great for studying enzyme kinetics and HTS applications.

STABILITY: One year from date of receipt when stored as directed.

REFERENCE(S):

Gee, K.R., *et al.*, *Anal Biochem*, 1999 Aug 15; **273(1)**:41-8.
Brown-Shimer, S., *et al.*, 1992 Jan 15; **52(2)**:478-82.

ASSAY PROTOCOL:

All samples and controls should be tested in duplicate.

- 1) Prepare **1X PTP Assay Buffer** with 1 mM DTT from **5X PTP Assay Buffer**. For example, add 200 µl **5X PTP Assay Buffer** and 2 µl 0.5M DTT to 798 µl distilled H₂O to make 1 ml **1X PTP Assay Buffer**.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**
Or you can Email us at: info@bpsbioscience.com
Please visit our website at: www.bpsbioscience.com

- 2) Prepare the master mixture: N wells × (24.5 µl **1X PTP Assay Buffer** (with DTT) + 0.5 µl 0.5 mM PTP Substrate).
- 3) Add 25 µl of master mixture to each well (Final concentration of the PTP substrate in a 50 µl reaction is 5 µM).
- 4) Prepare the inhibitor solution.

If the inhibitor compound is dissolved in DMSO, make a 100-fold higher concentration of the compound in DMSO than the highest concentration you want to test in the assay. Then make a 10 fold dilution in 1X assay buffer (at this step the compound concentration is 10-fold higher than the final concentration). If you want to run an IC50 or test lower concentrations of the compound, make a series of further dilutions using 1X assay buffer containing 10% DMSO, so the final concentration of DMSO will be 1% in all samples.

If the inhibitor compound is dissolved in water, make a solution of the compound in 1X assay buffer that is 10-fold higher than the final assay concentration.

- 5) Add 5 µl of the inhibitor solution to the well designed as “Test Sample”. Add 5 µl of the inhibitor buffer (without inhibitor) to the wells designed as “Blank”, and “Positive Control”.

| Component | Positive Control | Test Sample | Blank |
|----------------------------------|------------------|--------------|--------------|
| 1X assay buffer with DTT | 24.5 µl | 24.5 µl | 44.5 µl |
| Substrate | 0.5 µl | 0.5 µl | 0.5 µl |
| Test Inhibitor | – | 5 µl | – |
| Inhibitor Buffer (no inhibitor)* | 5 µl | – | 5 µl |
| PTP1B (2 pg/µl) | 20 µl | 20 µl | – |
| Total | 50 µl | 50 µl | 50 µl |

* Inhibitor buffer typically represents 1x PTP assay buffer with proper concentration of DMSO.

- 6) Thaw **Recombinant Human PTP1B (1-321)** on ice. Upon first thaw, briefly spin tube containing enzyme to recover full content of the tube. Aliquot **Recombinant Human PTP1B (1-321)** into single use aliquots. Store remaining undiluted enzyme in aliquots at -80°C. *Note: PTP1B enzyme is sensitive to freeze/thaw cycles. Do not re-use thawed aliquots or diluted enzyme.*
- 7) Dilute **Recombinant Human PTP1B (1-321)** in **1X PTP Assay Buffer** at 2 pg/µl (40 pg per reaction).

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.


To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com

- 8) Add 20 μ l diluted **Recombinant Human PTP1B (1-321)** solution to wells designated "Positive Control" and "Test Sample". Add 20 μ l 1X assay buffer to "Blank" wells.
- 9) Incubate at room temperature for 30 minutes. Measure the fluorescence intensity in a microtiter plate-reading fluorimeter capable of excitation at a wavelength 360 nm and detection of emission at a wavelength 460 nm. You can also measure the fluorescence intensity kinetically. "Blank" value is subtracted from all other values.

Example of Assay Results:



PTP1B enzyme activity, measured using the *Fluorogenic PTP1B (Catalytic Region) Assay Kit* (BPS Bioscience #79764). Fluorescence intensity was measured using a Tecan fluorescent microplate reader. *Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.com.*

Related Products

| <u>Product</u> | <u>Cat. #</u> | <u>Size</u> |
|--------------------------------------|---------------|-------------|
| Human PTP1B (1-321), GST-tag | 30010 | 20 μ g |
| PTP1B (1-321) Colorimetric Assay Kit | 30019 | 96 rxns. |
| Mouse PTP1B (1-321), GST-tag | 30012 | 20 μ g |
| Rat PTP1B (1-321), GST-tag | 30011 | 20 μ g |
| PTP1B (PTPN1) full length, GST-tag | 30009 | 20 μ g |
| 10x PTP1B Colorimetric Substrate | 79693 | 5 ml |

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401**, Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com