



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

Data Sheet
UBCH5c TR-FRET Assay Kit
Catalog # 79901
Size: 384 reactions

DESCRIPTION:

UBCH5c (UBE2D3) is an E2 conjugating ubiquitin protein that functions in the ubiquitination of the tumor-suppressor protein p53. It is a promising drug target in cancer immunotherapy. The *UBCH5c TR-FRET Assay Kit* is designed to measure ubiquitination activity in a homogeneous 384 reaction format. It utilizes biotin-labeled Ubiquitin and a terbium-labeled antibody recognizing the His-tagged UBCH5c protein to complete the TR-FRET pairing. This FRET-based assay requires no time-consuming washing steps, making it especially suitable for high throughput screening applications.

COMPONENTS:

Catalog #	Component	Amount	Storage	
80301	UBE1 (E1)	25 µg	-80°C	<i>Avoid freeze/ thaw cycles!</i>
80313	UBCH5c (E2)	5 µg	-80°C	
	Biotin-Ubiquitin	400 µl	-80°C	
	ATP (10 mM)	400 µl	-80°C	
	CBL assay buffer	2 x 10 ml	-80°C	
30017	Anti-His Tb-labeled donor	10 µl	-20°C	
	Dye-labeled acceptor	10 µl	-20°C	
	Methyl-ubiquitin (1 mM in 50 mM Tris, pH 7.4)	20 µl	-80 °C	
	White Corning microtiter plate	1	Room temp.	

MATERIALS OR INSTRUMENTS REQUIRED BUT NOT SUPPLIED:

Fluorescent microplate reader capable of measuring Time Resolved Fluorescence Resonance Energy Transfer (TR-FRET)
Adjustable micropipettor and sterile tips

APPLICATIONS: Great for screening small molecular inhibitors for drug discovery and HTS applications.

STABILITY: At least 6 months from date of receipt when stored as directed.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

REFERENCES: 1. Saville, M.K., *et al.*, *J. Biol. Chem.* 2004; **279(40)**: 42169-81.
2. Kim, J.-H., *et al.*, *BMB Rep.* 2015; **48(1)**: 25-29.

ASSAY PROTOCOL:

All samples and controls should be tested in triplicates.

- 1) Thaw **UBE1**, **UBCH5c**, **Biotin-Ubiquitin**, **CBL assay buffer**, and **ATP** on ice. Aliquot each protein, **CBL assay buffer**, and **ATP** into single-use aliquots and store at -80°C immediately. *Note: **UBE1**, **UBCH5c**, **Biotin-Ub**, **CBL assay buffer**, and **ATP** are sensitive to freeze/thaw cycles. Avoid multiple freeze-thaw cycles.*
- 2) Carefully calculate the amount of proteins needed. Prepare appropriate amounts of diluted proteins; dilute only the amount required for the assay. Do not store diluted proteins

Dilute the **UBE1** in **CBL assay buffer** at 40 ng/μl
Dilute the **UBCH5c** in **CBL assay buffer** at 8 ng/μl
Keep the diluted reagents on ice until use.

- 3) Prepare the master mixture using diluted reagents: N wells × (1 μl **Biotin-Ub** + 1.5 μl diluted **UBE1** + 1.5 μl diluted **UBCH5c**). Add 4 μl of master mixture to each well designated for the "Substrate Control," "Positive Control," "Test Inhibitor." For the wells labeled "Blank," add 1 μl **Biotin-Ub** + 1.5 μl diluted **UBE1** + 4.5 μl **CBL assay buffer**.
- 4) Dilute your test inhibitor (100x in DMSO) 1:20 in distilled water. Add 2 μl of diluted test inhibitor solution to each well designated "Test Inhibitor." For the "Positive Control," "Substrate Control," and "Blank," add 2 μl of the same solution without the test inhibitor (inhibitor buffer), in this instance, 5% DMSO in water. Final DMSO concentration in the assay should be ≤1%.

Note: To make an IC50 using the methyl-ubiquitin, prepare serial dilutions in the same inhibitor buffer as your test inhibitor (typically 5% DMSO in water)

- 5) Add 4 μl of assay buffer to the well designated "Substrate Control." Add 3 μl of assay buffer to the wells designated "Positive Control," and "Test Inhibitor."
- 6) Initiate the reaction by adding 1 μl of **ATP (10 mM)** to the wells labeled "Positive Control," "Test Inhibitor," and "Blank." Incubate the reaction at 30°C for four hours. Cover the plate with a plate sealer if necessary.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
 San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

	Blank	Substrate Control	Positive Control	Test Inhibitor
Biotin-Ub	1 μ l	1 μ l	1 μ l	1 μ l
UBE1 (40 ng/ μ l)	1.5 μ l	1.5 μ l	1.5 μ l	1.5 μ l
UBCH5c (8 ng/ μ l)	-	1.5 μ l	1.5 μ l	1.5 μ l
Test Inhibitor/Activator	-	-	-	2 μ l
5% DMSO in water (Inhibitor buffer)	2 μ l	2 μ l	2 μ l	-
CBL assay buffer	4.5 μ l	4 μ l	3 μ l	3 μ l
ATP (10 mM)	1 μ l	-	1 μ l	1 μ l
Total	10 μl	10 μl	10 μl	10 μl

- 7) Dilute **Tb-labeled donor** (1:400) and **Dye-labeled acceptor** (1:400) in one step using CBL Assay Buffer. Prepare only the amount required for the assay. Add 10 μ l of diluted donor/acceptor mixture into each well. Incubate at room temperature for one hour.
- 8) Read the fluorescent intensity in a microtiter-plate reader capable of measuring TR-FRET. Blank value is subtracted from all other values.

Instrument Settings

Reading Mode	Time Resolved
Excitation Wavelength	340 \pm 20 nm
Emission Wavelength	620 \pm 10 nm
Lag Time	60 μ s
Integration Time	500 μ s
Excitation Wavelength	340 \pm 20 nm
Emission Wavelength	665 \pm 10 nm
Lag Time	60 μ s
Integration Time	500 μ s

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com



6042 Cornerstone Court W, Ste B
San Diego, CA 92121
Tel: 1.858.202.1401
Fax: 1.858.481.8694
Email: info@bpsbioscience.com

CALCULATING RESULTS:

Two sequential measurements should be conducted. Tb-donor emission should be measured at 620 nm followed by dye-acceptor emission at 665 nm. Data analysis is performed using the TR-FRET ratio (665 nm emission/620 nm emission).

When percentage activity is calculated, the FRET value from the negative control (Blank or Substrate Control) can be set as zero percent activity and the FRET value from the positive control can be set as one hundred percent activity.

$$\% \text{ Activity} = \frac{\text{FRET}_s - \text{FRET}_{\text{neg}}}{\text{FRET}_p - \text{FRET}_{\text{neg}}} \times 100\%$$

Where FRET_s = Sample FRET, FRET_{neg} = negative control FRET, and FRET_p = Positive control FRET.

Example of Assay Results:




Figure 1: Titration of UBCH5c activity using the *UBCH5c TR-FRET Assay Kit*, BPS Bioscience #79901. *Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.com.*

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com




Figure 2: Inhibition of UCBH5c Assay FRET signal by Methylated Ubiquitin, measured using the *UBCH5c TR-FRET Assay Kit*, BPS Bioscience #79901. *Data shown is lot-specific. For lot-specific information, please contact BPS Bioscience, Inc. at info@bpsbioscience.com.*

RELATED PRODUCTS

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
UBE1 (UBA1), FLAG-tag	#80301	100 µg
Ubiquitin, Biotin Labeled	#11236	50 µg
UbcH5a, His-Tag (Human)	#80315	100 µg
UbcH5b, His-Tag (Human)	#80314	100 µg
UbcH5c, His-Tag (Human)	#80313	100 µg
UBA6 (UBE1L2), FLAG-tag	#80303	100 µg
CBL-B, GST-Tag (Human)	#80415	100 µg
CBL-B, His-Avi-Tag	#80414	100 µg
CBL-B, Biotin-labeled (Human)	#80412	50 µg
CBL-B (Y363F), Biotin-labeled (Human)	#80413	50 µg

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.202.1401** Fax **1.858.481.8694**

Or you can Email us at: info@bpsbioscience.com

Please visit our website at: www.bpsbioscience.com