



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Product Information

## Annexin V Conjugates

Cat. No.	Product Name	Unit Size
29045-100uL	Annexin V, R-PE Conjugate	100 uL (20 tests)
29045-500uL		500 uL (100 tests)
29057-100uL	Annexin V, APC Conjugate	100 uL (20 tests)
29057-500uL		500 uL (100 tests)

### Storage and Handling

Store at 4°C, protected from light. Storage at -20°C is not recommended for R-PE/APC conjugates. Product is stable for at least 6 months from date of receipt when stored as recommended.

### Product Description

Fluorescent conjugates of Annexin V can be used to label apoptotic cells. The human anticoagulant Annexin V is a 35-36 kDa, Ca<sup>2+</sup>-dependent phospholipid-binding protein with high affinity for phosphatidylserine (PS). In normal viable cells, PS is located on the inner leaflet of the cytoplasmic membrane. However, in apoptotic cells, PS is translocated from the inner to the outer leaflet of the plasma membrane, where it is available for binding to fluorescently labeled Annexin V.

The fluorescent proteins R-phycoerythrin (R-PE) and allophycocyanin (APC) belong to the family of phycobiliproteins, which are derived from cyanobacteria and eukaryotic algae. R-PE is a 240 kDa protein with high extinction coefficient (~1,960,000 M<sup>-1</sup>cm<sup>-1</sup>) and high fluorescence quantum yield (~0.82). APC is a 105 kDa protein with high extinction coefficient (240,000 M<sup>-1</sup>cm<sup>-1</sup>) and high quantum yield (0.68) in the far-red region. The APC in Annexin V-APC conjugates are chemically cross-linked. Annexin V R-PE and APC conjugates are designed for detection by flow cytometry.

Biotium offers a variety of Annexin V conjugates including those labeled with our outstanding series of CF® dyes for fluorescence microscopy or flow cytometry. Please visit [www.biotium.com](http://www.biotium.com) for details.

### Staining Protocols

We recommend using our Annexin V Binding Buffer (Catalog no. 99902) with Annexin V conjugates. Staining can be performed in other calcium-containing buffers such as HEPES-buffered saline containing 2.5 mM CaCl<sub>2</sub>. The optimal staining concentration for each conjugate should be determined empirically. The protocols provided below are intended to serve as general guidelines.

Note: Annexin V cannot be used to stain fixed cells or tissues. After staining with Annexin V and washing, cells can be fixed with 2% formaldehyde. Annexin V staining is calcium dependent, therefore 2.5 mM CaCl<sub>2</sub> should be included in all buffers used for washing and fixation. Annexin V binds to a phospholipid in the plasma membrane, therefore staining is not compatible with alcohol-based fixation or detergent permeabilization.

### Suspension cells for flow cytometry

1. Induce apoptosis in cells by desired method. Include a control sample of untreated cells.
2. Dilute 5X Annexin V Binding Buffer (catalog number 99902) 1:5 in distilled water to obtain 1X Binding Buffer. HEPES-buffered saline containing 2.5 mM CaCl<sub>2</sub> can be used in place of 1X Binding Buffer.
3. Wash cells with PBS once and resuspend cells at 2-3x10<sup>6</sup> cells/mL in 1X Binding Buffer.
4. Aliquot 100 uL cells per tube.
5. Add 5 uL Annexin V R-PE or APC conjugate to each tube and mix well.
6. Incubate at room temperature for 15 minutes, protected from light.
7. Add 400 uL 1X Binding Buffer to each tube. Analyze the cells by flow cytometry within 1 hour of staining.

### Adherent cells for flow cytometry

1. Induce apoptosis in cells by desired method. Include a control sample of untreated cells.
2. Dilute 5X Binding Buffer 1:5 in distilled water to obtain 1X Binding Buffer. HEPES-buffered saline containing 2.5 mM CaCl<sub>2</sub> can be used in place of 1X Binding Buffer.
3. Wash cells with PBS twice and detach cells from cell culture plate or well by trypsin or cell dissociating buffer.
4. Pellet cells and discard supernatant. Resuspend cells at 2-3x10<sup>6</sup> cells/mL in 1X Binding Buffer.
5. Aliquot 100 uL cells per tube.
6. Add 5 uL Annexin V R-PE or APC conjugate to each tube and mix well.
7. Incubate at room temperature for 15 minutes, protected from light.
8. Add 400 uL 1X Binding Buffer to each tube. Analyze the cells by flow cytometry within 1 hour of staining.

### Related Products

Cat. No.	Product Name
99902	5X Annexin V Binding Buffer
10405	NucView™405 Caspase-3 Substrate, 1 mM in DMSO
10402	NucView™488 Caspase-3 Substrate, 1 mM in DMSO
10406	NucView™530 Caspase-3 Substrate, 1 mM in DMSO
30067	Dual Apoptosis Assay Kit with NucView™488 Caspase-3 Substrate & CF®™594 Annexin V
30062	NucView™488 and MitoView™ 633 Apoptosis Kit
30072	NucView™488 and RedDot™2 Apoptosis and Necrosis Kit
30060	CF® 488A Annexin V and 7-AAD Apoptosis Kit
30061	CF® 488A Annexin V and PI Apoptosis Kit
32010	Live-or-Dye NuCF® ix™ Red Staining Kit
32002-32009	Live-or-Dye™ Fixable Viability Staining Kits

A full selection of CF® dye labeled products including secondary antibodies, streptavidin and anti-biotin antibodies, antibody labeling kits, and other bioconjugates such as phalloidins, lectins, and a-bungarotoxin are also available. Please visit the Biotium website at [www.biotium.com](http://www.biotium.com) for details.

Biotium products are for research use only, and are not intended for food, drug, or household use.