



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Product Information

## Apoptotic, Necrotic and Healthy Cells Quantitation Kit Plus

**Catalog Number:** 30066

**Unit Size:** 50 assays

### Kit Contents

Component	Cat. No.	Size
CF@488A Annexin V in TE/0.1% BSA/0.1% NaN <sub>3</sub>	99967	250 uL
Ethidium Homodimer III (EthD-III) in PBS	99968	250 uL
Hoechst 33342 in dH <sub>2</sub> O	30018C	250 uL
5X Annexin V Binding Buffer	99902	15 mL

### Storage and Handling

Store the kit at 4°C. Do not freeze. Protect CF@488A Annexin V, Ethidium Homodimer III and Hoechst 33342 from light. When stored as directed, the kit is stable for at least 6 months from the date it is received.

### Spectral Properties

Hoechst 33342 Ex/Em: 352/462 nm (with DNA)

CF@488A Ex/Em: 490/515 nm

Ethidium Homodimer III (EthD-III) Ex/Em: 532/625 nm\* (with DNA)

\*Ethidium Homodimer III also has a strong UV absorbance peak at 279 nm

### Product Description

The Apoptosis, Necrosis and Healthy Cell Quantitation Kit Plus provides a convenient assay for detecting apoptotic cells (green) and necrotic cells (red) within the same cell population by flow cytometry or fluorescence microscopy. The nuclear stain Hoechst is included for staining the entire cell population (live and dead) blue.

Apoptosis and necrosis are two major processes by which cells die. Apoptosis is an active, genetically regulated disassembly of the cell from within. During apoptosis, phosphatidylserine (PS) is translocated from the inner to the outer leaflet of the plasma membrane for phagocytic cell recognition. The human anticoagulant Annexin V is a 35 kDa, calcium-dependent phospholipid-binding protein with high affinity for PS. Annexin V labeled with CF@488A stains apoptotic cells green by binding to PS exposed on the cell surface. CF@488A is spectrally similar to fluorescein (FITC), with much brighter and more photostable fluorescence.

Necrosis normally results from a severe cellular insult. Both internal organelle and plasma membrane integrity are lost, resulting in the spilling of cell contents into the surrounding environment. Ethidium Homodimer III (EthD-III) is a highly positively charged nucleic acid probe, which is impermeant to live cells and early apoptotic cells, but stains necrotic cells and late apoptotic cells with red fluorescence.

EthD-III is a superior alternative to other red dead cell stains such as propidium iodide (PI) or Ethidium Homodimer I (EthD-I) due to its significantly higher affinity for DNA and higher fluorescence quantum yield, resulting in a brighter signal.

Membrane-permeant blue fluorescent nuclear dye Hoechst 33342 is included in the kit for staining the entire cell population.

### Assay Protocols

**Note:** This assay must be used on unfixed cells. Both Annexin V and EthD-III rely upon the presence of intact membranes in healthy cells to accurately distinguish healthy cells from apoptotic or necrotic cells. The dyes cannot be used for live/dead discrimination in fixed cells or tissues. The dyes cannot withstand fixation after staining.

**Note:** We recommend including single stain controls (i.e., staining cells with only CF488A Annexin V, EthD-III, or Hoechst 33342).

### Suspension cell staining

1. Prepare 1X Binding Buffer by diluting 5X Annexin V Binding Buffer with dH<sub>2</sub>O (mix 1 volume of 5X buffer with 4 volumes of dH<sub>2</sub>O).
2. Wash cells with PBS and resuspend cells at 2-3x10<sup>6</sup> cells/mL in 1X Binding Buffer.
3. Pipet 100 uL cell suspension into a microcentrifuge tube.
4. Add 5 uL of CF@488A Annexin V, 5 uL of EthD-III, and 5 uL Hoechst 33342 to each tube.  
**Note:** Reducing the concentration of EthD-III may result in better signal-to-noise for some cell types.
5. Incubate at room temperature for 15 minutes in the dark.
6. For flow cytometry analysis, add 400 uL 1X Binding Buffer to each tube and measure fluorescence in Pacific Blue®, FITC, and PE channels within 1 hour of staining.
7. For fluorescence microscopy analysis, wash cells with 1X Binding Buffer, resuspend cells in 1X Binding Buffer, and observe fluorescence using DAPI, FITC, and Texas Red® or Cy@3 filter sets.

### Adherent cell staining for fluorescence microscopy

1. Prepare 1X Binding Buffer by diluting 5X Annexin V Binding Buffer with dH<sub>2</sub>O (mix 1 volume of 5X buffer with 4 volumes of dH<sub>2</sub>O).
2. Wash cells twice with PBS.
3. Prepare staining solution by adding 5 uL of CF@488 Annexin V, 5 uL EthD-III, and 5 uL Hoechst 33342 per 100 uL 1X Binding Buffer. Prepare enough staining solution to cover cells.  
**Note:** Reducing the concentration of EthD-III may result in better signal-to-noise for some cell types.
4. Incubate samples 15 minutes at room temperature, protected from light.
5. Wash cells with 1X Binding Buffer 1-2 times.
6. Cover cells with 1X Binding Buffer and observe fluorescence using DAPI, FITC, and Texas Red® or Cy@3 filter sets.

### Staining adherent cells for flow cytometry

1. Detach cells from cell culture plate or well using trypsin or other cell dissociation method. Pellet cells and discard supernatant.
2. Follow protocol for suspension cell staining.

## Related Products

Cat. No.	Product
99902	5X Annexin V Binding Buffer
29003... 29085	Annexin V CF® Dye Conjugates
29006... 29082	Annexin V Near IR CF® Dye Conjugates
29003R... 29085R	Annexin V CF® Dye Conjugates, Azide-Free, Lyophilized
30072	NucView® 488 and RedDot™2 Apoptosis and Necrosis Kit
10403	NucView® 488 Caspase-3 Substrate, 1 mM in PBS
10407	NucView® 405 Caspase-3 Substrate, 1 mM in PBS
10408	NucView® 530 Caspase-3 Substrate, 1 mM in PBS
30062	NucView® 488 and MitoView™ 633 Apoptosis Kit
40081-40082	NucSpot® Live Cell Nuclear Stains
40083... 41040	NucSpot® Nuclear Stains
10059	MTT
10060	XTT
30063... 30074	TUNEL Assay Apoptosis Detection Kit
70055	MitoView™ 633
32002... 32018	Live-or-Dye™ Fixable Viability Staining Kits
32010	Live-or-Dye™ NucFix Red Staining Kit
30020	ATP-Glo™ Bioluminometric Cell Viability Assay
30068	ViaFluor® SE 405 Cell Proliferation Kit
30086	ViaFluor® SE 488 Cell Proliferation Kit
30050	ViaFluor® SE CFSE Cell Proliferation Kit

Please visit our website at [www.biotium.com](http://www.biotium.com) for information on our life science research products, including environmentally friendly EvaGreen® qPCR master mixes, fluorescent CF® Dye antibody conjugates and reactive dyes, apoptosis reagents, fluorescent probes, and kits for cell biology research.

Texas Red and Pacific Blue are registered trademarks of Thermo Fisher Scientific; Cy is a registered trademark of Cytiva.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.