



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# PRODUCT INFORMATION



## Cysteinyl Leukotriene Affinity Sorbent

Item No. 400396

### Overview and Properties

<b>Contents:</b>	This vial contains 5 ml of cysteinyl leukotriene affinity sorbent.
<b>Storage:</b>	4°C (as supplied)
<b>Stability:</b>	As supplied, 1 year from the QC date provided on the Certificate of Analysis, when stored properly

### Procedures

This vial contains 5 ml of cysteinyl leukotriene (cysLT) affinity sorbent (mouse monoclonal anti-LTC<sub>4</sub>/LTD<sub>4</sub>/LTE<sub>4</sub> covalently bound to Sepharose 4B) supplied as a 50:50 mixture in in Eicosanoid Affinity Column Buffer. This sorbent is stable for at least two years if stored at 4°C. Do NOT freeze. The sorbent has a binding capacity of 10 ng CysLT per ml of sorbent determined using LTE<sub>4</sub>.

#### Sorbent Preparation

Wash the sorbent with 2 to 4 volumes of Eicosanoid Affinity Column Buffer (Item No. 400220). Briefly centrifuge at 500 x g to sediment the sorbent. Repeat wash. Gently mix the sorbent prior to use by vortexing or inversion to completely resuspend the sorbent in the buffer. Best results are obtained when a moderate excess of the sorbent is used to bind the CysLTs from biological samples. We recommend using 50 µl-200 µl of the resuspended sorbent per 1 ml of sample.

#### Eicosanoid Affinity Column Buffer Preparation

Prepare a 0.1 M phosphate buffer solution by combining 13.3 g K<sub>2</sub>HPO<sub>4</sub>, 3.22 g KH<sub>2</sub>PO<sub>4</sub>, 0.5 g sodium azide, and 29.2 g sodium chloride. Dilute to a total volume of 1.0 liter with UltraPure water. The pH of this buffer will be 7.0. This buffer may also be purchased as a 5X concentrated buffer (Item No. 400220).

#### Sample Preparation

This sorbent has been validated with the LTE<sub>4</sub> and CysLT Express EIA Kits (Item No. 520411 and 10009291, respectively) for human urinary cysLT extraction. Purification of other sample matrices have not been fully validated.

#### Urine

When known amounts of LTE<sub>4</sub> are spiked into human urine samples it was found that 200 µl of sorbent is sufficient to extract between 20 pg and 2 ng LTE<sub>4</sub> from 1 ml of normal human urine.

#### Purification Protocol

1. Aliquot 1 ml of sample to a 1.5 ml microfuge or similar tube. Larger samples can be used. It is important to also include a blank extract (sorbent incubated in buffer alone). We recommend using plastic conical centrifuge tubes to facilitate collection of the sorbent by centrifugation. When using larger volumes of sample, the amount of sorbent and incubation time may need to be increased.
2. Gently mix the sorbent prior to use to resuspend. Add 50-200 µl of resuspended sorbent mixture to each 1 ml sample (or blank) and mix gently for 30-60 minutes. The amount of sorbent needed will depend on the concentration of CysLT in the sample.
3. Briefly centrifuge the samples at 500 x g to sediment the sorbent.
4. Carefully remove the supernatant by aspirating with a pipette. Care must be taken to retain all the sorbent as this contains the bound CysLT.

**WARNING**  
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

**SAFETY DATA**  
This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

**WARRANTY AND LIMITATION OF REMEDY**  
Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 07/06/2017

**CAYMAN CHEMICAL**  
1180 EAST ELLSWORTH RD  
ANN ARBOR, MI 48108 · USA  
PHONE: [800] 364-9897  
[734] 971-3335  
FAX: [734] 971-3640  
CUSTSERV@CAYMANCHEM.COM  
WWW.CAYMANCHEM.COM

# PRODUCT INFORMATION



5. Wash the sorbent once with 1 ml Column Buffer, centrifuge, and carefully remove and discard the supernatant.
6. Resuspend the sorbent in 0.5 ml ice-cold methanol. Vortex briefly.
7. Centrifuge to sediment the sorbent and carefully remove the methanol and transfer to a clean tube. **DO NOT DISCARD** the methanol as this contains the eluted CysLT. Repeat Steps 6 and 7.
8. Combine the methanol eluates and evaporate to dryness under nitrogen or vacuum. If the analysis cannot be performed at once, store the methanol eluates at  $-80^{\circ}\text{C}$ ; they will be stable for at least six months.
9. Immediately dissolve the CysLT in the buffer or solvent appropriate for your application. If you are assaying the sample with one of our EIA kits (Cysteinyl Leukotriene EIA Kit - Item No. 500390; Cysteinyl Leukotriene Express EIA Kit - Item No. 10009291; Leukotriene E4 EIA Kit - Item No. 520411), dissolve the sample in the buffer or solvent recommended in the instructions provided with the kit. The amount of buffer depends on the original volume and the expected concentration of CysLT in the sample. A dilute sample may be concentrated by dissolving the residue in a smaller volume of buffer than the original sample volume. *NOTE: Remember to also subtract the value of the blank (i.e., buffer extracted sample) from all sample values in the final analysis.*

CAYMAN CHEMICAL  
1180 EAST ELLSWORTH RD  
ANN ARBOR, MI 48108 · USA  
PHONE: [800] 364-9897  
[734] 971-3335  
FAX: [734] 971-3640  
CUSTSERV@CAYMANCHEM.COM  
WWW.CAYMANCHEM.COM