



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ASSAY NAME: SARS CoV2 N1P (Omicron)

Quantity: 100 x 20µL PCR reactions

2-plex assay: SARS CoV2 N1P and human RPP30 mRNA

Detects: Omicron variants: BA.1, BA.2, BA.4 and BA.5 through JN.1 Targets: N gene: Δ31-33

SKU#'s:

PNP-N1P-D-BR (Bio-Rad with control assay)

PNP-N1P-N-BR (Bio-Rad without control assay)

PNP-N1P-D-QS (QuantStudio with control assay)

PNP-N1P-N-QS (QuantStudio without control assay)

PNP-N1P-D-MIC (BMS MIC with control assay)

PNP-N1P-N-MIC (BMS MIC without control assay)

(RUO). Research Use Only. Not for use in Diagnostic Procedures.

SCOPE OF THIS DOCUMENT

The oligonucleotide recipes are optimized for each instrument (BioRad, QuantStudio, MIC). The pre-validation data presented in this document were performed using PNP-N1P-D-BR-100 on a BioRad CFX96. The performance of the other SKUs on their corresponding instrument should be similar. Contact PCRassays.com if need to use a different qPCR instrument.

CONTENTS

The primers and probes in the SARS CoV2 N1P kit are provided in Tube 1 as a 20X concentrated working solution. The same mix also contains primers/probe targeting human RPP30RNA splice junction as an RT-qPCR positive control assay for human samples. The probes are designed as TaqMan⁵ cleavage mechanism and thus the reaction requires a DNA polymerase with 5'-exonuclease activity (we recommend InhibiTaq Standard RT-qPCR Master Mix).

Table of Dyes used in this kit:

Pathogen/Target	Dyes	Quencher	Refs.
RPP30-RNA control	HEX	BHQ-1	3, 4
SARS CoV2 N1P	FAM	BHQ-2	1, 2

Tube 2 “double positive control” contains a mixture of synthetic 500 bp DNA constructs containing the amplicon regions of all targets and human RPP30 is provided as a positive extraction control. The concentration of each DNA construct is approximately 5,000 copies/µL. The DNA constructs are for validation purposes only and **Tube 2 should NOT be added to wells for specimen unknowns.**

Note: molecular biology grade water should be used to prepare the PCR reactions, which is NOT included in this kit.

Assay contents:

Tube 1: 20X Primer/Probe mix for SARS CoV2 N1P. If you order SKU#: PNP-N1P-D, then primers and probes for hRPP30 intron are also included.

Tube 2: (optional if ordered) 5000 copies/µl Positive controls of synthetic 500 bp DNA fragments of SARS CoV2 N1P and hRPP30.

Tube 3: (optional if ordered) Spike-in control. 1.0E6 copies/µL of synthetic 500 BP regions of SARS CoV2 N1P and human RPP30 intron.

Tube 4: (optional if ordered) InhibiTaq Standard RT-qPCR enzyme Mastermix (enough for 100 rxns. with 20 µL total volume). This is a custom formulation from Fortis Life Sciences to the specifications of PCRassays.com.



KIT HANDLING AND CONTAMINATION

The SARS CoV2 N1P assay is shipped at ambient temperature, and should be stored at -20 °C. The kit should be kept on ice once thawed. Do not subject the enzyme to multiple freeze-thaw cycles.

Any contamination should be avoided by using appropriate personal protective equipment (PPE), powder free gloves, aerosol barrier pipette tips, and a clean hood.

EXPERIMENTAL

Set up your reaction (20 µL) as follows on ice:

Component	Volume (µL)
InhibiTaq Standard RT-qPCR enzyme mastermix (2X)	10
Primer/Probe mix (20X)	1
Sample	2
Water	7

Notes: To improve assay sensitivity, up to 9 µL of sample can be added (water volume adjusted accordingly) for a total reaction volume of 20 µL. For positive control rxns., add 2 µL of the solution from Tube 2 (i.e., the “sample”).

An RT-PCR protocol was used at PCRassays.com for verification on a Bio-Rad CFX96™ Real-Time System, with the following program:

Step	Thermocycling Protocol:
1	Incubate @ 50 °C for 5 minutes
2	Incubate @ 94 °C for 3 minutes
3	Incubate @ 94 °C for 5 seconds
4	Incubate @ 60 °C for 15 seconds
5	Plate Read
6	Go to Step 3, repeat 44x more
7	(optional) Incubate @60 °C for 3 minutes

For QuantStudio instruments, we recommend a Step 3 cycle time of 22 seconds at 60 °C.

RESULT INTERPRETATION

After running the qPCR reaction, perform a regression analysis on the data to determine the quantification cycle, Cq. (Cq is preferred over Ct). Each fluorescence channel with a Cq < 37 cycles and final RFU > “threshold” is considered “positive” or “+” in the Table below. The “threshold” is 2.0 on the BMS MIC, 200 on BioRad instruments and 200,000 on QuantStudio 5, 6, 7, 12K instruments.

SARS CoV2 N1P FAM™	hRPP30 HEX™	Recommended Interpretation
–	–	The PCR reaction failed. Please repeat the experiment.
–	+	The sample does not contain bacterial DNA of interest. The sample contains human RPP30 DNA.
+	–	The sample contains SARS CoV2 N1P DNA. The sample may not contain human RPP30 DNA.
+	+	The sample contains SARS CoV2 N1P DNA and human RPP30 DNA.

VERIFICATION EXPERIMENTS

Experiments were performed in triplicate using the experimental procedure given above, but with different samples added to each reaction. The samples used for the validation experiments contained 1×10⁵ copies/reaction synthetic viral RNA obtained from Twist Biosciences as follows:

SARS-CoV-2 Delta variant (Twist® Standard RNA #18)
SARS-CoV-2 Omicron BA.1 variant (Twist® Standard RNA #48)
SARS-CoV-2 Omicron BA.2 variant (Twist® Standard RNA #50)

The samples also contained human brain RNA (1500 copies) from Roche and human genomic DNA (3100 copies) from Clontech. The RPP30 control primers and probe specifically reverse transcribe and amplify the human RPP30 mRNA and not genomic DNA (See DNAS Product insert about RPP30 RNA control for more information). The presence of the human genomic DNA in the reaction appears to have no effect on the amplification of SARS-CoV-2 RNA with this kit (data not shown). The results of these experiments are shown in **Figure 1**.




Figure 1: Verification experiments with single or double target(s) (given in text boxes for each panel). All sets of probes and primers are present in every reaction, but positive signal is only observed when the target(s) is present, indicating that the amplification is specific. The **FAM** probe detects SARS CoV2 for both Omicron variants BA.1 and BA.2. The **HEX** probe detects human RPP30 mRNA.




Figure 2: Serial dilution experiments show LOD <20 molecules for each target. For the bottom panel (BA.2), 1 out of the 3 NTC reactions showed an amplification with Cq=37.4 (the other 2 NTC reactions were flat). We think there may have been a single molecule contamination in that reaction since we have not observed that in any of the other many NTC reactions run for this assay.

Conclusion: The data in **Figure 1** indicates that the SARS CoV2 N1P primers and probe are compatible with the human RPP30 mRNA positive control primers and probe in the human genomic DNA and RNA.

Limit of detection (LOD) was estimated by performing serial dilution experiments in triplicate (**Figure 2**). For dilution series only one SARS CoV-2 template RNA was added (*i.e.*, Omicron BA.1, or Omicron BA.2). The results show a limit of detection (LOD) <20 copies/reaction for both templates.

CONTACT US

For assistance, please contact DNA Software using the link:
<https://dnasoft.jira.com/servicedesk/customer/portals>

Address: Michigan Life Science and Innovation Center,
46701 Commerce Center Dr, Plymouth, MI 48170

Phone: (734) 222-9080

NOTES

- ¹ FAM™ (Carboxyfluorescein), a trademark of Life Technologies Corporation.
- ² BHQ-2™ (Black Hole Quencher) is a trademark of Biosearch Technologies, Inc.
- ³ HEX™ (Hexachloro-fluorescein) is a trademark of Applera Corp.
- ⁴ BHQ-1™ (Black Hole Quencher) is a trademark of Biosearch Technologies, Inc.
- ⁵ “TaqMan” is a trademark of Roche Molecular Systems, Inc.

