



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®]
www.cedarlanelabs.com



Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Elite PAP Pen

CLSG80125
CLSG90451

Intended Use: To provide a hydrophobic barrier when a circle is drawn around a specimen on the glass slide.

Introduction: Elite PAP pen provides a thin, film like water repellent barrier when a circle is drawn around the specimen on the glass slide. This barrier creates the proper surface tension to hold antibody solution within the target area on the slide. This avoids the repeated wipeout around the specimen during the staining procedure. Since the reagent is guarded by the hydrophobic circle around the tissue, it minimizes the undesirable dilution of the reagents with the wash buffer. This pen can be used for PAP, ABC, immunofluorescence method, frozen sections and in situ hybridization procedures. Two specimens separated by two circles can be applied to the same slide. The Mini PAP pen (**CLSG90451**) is especially useful when working with multiple sections on a single slide because its tip is finer than the PAP pen (**CLSG80125**).

Properties:

1. Heat resistance up to 120°C.
2. Ready to be used directly on the wet slides.
3. Applicable to microwave staining and in situ hybridization.
4. Circle made by Elite PAP pen is insoluble in alcohol, acetone.
5. Withstand rehydration and dehydration steps in alcohol.
6. If desired, circle can be removed by xylene after the staining procedure is completed.

Preparation before use: Before use, hold the pen point up and press it with fingertip to take blocked air out from the pen barrel to prevent overflow.

Procedure: To start, press the tip of the pen straight down until fluid starts soaking the tip and then draw the circle around the specimen.

Storage: Always keep the pen tightly capped and store horizontally. Press down on the tip when pen starts to dry out. Store pen at room temperature.

Laboratory Reagent For Research Use Only

JK 11/20/15

Visit our website for your local distributor.

CEDARLANE[®]



www.cedarlanelabs.com

An ISO 9001:2000 and ISO 13485:2003
registered company.

In CANADA: **Toll Free: 1-800-268-5058**

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: **Toll Free: 1-800-721-1644**

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com