



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Technically
Speaking

CEDARLANE[®]
www.cedarlanelabs.com



Conveniently Delivering You Today's Innovations
for the Science of Tomorrow™

Streptavidin Coated 96-well Microtiter Plates

CAT#: CLSTAN-5P

INTRODUCTION:

Streptavidin is a 52.8 KDa protein originally isolated from the bacterium *Streptomyces avidinii* that naturally forms a homo-tetramer. Each of the four subunits of streptavidin contain a single binding site for biotin. Streptavidin has an extremely high affinity for biotin, having a dissociation constant (Kd) on the order of ~10–14 mol/L; this makes it one of the strongest non-covalent interactions known in nature. Additionally, Streptavidin has no carbohydrate groups and a near neutral isoelectric point, resulting in low nonspecific interactions. Taken together, the high affinity and low background interactions of streptavidin and biotin make it an ideal system for use in biological assays.

Cedarlane's Streptavidin Coated Microtiter Plates are coated with highly purified streptavidin. After coating, the plates are blocked to ensure low non-specific binding and stability as well as allowing for direct use. The plates are useful for capturing biotinylated proteins, peptides and DNA and can be used in various applications including ELISA and measure of protein-protein interactions.

STORAGE/STABILITY:

The microplates, if unopened, are stable refrigerated until the expiration date printed on the label. If opened, store in closed pouch with desiccant and use within 2- 4 weeks.

SPECIFICATIONS:

Components

Five (5) individually pouched 96-well plates, configured in 12 removable 8-well strips.

Coating

Purified Streptavidin is coated using 100 µl/well. Nominal binding capacity is ~10 pmol biotin/well. The strips are blocked with 200 µl/well to ensure low non-specific binding and long-term stability.

Sensitivity

Biotinylated IgG was detected at concentrations significantly above background in an ELISA format using anti-rabbit or mouse IgG-HRP (at 0.1 µg/ml) as detector and TMB as substrate (see Figure 1 below), as follows:

Rabbit IgG – 0.156 ng/well

Mouse IgG – 0.234 ng/well

Continued Overleaf....

Visit our website for your local distributor.

CEDARLANE[®]
www.cedarlanelabs.com



An ISO 9001:2015 and ISO 13485:2016
registered company.

In CANADA: Toll Free: 1-800-268-5058

4410 Paletta Court, Burlington, ON L7L 5R2 ph: (289) 288-0001, fax: (289) 288-0020
e-mail: general@cedarlanelabs.com

In the USA: Toll Free: 1-800-721-1644

1210 Turrentine Street, Burlington, NC 27215 ph: (336) 513-5135, fax: (336) 513-5138
e-mail: service@cedarlanelabs.com




Figure 1. Capture of biotinylated rabbit and mouse IgG.

APPLICATION:

Example ELISA Procedure:

A. Materials Required

- Wash Buffer / Assay buffer: 1x phosphate-buffered saline, 0.05% Tween[®]-20, pH 7.4 or alternative buffer (assay dependent).
- Biotinylated capture antibody to antigen of interest diluted in wash buffer.
- Antigen adjusted to appropriate concentration with wash buffer.
- Enzyme-labeled (HRP or ALP) primary antibody, or primary antibody + enzyme-labeled secondary antibody diluted in wash buffer.
- Appropriate enzyme substrate (TMB or Phosphatase); stop solution (1M H₂SO₄).

B. Method

1. Add 100μL of the biotinylated capture antibody to each well and incubate for 1 hour at 37°C.
2. Wash each well five times with 300μL of Wash Buffer. Make a serial dilution of the antigen and add 100μL to each well. Incubate plate for 1 hour at 37°C.
3. Wash each well five times with 300μL of Wash Buffer. Add 100μL of the primary antibody (unconjugated) or primary antibody-HRP conjugate to each well and incubate plate for 1 hour at 37°C.
4. Wash each well five times with 300μL of Wash Buffer.
5. Proceed to TMB step in the case of antibody-HRP conjugate or add 100μL of the HRP-labeled secondary antibody to each well. Incubate plate for 1 hour at 37°C. Wash each well five times with 300μL of Wash Buffer.
6. Add 100μL of TMB and incubate at RT x 15-30 min.
7. Add 100μL of Stop Solution, read absorbance at 450nm and process data.

For Research Use Only

® is a registered trademark of Cedarlane Corporation