Duplica^{RealTime} *Neisseria gonorrhoeae* 2nd Generation Detection Kit

REV. EBR018032_IFU_REV.02D_ENITA

REF EBR018032 -32 tests

Instructions For Use

INTENDED USE

Duplica^{RealTime} *Neisseria gonorrhoeae* 2nd Generation Detection Kit is an *in vitro* nucleic acid amplification test for the detection of *Neisseria gonorrhoeae* DNA extracted from clinical specimens or from BD ProbeTec ET™ lised samples.

INTRODUCTION

Neisseria gonorrhoeae is responsible for the Sexually Transmitted Disease (STD) Gonorrhoea. Gonorrhoea is a serious public health problem and is the most commonly reported infectious disease. *Neisseria gonorrhoeae* is a species of Gram-negative kidney bacteria most frequently spread during sexual contact. However, it can also be transmitted from the mother's genital tract to the newborn during birth. In women, the cervix is the most common site of Gonorrhoea, resulting in endocervicitis and urethritis, which can be complicated by pelvic inflammatory disease (PID). In men, Gonorrhoea causes anterior urethritis.

PRINCIPLE OF THE TEST

Duplica^{RealTime} *Neisseria gonorrhoeae* 2nd Generation Detection Kit is an assay for the DNA-based detection of *Neisseria gonorrhoeae*. The kit was designed to identify the sequences of **Cryptic plasmid pJD1** and **16S rRNA**. The reagents for the amplification are ready-to-use and provided in 3 separate tubes:

- **AMPLIFICATION MIX:** with Hot Start Taq DNA polymerase, nucleotides, MgCl₂ and buffer.
- **OLIGO MIX:** with primers and fluorogenic probes.
- **INTERNAL CONTROL:** to check the successful extraction and/or amplification process.

Duplica^{RealTime} *Neisseria gonorrhoeae* 2nd Generation Detection Kit is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulation of PCR amplification products by fluorescent DNA probes. The probe designed to detect the target carries the fluorophore FAM (6-carboxy-fluorescein) at the 5' end, while the probe detecting the internal control (IC), is labelled with the fluorophore HEX (hexa-chloro-fluorescein). Both probes have a non-fluorescent quencher at the 3' end. If excited, the whole probe doesn't emit fluorescence, since the proximity of the quencher to the reporter prevents the emission of the fluorescence from the reporter (quenching effect).

REAGENTS PROVIDED

Each kit contains enough reagents to perform 32 tests when used in 4 analytical sessions with 6 **samples** and 2 **controls** each or 5 **samples** and 3 **controls** each.

Kit Components

Reagent	Color Code	Storage (range, °C)	Volume (µl)	Quantity (tubes)
Oligo Mix (OM)*	Green Cap	-22÷-18	400	1
Amplification Mix (AM)	Blue Cap	-22÷-18	400	1
Control 1 (C1, Positive C.)	Red Cap	-22÷-18	>50	1
Control 2 (C2, Internal C.)**	Yellow Cap	-22÷-18	400	1
Reaction Blank (BM)***	White Cap	-22÷-18	1500	1

***protect the tube from direct light**

**** The Internal Control (C2) can be used as an Extraction/Amplification control (EAC) by adding it directly to the sample/lysis buffer mix during the DNA extraction procedure (for details see the paragraph "a) DNA purification" and refer to the specific protocol of the extraction kit) or as an Amplification Control by adding it directly to the PCR Mix (for details see the paragraph "c) PCR mix preparation")**

***** Used both in the extraction procedure as an extraction negative control (ENC, see paragraph "a) DNA purification"), and in amplification as reaction blank (see the paragraph "c) Preparation of the PCR mix").**

STORAGE AND HANDLING

All reagents must be stored at **-22÷-18°C** and can be used until the expiry date printed on the labels. Do not freeze and thaw the products more than six times.

MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to the specific handbook's section)
- Optical tubes or microplate for Real Time PCR
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 µl, 20-200 µl and 100-1000 µl)
- Disposable RNase/DNase-free tips with aerosol barriers
- Tube racks
- Desktop centrifuge
- PCR box
- Refrigerator
- Deep-freezer
- Thermal cycler for Real Time PCR

The kit has been developed to be used on SmartCycler® (Cepheid) and Rotor-Gene® Q (Qiagen) thermal cyclers. Other makes and models should be fully tested and evaluated for optimal performance by the user before reporting results.

The equipment should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure accurate PCR cycling and optimal performance.

PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood.
- If required, CLONIT offers the necessary technical support for the correct use of the kit.
- Carefully read this Instruction for Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipment. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite
- Use powder-free gloves. Do not leave fingerprints on optical section of Smart Tubes.
- Materials containing or potentially-containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood.
- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use.
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet (MSDS, available on request).
- The kit reagents, individual protective equipment, used materials, biological samples and test residuals must be disposed in accordance with local regulations.
- Patient Drug treatment may interfere with the final result of the molecular biology analysis.

OPERATING PROCEDURE

a) DNA purification

1. DNA extraction for clinical specimens

- Collect the first morning urine into a clean, sealed container with no preservative.
- The specimens should be stored between 2-8°C and ideally tested immediately or within 24 hours. Alternatively, the specimens should be frozen immediately upon receipt and stored up to ten days at -22±-18°C or colder. A single freeze thaw cycle should not affect the results.
CLONIT recommends for Automatic extraction Duplica®Prep Body Fluid kit (ref. EDI004200) with the Duplica®Prep instrument and GXT NA Extraction kit (BKR120802) with the GenoXtract® system (ref. DS083101)
- For Manual Extraction is recommended to use Bact Extra Pure Kit (ref. EDR004050).

2. DNA extraction for BD ProbeTec™ samples

- From previously lysed BD ProbeTec™ samples, please either use the DuplicaPrep Body Fluid kit (ref. EDI004200) according to the plasma protocol or Bact Extra Pure Kit (ref. EDR004050). Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

For **Duplica®Prep Body Fluid kit** and **Bact Extra Pure kit**, use **200µl** of the biological sample. If using the Internal Control (C2) as an Extraction/Amplification Control add **10µl** of the **C2** during the DNA purification procedure directly into the sample/lysis buffer mix. Elute in **50µl** of elution buffer.

For GXT NA Extraction kit, use **250µl** of the biological sample. If using the Internal Control (C2) as an Extraction/Amplification Control add **10µl** of the **C2** during the DNA purification procedure directly to the sample together with Proteinase K (not supplied). The use of polyA (not supplied) is optional. Elute in **50µl** of elution buffer.

For all extraction methods, if using the C2 as an Extraction/Amplification Control, prepare an **Extraction Negative Control (ENC)** by using the supplied Reaction Blank (**BM**), to which the **C2** should be added similarly to all samples to be extracted. In all cases, avoid mixing the C2 directly with the sample prior to the lysis step. The Positive Control (**C1**) does not require an extraction step and must be used directly during the amplification phase.

b) Thermal cycler Setup

Important points before starting: Refer to the specific handbook of the equipment used to set the thermal profile indicated in the **Thermal Profile Table**. We recommend to switch on the instrument and to set the thermal profile before preparing the reaction mix.

*N.B.: before starting the run is recommended to save the file as "**Neisseria gonorrhoeae Test**"; in this way it is possible to save the Thermal Profile and settings and recall them in subsequent runs.*

SmartCycler® platform

- Open the SmartCycler® software and click "Define Protocols".
- Click "New Protocol" (bottom left).
- Enter Name "**Neisseria gonorrhoeae Test**".
- When ready click "Save Protocol".

Rotor-Gene® Q platform

- Start the software and on the box *New Run* select *Advanced*
- Select a *new template* in *Empty Run* or a pre-existing one
- Select the Rotor Type of your instrument and then *Next*
- Type 25 µl in the reaction volume and then *Next*
- Select *Edit Profile* and set up the correct Thermal Profile as indicated in the table below
- Select *Gain Optimisation* and then flag the option *Perform Optimisation before 1st acquisition*
- On *Channel Settings* select the green/yellow fluorophores and tube position "1" to perform the optimization. Then close the window and select *Next* and *Start Run*

Thermal Profile Table (common for all the supported platforms)

TIME	TEMPERATURE	CYCLES
10 min	95°C	1
15 sec	95°C	45
60 sec	62°C	Fluorescence Acquisition

c) PCR mix preparation

The total reaction volume is **25 µl**.

For each experiment prepare a PCR mix for the required number of controls and **n+1** samples. If using the C2 as an Extraction/Amplification Control (EAC), these controls include 1 Extraction Negative Control (**ENC**) and 1 Positive Control (C1). CLONIT also recommends to include 1 Reaction Blank (**BM**) in every session to evaluate the presence of any contamination; in this case use directly 5µl of Reaction Blank (**BM**) provided, considering it as an extracted sample.

The mix must be prepared by mixing the reagents as indicated in the table:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Extracted DNA with EAC	5

The PCR mix has to be freshly prepared every time

After its preparation, aliquot **20µl of Master Mix** in the tubes or in the microplates well for PCR then add in each tube/well **5µl** from the **extracted DNA (or ENC)** or from **amplification controls (C+ and BM)**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermal cycler.

Alternatively, the Internal Control (C2) can also be used as a control in the amplification phase only. In this case, extract the samples according to the specific protocol of the kit in use, as described in the paragraph "a) DNA purification", excluding the step of the addition of the C2. The internal control (C2) should be added directly to the PCR MIX, prior to the aliquotation. The reagents of the PCR mix have to be mixed as indicated in the table below:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Internal Control (C2)	1
Extracted DNA without EAC	4

After its preparation, aliquot **21µl of Master Mix** in the tubes or in the microplates well for PCR then add in each tube/well **4µl** from the **extracted DNA** or from **amplification controls (C+ and BM)**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermal cycler.

d) ANALYSIS and INTERPRETATION of RESULTS

Important points before starting: For a detailed description on how to analyse data, refer to *System User's manuals*. **Always visually inspect the amplification plot for each sample tested versus C_T values obtained with the software.**

Results interpretation

Refer to the instrument-specific user guide to visualize the amplification plots for the entire plate/rotor. Detailed analysis of raw data depends on the real-time PCR instrument used. Baseline noise levels should either be set automatically or at predefined cycles.

The fluorescence in each channel indicates the hybridisation of the target specific probes:

- **Channel 1** for **FAM/Green= Target probe**
- **Channel 2** for **HEX/Cy3/Yellow= Internal Control probe.**

- If a sample shows a fluorescence in **FAM/Green (C_T>0)**, the sample is surely positive and the signal detected by **HEX/Cy3/Yellow** fluorophore (C_T≥0) is not relevant.
- When no signal at fluorophore **FAM/Green (C_T=0)** is detected, to confirm the negative result, the completed amplification of internal control and therefore the appearance of a fluorescent signal at **HEX/Cy3/Yellow (C_T>0)** fluorophore level, must be verified. Only in this case we can state that the sample is definitely negative.
- Condition in which no signal is detected indicates no DNA extraction or PCR inhibition. The sample must be repeated; a dilution 1:10 of the target DNA or complete DNA extraction is suggested.

In case the Internal Control (C2) has been used as an Extraction/Amplification Control (EAC), please be aware that:

- The Reaction Blank (**BM**) must be negative in both the **FAM/Green** and the **HEX/Cy3/Yellow** channel.
- The Positive Control (**C1**) must be positive in the **FAM/Green** channel.
- The Extraction Negative Control (**ENC**) must be negative in the **FAM/Green** channel, but positive in the **HEX/Cy3/Yellow** channel.

In case the Internal Control (C2) has been used only during amplification, adding it to the reaction mix, the signal of the HEX/Cy3/Yellow channel will also be present in the Positive Control (C1) and in the Reaction Blank (BM). For the Positive Control (C1), the **HEX/Cy3/Yellow** signal, due to competition with **FAM/Green** probe, may arise later or be absent. In this case the **HEX/Cy3/Yellow** signal is not relevant and the run is valid anyway.

If all these conditions have been met, the run is valid and it's possible to analyse the data; otherwise the run is not valid. It's responsibility of the user to validate the run.

Results Interpretation Table

FAM/Green	HEX/ Cy3/Yellow	Results
C _T > 0	Not relevant (C _T ≥ 0)	POSITIVE
C _T = 0	C _T > 0	NEGATIVE
C _T = 0	C _T = 0	INHIBITION

TROUBLESHOOTING

Problem 1: Weak or no signal of Positive Control, C1.

1. The PCR conditions didn't comply with the instructions. All sample results are IVALID:
 - Check the amplification protocol and select the fluorescence channel reported in the manual.
2. Deterioration of dyes and/or primers. The reagents storage conditions didn't comply with the instructions:
 - Check storage conditions.

Problem 2: Weak or no signal of the Internal Control, C2, in unknown samples AND Reaction Blank BM.

1. The PCR was inhibited:
 - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.
2. Pipetting error due to omitted reagents or samples:
 - Repeat the analysis starting from the PCR.
3. Deterioration of dyes and/or primers. The reagents storage conditions didn't comply with the instructions:
 - Check storage conditions.
4. Very low starting amount and/or low purity of genomic DNA. Improper DNA extraction:
 - Repeat the DNA purification.
5. Wrong channel/filter was chosen. The PCR conditions didn't comply with the instructions:
 - Check the PCR conditions and select the fluorescence channels reported in the protocol for the Unknown Sample detection.

Problem 3: FAM signal in Reaction Blank BM or in Extraction Negative Control (ENC).

1. Contamination during DNA Extraction and/or PCR preparation procedure. All samples results are INVALID:
 - Decontaminate the working area and all instruments.
 - Pipette the positive control (C1) at last.
 - Repeat the DNA Extraction and PCR preparation using a new set of reagents.

Problem 4: No FAM and HEX signals in unknown sample.

1. The PCR was inhibited:
 - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.

Problem 5: Wide Fluctuations in fluorescence values.

1. The PCR Master Mix is not well prepared:
 - Carefully repeat the PCR preparation procedure.
2. Air bubbles trapped in the PCR tubes:
 - Check the presence of air bubbles before starting a new run.

Problem 6: Absence of any fluorescent signal.

1. Verify the performance of the thermal cycler:
 - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions didn't comply with the instructions:
 - Check storage conditions.
 - Check the expiry date of the kit.

Problem 7: The thermal cycler gives an error message.

1. Refer to the real-time PCR instrument user manual or contact the local technical support of the real-time PCR instrument company.

Problem 8: The kit reagents left out of the storage range temperature.

1. These reagents must be stored **as indicated**. for a proper execution of the test. The performance of the product is not guaranteed if the reagents have not been properly stored.

Duplica^{RealTime} *Neisseria gonorrhoeae* 2nd Generation Detection Kit

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CLONIT S.r.l.