Duplica^{RealTime} Advanced *C. Difficile* A/B Kit

REV: EBR026032_IFU_REV.06D_ENITA

REF: EBR026032 -32 tests

Instructions For Use

INTENDED USE/PURPOSE

Duplica^{RealTime} Advanced *Clostridium Difficile* A/B Kit is an *in vitro* nucleic acid amplification test for the detection of the DNA of *Clostridium Difficile* toxins A and B in stool samples. The device is intended for detecting the presence of an infectious agent and for determining infective disease status.

INTRODUCTION

Clostridium difficile is regarded as an emerging pathogen and in recent years has been recognized as an increasingly important cause of nosocomial disease. It is a Gram positive spore-forming anaerobe and is the causative agent of several antibiotic-associated diarrheal diseases that are induced by treatment with antibiotics or by disruption of the normal gastrointestinal flora. These syndromes are collectively known as *C. difficile* infections or CDI, and include pseudomembranous colitis, toxic megacolon, and may lead to chronic infection, shock and death. Its spore-forming ability and resultant persistence has allowed *C. difficile* to become common in the hospital environment and it is now the leading cause of infectious diarrhoea in hospitals worldwide.

The major virulence factors of *C. difficile* are thought to be toxin A and toxin B. These toxins are encoded by the *tcdA* and *tcdB* genes, respectively, which are located within a 19.6 kb pathogenicity locus, the PaLoc.

Both toxin A and toxin B are proinflammatory and cytotoxic, causing disruption of the actin cytoskeleton and impairment of tight junctions in human intestinal epithelial cells, with resulting fluid accumulation and extensive damage to the large intestine.

PRINCIPLE OF THE TEST

Duplica^{RealTime} Advanced *Clostridium Difficile* A/B Kit is an assay for the DNA-based detection of the *Clostridium Difficile*. The reagents for the amplification are ready-to-use and provided in 3 separate tubes:

- **AMPLIFICATION MIX:** with Hot Start Taq DNA polymerase, nucleotides, MgCl₂ and buffer.
- **OLIGO MIX:** with primers and fluorogenic probes.
- **INTERNAL CONTROL OF AMPLIFICATION:** to check the successful of amplification reaction.

Duplica^{RealTime} Advanced *Clostridium Difficile* A/B Kit is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulated PCR amplification products by fluorescent DNA probes. The probe, designed to detect the target, carries the fluorophore FAM (6-carboxy-fluorescein) at the 5' end, while the probe specific for the internal control (IC), is labelled with the fluorophore HEX (hexa-chloro-fluorescein). Both the probes have a non-fluorescent quencher at the 3' end. If excited, the whole probe does not emit fluorescence since the proximity of the quencher to the reporter prevents the emission of the fluorescence from the reporter (quenching effect).

REAGENTS PROVIDED

Each kit contains enough reagents to perform 32 tests when used in 4 analytical sessions with 6 samples, 1 **Positive Control of Amplification (Control 1, C1)**, and 1 **Reaction Blank (BM)** each.

Kit Components

Reagent	Color Code	Storage (range, °C)	Volume (µl)	Quantity (tubes)
Oligo Mix (OM)*	Green Cap	-22÷-18	400	1
Amplification Mix (AM)	Blue Cap	-22÷-18	400	1
Control 1 (C1, Positive C.)	Red Cap	-22÷-18	>50	1
Control 2 (C2, Internal C.)	Yellow Cap	-22÷-18	>50	1
Reaction Blank (BM)	Neutral Cap	-22÷-18	>50	1

***protect the tube from direct light**

STORAGE AND HANDLING

All reagents must be stored at **-22÷-18°C** and can be used until the expiry date printed on the labels. Do not freeze and thaw the products more than six times.

MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to the specific handbook's section)
- Optical tubes or microplate for Real Time PCR
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 µl, 20-200 µl and 100-1000 µl)
- Disposable RNase/DNase-free tips with aerosol barriers
- Tube racks
- Vortex
- Desktop centrifuge
- PCR box
- Refrigerator
- Deep freezer
- Thermal cycler for Real Time PCR

The kit has been developed for Real time PCR platforms for use on RotorGene® Q (Qiagen), DX®/CFX (Biorad), AB7500 (Life Technologies) and q16/q32 (Novacyt) thermocyclers.

Other makes and models should be fully tested and evaluated for optimal performance by the user before reporting results.

The equipment should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure accurate PCR cycling and optimal performance.

PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood.
- If required, Clonit offers the necessary technical support for the correct use of the kit.
- Carefully read this Instruction for Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipment. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite
- Use powder-free gloves. Do not leave fingerprints on optical section of PCR Tubes.
- Materials containing or potentially containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood.

- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use.
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet (MSDS, available on request).
- The kit reagents, individual protective equipment, used materials, biological samples and test residuals must be disposed in accordance with local regulations.
- Patient Drug treatment may interfere with the final result of the molecular biology analysis.

OPERATING PROCEDURE

a) Specimens handling

- Collect the stool specimens into a clean, sealed container with no preservative.
- The specimens should be stored between +2°C and +8°C and ideally tested immediately or within 24 hours. Alternatively, the specimens should be frozen immediately upon receipt and stored up to ten days at -20°C or colder. A single freeze thaw cycle should not affect the results.

b) Specimen Preparation

- Mix stool as completely as possible.
- Clonit recommends Copan Fecal Swab for sample preparation. Open the swab's pouch; remove the tube and the swab and collect the specimen.
- For solid stools: Collect a small 3-4 mm diameter portion of mixed stool.
- Unscrew the cap aseptically and insert the swab into the tube; break the applicator at the point indicated by the colored line; replace the cap on the tube and screw the cap tightly.
- For liquid or semi-liquid stool: twist three times the swab into the specimen (in case of semi-liquid stool do not collect more than 3-4 mm diameter of mixed stool).
- Unscrew the cap aseptically and insert the swab into the tube; break the applicator at the point indicated by the colored line; replace the cap on the tube and screw the cap tightly.

c) DNA purification

- Vortex the tube containing the swab for 5-10 seconds until the solution becomes homogeneous.
- Collect 200 µl of the transport media and process with the Duplic α ®Prep equipment.
- Samples should be purified on the same day of the collection

For DNA purification Clonit recommends to use:

-Duplic α Prep NA Body Fluid kit (ref. EDI004200) with Duplic α PREP Automatic Extractor (ref. EDI001). Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

d) Thermalcycler Setup

Important points before starting: Refer to the specific handbook of the equipment used to set the thermal profile indicated in the **Thermal Profile Table**. We recommend to switch on the instrument and to set the thermal profile before preparing the reaction mix.

*N.B.: before starting the run is recommended to save the file as "**Advanced Clostridium Difficile A/B Test**"; this way it is possible to save the Thermal Profile and settings and recall them in subsequent runs.*

Rotor-Gene® Q platform

- Start the software and on the box *New Run* select *Advanced*
- Select a *new template* in *Empty Run* or a pre-existing one
- Select the Rotor Type of your instrument and then *Next*
- Type 25 µl in the reaction volume and then *Next*
- Select *Edit Profile* and set up the correct Thermal Profile as indicated in the table below
- Select *Gain Optimisation* and then flag the option *Perform Optimisation before 1st acquisition*
- On *Channel Settings* select the green/yellow fluorophores and tube position "1" to perform the optimization. Then close the window and select *Next* and *Start Run*

DX®/CFX platform

- Select *Create a new Experiment*
- In the *Protocol* section select *Create New* to set the thermal profile
- In *Plate Editor*: set FAM and as *Fluorophores* and define the sample name.

Applied Biosystems® 7500 platform

- Select *Create a new Experiment*
- In the *Experiment Properties* select *Quantitation-Standard curve* as experiment type and *Taqman® Reagents* as reagent type
- Select the Run Mode (*7500* or *7500 Fast*) and the Ramp Speed (*Standard*)
- In the *Setup* menu select *Run Method* to set the thermal profile
- In the *Setup* menu select *Plate Setup-Define Targets and Samples* to select Samples and assign VIC target to internal control (C2) and FAM to the positive control (C1). Leave the default value (*NFQ-MGB*) as quencher type
- In *Define Samples* insert the sample name
- Select the used wells in *Plate Setup-Assign Targets and Samples*
- Select *None* as passive reference

q16/q32 platform

- Select "New" and select the correct type of assay (the run conditions are automatically loaded and do not need to be entered manually).
- In "Samples" enter the sample name.
- The positive and negative controls are required to be placed in the specific well positions as indicated by the software.
- For the Q32 instrument ensure that positions A1, A8, D1, D8 will contain a tube to balance the heated lid.
- To start the run select the "Start Run" option, choose the file location to save the run file and the appropriate instrument.
- To start the run from a USB drive, when asked to select an instrument select "Start run from USB" and select the drive where the USB is located. The USB drive can now be placed in the instrument and the run will start automatically (confirmed by a sound in Q16 or by LEDs in Q32 changing from orange to blue).

Thermal Profile Table (common for all the supported platforms)

TIME	TEMPERATURE	CYCLES
5 min	95 °C	1
15 sec	95 °C	X 15 Fluorescent Acquisition off
60 sec	63 °C	
15 sec	95 °C	X 45 Fluorescent Acquisition on
60 sec	53 °C	

d) PCR mix preparation

The total reaction volume is **26 µl**.

For each experiment prepare a PCR mix for the positive control (**C1**), the **Reaction Blank (BM)** and **n+1** samples. The internal control of amplification (**C2**) should be added directly to the PCR MIX, prior the aliquotation. The reagents of the PCR mix have to be mixed as indicated in the table below:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Internal Control (C2)	1
Extracted DNA	5

The PCR mix has to be freshly prepared every time

After its preparation, aliquot **21 µl of Master Mix** in the tubes or in the microplates well for PCR then add

in each tube/well **5 µl** from the **extracted DNA** or **Control 1 DNA**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermalcycler.

e) ANALYSIS and INTERPRETATION of RESULTS

Important points before starting: For a detailed description on how to analyze data, refer to *System User's manuals*. **Always visually inspect the amplification plot for each sample tested versus C_T values obtained with the software.**

Results interpretation

Refer to the instrument-specific user guide to visualize the amplification plots for the entire plate/rotor. Detailed analysis of raw data depends on the real-time PCR instrument used. Baseline noise levels should either be set automatically or at predefined cycles. For q16/q32 instruments, the calculated **CT** (C_q) values can be found in the Results section. Thresholds for C_q determination are automatically applied by the software and cannot be changed.

- The fluorescence in each channel indicates the hybridisation of the specific probes:
 - o **Channel 1 for FAM/Green= Target probe**
 - o **Channel 2 for HEX/Cy3/Yellow= Internal Control probe.**
- If a sample shows a fluorescence in **FAM/Green (C_T>0)**, the sample is positive and the signal detected by **HEX/Cy3/Yellow** fluorophore (C_T≥0) is not relevant.
- When no signal at fluorophore **FAM/Green (C_T=0)** is detected, to confirm the negative result, the completed amplification of internal control and therefore the appearance of a fluorescent signal at **HEX/Cy3/Yellow (C_T>0)** fluorophore level, must be verified. Only in this case we can state that the sample is negative.
- Condition in which no signal is detected indicates PCR inhibition. The sample must be repeated and a dilution 1:10 of the target DNA is suggested.

Please be aware that:

- The Blank is positive in the **HEX/Cy3/Yellow** channel, but negative in the **FAM/Green** channel.
- The Positive Control is positive in both the **FAM/Green** and the **HEX/Cy3/Yellow** channels. The **HEX/Cy3/Yellow** signal, due to competition with **FAM/Green** probe, may arise later or be absent. In this case the **HEX/Cy3/Yellow** signal is not relevant and the run is valid.

If all these conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run.

Results Interpretation Table

FAM/Green	HEX/CY3/Yellow	Results
C _T > 0	Not relevant (C _T ≥ 0)	POSITIVE
C _T = 0	C _T > 0	NEGATIVE
C _T = 0	C _T = 0	INHIBITION

ANALYTICAL AND CLINICAL PERFORMANCE

a) ANALYTICAL SENSITIVITY AND SPECIFICITY

The **Analytical Sensitivity** of the test was determined using purified and quantified genomic DNA of known concentration. Test was carried out in four different days on twenty replicates; the analytical sensitivity obtained is less than 10 copies/reaction.

The **Analytical Specificity** of Duplica^{Real Time} *Advanced C. Difficile A/B* Kit is foremost granted by the selection of the primers and probes, as well as the selection of stringent reaction conditions. The primers and probes have been checked for possible homologies against the major genetic data banks. The analytical specificity was evaluated by experimentally testing the cross-reactivity of the Duplica^{Real Time} *Advanced C. Difficile A/B* Kit against potentially interfering microorganisms. None of the following

pathogens demonstrated reactivity: *Aspergillus*, *Salmonella*, *Listeria monocytogenes*, *Escherichia coli*, *Legionella ssp.*, *Adenovirus*, *Chlamydia trachomatis*, *Campylobacter jejuni*.

b) CLINICAL SENSITIVITY AND SPECIFICITY

The clinical sensitivity and specificity were tested on 47 biological samples comparing the kit with a reference test (Xpert C. difficile Method). The results demonstrated an overall sensitivity of 100% and specificity of 100%.

c) INTRA-ASSAY PRECISION

In order to demonstrate the reproducibility of the test within the same analytical session, 2 positive samples, 2 negative samples, 1 positive control and 1 negative control were tested in 4 replicates in the same run. Results show 100% of accuracy between expected and obtained results. CV% is below 10% in agreement with the acceptance criteria.

d) INTER-ASSAY PRECISION

In order to demonstrate the reproducibility of the test among different analytical sessions, 2 positive samples, 2 negative samples, 1 positive control and 1 negative control were tested as single replicates in 4 independent runs. Results show 100% of accuracy between expected and obtained results. CV% is below 10% in agreement with the acceptance criteria.

BIBLIOGRAPHY

"Toxin B is essential for Virulence of Clostridium difficile", Nature-07822 – D.Lyras et al., 2009

TROUBLESHOOTING

Problem 1: Weak or no signal of Positive Control, C1.

1. The PCR conditions did not comply with the instructions. All sample results are IVALID:
 - Check the amplification protocol and select the fluorescence channel reported in the manual.
2. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.

Problem 2: Weak or no signal of the Internal Control, C2 in unknown samples AND Reaction Blank BM.

1. The PCR was inhibited:
 - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.
2. Pipetting error due to omitted reagents or samples:
 - Repeat the analysis starting from the PCR.
3. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.
4. Very low starting amount and/or low purity of genomic DNA. Improper DNA extraction:
 - Repeat the DNA purification.
5. Wrong channel/filter was chosen. The PCR conditions did not comply with the instructions:
 - Check the PCR conditions and select the fluorescence channels reported in the protocol for the Unknown Sample detection.

Problem 3: FAM signal in Reaction Blank BM.

1. Contamination during PCR preparation procedure. All samples results are INVALID:
 - Decontaminate the working area and all instruments.
 - Pipette the controls C1 and C2 at last.
 - Repeat the PCR preparation using a new set of reagents.

Problem 4: No FAM and HEX signals in unknown sample.

1. The PCR was inhibited:

- Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.

Problem 5: Wide Fluctuations in fluorescence values.

1. The PCR Master Mix is not well prepared:
 - Carefully repeat the PCR preparation procedure.
2. Air bubbles trapped in the PCR tubes:
 - Check the presence of air bubbles before starting a new run.

Problem 6: Absence of any fluorescent signal.

1. Verify the performance of the thermal cycler:
 - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions did not comply with the instructions:
 - Check storage conditions.
 - Check the expiry date of the kit.

Problem 7: The thermal cycler gives an error message.

1. Refer to the real-time PCR instrument user manual or contact the local technical support of the real-time PCR instrument company.

Problem 8: The kit reagents left out of the storage range temperature.

1. These reagents must be stored **as indicated** for a proper execution of the test. The performance of the product is not guaranteed if the reagents have not been properly stored.

This page left intentionally blank

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CLONIT S.r.l.