For Professional Use Only

Duplica^{RealTime} Advanced *Mycoplasma Genitalium* Kit

REV. EBR027032_IFU_REV.03D_ENITA

REF EBR027032 -32 tests

Instructions For Use

INTENDED USE

Duplica^{RealTime} Advanced *Mycoplasma Genitalium* Kit is an *in vitro* nucleic acid amplification test for the detection of *Mycoplasma genitalium* (MG) DNA, extracted from clinical specimens (such as urine and swab).

INTRODUCTION

Mycoplasma genitalium (MG) is a member of genital mycoplasmas, which is emerging as an important causative agent of sexually transmitted infections both in males and females. The advent of polymerase chain reaction and other molecular methods have made studies on *M. genitalium* more feasible, which is otherwise a difficult organism to isolate. Besides *Chlamydia trachomatis*, *M. genitalium* is now an important and established cause of non gonococcal urethritis (NGU) in men, more so in persistent and recurrent NGU. Multiple studies have also shown a positive association of *M. genitalium* with mucopurulent cervicitis and vaginal discharge in females as well. The evidence for *M. genitalium* pelvic inflammatory diseases and infertility are quite convincing and indicate that this organism has potential to cause ascending infection. Lack of clear association with *M. genitalium* has been reported for bacterial vaginosis and adverse pregnancy outcomes. Diagnosis of *M. genitalium* infections is performed exclusively using nucleic acid amplification tests (NAATs), owing to poor or slow growth of bacterium in culture.

PRINCIPLE OF THE TEST

Duplica^{RealTime} Advanced *Mycoplasma Genitalium* Kit is an assay for the DNA-based detection of the *Mycoplasma genitalium*. The kit is designed to amplify a portion of MgPa and 16S genes. The reagents for the amplification are ready-to-use and provided in 3 separate tubes:

- **AMPLIFICATION MIX:** with Hot Start Taq DNA polymerase, nucleotides, MgCl₂ and buffer.
- **OLIGO MIX:** with primers and fluorogenic probes.
- **INTERNAL CONTROL OF AMPLIFICATION:** to check the successful extraction and/or amplification process.

Duplica^{RealTime} Advanced *Mycoplasma Genitalium* kit is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulation of PCR amplification products by fluorescent DNA probes. The probe designed to detect the target carries the fluorophore FAM (6-carboxy-fluorescein) at the 5' end, while the probe detecting the internal control (IC), is labelled with the fluorophore HEX (hexa-chloro-fluorescein). Both the probes have a non fluorescent quencher at the 3' end. If excited, the whole probe does not emit fluorescence since the proximity of the quencher to the reporter prevents the emission of the fluorescence from the reporter (quenching effect).

REAGENTS PROVIDED

Each kit contains enough reagents to perform 32 tests when used in 4 analytical sessions with 6 **samples** and 2 **controls** each or 5 **samples** and 3 **controls** each.

Kit Components

Reagent	Color Code	Storage (range, °C)	Volume (μl)	Quantity (tubes)
Oligo Mix (OM)*	Green Cap	-22÷-18	400	1
Amplification Mix (AM)	Blue Cap	-22÷-18	400	1
Control 1 (C1, Positive C.)	Red Cap	-22÷-18	>50	1
Control 2 (C2, Internal C.)**	Yellow Cap	-22÷-18	400	1
Reaction Blank (BM)***	White Cap	-22÷-18	1500	1

***protect the tube from direct light**

**** The Internal Control (C2) can be used as an Extraction/Amplification control (EAC) by adding it directly to the sample/lysis buffer mix during the DNA extraction procedure (for details see the paragraph "a) DNA purification" and refer to the specific protocol of the extraction kit) or as an Amplification Control by adding it directly to the PCR Mix (for details see the paragraph "c) PCR mix preparation")**

***** Used both in the extraction procedure as an extraction negative control (ENC, see paragraph "a) DNA purification"), and in amplification as reaction blank (see the paragraph "c) Preparation of the PCR mix").**

STORAGE AND HANDLING

All reagents must be stored at **-22÷-18°C** and can be used until the expiry date printed on the labels.

Do not freeze and thaw the products more than six times.

MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to the specific handbook's section)
- Optical tubes or microplate for Real Time PCR
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 μ l, 20-200 μ l and 100-1000 μ l)
- Disposable RNase/DNase-free tips with aerosol barriers
- Tube racks
- Vortex
- Desktop centrifuge
- PCR box
- Refrigerator
- Deep-freezer
- Thermal cycler for Real Time PCR

The kit has been developed to be used on SmartCycler® (Cepheid) Rotor-Gene® Q (Qiagen), DX®/CFX (Bio-Rad) and Applied Biosystems® 7500 (Life Technologies™) thermal cyclers.

Other makes and models should be fully tested and evaluated for optimal performance by the user before reporting results.

The equipments should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure accurate PCR cycling and optimal performance.

PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood.
- If required, Clonit offers the necessary technical support for the correct use of the kit.
- Carefully read this Instruction for Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipments. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite
- Use powder-free gloves. Do not leave fingerprints on optical section of Smart Tubes.
- Materials containing or potentially containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood.

- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use.
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet (MSDS, available on request).
- The kit reagents, individual protective equipments, used materials, biological samples and test residuals must be disposed in accordance with local regulations.
- Patient Drug treatment may interfere with the final result of the molecular biology analysis.

OPERATING PROCEDURE

a) Specimens handling

- Collect the first morning urine into a clean, sealed container with no preservative.
- The specimens should be stored between 2 ÷ 8°C and ideally tested immediately or within 24 hours. Alternatively, the specimens should be frozen immediately upon receipt and stored up to ten days at -20°C or colder. A single freeze thaw cycle should not affect the results.

b) DNA purification

- Vortex the tube containing the swab for 5-10 seconds until the solution becomes homogeneous or resuspend the urine. Collect 200 µL of the transport media or urine.
- Samples should be purified on the same day of the collection.

For genomic DNA purification Clonit recommends to use:

- Duplica®Prep NA Body Fluid kit (ref. EDI004200) with Duplica® PREP Automatic Extractor (ref. EDI001).
- GXT NA Extraction kit (BKR120802) with the GenoXtract® system (ref. DS083101)
- Bact Extra Pure Kit (ref. EDR004050) for manual extraction.

Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

For **Duplica®Prep Body Fluid kit** and **Bact Extra Pure kit**, use **200µl** of the biological sample. If using the Internal Control (C2) as an Extraction/Amplification Control add **10µl** of the **C2** during the DNA purification procedure directly into the sample/lysis buffer mix. Elute in **50µl** of elution buffer.

For GXT NA Extraction kit, use **250µl** of the biological sample. If using the Internal Control (C2) as an Extraction/Amplification Control add **10µl** of the **C2** during the DNA purification procedure directly to the sample together with Proteinase K (not supplied). The use of polyA (not supplied) is optional. Elute in **50µl** of elution buffer. For all extraction methods, if using the C2 as an Extraction/Amplification Control, prepare an **Extraction Negative Control (ENC)** by using the supplied Reaction Blank (**BM**), to which the **C2** should be added similarly to all samples to be extracted. In all cases, avoid mixing the C2 directly with the sample prior to the lysis step.

The Positive Control (**C1**) does not require an extraction step and must be used directly during the amplification phase

c) Thermal cycler Setup

Important points before starting: Refer to the specific handbook of the equipment used to set the thermal profile indicated in the **Thermal Profile Table**. We recommend to switch on the instrument and to set the thermal profile before preparing the reaction mix.

N.B.: before starting the run is recommended to save the file as "**Advanced Mycoplasma Genitalium Test**"; in this way is possible to save the Thermal Profile and settings and recall them in subsequent runs.

SmartCycler® platform

- Open the SmartCycler® software and click "Define Protocols".
- Click "New Protocol" (bottom left).
- Enter Name "**Advanced Mycoplasma Genitalium Test**".
- When ready click "Save Protocol".

Rotor-Gene® Q platform

- Start the software and on the box *New Run* select *Advanced*
- Select a *new template* in *Empty Run* or a pre-existing one
- Select the Rotor Type of your instrument and then *Next*
- Type 25 µl in the reaction volume and then *Next*
- Select *Edit Profile* and set up the correct Thermal Profile as indicated in the table below
- Select *Gain Optimisation* and then flag the option *Perform Optimisation before 1st acquisition*

- On *Channel Settings* select the green/yellow fluorophores and tube position "1" to perform the optimization. Then close the window and select *Next* and *Start Run*

DX®/CFX platform

- Select *Create a new Experiment*
- In the *Protocol* section select *Create New* to set the thermal profile
- In *Plate Editor*: set FAM and HEX as *Fluorophores* and define the sample name.

Applied Biosystems® 7500 platform

- Select *Create a new Experiment*
- In the *Experiment Properties* select *Quantitation-Standard curve* as experiment type and *Taqman® Reagents* as reagent type
- Select the Run Mode (*7500* or *7500 Fast*) and the Ramp Speed (*Standard*)
- In the *Setup* menu select *Run Method* to set the thermal profile
- In the *Setup* menu select *Plate Setup-Define Targets and Samples* to select Samples and assign VIC target to internal control (C2) and FAM to the positive control (C1). Leave the default value (*NFQ-MGB*) as quencher type
- In *Define Samples* insert the sample name
- Select the used wells in *Plate Setup-Assign Targets and Samples*
- Select *None* as passive reference

Thermal Profile Table (common for all the supported platforms)

TIME	TEMPERATURE	CYCLES
5 min	95 °C	1
15 sec	95 °C	X 10 Fluorescent Acquisition off
60 sec	62 °C	
15 sec	95 °C	X 45 Fluorescent Acquisition on
60 sec	52 °C	

d) PCR mix preparation

The total reaction volume is **25 µl**.

For each experiment prepare a PCR mix for the required number of controls and **n+1** samples. If using the C2 as an Extraction/Amplification Control (EAC), these controls include 1 Extraction Negative Control (**ENC**) and 1 Positive Control (C1). Clonit also recommends to include 1 Reaction Blank (**BM**) in every session to evaluate the presence of any contamination; in this case use directly 5µl of Reaction Blank (**BM**) provided, considering it as an extracted sample.

The mix must be prepared by mixing the reagents as indicated in the table:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Extracted DNA with EAC	5

The PCR mix has to be freshly prepared every time

After its preparation, aliquot **20µl of Master Mix** in the tubes or in the microplates well for PCR then add in each tube/well **5µl** from the **extracted DNA (or ENC)** or from **amplification controls (C+ and BM)**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermal cycler.

Alternatively, the Internal Control (C2) can also be used as a control in the amplification phase only. In this case, extract the samples according to the specific protocol of the kit in use, as described in the paragraph "a) DNA purification", excluding the step of the addition of the C2. The internal control (C2) should be added directly to the PCR MIX, prior to the aliquotation. The reagents of the PCR mix have to be mixed as indicated in the following table:

REAGENT	VOLUME (μ l)
Amplification mix	10
Oligo mix	10
Internal Control (C2)	1
Extracted DNA without EAC	4

After its preparation, aliquot **21 μ l of Master Mix** in the tubes or in the microplates well for PCR then add in each tube/well **4 μ l** from the **extracted DNA** or from **amplification controls (C+ and BM)**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermal cycler.

e) ANALYSIS and INTERPRETATION of RESULTS

Important points before starting: For a detailed description on how to analyze data, refer to *System User's manuals*. **Always visually inspect the amplification plot for each sample tested versus C_T values obtained with the software.**

Results interpretation

Refer to the instrument-specific user guide to visualize the amplification plots for the entire plate/rotor. Detailed analysis of raw data depends on the real-time PCR instrument used. Baseline noise levels should either be set automatically or at predefined cycles.

The fluorescence in each channel indicates the hybridisation of the specific probes:

- **Channel 1** for **FAM/Green= Target probe**
- **Channel 2** for **HEX/Cy3/Yellow= Internal Control probe.**
- If a sample shows a fluorescence in **FAM/Green ($C_T > 0$)**, the sample is positive and the signal detected by **HEX/Cy3/Yellow** fluorophore ($C_T \geq 0$) is not relevant.
- When no signal at fluorophore **FAM/Green ($C_T = 0$)** is detected, to confirm the negative result, the completed amplification of internal control and therefore the appearance of a fluorescent signal at **HEX/Cy3/Yellow ($C_T > 0$)** fluorophore level, must be verified. Only in this case we can state that the sample is negative.
- Condition in which no signal is detected indicates no DNA extraction or PCR inhibition. The sample must be repeated; a dilution 1:10 of the target DNA or complete DNA extraction is suggested.

In case the Internal Control (C2) has been used as an Extraction/Amplification Control (EAC), please be aware that:

- The Reaction Blank (**BM**) must be negative in both the **FAM/Green** and the **HEX/Cy3/Yellow** channel.
- The Positive Control (**C1**) must be positive in the **FAM/Green** channel.
- The Extraction Negative Control (**ENC**) must be negative in the **FAM/Green** channel, but positive in the **HEX/Cy3/Yellow** channel.

In case the Internal Control (C2) has been used only during amplification, adding it to the reaction mix, the signal of the HEX/Cy3/Yellow channel will also be present in the Positive Control (C1) and in the Reaction Blank (BM). For the Positive Control (C1), the **HEX/Cy3/Yellow** signal, due to competition with **FAM/Green** probe, may arise later or be absent. In this case the **HEX/Cy3/Yellow** signal is not relevant and the run is valid anyway.

If all these conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run.

Results Interpretation Table

FAM/Green	HEX/CY3/Yellow	Results
$C_T > 0$	Not relevant ($C_T \geq 0$)	POSITIVE
$C_T = 0$	$C_T > 0$	NEGATIVE
$C_T = 0$	$C_T = 0$	INHIBITION

TROUBLESHOOTING

Problem 1: Weak or no signal of Positive Control, C1.

1. The PCR conditions did not comply with the instructions. All sample results are INVALID:
 - Check the amplification protocol and select the fluorescence channel reported in the manual.
2. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.

Problem 2: Weak or no signal of the Internal Control, C2 in unknown samples AND Reaction Blank BM.

1. The PCR was inhibited:
 - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.
2. Pipetting error due to omitted reagents or samples:
 - Repeat the analysis starting from the PCR
3. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.
4. Very low starting amount and/or low purity of genomic DNA. Improper DNA extraction:
 - Repeat the DNA purification.
5. Wrong channel/filter was chosen. The PCR conditions did not comply with the instructions:
 - Check the PCR conditions and select the fluorescence channels reported in the protocol for the Unknown Sample detection.

Problem 3: FAM signal in Reaction Blank BM or in Extraction Negative Control (ENC).

Contamination during DNA Extraction and/or PCR preparation procedure. All samples results are INVALID:

1. Decontaminate the working area and all instruments.
2. Pipette the controls C1 at last.
3. Repeat the DNA Extraction and the PCR preparation using a new set of reagents.

Problem 4: No FAM and HEX signals in unknown sample.

1. The PCR was inhibited:
 - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.

Problem 5: Wide Fluctuations in fluorescence values.

1. The PCR Master Mix is not well prepared:
 - Carefully repeat the PCR preparation procedure.
2. Air bubbles trapped in the PCR tubes:
 - Check the presence of air bubbles before starting a new run.

Problem 6: Absence of any fluorescent signal.

1. Verify the performance of the thermal cycler:
 - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions did not comply with the instructions:
 - Check storage conditions.
 - Check the expiry date of the kit.

Problem 7: The thermal cycler gives an error message.

1. Refer to the real-time PCR instrument user manual or contact the local technical support of the real-time PCR instrument company.

Problem 8: The kit reagents left out of the storage range temperature.

1. These reagents must be stored **as indicated** for a proper execution of the test. The performance of the product is not guaranteed if the reagents have not been properly stored.

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CLONIT S.r.l.