Duplica^{RealTime} Mix & Match HFE S65C Kit

REV. EER046032_IFU_REV.02D_ENITA

REF EER046032-32 tests

Instructions For Use

INTENDED USE

Duplica^{RealTime} Mix & Match HFE S65C Kit is an *in vitro* nucleic acid amplification test for the detection of the aminoacidic substitution S65C (allelic variant g.26091185A>T) of HFE gene, in human genomic DNA extracted from peripheral whole blood samples, collected in EDTA.

Clonit **Duplica^{RealTime} Mix & Match kits** belong to the *HFE panel* (EER044032, EER045032 and EER046032) and share a common thermal profile.

INTRODUCTION

Hereditary haemochromatosis (HH) is an inherited autosomal recessive disorder of iron metabolism. Due to excessive intestinal absorption, iron accumulates in parenchymal cells of the liver, pancreas, heart and other organs with resulting damage to their structure and impairment of their function. It is one of the most common genetic diseases in Caucasians with a prevalence of nearly 1 in 300. Although the symptoms of the disease are often nonspecific, much of the organ damage is irreversible once it has occurred. Early detection and therapy are therefore very important as a part of preventive medicine. The discovery of the responsible gene HFE in 1996 enabled molecular analysis to be included in the diagnostic strategy for HH. A number of different HFE mutations have been reported so far. The majority of HH cases (52–96%) in European regions are associated with a homozygous **g.845G>A** mutation within **exon 4** of the HFE gene, which results in amino acid change at position **282** from **cysteine to tyrosine (C282Y)**. A second mutant allele **g.187C>G** detected with relatively high frequency occurs within the **exon 2** of the HFE gene where **aspartate** replaces **histidine** at aminoacidic position **63 (H63D)**. The contribution of this allele to iron overload is most relevant in the case of combined heterozygosity with C282Y allele (C282Y/H63D). The third common mutation of HFE is **g.193A>T** substitution in **exon 2 (serine** replaced by **cysteine** in the protein sequence, **S65C**) and was shown to be generally benign, although a C282Y/S65C genotype may confer a slight increase in disease risk, contributing to a mild disease phenotype.

PRINCIPLE OF THE TEST

Duplica^{RealTime} Mix & Match HFE S65C Kit is designed to identify the aminoacidic substitution S65C in the human HFE gene (point mutation 193A>T). Two reaction mixes are provided for the amplification:

- **AMPLIFICATION MIX**, containing Hot Start Taq DNA polymerase, nucleotides, MgCl₂ and buffer.
- **OLIGO MIX**, containing primers and fluorogenic probes.

Duplica^{RealTime} Mix & Match HFE S65C Kit is based on specific recognition and amplification of target sequences by PCR, and the simultaneous detection of the accumulation of PCR amplification products by fluorescent DNA probes. In particular, the probe designed to detect the Wild Type allele carries the fluorophore FAM (6-carboxy-fluorescein) at the 5' end, while the probe detecting the Mutated allele, is labelled with the fluorophore HEX (hexa-chloro-fluorescein). Both the probes have a non fluorescent black quencher at the 3' end. If excited, the whole probe does not emit fluorescence since the proximity of the quencher to the reporter prevents the emission of the fluorescence from the reporter (quenching effect).

REAGENTS PROVIDED

Each kit contains enough reagents to perform 32 tests when used in 4 analytical sessions with 5 **samples**, 1 **Wild Type control (Control 1, C1)**, 1 **Mutated control (Control 2, C2)** and one **Reaction Blank (BM)** each.

Kit Components

Reagent	Color Code	Storage (range, °C)	Volume (µl)	Quantity (tubes)
Oligo Mix (OM)*	Green Cap	-22÷-18	400	1
Amplification Mix (AM)	Blue Cap	-22÷-18	400	1
Control 1 (C1, Wild Type)	Red Cap	-22÷-18	>50	1
Control 2 (C2, Mutated)	Yellow Cap	-22÷-18	>50	1
Reaction Blank (BM)	White Cap	-22÷-18	>50	1

* **protect the tube from direct light**

STORAGE AND HANDLING

All reagents must be stored at **-22÷-18°C** and can be used until the expiry date printed on the labels. Do not freeze and thaw the products more than six times.

MATERIAL REQUIRED BUT NOT PROVIDED

- Extraction kit for DNA purification (refer to the specific handbook's section)
- Optical tubes or microplate for Real Time PCR
- Disposable powder-free gloves and laboratory coat
- Variable volumes pipettes (5-20 µl, 20-200 µl and 100-1000 µl)
- Disposable RNase/DNase-free tips with aerosol barriers
- Tube racks
- Desktop centrifuge
- PCR box
- Refrigerator
- Deep-freezer
- Thermal cycler for Real Time PCR

The kit has been optimized to be used on DX®/CFX (Bio-Rad), Rotor-Gene® Q (Qiagen) and Applied Biosystems® 7500 (Life Technologies™) thermal cyclers. Other makes and models should be fully tested and evaluated for optimal performance by the user before reporting results.

The equipments should be regularly maintained, in accordance with the manufacturer's instructions, and calibrated to ensure an optimal performance.

PRECAUTIONS AND WARNINGS

- In compliance with Good Laboratory Practice, define three separate laboratory's areas for: DNA extraction, PCR reaction mix preparation; manipulation of controls provided with the kit. Each area must have dedicated pipettes and laminar flow hood
- If required, Clonit offers the necessary technical support for the correct use of the kit
- Carefully read this Instruction for Use before using the kit
- Do not use the reagents after the expiry date
- Thaw and carefully mix the reagents of the kit before use
- Do not mix the reagents from different lots of the product
- Use calibrated and regularly checked pipettes and instrumentation only
- Use dedicated laboratory equipments. Change gloves frequently
- Periodically wipe the working area with 0,5% hypochlorite
- Use powder-free gloves. Do not leave fingerprints on optical caps. Do not write on caps as this may cause an interference with fluorescent detection
- Materials containing or potentially containing infectious agents must always be manipulated in a separated microbiological safety room under a Biohazard biological hood
- In case of damaged package, contact the technical support before using the kit
- Do not use the product when stored at temperatures other than those indicated on the labels or described in this Instructions For Use
- In case of spillage of the kit contents, please refer to the specific Material Safety Data Sheet (MSDS, available on request)
- The kit reagents, individual protective equipments, used materials, biological samples and test residuals must be

- disposed in accordance with local regulations
- Patient Drug treatment may interfere with the final result of the molecular biology analysis

OPERATING PROCEDURE

a) DNA purification

For genomic DNA purification Clonit recommends to use:

- Duplicα Blood DNA kit (ref. EDI002250) and Duplicα NA Body Fluid kit (ref. EDI004200) with Duplicα® PREP Automatic Extractor (ref. EDI001) for automated purification.

For manual extraction/purification:

- Fasset DNA Releaser (ref. EMR057050) for **fresh** peripheral blood (i.e. stored for up to 24 hours at 2÷8°C).
- Spin DNA Purification Kit (ref. EMR061050) for **frozen** blood or stored at 2÷8°C for more than 24 hours.

Other extraction reagents and methods should be fully tested and evaluated for optimal performance by the user before reporting results.

Attention! Use blood samples in EDTA anti-coagulant solution.

b) Thermalcycler Setup

Refer to the specific handbook of the equipment used to set the thermal profile indicated in the **Thermal Profile Table**. We recommend to switch on the instrument and to set the thermal profile before preparing the reaction mix.

DX®/CFX platform

- Select *Create a new Experiment*
- In the *Protocol* section select *Create New* to set the Thermal Profile
- In *Plate Editor*: set FAM (Wild Type genotype, Allele 1) and HEX (Mutated genotype, Allele 2) as *Fluorophores* and define the sample name

Applied Biosystems® 7500 platform

- Select *Create a new Experiment*
- In the *Experiment Properties* select *Quantitation-Standard curve* as experiment type and *Taqman® Reagents* as reagent type
- Select the Run Mode (*7500* or *7500 Fast*) and the Ramp Speed (*Standard*)
- In the *Setup* menu select *Run Method* to set the Thermal Profile
- In the *Setup* menu select *Plate Setup-Define Targets and Samples* to select Samples and assign VIC target to MUT (Mutated genotype) and FAM target to WT (Wild-Type genotype). Leave the default value (*NFQ-MGB*) as quencher type
- In *Define Samples* insert the sample name
- Select the used wells in *Plate Setup-Assign Targets and Samples*
- Select *None* as passive reference

Rotor-Gene® Q platform

- Start the software and on the box *New Run* select *Advanced*
- Select a *new template* in *Empty Run* or a pre-existing one
- Select the Rotor Type of your instrument and then *Next*
- Type 25 µl in the reaction volume and then *Next*
- Select *Edit Profile* and set up the correct Thermal Profile as indicated in the table below
- Select *Gain Optimisation* and then flag the option *Perform Optimisation before 1st acquisition*
- On *Channel Settings* select the green/yellow fluorophores and tube position "1" to perform the optimization. Then close the window and select *Next* and *Start Run*

Thermal Profile Table (common for all the supported platforms)

Time	Temperature	Cycles
10 min	95°C	1
15 sec	95°C	35
60 sec	60°C	Fluorescence Acquisition

c) Preparation of PCR mix

The total reaction volume is **25 µl**.

For each experiment prepare a PCR mix for the **2 controls (C1 and C2)**, **1 Reaction Blank (BM)** and **n+1 samples**. The reagents of the PCR mix have to be mixed as indicated in the table below:

REAGENT	VOLUME (µl)
Amplification mix	10
Oligo mix	10
Extracted DNA	5

The PCR mix has to be freshly prepared every time.

After its preparation, aliquot **20 µl of Master Mix** in the tubes or in the microplates well for PCR then add in each tube/well **5 µl** (100-250 ng/reaction) from the **extracted DNA** or **control DNA**, place in order the tubes/microplate in the instrument and start the program of amplification. At the end of the program remove the tubes/microplate from the thermal cycler.

d) ANALYSIS and INTERPRETATION of the RESULTS

Real Time PCR curves analysis

Refer to the instrument specific user guide to visualize the amplification plots for the entire plate/rotor. Detailed analysis of raw data depends on the Real Time PCR instrument used. Baseline noise levels should either be set automatically or at predefined cycles.

The fluorescence in each channel indicates the hybridisation of the allelic specific probes: **Channel 1** for **FAM/green= Wild Type allele probe** and **Channel 2** for **HEX/VIC/Yellow= Mutated allele probe**. If a sample shows a fluorescence in **FAM/Green**, the sample has the **Wild Type allele**. If a sample shows a fluorescence in **HEX/VIC/Yellow**, the sample has the **Mutated allele**.

Therefore, if only a **FAM/Green signal** is detected the sample is **Homozygous Wild Type**, whereas if only a **HEX/VIC/Yellow signal** is detected the sample is **Homozygous Mutated**. Finally, if both **FAM/Green** and **HEX/VIC/Yellow** are detected the sample is **Heterozygous**.

Condition in which no signal is detected indicates PCR inhibition. In this case, refer to troubleshooting.

Control 1 (Wild Type), Control 2 (Mutant) and reaction blank (BM) are provided in order to properly set the threshold line before samples analysis. After the run, the threshold line has to be set so that: **Control 1 is detected in FAM/Green and not detected in HEX/VIC/Yellow, whereas Control 2 is detected in HEX/VIC/Yellow and not detected in FAM/Green. The Reaction Blank must not be detected in any channel.**

If these three conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run.

Results Interpretation Table

FAM/Green	HEX/VIC/Yellow	Results
Detected	Not detected	Wild Type
Not detected	Detected	Mutated
Detected	Detected	Heterozygous
Not detected	Not detected	Inhibition


Scatter Plot analysis

Additionally, the **Scatter Plot Analysis** can be used to visualize the results for **Rotor-Gene®Q** and **DX®/CFX** platforms. In this analysis each sample is identified by a spot representing the end-point fluorescence values for the two dyes.

Rotor-Gene® Q platform

- In the *Analysis* menu select *Other*
- Select *Scatter Graph Analysis*, highlight both *Green* and *Yellow* channels and press *Show* to display the scatter plot. The scatter plot is divided in four areas corresponding to samples classified as Wild Type, Mutant, Heterozygous.


In the *Scatter Analysis graph* window that appears select with the mouse the area corresponding to a specific Genotype and name it as show in the figure below (*Green* → *Wild-type*; *Yellow* → *Mutated*; and *Green/Yellow* → *Heterozygous*):

DX®/CFX platform

- Select *Ct Determination Mode* and *Regression* in the *Settings* menu
- In the *View/Edit Plate* menu set the reaction blank as *NTC*
- Select *Quantitation*. In the *Ct* table in the lower right part of the screen no amplification must be detected (*Ct* value=N/A for both fluorophores) for the **blank/NTC** sample. **Control 1 (C1)** must be detected (*Ct* value≠N/A) in FAM and non detected (*Ct* value=N/A) in HEX, whereas **Control 2 (C2)** detected in HEX and not detected in FAM. If these conditions have been met, the run is valid and it is possible to analyse the data; otherwise the run is not valid. It is responsibility of the user to validate the run by checking that these conditions have been met
- **Unknown (Unk)** samples must be detected in at least one fluorescence channel. Condition in which no signal is detected (*Ct* value=N/A for both fluorophores) indicates PCR inhibition and **samples presenting this condition must be excluded from the following Allelic Discrimination analysis**. For troubleshooting refer to the point 1 and 2 of the Troubleshooting section.
- In the *Allelic Discrimination* menu select *RFU* and flag the *Normalize Data* in the *Display Mode* to display the relative Scatter Plot
- For both Fluorophores set manually the Threshold to 0,15

The Scatter Plot displays three main quadrants in which each spot represents the two fluorescent signals for a single sample/control:

The corresponding genotype of the samples is reported in a table displayed close to the Scatter Plot:

Bottom Right Quadrant	Wild Type Samples (WT- <i>Allele 1</i>)
Upper Left Quadrant	Mutated Samples (MUT- <i>Allele 2</i>)
Upper Right Quadrant	Heterozygote Samples (HET- <i>both Alleles</i>)

TROUBLESHOOTING

Problem 1: Weak or no signal of unknown samples.

1. The PCR was inhibited:
 - Make sure to use a recommended DNA purification method and carefully follow the manufacturer's instructions.
2. Pipetting error due to omitted reagents or samples:
 - Repeat the analysis starting from the PCR.
3. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.
4. Very low starting amount and/or low purity of genomic DNA. Improper DNA extraction:
 - Repeat the DNA purification.
5. Wrong channel/filter was chosen. The PCR conditions did not comply with the instructions:
 - Check the PCR conditions and select the fluorescence channels reported in the protocol for the Unknown Sample detection.

Problem 2: Weak or no signal of the Wild Type Control C1 and/or Mutated Control C2.

1. The PCR conditions did not comply with the instructions:
 - Check the amplification protocol and select the fluorescence channel reported in the manual.
2. Deterioration of dyes and/or primers. The reagents storage conditions did not comply with the instructions:
 - Check storage conditions.

Problem 3: Any signal with Reaction Blank BM.

1. Contamination during PCR preparation procedure. All samples results are INVALID:
 - Decontaminate the working area and all instruments.
 - Pipette the controls C1 and C2 at last.
 - Repeat the PCR preparation using a new set of reagents.

Problem 4: HEX signal in Wild Type Control C1 and/or FAM signal in Mutated Control C2.

1. Contamination during PCR preparation procedure. All samples results are INVALID:
 - Decontaminate the working area and all instruments.
 - Use tips with aerosol barriers only. Change tips between tubes/wells.
 - Repeat the PCR preparation using a new set of reagents.

Problem 5: Fluorescence intensity varies.

1. The PCR Master Mix is not well prepared:
 - Carefully repeat the PCR preparation procedure.
2. Air bubbles trapped in the PCR tubes:
 - Check the presence of air bubbles before starting a new run.

Problem 6: Absence of any fluorescent signal.

1. Verify the performance of the thermal cycler:
 - Calibrate the equipment.
2. Deterioration of dyes and/or primers. The storage conditions did not comply with the instructions:
 - Check storage conditions.
 - Check the expiry date of the kit.

Problem 7: The thermal cycler gives an error message.

1. Refer to the Real Time PCR instrument user manual or contact the local technical support of the Real Time PCR instrument company.

Problem 8: The kit reagents left out of the storage range temperature.

1. These reagents must be stored as indicated for a proper execution of the test. The performance of the product is not guaranteed if the reagents have not been properly stored.

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CLONIT S.r.l.