



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# CRP Control HumaFIA SR

Package Size 2 x 0.1 mL

[REF] 16090/21

[IVD]

## Intended Purpose

CRP Control HumaFIA SR is a control to check the accuracy and reliability of the CRP HumaFIA SR assay.

The control is designed to be used together with the CRP HumaFIA SR reagent on semi-automated HumaFIA analyzer.

For professional in vitro diagnostic use only.

## Reagents and Contents

[CTR]1	0.1 mL	Control Level 1 (green), liquid, ready-for-use
[CTR]2	0.1 mL	Control Level 2 (brown), liquid, ready-for-use

The [CTR] concentrations levels are lot specific, see target value sheet.

## Applicable system

[REF]	
16090	HumaFIA
16090/20	CRP HumaFIA SR

## Storage/Stability

Unopened at 2...8°C	Up to the stated expiration date
Opened vial at 2...8°C	5 days open vial. Do not freeze.

Store upright and protected from light.

## Procedure

Follow the procedure exactly as described.

### Procedural Notes

1. Do not use [CTR] after the expiration date.
2. Bring [CTR] to room temperature before use, and keep it at 18-25°C for 20 minutes.
3. Thoroughly mix the contents of the vial before each use by gently inverting for at least 10 times avoiding the formation of foam.
4. Treat the control in the same manner as patient specimen in the assay procedure. **Test should be completed at 22-24°C.**
5. Read the control results, refer to the target value sheet.
6. Close the cap and store the remaining volume at 2-8°C immediately after each use.

### Notes

1. This product is not intended to be used as standard material.

### Assay Procedure

Refer to the instruction for the use of CRP HumaFIA SR

[REF] 16090/20.

## Safety Notes

All patient specimens, calibrators and controls should be handled as potentially infectious. All materials of animal origin avoid many risks associated with the use of human serum (e.g. Hepatitis B and C, HIV). Nevertheless, all material of human or animal origin should still be treated as potentially infectious material.

For users in the European Union only: Please report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

P234 Keep only in original container.

P260 Do not breathe dust/fume/gas/mist/vapors/spray.

P262 Do not get in eyes, on skin, or on clothing.

P281 Use personal protective equipment as required.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

P401 Store in accordance with local/regional/national/international regulations.

P501 Dispose of contents/container in accordance with local/regional/ national/international regulations.



Beijing HumaDX Tech Co., Ltd.

Address: Room 1309, 3rd Floor, 1st Block, No. 55 Jiachuang 2nd

Road, BDA 101111, Beijing, China

E-mail: service@humadx.com.cn Tel: +86-10-80828658



HUMAN Gesellschaft für Biochemica und Diagnostica mbH

Address: Max-Planck-Ring 21, 65205, Wiesbaden, Germany

Tel: +49-6122-99880, E-Mail: human@human.de



Effective Date: Feb 5, 2024  
Version: HD/CE 21 202402 V001