SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

ImmunoCard[®] *H. pylori*

A Rapid Enzyme Immunoassay for the Detection of IgG to *Helicobacter pylori*
in Human Serum, Plasma and Whole Blood

REF Catalog No. 710030

IVD For *in vitro* diagnostic use only

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
Meridian
Bioscience, Inc.

ImmunoCard® *H. pylori*

A Rapid Enzyme Immunoassay for the Detection of IgG to *Helicobacter pylori*
in Human Serum, Plasma and Whole Blood

REF Catalog No. 710030

IVD For *in vitro* diagnostic use only

INTENDED USE

The ImmunoCard® *H. pylori* enzyme immunoassay (EIA) is an *in vitro* qualitative procedure for the detection of IgG to *Helicobacter pylori* in human serum, plasma, or whole blood. Test results are intended to aid in the diagnosis of *H. pylori* infection. The ImmunoCard® *H. pylori* assay may be performed in clinical laboratories and physician's offices.

EXPLANATION

The importance of *Helicobacter pylori* in gastrointestinal diseases has increased greatly since Marshall and Warren described the presence of campylobacter-like organisms in the antral mucosa of patients with histological evidence of antrum gastritis and peptic ulcers especially duodenal ulcers (15,3). Though a causative relationship has not yet been fully established, the strong correlation between the presence of *H. pylori* and histologically confirmed gastritis, peptic ulcer disease and gastric carcinoma can no longer be debated (3,11,20).

The ecological niche in humans appears to be restricted to the stomach and duodenum. Patients who harbor the organism are divided into two basic groups. The first are those who are said to be "colonized". These patients have the organism, yet have no signs of gastrointestinal disease. Those with symptoms and *H. pylori* are considered to be infected. The process by which a colonized individual becomes infected remains unclear. The process by which patients become colonized is also still under investigation (3,7,11,13,14).

The diagnostics strategies for the determination of *H. pylori* have been developed along two lines. The first involves direct detection of the organism and, to date, is invasive. The second, a noninvasive approach, involves the detection of antibodies made against *H. pylori* or the detection of the organism's urease (8,9,22).

The invasive route requires that a biopsy be taken from the upper gastrointestinal tract. The presence of *H. pylori* is then confirmed by direct microscopic examination or culturing of the organism from the biopsy material. This strategy has the advantage of being able to detect active infections while being highly specific with a very high positive predictive value. The difficulties associated with this approach is that there is a risk and discomfort to the patient. In addition, *H. pylori* tends to colonize in patches and may be missed totally by the biopsy. Culturing of *H. pylori* from biopsy material tends to be difficult and time consuming. These technical difficulties may lead to false negative results (1,2,3,19).

The urea breath test is a noninvasive method that detects *H. pylori* by exploiting its highly active urease (3,8,10).

The most common noninvasive approach to the detection of *H. pylori* is the serological identification of specific antibodies in infected patients. Through an indirect approach, the correlation between histological gastritis, the presence of *H. pylori*, and seropositivity is extremely strong. The disadvantages of these tests are that most require instruction to be interpreted and many have a lower specificity due to cross reactions from other organisms. Their major advantage is that they are rapid and offer a high sensitivity (5,18,21).

ImmunoCard *H. pylori* is an enzyme immunoassay in a self contained device which detects IgG antibodies to *H. pylori*. Results are interpreted by observing for a visual color change within seven minutes after addition of specimen.

BIOLOGICAL PRINCIPLES

The test detects IgG antibodies to *H. pylori* through the use of an enzyme immunoassay. Each test card consists of a plastic housing and a fiber membrane which provides the solid support for the EIA. There are two large sample ports and two small reaction ports. The left reaction port serves as a CONTROL and has human IgG antibodies immobilized on the fiber membrane. The right reaction port is the TEST port and contains immobilized *H. pylori* antigens.

Patient serum, plasma, or whole blood is added to each of the two sample ports. After allowing the sample to enter the device and migrate along the membrane to the reaction ports, anti-human IgG alkaline phosphatase conjugate is added to the reaction ports and allowed to enter the device. Wash Buffer and Substrate are then added to each of the reaction ports. Reaction ports are observed for the development of blue color.

The development of blue color in the TEST reaction port indicates a positive test result. The absence of blue color in the TEST reaction port denotes a negative test result. The CONTROL reaction port serves as the procedural control. Color development confirms that the reagents are performing and that the procedure was done correctly.

MATERIALS PROVIDED

1. Test Cards (30) — Individually foil pouched cards containing immobilized *H. pylori* antigens (TEST reaction port) and human IgG (CONTROL reaction port).
2. Positive Control (2.0ml) — Human Serum containing IgG antibodies to *H. pylori* in a buffer (pH 7.3) with 0.1% sodium azide.
3. Negative Control (2.0ml) — Nonreactive Human Serum in a buffer (pH 7.3) with 0.1% sodium azide.
4. Enzyme Conjugate (5.0ml) — Monoclonal anti-human IgG labeled with alkaline phosphatase in a buffer (pH 7.25) containing 0.1% sodium azide.
5. Wash Buffer (15ml) — Buffer (pH 7.2) containing 9.5% (weight/vol.) guanidine hydrochloride.
6. Substrate Reagent (15ml) — Buffer (pH 9.8) containing 5-bromo-4-chloro-3-indolyl phosphate and 0.1% sodium azide.
7. Whole Blood Sample Diluent (7.5ml) — Buffer (pH 7.0) containing sodium dodecyl sulfate and 0.1% sodium azide.
8. Transfer Pipettes (60).

MATERIALS NOT PROVIDED

1. Timer.
2. Test tubes (for whole blood only).
3. Capillary tubes (for whole blood only).

PRECAUTIONS

1. All reagents are for *in vitro* diagnostic use only.
2. Reagent concentrations, incubation times and temperature (22°-27°C) have been optimized for sensitivity and specificity. Best results are obtained by adhering to these specifications. Once the assay has been started, complete all subsequent steps without interruption.

- Patient specimens and controls may contain infectious agents. The positive and negative controls contain human sera which were screened for HBsAg and antibody to HIV and found to be negative. However, no test can offer complete assurance that human blood will not transmit HIV, hepatitis or other infectious agents. Therefore, patient specimens, controls, used Test Cards and other potentially hazardous material should be handled at Biosafety Level 2 as recommended in the CDC/NIH manual, "Biosafety in Microbiology and Biomedical Laboratories"; 1988.
 - All reagents and patient specimens should be at room temperature (22°-27°C) before use. Additionally, all reagents should be removed from the foam insert and mixed prior to use.**
 - Do not interchange reagents from different kit lot numbers or use expired reagents.
 - Hold reagent vials vertically to insure proper drop size and delivery. Do not allow the tips of bottles or pipets to touch either the sample or reaction ports.
 - Replace caps on correct vials.
 - The Substrate Reagent may be light sensitive and should not be exposed to excessive illumination. If this reagent exhibits a blue color, it should be discarded.
 - Use only one transfer pipet per specimen. Discard after use. Do not attempt to reuse.
 - Disregard any color reactions in the sample (lower) ports. Results are determined by color development in the reaction (upper) ports only.
- WARNING:** Some reagents in this kit contain sodium azide, which is a skin irritant. Avoid skin contact with the kit reagents. Disposal of reagents containing sodium azide into lead or copper plumbing can result in the formation of explosive metal azides. This can be avoided by flushing with a large volume of water during such disposal.

SHELF LIFE AND STORAGE


The expiration date is indicated on the kit label. Store at 2°-8°C and return the kit promptly to the refrigerator after each use.

SPECIMEN HANDLING

- Serum, plasma, or whole blood can be used in this assay. Sodium heparin, lithium heparin, sodium citrate and EDTA are acceptable anticoagulants. The performance of the ImmunoCard® *H. pylori* assay with plasma obtained using other anticoagulants is unknown. Specimens must be collected using established laboratory procedures.
- Serum and plasma can be stored in the refrigerator at 2°-8°C for up to five days. If testing is to be performed at a later date, the specimen should be frozen. Avoid repeated freezing and thawing. A cloudy serum/plasma specimen, or one containing particulate matter, may be centrifuged prior to testing.
- Whole blood with anticoagulant can be used without centrifugation. With lithium heparin, EDTA, or sodium citrate, whole blood specimens can be stored at 2°-8°C for up to 72 hours.
- For fingerstick whole blood, completely fill capillary tube and test immediately.

PROCEDURAL NOTES

- The ImmunoCard® Test Card format is diagrammed below:

- Batch processing of samples or controls is possible provided all procedural steps and incubation times are maintained for each Test Card.
- The **CONTROL** (left) side of each card provides a procedural control which detects problematic specimens and tests the reagent performance.

PROCEDURE

This test should be performed by qualified personnel per local regulatory requirements.

PLASMA, SERUM AND CONTROLS

- Remove the appropriate number of Test Cards from their envelopes. Label with appropriate identification. Use one card for each control or specimen to be tested. Test Cards must be maintained in a horizontal position during the performance of the assay procedure.
- Using a transfer pipet, dispense two drops of serum or plasma to both lower sample ports.
Optional Controls: The Positive Control is run by directly adding two drops to both lower sample ports of a card. The negative Control can be run in the same manner on a separate card.
- Wait one minute after the sample or control has been added to both lower ports.
- Add two drops of Enzyme Conjugate to both upper ports.
- As soon as the Enzyme Conjugate is fully absorbed, add two drops of Wash Buffer to both upper ports.
- As soon as the Wash Buffer is fully absorbed, add two drops of Substrate Solution to both upper ports.
- Wait five minutes.**
- Without delay, visually read upper wells only for the presence or absence of blue color in the same order that the Substrate Solution was added.**

WHOLE BLOOD AND CONTROLS

- Remove the appropriate number of Test Cards from their envelopes. Label with appropriate identification. Use one card for each specimen to be tested. Test Cards must be maintained in a horizontal position during the performance of the assay procedure.
- Add four drops of Whole Blood Sample Diluent to labeled glass test tube.
- Using a transfer pipet, add two drops of whole blood to the test tube.
NOTE: For fingerstick procedure, completely fill a capillary tube with whole blood. Add the entire contents to a test tube that contains four drops of Whole Blood Sample Diluent.
Optional Controls: The Positive and Negative Controls are tested by adding two drops to the test tube which contains four drops of Whole Blood Sample Diluent.
- Using a **new** transfer pipet, dispense two drops of the diluted whole blood or controls to both lower sample ports.
- Wait three minutes after the sample or control has been added to both lower ports.
CAUTION: Some whole blood specimens may clog the membrane. If the sample does not completely wet the test and control reaction ports, the test is invalid.
- Add two drops of Enzyme Conjugate to both upper ports.
- As soon as the Enzyme Conjugate is fully absorbed, add two drops of Wash Buffer to both upper ports.
- As soon as the Wash Buffer is fully absorbed, add two drops of Substrate Solution to both upper ports.
- Wait five minutes.
- Without delay, visually read upper wells only for the presence or absence of blue color in the same order that the Substrate Solution was added.

QUALITY CONTROL

The Positive and Negative Controls should be assayed upon receipt of the kit and must yield the following results:

1. Positive Control: Visually detectable blue color in both reaction (upper) ports (Figure 1).
 2. Negative Control: Visually detectable blue color in CONTROL (upper left) reaction port only. The TEST (upper right) reaction port should be colorless (Figure 2).
- The Procedural Control present in the upper left port of each Test Card tests the individual specimen and reagent performance. Failure of the Procedural Control to yield a blue color with any specimen or control reagent indicates an invalid test and the test should be repeated.

At the time of each use, kit components should be visually examined for obvious signs of microbial contamination, freezing, or leakage. It is suggested that the results of each quality control check be recorded in an appropriate log book to maintain high quality testing records.

The Positive Control stresses the test at the assay cut-off.

The Procedural Control present in the upper left portion of each Test Card verifies the reactivity of each reagent.

If the Procedural, Positive or Negative Control does not generate the appropriate and expected results, the test is invalid and the patient results should not be reported. If the expected reactions are not observed and the reagents are still within their expiration date, please contact Meridian Bioscience Technical Support Center at 1-800-343-3858.

INTERPRETATION OF RESULTS


Fig. 1 Reactive Result Fig. 2 Nonreactive Result Fig. 3 Invalid Result Fig. 4 Invalid Result

READ ONLY UPPER REACTION PORTS TO INTERPRET RESULTS

Positive Test Result: Visually detectable blue color in **BOTH** reaction ports. Occasionally a positive test reaction may show evidence of a gradation of blue color from the bottom of the port to the top. This should be considered a positive result. A positive test result indicates the presence of IgG to *H. pylori*.

Negative Test Result: Visually detectable blue color in **CONTROL** reaction port (upper left) only. **The TEST reaction port (upper right) should be colorless to faint grey.** Occasionally, the **TEST** reaction port (upper right) may show evidence of a hint of blue color in the right or left side of the port, with the rest of the port remaining colorless. This should be considered a negative test result. A negative test result indicates that IgG to *H. pylori* is absent or below the level of detection of the assay.

Invalid Test Result: No detectable color in **CONTROL** reaction port. Invalid test results may be due to an **ImmunoCard**® reagent problem, a procedural error, or clogging of the membrane by specimen: repeat test.

LIMITATIONS OF THE PROCEDURE

1. **ImmunoCard**® *H. pylori* test results should be used in conjunction with information from the patient clinical evaluation and other diagnostic procedures.
2. This test should only be used on patients who present symptoms suggestive of gastrointestinal disease. **ImmunoCard**® *H. pylori* is not intended for use with asymptomatic patients.
3. A Positive test result does not distinguish between an active infection and colonization by *H. pylori*.
4. **ImmunoCard**® *H. pylori* is a qualitative test that only indicates the presence of IgG antibody to *H. pylori*. This test does not indicate that gastrointestinal disease is present.
5. A Negative test result indicates the absence of IgG antibody to *H. pylori* or that the IgG level is below that which can be detected by this test.
6. **ImmunoCard**® *H. pylori* has not been tested in pediatric populations. Therefore, no performance characteristics have been established for this group.
7. **ImmunoCard**® *H. pylori* was evaluated for possible interferences from high levels of hemoglobin, lipid, and bilirubin. No interference was observed in serum specimens that were either positive or negative for antibodies to *H. pylori*. Whole blood specimen results were not evaluated for interferences by these compounds.
8. Whole blood specimens with hematocrit levels >61 volume % have not been evaluated and may clog the membrane.

EXPECTED VALUES

Studies on the epidemiology of *H. pylori* have shown that this organism is present worldwide (12,16). Gastritis caused by *H. pylori* has been shown to correlate with age, ethnic background, family size and socioeconomic class (6,10). While the epidemiological studies on *H. pylori* infection or colonization are ever expanding, the route by which individuals acquire this organism is, at present, unknown (3,7,11,13,14).

Studies indicate that the incidence of infection in the United States may increase 1-2% annually (10). It is not uncommon to see positivity rates of 50% in those who are 60 years or older. The frequency of *H. pylori* infection in patients diagnosed with duodenal ulcers, however, has been shown for patients in every age group to be approximately 80% (12). The positivity rates seen with *H. pylori* clearly depend on the sample population and how that population is defined.

The **ImmunoCard**® *H. pylori* test detects the presence of IgG antibodies to *H. pylori* in serum, plasma, and whole blood. Sera from 43.2% of a general symptomatic population were antibody positive by **ImmunoCard**® *H. pylori*. There was a 99.1% correlation of whole blood and serum results in a general population of symptomatic and asymptomatic individuals. Expected values for a given population should be determined for each laboratory. The rate of positivity may vary depending on geographic location, method of specimen collection, handling and transportation, test employed, and general health environment of the patient population under study.

PERFORMANCE CHARACTERISTICS

The **ImmunoCard**® *H. pylori* test was evaluated at two hospitals and one gastroenterology laboratory (4). A total of 382 serum specimens were obtained from patients with gastrointestinal complaints. All patients underwent endoscopy. The **ImmunoCard**® *H. pylori* test was compared to results obtained from the biopsied material using rapid biopsy urease tests.

ImmunoCard® <i>H. pylori</i>	Biopsy Urease Test	
	Positive	Negative
Positive	115	50
Negative	2	215

n= 382
Sensitivity = 98.2% (94.0-99.5)*
Specificity = 81.1% (76.0-85.4)
Agreement = 86.4% (82.5-89.5)

95% confidence intervals calculated by the normal method are indicated in parenthesis

Forty-seven biopsy urease test negative/ImmunoCard® *H. pylori* positive specimens were positive by a competitor's microtiter well EIA. Forty-six of these specimens were positive by histological staining and/or the pathology diagnosis was associated with *H. pylori* disease. The additional information indicates that 47 of the 50 relative false positive results were, in fact, true antibody positives that were not positive by biopsy urease test. The two urease test positive/ImmunoCard® *H. pylori* negative specimens were antibody negative by a competitor's microtiter well EIA, indicating that these relative false negative results were true antibody negative results.

ImmunoCard® *H. pylori* was compared to a competitor's microwell based EIA also for the detection of IgG to *H. pylori*.

ImmunoCard® <i>H. pylori</i>	Competitor EIA		
	Positive	Negative	I
Positive	161	4	0
Negative	5	211	1

n= 382
Relative Sensitivity = 97.0% (93.1-98.7)*
Relative Specificity = 98.1% (95.3-99.3)
Relative Agreement = 97.6% (95.6-98.7)

One EIA negative/ImmunoCard® *H. pylori* positive specimen was positive by the biopsy urease test. Three EIA negative/ImmunoCard® *H. pylori* positive specimens were negative by the biopsy urease test. All five ImmunoCard® *H. pylori* negative/EIA positive specimens were negative by the biopsy urease test. The one EIA indeterminate/ImmunoCard® *H. pylori* negative specimen was biopsy urease test negative.

A specimen matrices study evaluated whole blood collected in lithium heparin, EDTA and sodium citrate and serum collected from 108 symptomatic and asymptomatic individuals. The relative comparison of whole blood to serum was not traceable to biopsy results. There was a 99.1% (95% CI = ±1.8%) agreement between whole blood results, regardless of anticoagulant (30 positives) and serum results (31 positives).

A comparison of the results of capillary whole blood obtained by finger stick and serum from 77 individuals indicated 100% agreement; 15 positives and 62 negatives. The specimen results were not traceable to biopsy results.

REPRODUCIBILITY

Two clinical test sites ran a negative, low positive and high positive serum specimen ten times a day each of three days. The Positive and Negative Controls were run three times a day over this testing period. Three clinical test sites ran the Positive and Negative Controls and four negative and four positive serum specimens three times a day for three days. A total of 450 assays were performed at all the testing sites. The inter and intra-assay reproducibility of ImmunoCard® *H. pylori* with all specimens was 100% (4).

ASSAY SPECIFICITY

A. Bacterial Studies

Adsorption studies demonstrated that ImmunoCard® *H. pylori* does not suffer cross reactivity problems with the following clinical isolates (CI) or ATCC strains (4).

<i>Campylobacter coli</i>	CI	<i>Klebsiella pneumoniae</i>	ATCC 13883
<i>Campylobacter fetus</i>	CI	<i>Proteus vulgaris</i>	ATCC 8427
<i>Campylobacter hyointestinalis</i>	ATCC 29543	<i>Pseudomonas aeruginosa</i>	ATCC 9027
<i>Campylobacter jejuni</i>	ATCC 33292	<i>Pseudomonas fluorescence</i>	CI
<i>Campylobacter lari</i>	CI	<i>Salmonella sp</i>	CI
<i>Campylobacter rectus</i>	ATCC 35217	<i>Serratia liquefaciens</i>	ATCC 3551
<i>Candida albicans</i>	ATCC 10231	<i>Shigella flexneri</i>	CI
<i>Citrobacter freundii</i>	ATCC 8090	<i>Shigella sonnei</i>	CI
<i>Enterobacter aerogenes</i>	ATCC 13048	<i>Staphylococcus aureus</i> (Cowan I)	ATCC 12598
<i>Enterobacter cloacae</i>	ATCC 13407	<i>Streptococcus faecalis</i>	CI
<i>Escherichia coli</i>	ATCC 8739	<i>Wolinella succinogenes</i>	ATCC 29543
<i>Helicobacter cinaedi</i>	ATCC 35683	<i>Yersinia enterocolitica</i>	CI

B. Serum Studies

Sera from patients who were culture positive or had serologic evidence of the following diseases did not react when tested in ImmunoCard® *H. pylori* (4).

Anti-Nuclear Antibody	Cytomegalovirus	Rheumatoid Factor	Rubella
Lupus Erythematosus	Toxoplasma	Influenza A	Coccidioides
Histoplasma	Mononucleosis	Herpes Simplex Virus	Cryptococcus

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REFERENCES

1. Alpert, L.C., D.Y. Graham, D.J. Evans, Jr., H.H. Yoshimura, Sa.L. Hazell, D.G. Evans, and P.D. Klein. 1989. Diagnostic possibilities for *Campylobacter pylori* infection. *Eur. J. Gastroenterol. Hepatol.* 1:17-26.
2. Bartel, J.S., and E.D. Everett. 1990. Diagnosis of *Campylobacter pylori* infections: The "gold standard" and the alternatives. *Rev. Infect. Dis.* 12:S107-S114.
3. Buck, G.E. 1990. *Campylobacter pylori* and gastroduodenal disease. *Clin. Microbiol. Rev.* 3:1-12.
4. Data on file, Meridian Bioscience, Inc.
5. Evans, D.J. Jr., D.G. Evans, D.Y. Graham, and P.D. Klein. 1989. A sensitive and specific serologic test for detection of *Campylobacter pylori* infection. *Gastroenterology*. 96:1004-1008.
6. Fiedorek, S.C., H.M. Malaty, D.L. Evans, C.L. Pumphrey, H.B. Casteel, D.J. Evans, Jr., and D.Y. Graham. 1991. Factors influencing the epidemiology of *Helicobacter pylori* infection in children. *Pediatrics*. 88:578-582.
7. Graham, D.Y. 1991. I. *Helicobacter pylori*: Its epidemiology and its role in duodenal ulcer disease. *J. Gastroenterol. Hepatol.* 6:105-113.
8. Graham, D.Y., P.D. Klein, D.J. Evans, Jr., D.G. Evans, L.C. Alpert, A.R. Opekun, and T.W. Boutton. 1987. *Campylobacter pylori* detected noninvasively by the ¹³C-urea breath test. *Lancet* i: 1174-1177.
9. Graham, D.Y., P.D. Klein, A.R. Opekun, and T.W. Boutton. 1988. Effect of age on the frequency of active *Campylobacter pylori* infection diagnosed by the [¹³C] Urea breath test in normal subjects and patients with peptic ulcer disease. *J. Infect. Dis.* 157:777-780.
10. Graham, D.Y., H.M. Malaty, D.G. Evans, D.J. Evans, Jr., P.D. Klein, and E. Adam. 1991. Epidemiology of *Helicobacter pylori* in an asymptomatic population in the United States. Effect of age, race, and socioeconomic status. *Gastroenterology*. 100:1495-1501.
11. Gregson, D.B. and A.E. Simor. 1991. *Helicobacter pylori* and inflammatory disease of the stomach and duodenum. *Stomach and Bowel*, Sept. 1993, pp.29-30,35-37.
12. Loffeld, R.J.L.F., E. Stobberingh, J.P. Van Spreuwel, J.A. Fiendrig, and J.W. Arends. 1991. The prevalence of anti-*Helicobacter (Campylobacter) pylori* antibodies in patients and healthy blood donors. *J. Med. Microbiol.* 32:105-109.
13. Malaty, H.M., D.I. Graham, P.D. Klein, D.G. Evans, E. Adam, and D.J. Evans. 1991. Transmission of *Helicobacter pylori* infection. Studies in families of healthy individuals. *Scand. J. Gastroenterol.* 26:927-932.
14. Marshall, B.J., J.A. Armstrong, D.B. McGeachie and R.J. Glancy. 1985. Attempt to fulfill Koch's postulates for pyloric campylobacter. *Med. J. Aust.* 142:436-439.
15. Marshall, B.J. and J.R. Warren. 1984. Unidentified curved bacilli in the stomach of patients with gastritis and peptic ulceration. *Lancet* i:1311-1314.
16. Morris, A., G. Nicholson, G. Lloyd, D. Haines, A. Rogers, and D. Taylor. 1986. Seroprevalence of *Campylobacter pyloridis*. *N. Z. Med. J.* 99:657-659.
17. Murray, D.M. 1993. Clinical relevance of infection by *Helicobacter pylori*. *Clin. Microbiol. Lett.* 15:33-37.
18. Newell, D.G. 1987. Identification of the outer membrane proteins of *Campylobacter pyloridis* and antigenic cross-reactivity between *C. pyloridis* and *C. jejuni*. *J. Gen. Microbiol.* 133:163-170.
19. Nichols, L., M. Sughayer, P.C. DeGirolami, K. Balogh, D. Pleskow, K. Eichelberger, and M. Santos. 1991. Evaluation of diagnostic methods for *Helicobacter pylori* gastritis. *Am. J. Clin. Path.* 95:769-773.
20. Parsonnet, J., G.D. Friedman, D.P. Vandersteen, Y. Chang, J.H. Vogelman, N. Orentreich, and R.K. Sibley. 1991. *Helicobacter pylori* infection and the risk of gastric carcinoma. *N. Engl. J. Med.* 325:1127-1131.
21. Tailey, N.J., D.G. Newell, J.E. Ormand, H.A. Carpenter, W.R. Wilson, A.R. Zinsmeister, G.I. Perez-Perez, and M.J. Blaser. 1991. Serodiagnosis of *Helicobacter pylori*. Comparison of enzyme-linked immunosorbent assays. *J. Clin. Microbiol.* 29:1635-1639.
22. Westblom, T.U., E. Madan, J. Kemp, and M.A. Subik. 1988. Evaluation of a rapid urease test to detect *Campylobacter pylori* infection. *J. Clin. Microbiol.* 26:1393-1394.

USA/Corporate Office
 Manufactured by 3471 River Hills Drive
 Cincinnati, Ohio 45244
 Telephone: (513) 271-3700
 Orders/Customer Service: 1-800-543-1980
 Technical Support Center: 1-800-343-3858
 Information Fax: (513) 272-5432
 Ordering Fax: (513) 271-0124

INTERNATIONAL SYMBOL USAGE

You may see one or more of these symbols on the labeling/packaging of this product

	Manufactured by		For Performance Evaluation Only
	Authorized Representative		Temperature Limitations
	Catalog Number		Use By/Expiration Dating Information
	In vitro Diagnostics		Sufficient for "X" Tests
	Batch Code/Lot Information		CE Symbol

Authorized Representative
Meridian Bioscience Europe
 Via dell'Industria, 7
 20020 Villa Cortese (MI)
 Italy
 Tel.: +39 0331 433636
 Fax: +39 0331 433616
 e-mail: info@mdeur.com

Meridian Bioscience Europe s.a./n.v.
 Rue de l'Industrie 7
 1400 Nivelles
 Belgium
 Tel.: +32 (67) 895959
 Fax: +32 (67) 895958
 e-mail: info@mdeur.be

Meridian Bioscience Europe France
 455, Promenade des Anglais
 06299 Nice Cedex-3
 France
 Tel.: +33 (0) 493187210
 Fax: +33 (0) 493187211
 e-mail: mdeur.info@fnac.net

Meridian Bioscience Europe b.v.
 Halderheiweg 6
 5282 SN Boxtel
 Holland
 Tel.: +31 (411) 62 11 66
 Fax: +31 (411) 62 48 41
 e-mail: meridian@wxs.nl

OTHER MERIDIAN PRODUCTS

MERIDIAN PRODUCT

Premier™ Platinum HpSA
 Premier™ *H. pylori*
 ImmunoCard® *C. difficile*
 ImmunoCard® *Mycoplasma*

CATALOG NUMBER

601348
 606096
 706050
 709030

For technical assistance, please call our Technical Support Center at 1-800-343-3858 between the hours of 8 AM and 6 PM, Eastern Standard Time. To place an order, please call our Customer Service Department at 1-800-543-1980.

SN 11233

Printed in the U.S.A. on recycled paper
 Copyright © Meridian Bioscience, Inc.

Rev. 4/03