SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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Vivalytic MRSA/SA – Annex

Table 1 – Limit of Detection (per Test)

MRSA*	300 colony forming units
MSSA**	90 colony forming units

* tested with MRSA strain ATCC 1683 (SCCmec type IV a)

** tested with MSSA strain DSM 2569

Strains in Copan Liquid Amies Medium, supplemented with 10,000 Human Bronchial/Tracheal Epithelial Cells (HBTEC)

Table 2 – Inclusivity

MRSA-strains*	MRSA-SCCmec types
MRSA-strain BAA-41	Type II
MRSA-strain BAA-42	Type VI
MRSA-strain BAA-44	Type I
MRSA-strain BAA-2094***	Type V
MRSA-strain BAA-2313	Type XI
MRSA-strain ATCC 33592	Type III
MRSA-strain ATCC 43300***	Type II
MRSA-strain ATCC 1683	Type IV a
MSSA-strains**	
MSSA isolate 129	
MSSA strain DSM 1104	
MSSA strain DSM 2569	

MRSA strains were tested at

*300 CFU/test in Copan Liquid Amies Medium, supplemented with 10,000 Human Bronchial/Tracheal Epithelial Cells (HBTEC) and 33.33 % (v/v) respiratory swab matrix negative control (MC110). All strains were detected except for two. These were re-tested with

***900 CFU/test, verifying inclusivity for MRSA ATCC 43300 whereas a slightly decreased detection rate for MRSA strain BAA-2094 was detected.

MSSA strains were tested at

**75 CFU/test in Copan Liquid Amies Medium, supplemented with 10,000 Human Bronchial/Tracheal Epithelial Cells (HBTEC) and 33.33 % (v/v) respiratory swab matrix negative control (MC110). All strains were detected.

Table 3a– Exclusivity: Strains investigated by *in-vitro* analysis*

Staphylococcus hominis Z021

Staphylococcus epidermidis HER1292

Staphylococcus capitis isolate 48189

Staphylococcus epidermidis RV1625391

Staphylococcus warneri NTC 5955

Staphylococcus capitis isolate 644

Staphylococcus haemolyticus ATCC 29970

* tested at 10⁴–10⁵ CFUs/ml in Copan Liquid Amies Medium, supplemented with 10,000 Human Bronchial/Tracheal Epithelial Cells (HBTEC) and 33.33 % (v/v) respiratory swab matrix negative control (MC110)

Tested strains indicated no cross reactivity

Table 3b – Exclusivity: Strains investigated by *in-silico* analysis**Bacteria**

<i>Acinetobacter baumannii</i>	<i>Lactobacillus crispatus</i>
<i>Acinetobacter lwoffii</i>	<i>Legionella pneumophila</i>
<i>Aeromonas hydrophila</i>	<i>Listeria monocytogenes</i>
<i>Aspergillus fumigatus</i>	<i>Micrococcus luteus</i>
<i>Aspergillus niger</i>	<i>Moraxella catarrhalis</i>
<i>Bacillus cereus</i>	<i>Neisseria meningitidis</i>
<i>Bacteroides fragilis</i>	<i>Proteus vulgaris</i>
<i>Bordetella pertussis</i>	<i>Pseudomonas aeruginosa</i>
<i>Campylobacter coli</i>	<i>Salmonella enteritidis</i>
<i>Campylobacter jejuni</i>	(<i>Salmonella enterica</i> subsp. <i>enterica</i> – Serovar Enteritidis)
<i>Candida albicans</i>	<i>Serratia marcescens</i>
<i>Candida glabrata</i>	<i>Shigella sonnei</i>
<i>Citrobacter freundii</i>	<i>Staphylococcus capitis</i>
<i>Clostridium difficile</i>	<i>Staphylococcus epidermidis</i>
<i>Clostridium perfringens</i>	<i>Staphylococcus epidermidis</i>
<i>Corynebacterium diphtheriae</i>	<i>Staphylococcus haemolyticus</i>
<i>Corynebacterium jeikeium</i>	<i>Staphylococcus hominis</i>
<i>Corynebacterium striatum</i>	<i>Staphylococcus saprophyticus</i>
<i>Cryptococcus neoformans</i>	<i>Staphylococcus warneri</i>
<i>Enterobacter aerogenes</i>	<i>Streptococcus agalactiae</i>
<i>Enterobacter cloacae</i>	<i>Streptococcus mitis</i>
<i>Enterococcus faecalis</i>	<i>Streptococcus mutans</i>
<i>Enterococcus faecium</i>	<i>Streptococcus pneumoniae</i>
<i>Escherichia coli</i> (containing O157:H7)	<i>Streptococcus pyogenes</i>
<i>Haemophilus influenzae</i>	<i>Streptococcus salivarius</i> subsp. <i>thermophilus</i>
<i>Klebsiella oxytoca</i>	<i>Streptococcus uberis</i>
<i>Klebsiella pneumoniae</i>	<i>Vibrio parahaemolyticus</i>
	<i>Yersinia enterocolitica</i>

Viruses

<i>Adenoviridae</i> (containing Adenovirus 7A)	Influenza A
Coronavirus 229E	Influenza B
Cytomegalovirus	Measles morbillivirus
Enterovirus (containing Rhinovirus)	Mumps orthorubulavirus
Epstein Barr virus	Paramyxoviridae (Parainfluenza)
Human Metapneumovirus	Human orthopneumovirus (containing Respiratory Syncytial Virus)
	SARS-CoV-2

Tested strains indicated no cross reactivity.

Table 4 – Precision

c*	Number of total tests	Number of positive tests (MRSA)	Number of positive tests (MRSA/MRSA)	Proportion of positive/negative results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
blank	108	0	107	99.07	94.94 – 99.84	99.95 – 99.98
c95	108	106	0	98.15	93.5 – 99.49	93.47 – 99.77
c100	108	108	0	100	96.57 – 100	96.64 – 100

*blank = negative control; c95 = 95 % predicted positive agreement;
 c100 = 100 % predicted positive agreement

Table 5 – Interference substance	Final concentration of interfering substance
Human K2-EDTA blood	2 % v/v
Salbutamol-ratiopharm® inhalation solution	5 % v/v
Rhinomer nasal spray	10 % v/v
Nasivin® nasal spray	5 % v/v
Neosynephrin-POS® 5 % eyedrops	5 % v/v
Turixin 2 % nasal ointment Mupirocin (20 mg mupirocin per 1 g nasal ointment)	2.5 % w/v
TOBRAZID® (80 mg Tobramycin per 2 mL vial) injection-/infusion-solution	4 µg/mL (4 µg Tobramycin per mL)
Otri-Allergie nasal spray Fluticason	5 % v/v
MometaHEXAL® hayfever spray	10 % v/v
Polihexanide (Prontoderm Gel light)	2 % (w/v)

No interferences were observed at tested concentrations.

Table 6 – Diagnostic Sensitivity (PPA) and Specificity (NPA), 95%-width of CI*	
MRSA	MSSA
PPA 93.5 % (82.1 – 98.63 %)	PPA 85.96 % (74.21 – 93.74 %)
NPA 97.0 % (91.56 – 99.38 %)	NPA 95.12 % (83.47 – 99.4 %)

*Clinical performance data of Vivalytic MRSA/SA in comparison to the reference test system Xpert® SA Nasal Complete G3 (Cepheid) after discrepancy analysis. Retrospective clinical study with in total 146 valid nasal and oropharyngeal swab samples tested at one study site. In case of an insufficient number of positive swab specimens, samples were spiked (cultured bacterial cells) in a concentration of 10^4 – 10^8 CFU/ml for MRSA. In total, 20 samples for MRSA were spiked. The discordant test results between the two test systems, were further analysed with MRSA selective agar.

Table 7 – Document History	
Revision 01	Initial document
Revision 02	Changed test volume from 300µl to 600µl, update of performance data