SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

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Intended Use

The Bosch Vivalytic MG, MH, UP/UU test is an automated qualitative *in vitro* diagnostic test based on real-time polymerase chain reaction (PCR) for the detection of nucleic acids from *Mycoplasma genitalium*, *Mycoplasma hominis* and *Ureaplasma parvum/urealyticum* from human urine samples or human vaginal, cervical, urethral, rectal swab samples to aid in the diagnosis of sexually transmitted infections of symptomatic and asymptomatic individuals.

Results should not be used as the sole basis for diagnosis, treatment or other patient management decisions. Positive results do not exclude co-infection with other pathogens. The agent detected may not be the definite cause of disease. Negative results do not exclude a Mycoplasma- and Ureaplasma related infection or another infection. Results must be clinically correlated with patient history, clinical observations and epidemiological information. Other diagnostic information is necessary to determine patient infection status.

Intended for use with a Vivalytic *one* analyser by healthcare professionals only in both near patient testing such as doctor's offices, emergency departments and nursing homes as well as laboratory settings such as hospital laboratories and reference laboratories.

Safety Information

The present Instructions for Use contains test-specific information only. For additional warnings and instructions refer to the Instructions for Use provided with your Vivalytic *one* analyser (chapter device safety information).

Only use Vivalytic cartridges and accessories approved for the Vivalytic *one* analyser. Take care to avoid any contamination when handling patient samples and cartridges. When sample was spilled on the cartridge, do not use the cartridge and dispose it.

For *in vitro* diagnostic use by trained healthcare professionals.

WARNING

- Always follow good laboratory practice to ensure the proper performance of this test.
- Make sure to wear appropriate personal protective equipment.
- Do not use a cartridge if the sealed pouch or the cartridge itself is visibly damaged.
- Do not touch or scratch the detection area of the cartridge.
- Do not reuse a cartridge.
- Do not use expired cartridges. The expiration date can be found on the packaging and the cartridge label.
- Do not wait longer than 15 minutes after opening the cartridge pouch to begin the test. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.
- Do not shake a cartridge that contains a sample.
- Do not turn the cartridge upside down.
- Place the cartridge on a clean and flat surface only.
- Do not use sample types, media and volumes that are not approved for the test.
- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions. Handle potentially infectious patient samples and cartridges according to national laboratory standards and dispose samples and cartridges according to regional and laboratory standards.
- Be compliant with the national safety regulations and practices.

Note: Further information can be found in the safety data sheet (SDS) of the product. Please contact the customer support of your local distributor.

Vivalytic MG, MH, UP/UU – Instructions for Use

Additional Equipment & Consumables Required but not Provided

- Bosch Vivalytic *one* analyser (reference number F 09G 300 115)
- Pipettor (100–1000 µl)
- Sterile pipette tips 100 –1000 µl
- Swab collection kits
 - Regular Flocked Swab FLOQSwabs® 552C (COPAN Italia s.p.a.)
 - Minitip Flocked Swab FLOQSwabs® 551C (COPAN Italia s.p.a.)
- Transport medium
 - eNAT® 608C (COPAN Italia s.p.a.)
- Suitable protective clothing

Test Principle

Vivalytic MG, MH, UP/UU is a qualitative real-time PCR based test.

Storage und Usage Conditions

Product is stable until the expiry date if stored at +15 °C to +25 °C. Storage and usage conditions can be taken either from the cartridge, pouch, or box label. Cartridge has to be used at +15 °C to +25 °C, relative humidity < 65 %, and within 15 minutes upon pouch opening. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.

Reagents

All reagents necessary for the sample processing are integrated into the cartridge. The processing includes cell lysis, nucleic acid extraction, DNA amplification, and detection.

Reagents within the cartridge include PCR bead, binding buffer, washing buffer, and elution buffer. The PCR bead contains the DNA polymerase, primers and probes. Binding buffer facilitates binding of nucleic acids during the purification process. Washing buffer is a formulation of different salts and solvents to remove impurities e. g. proteins during the extraction process. Elution buffer is a low-salt buffer and contains the purified nucleic acids at the end of the extraction process.

Sample Types / Medium

The test is intended for use with urine in eNAT® transport medium (eNAT® 608C, COPAN Italia s.p.a.) or vaginal, cervical, and rectal swab samples (Regular Flocked Swab FLOQSwabs® 552C, COPAN Italia s.p.a.) or urethral swab samples (Minitip Flocked Swab FLOQSwabs® 551C, COPAN Italia s.p.a.) in 1 ml eNAT® transport medium (eNAT® 608C, COPAN Italia s.p.a.). Collect and store samples as indicated in the manufacturer's instructions.

Sample Preparation

Vaginal, Cervical, Rectal and Urethral Swabs

Use a COPAN Regular Flocked Swab (FLOQSwabs® 552C) for collection of vaginal, cervical, rectal patient samples or a COPAN Minitip Flocked Swab (FLOQSwabs® 551C) for collection of urethral patient samples. Each swab is transferred into 1 ml eNAT® transport medium.

Break the swab, close the tube, and shake the sample tube containing the swab sample and eNAT® medium for homogenization (avoid foaming).

Fill in 300 µl of homogenized patient sample in the sample input of the cartridge and close the cartridge lid. Insert the cartridge in the Vivalytic *one* analyser.

Do not use viscous samples that are difficult to pipette.

Urine

Fill in 1 ml of urine into the 1 ml eNAT® transport medium to achieve a mixing ratio of 1:1 of urine in eNAT® medium.

Shake the sample tube containing the urine and eNAT® medium for homogenization (avoid foaming) and fill in 300 µl of homogenized patient sample in the sample input of the cartridge and close the cartridge lid. Insert the cartridge in the Vivalytic *one* analyser.

Do not use viscous samples that are difficult to pipette.

Test Result

After automatic processing of the sample with the Valytic *one* analyser the test result is shown on the screen of the Valytic *one* analyser. Time-to-result is about one hour.

The sample is classified either as valid positive (if at least one of the targets *M. genitalium*, *M. hominis*, or *U. parvum/urealyticum* is detected), valid negative (if no target and only Human Control is detected) or invalid (if neither target nor Human Control is detected). In case of a positive detection of one or more targets, the test is considered valid even if the Human Control is negative.

U. parvum and *U. urealyticum* are both detected in one channel. A positive result for UP/UU indicates the presence of either one or both species.

Detection of the human cell based whole process control (Human Control) shows a successful DNA purification and amplification procedure and excludes an inhibition of the PCR reaction.

Interpretation of results is shown in the table below.

Interpretation of Results					
<i>M. genitalium</i>	<i>M. hominis</i>	<i>U. parvum/urealyticum</i>	Human Control	Validity	Result
+	+	+	+/-	valid	Sample is considered positive for <i>M. genitalium</i> , <i>M. hominis</i> and <i>U. parvum/urealyticum</i>
-	+	+	+/-	valid	Sample is considered positive for <i>M. hominis</i> and <i>U. parvum/urealyticum</i>
+	-	+	+/-	valid	Sample is considered positive for <i>M. genitalium</i> and <i>U. parvum/urealyticum</i>
+	+	-	+/-	valid	Sample is considered positive for <i>M. genitalium</i> and <i>M. hominis</i>
+	-	-	+/-	invalid	Sample is considered positive for <i>M. genitalium</i>
-	+	-	+/-	valid	Sample is considered positive for <i>M. hominis</i>
-	-	+	+/-	valid	Sample is considered positive for <i>U. parvum/urealyticum</i>
-	-	-	+	valid	Sample is considered negative for <i>M. genitalium</i> , <i>M. hominis</i> and <i>U. parvum/urealyticum</i>
-	-	-	-	invalid	Not evaluable. ¹

¹ Retesting is recommended.

PCR – Curve and C_q Value

Real-time PCR curves (software-modified) are shown and classified as positive or negative by the software. In case of positive curves, the respective C_q value is displayed. Inconclusive results are marked by the software (Δ). Retesting is advised.

Invalid or Failed Tests

A test is considered invalid if neither target(s) nor Human Control is detected. Possible reasons for an invalid run might be poor sample quality due to a partial or complete absence of human cellular material in the sample. Results are displayed for an invalid test but are not allowed to be used for diagnostic interpretation.

Vivalytic MG, MH, UP/UU – Instructions for Use

Pay attention to use the correct sample type, sample collection kit and storage conditions of sample and cartridges prior to the test run. If required, repeat the analysis with a new sample. In case of a failed test, first check for correct operating conditions of the Vivalytic *one* analyser (refer to Vivalytic *one* analyser's Instructions for Use). Restart the Vivalytic *one* analyser. If the problem persists, contact the customer support of your local distributor.

Test Termination

As soon as a valid, positive result for one of the targets *M. genitalium*, *M. hominis*, or *U. parvum/urealyticum* is shown on the screen, the user has the option to finish the test. Only finish the test if no further targets need to be analyzed. Pending results will not be displayed.

Test Report

In the printed test report, pathogens, results, control and information on user, patient and Vivalytic *one* analyser are listed with a signature field.

Quality Control

If required by your local or laboratory standards, quality control testing must be performed. The integrated control (Human Control) covers the processes of sample collection, nucleic acid purification, amplification, and detection. In case a reference testing method for quality control is needed the following reference materials can be used:

- AMPLIRUN® MYCOPLASMA GENITALIUM DNA CONTROL, Ref. MBC085 (Vircell)
- AMPLIRUN® MYCOPLASMA HOMINIS DNA CONTROL, Ref. MBC084 (Vircell)
- AMPLIRUN® UREAPLASMA DNA CONTROL, Ref. MBC133, Ref. MBC112 (Vircell)
- AMPLIRUN® TOTAL CT/NG/TV/MGE CONTROL (SWAB), Ref. MBTC024-R (Vircell)

Please follow the manufacturer's instructions.

In case of unexpected results, repeat the analysis with another sample. If the result of a quality negative control sample, e.g. pure eNAT® medium remains positive, the Vivalytic *one* analyser or its environment might be contaminated. Stop using the Vivalytic *one* analyser and clean the device as described in the Instructions for Use of the Vivalytic *one* analyser. In case of repeated false results for quality control samples please contact the customer support of your local distributor.

Notice to Users in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Limitations

The results of the Vivalytic MG, MH, UP/UU test must be interpreted by a trained healthcare professional only. The results of the Vivalytic MG, MH, UP/ UU test must not be used as the sole parameter for diagnosis.

- A negative result does not exclude pathogens being present in the sample at a level below assay sensitivity or other pathogens being present not covered by this assay.
- There is a risk of false negative or false positive results due to improperly collected, transported or handled samples.
- In borderline cases atypical PCR characteristics (e.g. flat curve with low or high C_q value) can occur. In case of atypical characteristics results are not allowed to be used for diagnostic interpretation. Retesting is advised.
- Inconclusive results are marked by the software. Retesting is advised.
- A positive result does not necessarily mean that viable bacterial cells are present.
- Do not use the Vivalytic MG, MH, UP/UU test for neonates.
- Vivalytic MG, MH, UP/UU is a qualitative real-time PCR test and does not provide a quantitative result.

Analytical Performance Evaluation

Analytical Sensitivity (Limit of Detection)

The limit of detection (LoD) of the Vivalytic MG, MH, UP/UU test was determined as the concentration at which 95 % of the testings are expected to yield a positive result using a Probit approach.

The LoD for each target pathogen was assessed by analyzing serial dilutions of samples prepared from Amplirun® Total Controls (Vircell) with either swab or urine background matrix. Each dilution was tested in 20 replicates with two independent production LOTS of Vivalytic cartridges distributed on three days. [Table 1](#) shows the LoD for each target.

Inclusivity

The amplicon sequences corresponding to *M. genitalium*, *M. hominis*, *U. urealyticum*, and *U. parvum* were analyzed in-silico (BLAST alignment) for sequence inclusivity.

Inclusivity analysis using MegaBLAST identified 338, 34, 3, and 10 sequences for *M. genitalium*, *M. hominis*, *U. urealyticum*, and respectively, harboring 100 % amplicon coverage with >96 % sequence identity for all targets evaluated ([Table 2](#)).

Exclusivity / Analytical Specificity

To exclude cross-reactivity (exclusivity), specificity of primers and probes was ensured by performing an *in silico* analysis for possible cross-reactions (BLAST alignment). In addition, the corresponding target region of *M. genitalium*, *M. hominis*, *U. parvum*, and *U. urealyticum* was analyzed *in silico* against the genomic sequence of various other pathogens. No cross-reactivity for the pathogens listed in [Table 3](#) could be predicted.

Precision

Repeatability and reproducibility testing of the Vivalytic MG, MH, UP/UU test was performed at 3 test sites. For each site, 3 different concentrations of each target pathogen were tested by the same operator in 2 replicates on 2 days, with the same set of Vivalytic *one* analysers with 3 LOTS resulting in 324 observations per target pathogen. The summarized results are highlighted in [Table 4](#).

Interferences

Interferences were evaluated for endogenous and exogenous substances that are potentially present in the patient sample. Refer to [Table 5](#) for substances that have the potential to interfere with the test.

Clinical Performance Evaluation


Diagnostic Sensitivity and Specificity

The clinical performance of the Vivalytic MG, MH, UP/UU test was assessed in a retrospective clinical study with urine and vaginal, cervical, urethral, and rectal specimens whereas swab samples were collected in eNAT® medium (COPAN Italia s.p.a.) and urine samples were diluted in eNAT® medium (COPAN Italia s.p.a.), respectively.

All clinical samples were tested using the Allplex™ STI Essential Assay (Seegene) as the reference method. Sensitivity or Positive Percent Agreement (PPA) was calculated as $100 \% \times TP / (TP+FN)$. Specificity or Negative Percent Agreement was calculated as $100 \% \times TN / (TN+FP)$. The results of the clinical performance evaluation are shown in [Table 6](#).

Vivalytic MG, MH, UP/UU – Instructions for Use

Symbols

	Manufacturer		Serial number
	Date of manufacture		Temperature limit
	Expiry date		Do not use if package is damaged
	Lot number		For single use only
	Reference number		Consult Instructions for Use
	Contains <n> tests		<i>in vitro</i> diagnostic medical device
	CE mark		

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Table 2 – Inclusivity

Query Sequence Name(s)	Number of matches with 100 % amplicon coverage	Number of matches with ≥ 90 % amplicon coverage	% Query Coverage (Range)	% Identity (Range)
<i>Mycoplasma genitalium</i>	338	338	100	96.2 – 100
<i>Mycoplasma hominis</i>	34	36	98.1 – 100	96.2 – 100
<i>Ureaplasma urealyticum</i>	3	3	100	100
<i>Ureaplasma parvum</i>	10	10	100	100

Table 3 – Exclusivity	
Pathogens (virus, bacteria, fungi) in Alphabetical Order	
<i>Acinetobacter baumannii</i>	<i>Micrococcus luteus</i>
<i>Acinetobacter lwoffii</i>	<i>Mobiluncus curtisii</i> subsp. <i>curtisii</i>
<i>Atopobium vaginae</i>	<i>Mobiluncus mulieris</i>
<i>Bacteroides fragilis</i>	<i>Mycoplasma pneumoniae</i>
<i>Candida albicans</i>	<i>Neisseria cinerea</i>
<i>Candida glabrata</i>	<i>Neisseria flavescens</i>
<i>Candida krusei</i> (<i>Issatchenkia orientalis</i>) (<i>Pichia kudriavzevii</i>)	<i>Neisseria gonorrhoeae</i>
<i>Candida parapsilosis</i>	<i>Neisseria lactamica</i>
<i>Candida tropicalis</i>	<i>Neisseria meningitidis</i> Serogroup 29E
<i>Chlamydia trachomatis</i>	<i>Neisseria meningitidis</i> Serogroup A (A, B, C, W135, X, Y, Z)
<i>Chlamydophila pneumoniae</i>	<i>Neisseria polysaccharea</i>
<i>Citrobacter freundii</i>	<i>Neisseria sicca</i>
<i>Corynebacterium amycolatum</i>	<i>Neisseria subflava</i>
<i>Corynebacterium aurimucosum</i>	<i>Peptostreptococcus anaerobius</i>
<i>Corynebacterium glucuronolyticum</i>	<i>Prevotella intermedia</i>
<i>Corynebacterium urealyticum</i>	<i>Proteus mirabilis</i>
Cytomegalovirus (CMV/HHV-5)	<i>Proteus vulgaris</i>
<i>Enterobacter aerogenes</i> (<i>Klebsiella aerogenes</i>)	<i>Pseudomonas aeruginosa</i>
<i>Enterobacter cloacae</i>	<i>Serratia marcescens</i>
<i>Enterococcus faecalis</i>	<i>Staphylococcus aureus</i>
<i>Enterococcus faecium</i>	<i>Staphylococcus aureus</i> MRSA
Epstein-Barr virus (EBV/HHV-4)	<i>Staphylococcus aureus</i> MecA+
<i>Escherichia coli</i>	<i>Staphylococcus capitis</i>
<i>Fusobacterium nucleatum</i>	<i>Staphylococcus epidermidis</i>
<i>Gardnerella vaginalis</i>	<i>Staphylococcus lugdunensis</i>
<i>Haemophilus influenzae</i>	<i>Staphylococcus sciuri</i>
<i>Haemophilus ducreyi</i>	<i>Staphylococcus saprophyticus</i>
Human Herpes Virus type 6 (HBLV/HHV-6)	<i>Stenotrophomonas maltophilia</i>
Human Papilloma Virus (HPV 6, 11, 16, 18)	<i>Streptococcus agalactiae</i> (group B strep)
<i>Klebsiella oxytoca</i>	<i>Streptococcus intermedius</i>
<i>Klebsiella pneumoniae</i>	<i>Streptococcus parasanguinis</i>
<i>Lactobacillus acidophilus</i>	<i>Streptococcus pneumoniae</i>
<i>Lactobacillus crispatus</i>	<i>Streptococcus pyogenes</i>
<i>Lactobacillus jensenii</i>	<i>Treponema pallidum</i>
<i>Lactobacillus rhamnosus</i>	<i>Trichomonas vaginalis</i>
<i>Listeria monocytogenes</i>	Varicella Zoster Virus VZV/HHV-3
	<i>Vellionella parvula</i>

Table 4 – Precision

Overall results <i>M. genitalium</i> and <i>U. parvum</i> stratified by applied concentration in swab matrix											
c*	Number of total tests	Number of positive tests			Proportion of positive results (%)			95 % Wilson-Score confidence interval (%)			
		MG	UP	MG	UP	MG	UP	MG	UP	MG	UP
c10	108	26	49	24.1	45.4	17.0–32.9	36.3–54.8	16.4–33.3	35.8–55.2		
c95	108	104	105	96.3	97.2	90.9–98.6	92.2–99.1	90.8–99.0	92.1–99.4		
c100	108	108	108	100	100	96.6–100	96.6–100	96.6–100	96.6–100		

Overall results <i>M. hominis</i> and <i>U. urealyticum</i> stratified by applied concentration in swab matrix											
c*	Number of total tests	Number of positive tests			Proportion of positive results (%)			95 % Wilson-Score confidence interval (%)			
		MH	UU	MH	UU	MH	UU	MH	UU	MH	UU
c10	108	17	38	15.7	35.2	10.1–23.8	26.8–44.6	9.5–24.0	26.2–45.0		
c95	108	103	108	95.4	100	89.6–98.0	96.6–100	89.5–98.5	96.6–100		
c100	108	107	107	99.1	99.1	94.9–99.8	94.9–99.8	95.0–100	95.0–100		

*c10 = 10 % predicted positive agreement; c95 = 95 % predicted positive agreement; c100 = 100 % predicted positive agreement

Table 5 – Tested Substances for Interference

No interference detected

For swab samples	For urine samples
Aciclovir cream (2 %) 50 mg/g	Acidity (pH 4.0)
Clotrimazol cream (1 %)	Albumin (direct) 3 mg/ml
FEMSOFT powder (0.1 %)	Alkalinity (pH 9.0)
KY Jelly (2 %)	Antibiotics (Amikacin) 30 µg/ml
Multi-Gyn IntiSkin (2 %)	Antibiotics (Gentamycin) 30 µg/ml
Vagisan washing lotion (2 %)	Bilirubin 0.2 mg/ml
Whole blood (2 %)	Glucose 10 mg/ml
Xylocain gel (2 %)	

Possibly interfering substance

Clotrimazol cream (2 %) (swab samples)

Table 6 – Diagnostic Sensitivity (PPA) and Specificity (NPA)¹

Swab samples ²	Urine Samples
<i>M. genitalium</i> – swab 97.8 % (PPA; n=45) 98.2 % (NPA; n=55)	<i>M. genitalium</i> – urine 97.7 % (PPA; n=44) 100 % (NPA; n=85)
<i>M. hominis</i> – swab 90.0 % (PPA; n=40) 100 % (NPA; n=60)	<i>M. hominis</i> – urine 100 % (PPA; n=44) 98.8 % (NPA; n=85)
<i>U. parvum/U. urealyticum</i> – swab 96.1 % (PPA; n=51) 95.9 % (NPA; n=49)	<i>U. parvum/U. urealyticum</i> – urine 90.2 % (PPA; n=51) 97.4 % (NPA; n=78)

¹ Clinical performance data of Vivalytic MG, MH, UP/UU in comparison to the reference test system Allplex™ STI Essential Assay (Seegene). Retrospective clinical study with in total 239 samples tested at one study site (132 urine and 107 swab specimens). In case of an insufficient number of either positive swabs or urine specimens, samples were spiked with the respective target (inactivated bacterial cells) in a concentration of approximately 2–3×LoD.

² Swab samples originated from vaginal, cervical, urethral, or rectal swab samples.

Table 7 – Document History

Revision 01	Initial document
Revision 02	Translation in 22 European languages

For more information see
www.bosch-vivalytic.com