# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation

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
[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

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## Package Contents

15 Vivalytic SARS-CoV-2 DT, Flu A/B, RSV test cartridges for the detection of SARS-CoV-2, Influenza A virus (Flu A), Influenza B virus (Flu B) and Respiratory Syncytial Virus (RSV) specific gene sequences.

## Safety Information

These Instructions for Use contain test-specific information only. For additional warnings and instructions refer to the Instructions for Use provided with your Vivalytic *one* analyser (chapter device safety information). Only use Vivalytic cartridges and accessories approved for the Vivalytic *one* analyser. Make sure to insert samples properly. Remove spilled samples with a cloth soaked with the appropriate solution. You can use 70 % ethanol. If you need to disinfect the cartridge by spraying (for example after heavy spillage of sample) do not use the cartridge afterwards since it might be damaged by disinfectant.

### WARNING

- Do not use a cartridge if the sealed pouch or the cartridge itself is visibly damaged.
- Do not touch or scratch the detection area of the cartridge.
- Do not reuse a cartridge.
- Do not use expired cartridges. The expiration date can be found on the packaging and the cartridge label.
- Do not wait longer than 15 min after opening the cartridge pouch to begin the test.
- Do not use sample types, media and volumes that are not approved for the test.
- Do not shake a cartridge that contains a sample.
- Do not turn the cartridge upside down.
- Always follow good laboratory practice to ensure the proper performance of this test.
- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions. Handle potentially infectious patient samples and cartridges according to national laboratory standards and dispose samples and cartridges according to regional and laboratory standards. Make sure to wear appropriate personal protective equipment (PPE).
- Place the cartridge on a clean and flat surface only.
- Be compliant with the national safety regulations and practices.

### Additional Equipment required (required but not provided)

- Bosch Vivalytic *one* analyser (reference number F 09G 300 115)
- Pipette (1000 µl)
- Swab collection kit + transport buffer (FLOQSwabs® regular Ref. 552C, FLOQSwabs® flexible Ref. 553C, eNAT® transport medium Ref. 608C, COPAN Italia S.P.A.)

### Intended Use

The Bosch Vivalytic SARS-CoV-2 DT, Flu A/B, RSV is an automated qualitative *in vitro* diagnostic test based on real time polymerase chain reaction (PCR) for the detection of nucleic acids from SARS-Coronavirus-2, Influenza A virus (Flu A), Influenza B virus (Flu B) and Human Respiratory Syncytial virus (RSV) from human nasopharyngeal, oropharyngeal or "lolly"-saliva swabs to aid in the diagnosis of respiratory tract infections of symptomatic individuals. Results should not be used as the sole basis for diagnosis, treatment or other patient management decisions. Positive results do not exclude co-infection with other pathogens. The agent(s) detected may not be the definite cause of disease. Negative results do not exclude a SARS-CoV-2, Flu A/B or RSV infection or another respiratory tract infection. Results must be clinically correlated with patient history, clinical observations and epidemiological information. Other diagnostic information is necessary to determine patient infection status.

Intended for use with a Vivalytic *one* analyser by healthcare professionals only in both near patient testing such as doctor's offices, emergency departments and nursing homes as well as laboratory settings such as hospital laboratories and reference laboratories.

## Test Principle

Vivalytic SARS-CoV-2 DT, Flu A/B, RSV is a qualitative real-time PCR based test utilizing reverse transcription (RT).

## Storage und Usage Conditions

Product is stable until the expiry date if stored at +15 °C to +25 °C. Storage and usage conditions can be taken either from the cartridge, pouch, or box label. Cartridge has to be used at +15 °C to +25 °C, relative humidity < 65 %, within 15 min upon pouch opening. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.

## Reagents

All reagents necessary for sample processing are integrated into the cartridge. The processing includes cell lysis, nucleic acid extraction, reverse transcription, DNA amplification, and detection.

Reagents are PCR bead, binding buffer, washing buffer and elution buffer. The PCR bead contains the Reverse Transcriptase and Taq Polymerase enzymes. Binding buffer facilitates binding of nucleic acids during the purification process. Washing buffer is a formulation of different salts and solvents to remove impurities e. g. proteins during the extraction process. Elution buffer is a low-salt buffer and contains the purified nucleic acids at the end of the extraction process.

## Sample Type/Medium

The test is intended for use with nasopharyngeal swab samples in eNAT® medium (FLOQSwabs® flexible Ref. 553C, eNAT® transport medium Ref. 608C, COPAN Italia S.P.A.), oropharyngeal swab samples or “lolly”-saliva swab samples in eNAT® medium (FLOQSwabs® regular Ref. 552C, eNAT® transport medium Ref. 608C, COPAN Italia S.P.A.). Collect and store samples as indicated in the manufacturer’s instructions.

## Sample Preparation

### Nasopharyngeal/Oropharyngeal Swabs\*

Avoid eating, drinking, smoking, and chewing gum 30 minutes before oropharyngeal swab sample collection.

Use a COPAN FLOQSwabs® flexible swab for collection of nasopharyngeal swab samples and a COPAN FLOQSwabs® regular swab for the collection of oropharyngeal swab samples. Each swab is transferred into 1 ml eNAT® transport medium. Break the swab and close the tube. Shake the sample tube containing the swab sample in eNAT® medium for homogenization (avoid foaming). Fill 300 µl of homogenized patient sample in the sample input of the cartridge and close the cartridge lid. Insert the cartridge into the Vivalytic *one* analyser.

### Saliva “Lolly” Swabs\*

Avoid eating, drinking, smoking, and chewing gum 30 minutes before sample collection. Use a COPAN FLOQSwab® regular swab in dry tube. Carefully remove the swab out of the dry tube. Avoid touching the outer surface of the tube. Slightly press the swab against the palate using the tongue. Rotate the swab for at least 30 seconds in this position and put the swab back into the dry tube. Transport and process the saliva samples in a dry environment (relative humidity < 65 %) and at room temperature (15–25 °C) within 4 h to avoid sample degradation. Transfer the saliva swab samples into a 1 ml COPAN eNAT® tube. Rotate the swab in the eNAT® medium at the inner surface of the tube for at least 15 seconds to ensure the sample is released into the eNAT® medium. Dispose the swab afterwards. Shake the tube thoroughly but avoid foaming.

Fill 300 µl of homogenized patient sample in the sample input of the cartridge and close the cartridge lid. Insert the cartridge into the Vivalytic *one* analyser.

\*Do not use viscous samples that are difficult to pipette.

## Test Result

The test result (list of positively/negatively detected targets) is shown on the screen. The sample is classified either as positive, negative, or invalid. For SARS-CoV-2, Flu A, Flu B and RSV one target is detected in each case. RSV A and RSV B cannot be differentiated from each other.

Detection of the human cell based whole process control (human control) in negative human samples shows a successful extraction procedure and excludes an inhibition of the PCR reaction. In case of a positive pathogen detection, the test is considered valid. In this case the human control can but does not have to be detected.

The PCR control shows a successful PCR reaction. If either the pathogen target(s) or the human control is detected, the PCR control can, but does not have to be detected for a valid result. If only the PCR control is positive, the test is classified as invalid. If a pathogen has been classified as invalid, retesting should be considered (see table below).

The time to result is less than 1 hour. For high titer specimens results are available earlier and the test run can be terminated.

Target/s (Flu A/B, RSV) <sup>1</sup>	Target SARS-CoV-2	Human Control 1 <sup>2</sup>	Human Control 2 <sup>3</sup>	PCR Control <sup>4</sup>	Result
+	-	+/-	+	+/-	Sample is considered positive for Flu A/B and/or RSV Sample is considered negative for SARS-CoV-2
+	-	+/-	-	+	Sample is considered positive for Flu A/B and/or RSV Sample is considered negative for SARS-CoV-2
-	+	+	+/-	+/-	Sample is considered negative for Flu A/B and RSV Sample is considered positive for SARS-CoV-2
+	+	+/-	+/-	+/-	Sample is considered positive for Flu A/B and/or RSV and SARS-CoV-2
-	-	+	+	+/-	Sample is considered negative for Flu A/B, RSV and SARS-CoV-2
-	-	+	-	+	Sample is considered negative for Flu A/B, RSV and SARS-CoV-2
-	+	-	+/-	+/-	Sample is not evaluable for Flu A/B and RSV (invalid) <sup>5</sup> Sample is considered positive for SARS-CoV-2
-	-	-	+	+/-	Sample is not evaluable for Flu A/B and RSV (invalid) <sup>5</sup> Sample is considered negative for SARS-CoV-2
+	-	+/-	-	-	Sample is considered positive for Flu A/B and/or RSV Sample is not evaluable for SARS-CoV-2 (invalid) <sup>5</sup>
-	-	+	-	-	Sample is considered negative for Flu A/B and RSV Sample is not evaluable for SARS-CoV-2 (invalid) <sup>5</sup>
-	-	-	-	+/-	Invalid <sup>5</sup>

<sup>1</sup> Individual results for each target is given on the user interface (analyzer screen).

<sup>2</sup> Whole process control including performance of PCR for Flu A/B and RSV targets.

<sup>3</sup> Whole process control including performance of PCR for SARS-CoV-2 target.

<sup>4</sup> PCR control for SARS-CoV-2 target.

<sup>5</sup> Retesting is recommended.

## Vivalytic SARS-CoV-2 DT, Flu A/B, RSV – Instructions for use

### PCR – Curve and C<sub>q</sub> Value

Real-time PCR curves (software-modified) are shown and classified as positive or negative by the software. In case of positive curves, the respective C<sub>q</sub> value is displayed. Based on the classification of targets as either positive or negative the result is defined. Inconclusive results are marked by the software (△). Retesting is advised.

### Invalid or Failed Tests

A test is rated as invalid if no target RNA and no human control is detected. Possible reasons for an invalid run might be poor sample quality due to a partial or complete absence of human cellular material in the sample. Results are displayed for an invalid test but are not allowed to be used for diagnostic interpretation. Pay attention to use the correct sample type, sample collection and storage of the sample and cartridges prior to the test run. If required, repeat the analysis with a new sample. In case of a failed test, first check for correct operating conditions of the Vivalytic *one* analyser (refer to Vivalytic *one* analyser's Instructions for Use). Restart the Vivalytic *one* analyser. If the problem persists, contact the customer support of your local distributor.

### Test Report

In the printed test report, results, pathogens, tested targets, controls and information on user, patient and Vivalytic *one* analyser are listed with a signature field.

### Test Termination

As soon as a valid, positive result is shown on the screen, the user has the option to finish the test. Only finish the test if no further targets need to be analyzed. Pending results will not be displayed.

### Quality Control

If required by your local or laboratory standards, quality control testing has to be performed. You can either use pre-characterized patient samples that were investigated by a reference testing method or purchase quality control materials (e. g. AMPLIRUN® SARS-CoV-2 RNA CONTROL, Ref. MBC137-R, AMPLIRUN® Total SARS-CoV-2, Ref. MBTC030, AMPLIRUN Total SARS-CoV-2/FluA/FLuB/RSV Control (Swab), Ref. MBTC031, Vircell S.L. Spain). RNA controls cover the processes of nucleic acid purification, reverse transcription amplification, and detection, but do not cover sample collection and sample lysis. Please be aware that parts of the listed reference materials are for research use only (RUO). Please follow the manufacturer's instructions. In case of unexpected results, repeat the analysis with another sample. If the result of a negative quality control sample e. g. pure eNAT® medium remains positive, the Vivalytic *one* analyser or its environment might be contaminated. Stop using the Vivalytic *one* analyser and clean the device as described in the Instructions for Use of the Vivalytic *one* analyser. In case of repeated false results for quality control samples please contact the customer support of your local distributor.

### Notice to Users in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Limitations

The results of the Vivalytic SARS-CoV-2 DT, FluA/B, RSV test must be interpreted by a trained healthcare professional only. The results of the Vivalytic SARS-CoV-2 DT, FluA/B, RSV test must not be used as the sole parameter for diagnosis.

- A negative result does not exclude pathogens being present in the sample at a level below assay sensitivity or a pathogen that is not covered by this assay.
- There is a risk of false negative or false positive results due to improperly collected, transported, or handled samples.
- In borderline cases atypical PCR characteristics (e. g. flat curve with low or high  $C_q$ -value) can occur. In case of atypical characteristics results are not allowed to be used for diagnostic interpretation. Inconclusive results are marked by the software ( $\Delta$ ). Retesting is advised.
- A positive result does not necessarily mean that viable viral particles are present.
- Vivalytic SARS-CoV-2 DT, Flu A/B, RSV is a qualitative real-time PCR test and does not provide a quantitative result.

## Analytical Performance evaluation

### Sensitivity (Limit of Detection, 95 % Detection rate)

The limit of detection of the Vivalytic SARS-CoV-2 DT, Flu A/B, RSV test was determined as the concentration at which 95 % of the testing is expected to yield a positive result using a hit rate approach (Flu A, Flu B and RSV) or a Probit approach (SARS-CoV-2) ([Table 1](#)).

### Inclusivity and Exclusivity

Specificity was ensured by the selection of primers and probes and their *in-silico* analysis for possible cross-reactions based on publicly available nucleic acid sequences derived from the NCBI database.

To evaluate inclusivity, reference material of several virus strains was tested with Vivalytic SARS-CoV-2 DT, Flu A/B, RSV cartridges ([Table 2](#)).

To exclude cross-reactivity (exclusivity), an *in-silico* analysis (BLAST alignment) of the target regions was conducted against the genomic sequences of various other pathogens representing common respiratory pathogens or closely related species. Additionally, various respiratory organisms were tested with Vivalytic SARS-CoV-2 DT, Flu A/B, RSV cartridges. There was no evidence of an interference ([Table 3](#)).

### Precision

The reproducibility and repeatability of the Vivalytic SARS-CoV-2 DT, Flu A/B, RSV test for SARS-CoV-2 was established using a panel with 3 different concentrations. At 3 test sites, each mix was tested on the same set of Vivalytic *one* analysers by the same operator in 2 replicates on 2 days, respectively. Vivalytic cartridges were drawn from 3 LOTs and distributed in the overall testing setup yielding in a total of 324 observations per target pathogen. The obtained positivity rates for the different combinations were correlated to the expected positivity rates ([Table 4](#)).

The reproducibility and repeatability of the Vivalytic SARS-CoV-2 DT, Flu A/B, RSV test for Flu A/B and RSV was performed at one test site with a total of 4 mixes tested on the same set of Vivalytic *one* analysers by the same operator in 2 replicates on 3 days with 3 different LOTs. The obtained positivity rates for the different combinations were correlated to the expected positivity rates ([Table 5](#)).

### Interferences

Interferences were evaluated for endogenous and exogenous substances that are potentially present in the patient sample. There was no evidence of an interference ([Table 6](#)).

## Vivalytic SARS-CoV-2 DT, Flu A/B, RSV – Instructions for use

### Clinical Performance evaluation

#### Sensitivity and Specificity

The clinical performance of the Vivalytic SARS-CoV-2 DT, Flu A/B, RSV test was assessed in a retrospective clinical study with nasopharyngeal, oropharyngeal, and lolly-saliva swab samples in eNAT® medium (COPAN Italia S.P.A.). Sensitivity or Positive Percent Agreement (PPA) was calculated as  $100 \% \times TP / (TP+FN)$ . Specificity or Negative Percent Agreement was calculated as  $100 \% \times TN / (TN+FP)$ .

For SARS-CoV-2 sensitivity and specificity results derived from oropharyngeal and nasopharyngeal samples (positive and negative samples in eNAT®, COPAN Italia S.P.A.) and lolly-saliva swabs (negative samples in eNAT®). Samples were collected in a clinical setting and compared with those of a reference method ([Table 7](#)). In total, 666 samples were analysed. Different matrix types (nasopharyngeal, oropharyngeal and lolly-saliva) were shown to be equivalent.

For FluA, FluB, and RSV sensitivity and specificity results derived from oropharyngeal, nasopharyngeal and lolly-swab samples (positive and negative samples in eNAT®, COPAN Italia S.P.A) tested with two different reference methods. In total, 617 samples were analysed. The results of the clinical performance evaluation are shown in [Table 8](#).

#### Symbols

Manufacturer

Serial number

Date of manufacture

Temperature limit

Expiry date

Do not use if package is damaged

Lot number

For single use only

Reference number

Consult Instructions for Use

Contains <n> tests

*in vitro* diagnostic medical device

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\* determined using a spiking approach in artificial simulated matrix (Respiratory Swab Matrix, Vircell) in eNAT® (COPAN) plus 1000 hPBTEC (human primary bronchial/tracheal epithelial cells).

\*\* determined using Zeptomatrix material for a spiking approach in Lolly-saliva samples in eNAT® (COPAN).

SARS-CoV-2 B.1.1.7 (Vircell AMPLIRUN® RNA Control)*	Influenza B/Massachusetts/2/2012 (Yamagata-Line) (50 TCID <sub>50</sub> /ml, Zeptomatrix***)
SARS-CoV-2 B.1.351 (Vircell AMPLIRUN® RNA Control)*	Influenza B/Yamagata/16/88 (5 TCID <sub>50</sub> /ml, Zeptomatrix***)
SARS-CoV-2 P.1 (Vircell AMPLIRUN® RNA Control)*	RSV A (Isolate:2006 Isolate) (25 TCID <sub>50</sub> /ml, Zeptomatrix***)
SARS-CoV-2 B.1.617.2 (Vircell AMPLIRUN® RNA Control)*	RSV B (Strain: CH93(18) -18) (0,25 TCID <sub>50</sub> /ml, Zeptomatrix***)
SARS-CoV-2 (Vircell AMPLIRUN® Total Control)*	Influenza A A/Brisbane/59/2007 (H1N1)**
SARS-CoV-2 B.1.1.7 + Q.x**	Influenza A A/Texas/50/2012 (H3N2)
SARS-CoV-2 B.1.351 + SARS-CoV-2 B.1.351.2 + B.1.351.3**	Influenza A A/Victoria/361/2011 (H3N2)**
SARS-CoV-2 B.1.429 + B.1.427**	Influenza A A/Perth/16/2009
SARS-CoV-2 B.1.525**	Influenza A A/Brisbane/10/2007 (H3N2)**
SARS-CoV-2 B.1.526**	Influenza A A/Kansas/14/2017 (H3N2)**
SARS-CoV-2 B.1.617.1**	Influenza A A/Michigan/45/2015 (H1N1)**
SARS-CoV-2 B.1.617.2 + AY.x**	Influenza A A/Singapore/INFIMH-16-0019/2016
SARS-CoV-2 B.1.617.3**	Influenza A A/Cambodia/e0826360/2020 (H3N2)**
SARS-CoV-2 B.1.621 + B.1.621.1**	Influenza A A/Wisconsin/588/2019 (H1N1)**
SARS-CoV-2 C.37 + C.37.1**	Influenza A A/Hawaii/70/2019 (H1N1)**
SARS-CoV-2 P.1 + P.1.x**	Influenza B B/Massachusetts/2/2012 (Yamagata-Line)**
SARS-CoV-2 P.2**	Influenza B B/Wisconsin/1/2010 (Yamagata-Line)**

## Vivalytic SARS-CoV-2 DT, Flu A/B, RSV – Annex

SARS-CoV-2 P.3**	Influenza B B/Phuket/3073/2013 (Yamagata-Line)**
Influenza A/Texas/50/2012 (H3N2) (500 TCID <sub>50</sub> /ml, Zeptomatrix***)	Influenza B B/Washington/02/2019 (Victoria-Line)**
Influenza A/Perth/16/2009 (100 TCID <sub>50</sub> /ml, Zeptomatrix***)	B/Colorado/06/2017 (Victoria-Line)**
Influenza A/Brisbane/10/2007 H3 Culture fluid (5 TCID <sub>50</sub> /ml, Zeptomatrix***)	Influenza B B/Brisbane/60/2008 (Victoria-Line)**
Influenza A/Hong Kong/4801/2014 (H3N2) (0,5 TCID <sub>50</sub> /ml, Zeptomatrix***)	RSV A**
Influenza A/California/07/09 (50 TCID <sub>50</sub> /ml, Zeptomatrix***)	RSV B**

\* Experimentally verified, tested at 492 viral copies.

\*\* Strains validated for inclusivity via BLAST alignment.

\*\*\* Experimentally verified using a spiking approach in eNAT® (COPAN) plus 1000 hPBTEC (human primary bronchial/tracheal epithelial cells).

<i>Acinetobacter baumannii</i>	Human orthopneumovirus (HRSV)2
<i>Aspergillus terreus</i>	<i>Klebsiella pneumoniae</i>
<i>Bacillus anthracis</i>	<i>Lactobacillus (L. plantarum, L. sp.)</i>
Bocaparvovirus	<i>Legionella longbeachae</i>
<i>Bordetella parapertussis</i>	<i>Legionella pneumophila</i>
<i>Bordetella pertussis</i>	<i>Leptospirosis</i>
<i>Candida albicans</i>	Human adenovirus (Type 1, 4, 7, 31, 34, 37)
<i>Chlamydia psittaci</i>	Measles morbillivirus
<i>Chlamydia pneumoniae</i>	MERS-coronavirus
<i>Clostridium perfringens</i>	<i>Moraxella catarrhalis</i>
<i>Corynebacterium diphtheriae</i>	Mumps orthorubulavirus
<i>Coxiella burnetii</i>	<i>Mycobacterium tuberculosis</i>
Cytomegalovirus	<i>Mycoplasma pneumoniae</i>
Enterovirus (Echovirus, Enterovirus Typ 71, Human Coxsackie Virus, Human Rhinovirus)	<i>Neisseria elongate</i>
<i>Enterobacter cloacae</i>	<i>Neisseria meningitides</i>
Epstein-Barr Virus	Parechovirus
<i>Escherichia coli</i>	<i>Pneumocystis jirovecii</i>
<i>Haemophilus influenzae</i>	<i>Proteus vulgaris</i>
Human coronavirus (229E, HKU1, NL63, OC43)	<i>Pseudomonas aeruginosa</i>
Human cytomegalovirus	SARS-coronavirus (without SARS-CoV-2)
Human metapneumovirus (hMPV)	<i>Serratia marcescens</i>
Human respirovirus 1 (HPIV-1)	Herpes simplex virus
Human orthorubulavirus 2 (HPIV-2)	<i>Staphylococcus (Staphylococcus epidermidis, Staphylococcus aureus)</i>
Human respirovirus 3 (HPIV-3)	<i>Streptococcus (Streptococcus pneumonia, Streptococcus pyogenes, Streptococcus salivarius)</i>
Human orthorubulavirus 4 (HPIV-4)2	Varicella zoster virus

Table 4 – Precision SARS-CoV-2

c*	Number of total tests	Number of positive tests	Proportion of positive results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
c70	108	35	32.4	24.3–41.7	23.7–42.1
c98	108	108	100	96.6–100	96.6–100
c100	108	108	100	96.6–100	96.6–100

\* c70 = 70 % predicted positive agreement,  
 c98 = 98 % predicted positive agreement,  
 c100 = 100 % predicted positive agreement

Table 5 – Precision Flu A/B, RSV

Pathogen	Applied concentration	Positive detections out of total number of runs	Percentage positive detections (%)	95 % Wilson-Score confidence interval (%)
Influenza A	Blank*	0/18	0	n. a.
	1500 cp/test	54/54	100.0	93.4–100.0
Influenza B	Blank*	0/36	0	n. a.
	1500 cp/test	36/36	100.0	90.4–100.0
RSV	Blank*	0/18	0	n. a.
	1500 cp/test	18/18	100.0	82.4–100.0
	5000 cp/test	36/36	100.0	90.4–100.0

\*Blank = negative control

Table 6 – Tested Substances for Interference

Tested substances	Final concentration
Nasivin® Spray (0.5 mg/ml Oxymetazolin- Hydrochlorid)	5.0 %
Rhinomer (9 g/l NaCl)	10.0 %
Hay fever spray	5.0 %
Turixin®	2.5 %
TOBRAZID®	4.0 µg/ml
Simulated nasal matrix [(2.5 % (w/v) porcine mucin, 1.0 % (v/v) human whole blood in 0.85 % sodium chloride (NaCl) in 1x PBS with 15 % glycerol]	1.0 %
Salbutamol-ratiopharm®	5.0 %
Neosynephrin-POS® 10% eye drops	5.0 %
Otri-Allergie Spray Fluticason	5.0 %
Human K2-EDTA blood*	1.0 %

Interferences were experimentally verified for SARS-CoV-2 with 1476 viral copies/test using a spiking approach in artificial simulated matrix (Respiratory Swab Matrix, Vircell) in eNAT® (COPAN) plus 1000 hPBTEC (human primary bronchial/tracheal epithelial cells) and for FluA H5N1, Flu B, RSV A with 1500 viral copies/test for defined pathogen using a spiking approach in eNAT® (COPAN) plus 1000 hPBTEC (human primary bronchial/tracheal epithelial cells).

\*Not tested for SARS-CoV-2.

**Table 7 – Clinical Sensitivity (PPA) and Specificity (NPA) for SARS-CoV-2 samples in eNAT® (95 % confidence interval)**

**SARS-CoV-2**

PPA 97.5 % (93.80 – 99.32 %)\*

NPA 98.4 % (96.87 – 99.31 %)\*\*

Clinical Performance data gained from 504 negative samples (442 oropharyngeal, 8 nasopharyngeal, 54 lolly-saliva) and 162 positive samples (93 oropharyngeal, 69 nasopharyngeal) in comparison to reference test system Vivalytic SARS-CoV-2 (Bosch Healthcare Solutions) after discrepancy analysis. Retrospective study with in total 666 swab samples tested at one study site. .

\* 4 samples were tested as false negative

\*\* 8 samples were reported as false positive

**Table 8a – Clinical Sensitivity (PPA) and Specificity (NPA) for Flu A, Flu B and RSV lolly-saliva samples in eNAT® (95 % confidence interval)**

**Flu A**

PPA 51.1 % (40.44– 61.66 %)\*

NPA 99.6 % (97.85 – 99.99 %)\*\*

**Flu B**

PPA 81.1 % (74.93 – 86.35 %)\*\*

NPA 99.4 % (96.94 – 99.99 %)\*\*

**RSV A/B**

PPA 84.1 % (72.74 – 92.12 %)\*\*

NPA 100.0 % (98.83 – 100.00 %)

Clinical performance data of Vivalytic SARS-CoV-2 DT, Flu A/B, RSV in comparison to the RIDA® GENE Flu & RSV (R-Biopharm AG). Retrospective clinical study with in total 427 valid results for Flu A, Flu B, and RSV each, generated at one study site. For Flu A, Flu B and RSV, all positive samples were created by spiking of negative patient swab samples in COPAN eNAT® medium with respective reference material (viral cells) in different concentrations between 1000-64000 copies/test. For Flu A H1N1 and H3N2 strains were used for spiking.

\* 45 samples were tested as false negative.

\*\* 1 sample were reported as false positive.

\*\* 37 samples were reported as false negative.

\*\* 10 samples were reported as false negative.

**Table 8b – Clinical Sensitivity (PPA) and Specificity (NPA) for Flu A, Flu B and RSV nasopharyngeal and oropharyngeal swab samples in eNAT® (95 % confidence interval, retrospective clinical study)**

**Flu A**

PPA 78.6 % (49.20 – 95.34 %)\*

NPA 91.3 % (71.96 – 98.93 %)\*\*

**Flu B**

PPA 100.0 % (54.07 – 100.00 %)

NPA 100.0% (88.78 – 100.00 %)

**RSV A/B**

PPA 47.1 % (22.98 – 72.19 %)\*\*

NPA 100.0 % (83.16 – 100.00 %)

Clinical performance data of Vivalytic SARS-CoV-2 DT, Flu A/B, RSV in comparison to the Allplex™ RV Essential Assay (Seegene). Retrospective clinical study with in total 37 valid results for Flu A, Flu B, and RSV each, generated at one study site.

\* 3 samples were tested as false negative.

\*\* 2 samples were reported as false positive.

\*\* 9 samples were reported as false negative.

**Table 8c – Clinical Sensitivity (PPA) and Specificity (NPA) for Flu A and RSV oropharyngeal swab samples in eNAT® (95 % confidence interval, prospective clinical study)**

**Flu A**

PPA 93.3 % (86.05– 97.51 %)\*

NPA 100.0 % (97.60 – 100.00 %)

**RSV A/B**

PPA 97.0 % (92.23 – 99.18 %)\*\*

NPA 100.0 % (96.64 – 100.00 %)

Clinical performance data of Vivalytic SARS-CoV-2 DT, Flu A/B, RSV in comparison to the Allplex™ RV Master Assay (Seegene). Prospective clinical study with in total 242 valid results tested for Flu A and RSV each, gained from 125 patient samples tested with two different cartridge lots and generated at one study site.

\* 6 samples were tested as false negative.

\*\* 4 samples were reported as false negativ.

**Table 9 – Document History**

Revision 01	Initial document
Revision 02	Update of Performance data
Revision 03	Precision of result matrix, addition of further performance data, Correction of header of Table 1
Revision 04	Added missing paragraph of intended use in German version

For more information see  
[www.bosch-vivalytic.com](http://www.bosch-vivalytic.com)