SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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
mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

BOSCH

vivalytic

C. difficile

PCR test for detection of
Clostridioides difficile (tcdA/tcdB)

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Introduction

Clostridioides difficile (*C. difficile*) is a gram-positive, spore-forming anaerobic bacillus consisting of toxigenic and non-toxigenic strains.

C. difficile is one of the most common pathogens of nosocomial diarrhea, especially in patients with antibiotic-associated diarrhea. In recent years, the number and severity of cases has increased in Europe and worldwide.¹

The symptoms of *C. difficile* infection (CDI) ranges from mild diarrhea to severe life-threatening pseudomembranous colitis, although colonization of *C. difficile* does not necessarily lead to a symptomatic infection. The disruption of the balance of gut microorganisms e.g. due to antibiologically treatment may cause toxigenic *C. difficile* to establish, dominate colonization and might be the start of CDI.²

Factors involved in the pathogenesis of *C. difficile*, are the genes *tcdA* and *tcdB* that code for toxin A (enterotoxin) and toxin B (cytotoxin) which are located in a 19.6 kb chromosomal region called the pathogenicity locus (PaLoc). Most pathogenic strains are toxin A and B positive (A+B+), but some pathogenic variant isolates may be toxin A negative and B positive (A-B+). They are the cause of *C. difficile*-associated diseases.³

Some strains of toxigenic *C. difficile* also produce a toxin called *C. difficile* transferase (CDT), or binary toxin.⁴

Package Contents

15 Vivalytic *C. difficile* test cartridges for the qualitative detection of *Clostridioides difficile* specific *tcdA/tcdB* toxin gene sequences.

Intended Use

The Vivalytic *C. difficile* test is an automated qualitative *in vitro* diagnostic test based on real-time polymerase chain reaction (PCR) for the detection of nucleic acids from toxigenic *C. difficile* (toxin genes *tcdA/tcdB*) from liquid or soft human stool swabs to aid in the diagnosis of acute gastrointestinal infections of symptomatic individuals.

Results should not be used as the sole basis for diagnosis, treatment, or other patient management decisions. Positive results do not exclude co-infection with other pathogens. The agent(s) detected may not be the definite cause of disease. Negative results do not exclude a *C. difficile* infection or another gastrointestinal infection. Results must be clinically correlated with patient history, clinical observations, and epidemiological information. Other diagnostic information is necessary to determine patient infection status. Intended for use with a Vivalytic *one* analyser by healthcare professionals only in laboratory settings such as hospital laboratories and reference laboratories.

Safety Information

These Instructions for Use contain test-specific information only. For additional warnings and instructions refer to the Instructions for Use provided with your Vivalytic *one* analyser (chapter device safety information). Only use Vivalytic cartridges and accessories approved for the Vivalytic *one* analyser. Take care to avoid any contamination when handling patient samples and cartridges. When sample was spilled on the cartridge, do not use the cartridge, and dispose it.

For *in vitro* diagnostic use by trained healthcare professional.

**WARNING**

- Always follow good laboratory practice to ensure the proper performance of this test.
- Make sure to wear appropriate personal protective equipment (PPE).
- Do not use a cartridge if the sealed pouch or the cartridge itself is visibly damaged.
- Do not touch or scratch the detection area of the cartridge.
- Do not reuse a cartridge.
- Do not use expired cartridges. The expiration date can be found on the packaging and the cartridge label.
- Do not wait longer than 15 minutes after opening the cartridge pouch to begin the test. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.
- Do not shake a cartridge that contains a sample.
- Do not turn the cartridge upside down.
- Place the cartridge on a clean and flat surface only.
- Do not use sample types, media and volumes that are not approved for the test.
- Always follow good laboratory practice to ensure the proper performance of this test.
- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions. Handle potentially infectious patient samples and cartridges according to national laboratory standards and dispose samples and cartridges according to regional and laboratory standards.
- Be compliant with the national safety regulations and practices.

Note: Further information can be found in the safety data sheet (SDS) of the product. Please contact the customer support of your local distributor.

Additional Equipment & Consumables Required but not Provided

- Bosch Vivalytic *one* analyser (reference number F 09G 300 115)
- Pipettor (100–1000 µl)
- Sterile filter pipette tips 100 – 1000 µl
- Swab collection kits
 - Regular Flocked Swab FLOQSwabs® 552C (COPAN Italia s.p.a.)
- Transport medium
 - eNAT® 606C (COPAN Italia s.p.a.)
- Suitable protective clothing

Test Principle

Vivalytic C. difficile is a qualitative real-time PCR based test.

Storage und Usage Conditions

Product is stable until the expiry date if stored at +15 °C to +25 °C. Storage and usage conditions can be taken either from the cartridge, pouch, or box label. Cartridge has to be used at +15 °C to +25 °C, relative humidity <65 %, within 15 min upon pouch opening. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.

Reagents

All reagents necessary for the sample processing are integrated into the cartridge. The processing includes cell lysis, nucleic acid extraction, DNA amplification and detection.

Reagents are PCR bead, binding buffer, washing buffer, and elution buffer. The PCR bead contains the DNA polymerase, primers, and probes. Binding buffer facilitates binding of nucleic acids during the purification process. Washing buffer is a formulation of different salts and solvents to remove impurities e. g. proteins during the extraction process. Elution Buffer is a low-salt buffer and contains the purified nucleic acids at the end of the extraction process.

Sample Types / Medium

The test is intended for use with stool swab samples in eNAT® medium (Regular Flocked Swab FLOQSwabs® 552C, eNAT® transport medium Ref. 606C, COPAN Italia s.p.a.).

Collect and store samples as indicated in the manufacturer's instructions.

In case the sample is not processed immediately after sample collection, nucleic acids will be preserved in eNAT® transport medium for up to 4 weeks at room temperature and at 4°C and up to 6 months at -20° C to -80° C (see Instructions for Use eNAT® transport medium⁵).

Sample Preparation

Use a Regular Flocked Swab FLOQSwabs® (COPAN Italia s.p.a.) to collect a small amount of stool by inserting the tip of the flocked swab into the stool sample and rotate it. Bloody, slimy, or watery area of stools should be selected and collected. After collection, examine the swab to make sure there is fecal material visible on the tip. In case it is not, again insert the flocked swab into the stool sample and rotate taking care that all the area of the swab tip is in contact with the sample. After collection transfer the swab into the 2 ml Copan eNAT® tube with eNAT® transport medium. Hold the swab shaft between thumb and finger, mash and mix the stool specimen against the side of the tube to evenly disperse and suspend the specimen in the preservation medium. Hold the swab shaft close to the rim of the tube, bend it at a 180 degrees angle to break it off at the marked breakpoint. Discard the broken upper part of the swab shaft and tighten the cap. Shake the sample tube containing the swab sample and eNAT® medium (COPAN Italia s.p.a.) thoroughly for homogenization and use a pipettor in order to fill 300 µl of homogenized patient sample in the sample input of the cartridge. Do not use viscous samples that are difficult to pipette.

Test Result

After automatic processing of the sample with the Vivalytic *one* analyser the test result is shown on the screen of the Vivalytic *one* analyser. The time to result is about 50 min. For high titer specimens results are available after less than 35 min and the test run can be terminated earlier (see Chapter Test Termination).

The sample is classified either as *C. difficile* positive, *C. difficile* negative or invalid. In case of a positive detection of *C. difficile*, the test is considered valid even if the Human Control is negative.

Detection of the human cell based whole process control (Human Control) in negative samples shows a successful extraction procedure and excludes an inhibition of the PCR reaction. Interpretation of results is shown in the table below.

<i>C. difficile</i> (tcdA / tcdB)	Human Control	Validity	Result
+	+/-	valid	Sample is considered positive for <i>C. difficile</i> (tcdA/tcdB toxin gene detected).
-	+	valid	Sample is considered negative for <i>C. difficile</i> (tcdA/tcdB toxin gene not detected)
-	-	invalid	Not evaluable. ¹

¹ Retesting is recommended.

PCR – Curve and C_q Value

Real-time PCR curves (software-modified) are shown and classified as positive or negative by the software. In case of positive curves, the respective C_q value is displayed. Inconclusive results are marked by the software (△). Retesting is advised.

Vivalytic *C. difficile* – Instructions for Use

Invalid or Failed Tests

A test is rated as invalid if neither target DNA nor Human Control is detected. Possible reasons for an invalid run might be poor sample quality due to a partial or complete absence of human cellular material in the sample. Results are displayed for an invalid test but are not allowed to be used for diagnostic interpretation. Pay attention to use the correct sample type, sample collection and storage of the sample and cartridges prior to the test run. If required, repeat the analysis with a new sample.

In case of a failed test, first check for correct operating conditions of the Vivalytic *one* analyser (refer to Vivalytic *one* analyser's Instructions for Use). Restart the Vivalytic *one* analyser. If the problem persists, contact the customer support of your local distributor.

Test Termination

As soon as a valid, positive *C. difficile* result is shown on the screen, the user has the option to finish the test.

Test Report

In the printed test report, pathogen, results, control and information on user, patient and Vivalytic *one* analyser are listed with a signature field.

Quality Control

If required by your local or laboratory standards, quality control testing must be performed.

In case a reference testing method for quality control is needed the following reference materials can be used:

CLOSTRIDIUM DIFFICILE RT027 100K TOTAL CONTROL (NEUTRAL), Vircell S.L. Spain).

Please follow the manufacturer's instructions.

In case of unexpected results, repeat the analysis with another sample. If the result of a negative quality control sample e. g. pure eNAT® medium remains positive, the Vivalytic *one* analyser or its environment might be contaminated. Stop using the Vivalytic *one* analyser and clean the device as described in the Instructions for Use of the Vivalytic *one* analyser.

In case of repeated false results for quality control samples please contact the customer support of your local distributor.

Notice to Users in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Limitations

The results of the Vivalytic *C. difficile* test must be interpreted by a trained healthcare professional only. The results of the Vivalytic *C. difficile* test must not be used as the sole parameter for diagnosis.

- A negative result does not exclude pathogens being present in the sample at a level below assay sensitivity or other pathogens being present not covered by this assay.
- There is a risk of false negative or false positive results due to improperly collected, transported, or handled samples.
- In borderline cases atypical PCR characteristics (e. g. flat curve with low or high C_q -value) can occur. In case of atypical characteristics results are not allowed to be used for diagnostic interpretation. Inconclusive results are marked by the software. Retesting is advised.
- Vivalytic *C. difficile* is a qualitative real-time PCR test and does not provide a quantitative result.
- A positive result does not necessarily mean that viable pathogens are present.
- The test detects *C. difficile* carrying *tcdA/tcdB* genes, but does not differentiate between toxin producing and non-toxin producing strains
- A negative result does not preclude *C. difficile* infection. It is recommended that negative tested samples are interpreted in the context with additional laboratory data
- An excess amount of stool may have inhibitory effects on the assay performance.

Analytical Performance Evaluation

Analytical Sensitivity (Limit of Detection)

The limit of detection of the Vivalytic C. difficile test was determined as the lowest concentration of analyte that can be consistently detected ($\geq 95\%$ of samples tested under routine laboratory conditions using a defined type of sample).

Inclusivity

To evaluate inclusivity, an *in silico* analysis (BLAST alignment) of the genomic sequence of various relevant C. difficile strains against the sequence of the PCR primers and hydrolysis probe used in the Vivalytic C. difficile test for amplification and detection of the respective pathogens was performed. Inclusivity could be shown for strains listed in [Table 2](#).

Exclusivity / Analytical Specificity

To exclude cross-reactivity (exclusivity), an *in silico* analysis (BLAST alignment) of the target region of C. difficile against the genomic sequence of various other pathogens representing common gastrointestinal pathogens or closely related species was conducted. There was no evidence of an interference ([Table 3](#)).

Reproducibility

The reproducibility of the Vivalytic C. difficile test was established using a panel with 3 different concentrations of C. difficile. At 3 test sites, each mix was tested on the same set of Vivalytic instruments by the same operator with 3 LOTs in 4 replicates on 3 days, respectively, yielding in a total of 324 observations per target pathogen. The obtained positivity rates for the different combinations were correlated to the expected positivity rate ([Table 4a](#)).

Repeatability

The repeatability of the Vivalytic C. difficile test was established using a panel with 1 concentration (3x c95) of C. difficile. At 1 test site, the mix was tested on the same set of Vivalytic instruments by the same operator with 3 LOTs in 20 replicates on 1 day, respectively, yielding in a total of 60 observations per target pathogen. The obtained positivity rates for the different combinations were correlated to the expected positivity rate ([Table 4b](#)).

Interferences

Interferences were evaluated for endogenous and exogenous substances, that are potentially present in the patient sample. Refer to [Table 5](#) for substances that have the potential to interfere with the test.

Clinical Performance Evaluation

Sensitivity and specificity results derived from native liquid and soft human stool samples. Samples were collected in a clinical setting and compared with results of reference methods.

Samples for testing with Vivalytic C. difficile cartridges were freshly used or frozen for storage and prepared as describe above in eNAT® (COPAN Italia S.p.A.).

In case of reference testing samples were prepared according to recommendations of used reference methods. In total, 124 samples were analysed. Sensitivity or Positive Percent Agreement (PPA) was calculated as $100\% \times TP / (TP+FN)$. Specificity or Negative Percent Agreement was calculated as $100\% \times TN / (TN+FP)$. The results of the clinical performance evaluation are shown in [Table 6](#).

Technical Support

In case of questions please contact your local distributor or send a message to the customer service of Vivalytic (email: support-BHCS@bosch-vivalytic.com).

Vivalytic C. difficile – Instructions for Use

References

- ¹Kuijper EJ, Coignard B, Tüll P; ESCMID Study Group for Clostridium difficile; EU Member States; European Centre for Disease Prevention and Control. Emergence of Clostridium difficile-associated disease in North America and Europe. Clin Microbiol Infect. 2006 Oct;12 Suppl 6:2-18.
- ²Czepiel J, Drózdź M, Pituch H, Kuijper EJ, Perucki W, Mielimonka A, Goldman S, Wultańska D, Garlicki A, Biesiada G. Clostridium difficile infection: review. Eur J Clin Microbiol Infect Dis. 2019 Jul;38(7):1211-1221. doi: 10.1007/s10096-019-03539-6. Epub 2019 Apr 3. PMID: 30945014; PMCID: PMC6570665.
- ³Kordus SL, Thomas AK, Lacy DB. Clostridioides difficile toxins: mechanisms of action and antitoxin therapeutics. Nat Rev Microbiol. 2022 May; 20(5): 285-298.
- ⁴Gerding DN, Johnson S, Rupnik M, Aktories K. Clostridium difficile binary toxin CDT: mechanism, epidemiology, and potential clinical importance. Gut Microbes. 2014 Jan-Feb; 5(1): 15-27.
- ⁵COPAN Instructions for Use of Copan eNAT® Collection and Preservation System; Version EIFU002R01 Date 2021.06.

Symbols

Manufacturer

Serial number

Date of manufacture

Temperature limit

Expiry date

Do not use if package is damaged

Lot number

For single use only

Reference number

Consult Instructions for Use

Contains <n> tests

in vitro diagnostic medical device

CE mark

Pipette the indicated sample volume in the sample input of the cartridge as marked by the black triangle.

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Références

- ¹Kuijper EJ, Coignard B, Tüll P; ESCMID Study Group for Clostridium difficile; EU Member States; European Centre for Disease Prevention and Control. Emergence of Clostridium difficile-associated disease in North America and Europe. Clin Microbiol Infect. 2006 Oct;12 Suppl 6: 2-18.
- ²Czepiel J, Drózdź M, Pituch H, Kuijper EJ, Perucki W, Mielimonka A, Goldman S, Wultańska D, Garlicki A, Biesiada G. Clostridium difficile infection: review. Eur J Clin Microbiol Infect Dis. 2019 Jul; 38(7): 1211-1221. doi: 10.1007/s10096-019-03539-6. Epub 2019 Apr 3. PMID: 30945014; PMCID: PMC6570665.
- ³Kordus SL, Thomas AK, Lacy DB. Clostridioides difficile toxins: mechanisms of action and antitoxin therapeutics. Nat Rev Microbiol. 2022 May; 20(5): 285-298.
- ⁴Gerding DN, Johnson S, Rupnik M, Aktories K. Clostridium difficile binary toxin CDT: mechanism, epidemiology, and potential clinical importance. Gut Microbes. 2014 Jan-Feb; 5(1): 15-27.
- ⁵COPAN Instructions for Use of Copan eNAT® Collection and Preservation System; Version EIFU002R01 Date 2021.06.

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- ¹Kuijper EJ, Coignard B, Tüll P; ESCMID Study Group for Clostridium difficile; EU Member States; European Centre for Disease Prevention and Control. Emergence of Clostridium difficile-associated disease in North America and Europe. Clin Microbiol Infect. 2006 Oct;12 Suppl 6: 2-18.
- ²Czepiel J, Drózd M, Pituch H, Kuijper EJ, Perucki W, Mielimonka A, Goldman S, Wultańska D, Garlicki A, Biesiada G. Clostridium difficile infection: review. Eur J Clin Microbiol Infect Dis. 2019 Jul; 38(7): 1211-1221. doi: 10.1007/s10096-019-03539-6. Epub 2019 Apr 3. PMID: 30945014; PMCID: PMC6570665.
- ³Kordus SL, Thomas AK, Lacy DB. Clostridioides difficile toxins: mechanisms of action and antitoxin therapeutics. Nat Rev Microbiol. 2022 May; 20(5): 285-298.
- ⁴Gerding DN, Johnson S, Rupnik M, Aktories K. Clostridium difficile binary toxin CDT: mechanism, epidemiology, and potential clinical importance. Gut Microbes. 2014 Jan-Feb; 5(1): 15-27.
- ⁵COPAN Instructions for Use of Copan eNAT® Collection and Preservation System; Version EIFU002R01 Date 2021.06.

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Bibliografie

- ¹Kuijper EJ, Coignard B, Tüll P; ESCMID Study Group for Clostridium difficile; EU Member States; European Centre for Disease Prevention and Control. Emergence of Clostridium difficile-associated disease in North America and Europe. Clin Microbiol Infect. 2006 Oct; 12 Suppl 6: 2-18.
- ²Czepiel J, Drózd M, Pituch H, Kuijper EJ, Perucki W, Mielimonka A, Goldman S, Wultańska D, Garlicki A, Biesiada G. Clostridium difficile infection: review. Eur J Clin Microbiol Infect Dis. 2019 Jul; 38(7): 1211-1221. doi: 10.1007/s10096-019-03539-6. Epub 2019 Apr 3. PMID: 30945014; PMCID: PMC6570665.
- ³Kordus SL, Thomas AK, Lacy DB. Clostridioides difficile toxins: mechanisms of action and antitoxin therapeutics. Nat Rev Microbiol. 2022 May; 20(5): 285-298.
- ⁴Gerding DN, Johnson S, Rupnik M, Aktories K. Clostridium difficile binary toxin CDT: mechanism, epidemiology, and potential clinical importance. Gut Microbes. 2014 Jan-Feb; 5(1): 15-27.
- ⁵COPAN Instructions for Use of Copan eNAT® Collection and Preservation System; Version EIFU002R01 Date 2021.06.

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Reference

- ¹Kuijper EJ, Coignard B, Tüll P; ESCMID Study Group for Clostridium difficile; EU Member States; European Centre for Disease Prevention and Control. Emergence of Clostridium difficile-associated disease in North America and Europe. Clin Microbiol Infect. 2006 Oct;12 Suppl 6: 2-18.
- ²Czepiel J, Drózdź M, Pituch H, Kuijper EJ, Perucki W, Mielimonka A, Goldman S, Wultańska D, Garlicki A, Biesiada G. Clostridium difficile infection: review. Eur J Clin Microbiol Infect Dis. 2019 Jul;38(7): 1211-1221. doi: 10.1007/s10096-019-03539-6. Epub 2019 Apr 3. PMID: 30945014; PMCID: PMC6570665.
- ³Kordus SL, Thomas AK, Lacy DB. Clostridioides difficile toxins: mechanisms of action and antitoxin therapeutics. Nat Rev Microbiol. 2022 May; 20(5): 285-298.
- ⁴Gerding DN, Johnson S, Rupnik M, Aktories K. Clostridium difficile binary toxin CDT: mechanism, epidemiology, and potential clinical importance. Gut Microbes. 2014 Jan-Feb; 5(1): 15-27.
- ⁵COPAN Instructions for Use of Copan eNAT® Collection and Preservation System; Version EIFU002R01 Date 2021.06.

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Referencias

- ¹Kuijper EJ, Coignard B, Tüll P; ESCMID Study Group for Clostridium difficile; EU Member States; European Centre for Disease Prevention and Control. Emergence of Clostridium difficile-associated disease in North America and Europe. Clin Microbiol Infect. 2006 Oct;12 Suppl 6: 2-18.
- ²Czepiel J, Drózdź M, Pituch H, Kuijper EJ, Perucki W, Mielimonka A, Goldman S, Wultańska D, Garlicki A, Biesiada G. Clostridium difficile infection: review. Eur J Clin Microbiol Infect Dis. 2019 Jul; 38(7): 1211-1221. doi: 10.1007/s10096-019-03539-6. Epub 2019 Apr 3. PMID: 30945014; PMCID: PMC6570665.
- ³Kordus SL, Thomas AK, Lacy DB. Clostridioides difficile toxins: mechanisms of action and antitoxin therapeutics. Nat Rev Microbiol. 2022 May; 20(5): 285-298.
- ⁴Gerding DN, Johnson S, Rupnik M, Aktories K. Clostridium difficile binary toxin CDT: mechanism, epidemiology, and potential clinical importance. Gut Microbes. 2014 Jan-Feb; 5(1): 15-27.
- ⁵COPAN Instrucciones de uso del Copan eNAT® Collection and Preservation System; Version EIFU002R01 Date 2021.06.

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Guanidinium chloride; guanidine hydrochloride

Hazard statements

H302	Harmful if swallowed.
H315	Causes skin irritation.
H319	Causes serious eye irritation.

Precautionary statements

P264	Wash hands thoroughly after handling.
P280	Wear protective gloves and eye/face protection.
P302+P352	IF ON SKIN: Wash with plenty of water and soap.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.
P501	Dispose of contents/container in accordance with local/ regional/ national/ international regulation

Table 1 – Limit of Detection (LoD)

<i>C. difficile</i> (<i>tcdA/tcdB</i>)*	10000 copies /mL
---	------------------

*Determined by using Vircell material (surrogate strain R20291 Ribotype 027 (toxin A + B positive) for a spiking approach in eNAT® (COPAN) supplemented with approximately 1000 Human Bronchial/Tracheal Epithelial Cells (hPBTEC) and 6% (v/v) Stool Matrix Negative Control (artificial simulated matrix).

Table 2 – Inclusivity*Clostridioides difficile* toxin A gene*Clostridioides difficile* toxin B gene

Strains validated for inclusivity via BLAST alignment

Table 3 – Exclusivity

<i>Abiotrophia defectiva</i>	<i>Enterococcus faecium</i>
<i>Acinetobacter baumannii</i>	<i>Enterococcus gallinarum</i>
<i>Acinetobacter lwoffii</i>	<i>Enterococcus hirae</i>
Adenovirus	<i>Enterococcus raffinosus</i>
<i>Aeromonas caviae</i>	Enterovirus
<i>Aeromonas hydrophila</i>	<i>Escherichia blattae</i>
<i>Alcaligenes faecalis</i>	<i>Escherichia coli</i> (Shiga toxin producing) (STEC)
<i>Anaerococcus tetradius</i>	<i>Escherichia coli</i> (Enteroaggregative) (EAEC)
<i>Arcobacter butzleri</i>	<i>Escherichia coli</i> (Enteroinvasive) (EIEC)
<i>Arcobacter cryaerophilus</i>	<i>Escherichia coli</i> (Enteropathogenic) (EPEC)
<i>Ascaris lumbricoides</i>	<i>Escherichia coli</i> (Enterotoxigenic) (ETEC)
Astrovirus	<i>Escherichia coli</i> (non-pathogenic) (K-12 MG1655)

Table 3 – Exclusivity

<i>Bacillus cereus</i>	<i>Escherichia coli</i> (Enterohemorrhagic) (EHEC)
<i>Bacillus subtilis</i>	<i>Escherichia fergusonii</i>
<i>Bacteroides caccae</i>	<i>Escherichia hermannii</i>
<i>Bacteroides fragilis</i>	<i>Escherichia vulneris</i>
<i>Bacteroides merdae</i>	<i>Eubacterium rectale</i>
<i>Bacteroides stercoris</i>	<i>Flavonifractor plautii</i>
<i>Bacteroides thetaiotaomicron</i>	<i>Fusobacterium varium</i>
<i>Bifidobacterium adolescentis</i>	<i>Gardnerella vaginalis</i>
<i>Bifidobacterium bifidum</i>	<i>Gemella morbillorum</i>
<i>Bifidobacterium longum</i>	<i>Giardia lamblia</i>
<i>Blastocystis hominis</i>	<i>Giardia muris</i>
<i>Campylobacter coli</i>	<i>Grimontia hollisae</i>
<i>Campylobacter concisus</i>	<i>Hafnia alvei</i>
<i>Campylobacter curvus</i>	<i>Helicobacter fennelliae</i>
<i>Campylobacter fetus</i>	<i>Helicobacter pylori</i>
<i>Campylobacter gracilis</i>	<i>Helicobacter cinaedi</i>
<i>Campylobacter helveticus</i>	<i>Helicobacter hepaticus</i>
<i>Campylobacter hominis</i>	Hepatitis A virus
<i>Campylobacter hyointestinalis</i>	Herpes Simplex Virus 2
<i>Campylobacter jejuni</i>	<i>Hymenolepis nana</i>
<i>Campylobacter lari</i>	<i>Klebsiella oxytoca</i>
<i>Campylobacter mucosalis</i>	<i>Klebsiella ozaenae</i>
<i>Campylobacter rectus</i>	<i>Klebsiella pneumoniae</i>
<i>Campylobacter showae</i>	<i>Lactobacillus acidophilus</i>
<i>Campylobacter sputorum</i>	<i>Lactobacillus reuteri</i>
<i>Campylobacter upsaliensis</i>	<i>Lactococcus lactis</i>
<i>Candida albicans</i>	<i>Leminorella grimontii</i>
<i>Candida catenulata</i>	<i>Listeria grayi</i>
<i>Cedecea davisae</i>	<i>Listeria innocua</i>
<i>Chilomastix mesnili</i>	<i>Listeria monocytogenes</i>
<i>Chlamydia trachomatis</i>	<i>Megamonas hypermegale</i>
<i>Citrobacter amalonaticus</i>	<i>Morganella morganii</i>
<i>Citrobacter freundii</i>	Norovirus GIV
<i>Citrobacter koseri</i>	<i>Pentatrichomonas hominis</i>
<i>Citrobacter sedlakii</i>	<i>Peptoniphilus asaccharolyticus</i>
<i>Clostridium beijerinckii</i>	<i>Peptostreptococcus anaerobius</i>
<i>Clostridium bifermentans</i>	<i>Plesiomonas shigelloides</i>
<i>Clostridium boltea</i>	<i>Porphyromonas asaccharolytica</i>
<i>Clostridium botulinum</i>	<i>Prevotella melaninogenica</i>
<i>Clostridium butyricum</i>	<i>Proteus mirabilis</i>
<i>Clostridium chauvoei</i>	<i>Proteus penneri</i>
<i>Clostridium difficile</i> (non-toxigenic)	<i>Proteus vulgaris</i>
<i>Clostridium fallax</i>	<i>Providencia alcalifaciens</i>
<i>Clostridium haemolyticum</i>	<i>Providencia rettgeri</i>
<i>Clostridium histolyticum</i>	<i>Providencia stuartii</i>
<i>Clostridium innocuum</i>	<i>Pseudomonas aeruginosa</i>

<i>Clostridium methylpentosum</i>	<i>Pseudomonas fluorescens</i>
<i>Clostridium nexile</i>	<i>Pseudomonas putida</i>
<i>Clostridium novyi</i>	Rotavirus B
<i>Clostridium orbiscindens</i>	Rotavirus C
<i>Clostridium paraputrificum</i>	Rotavirus D
<i>Clostridium perfringens</i>	Rotavirus F
<i>Clostridium ramosum</i>	Rotavirus G
<i>Clostridium scindens</i>	Rotavirus H
<i>Clostridium septicum</i>	Rotavirus I
<i>Clostridium sordellii</i>	<i>Ruminococcus bromii</i>
<i>Clostridium sphenoides</i>	<i>Saccharomyces boulardii</i>
<i>Clostridium spiroforme</i>	<i>Salmonella enterica</i>
<i>Clostridium sporogenes</i>	<i>Salmonella enterica</i> spp. <i>enterica</i> serovar <i>Typhimurium</i>
<i>Clostridium symbiosum</i>	<i>Salmonella enteritidis</i>
<i>Clostridium tertium</i>	<i>Salmonella subterranea</i>
<i>Clostridium tetani</i>	Sapovirus
<i>Collinsella aerofaciens</i>	<i>Serratia liquefaciens</i>
<i>Corynebacterium genitalium</i>	<i>Serratia marcescens</i>
Coxsackie virus	<i>Shigella boydii</i>
<i>Cryptosporidium canis</i>	<i>Shigella dysenteriae</i>
<i>Cryptosporidium cuniculus</i>	<i>Shigella flexneri</i>
<i>Cryptosporidium felis</i>	<i>Shigella sonnei</i>
<i>Cryptosporidium fetus</i>	<i>Staphylococcus aureus</i>
<i>Cryptosporidium hominis</i>	<i>Staphylococcus aureus</i> subsp. <i>aureus</i>
<i>Cryptosporidium meleagridis</i>	<i>Staphylococcus epidermidis</i>
<i>Cryptosporidium muris</i>	<i>Stenotrophomonas maltophilia</i>
	<i>Streptococcus agalactiae</i>
Cytomegalovirus	<i>Streptococcus dysgalactiae</i>
<i>Desulfovibrio piger</i>	<i>Streptococcus dysgalactiae</i> subsp. <i>dysgalactiae</i>
<i>Dientamoeba fragilis</i>	<i>Streptococcus intermedius</i>
<i>Diphyllobothrium latum</i>	<i>Streptococcus pyogenes</i>
<i>Dysgonomonas capnocytophagoides</i>	<i>Streptococcus salivarius</i>
Echovirus	<i>Streptococcus uberis</i>
<i>Edwardsiella tarda</i>	<i>Taenia saginata</i>
<i>Eggerthella lenta</i>	<i>Trabulsiella guamensis</i>
<i>Encephalitozoon intestinalis</i>	<i>Veillonella parvula</i>
<i>Endolimax nana</i>	<i>Vibrio cholerae</i>
<i>Entamoeba coli</i>	<i>Vibrio mimicus</i>
<i>Entamoeba hartmanni</i>	<i>Vibrio parahaemolyticus</i>
<i>Entamoeba histolytica</i>	<i>Vibrio vulnificus</i>
<i>Entamoeba moshkovskii</i>	<i>Yersinia bercovieri</i>
<i>Entamoeba polecki</i>	<i>Yersinia enterocolitica</i>
<i>Enterobacter aerogenes</i>	<i>Yersinia enterocolitica</i> subsp. <i>enterocolitica</i>
<i>Enterobacter cloacae</i>	<i>Yersinia frederiksenii</i>

Table 3 – Exclusivity

<i>Enterococcus casseliflavus</i>	<i>Yersinia intermedia</i>
<i>Enterococcus cecorum</i>	<i>Yersinia mollaretii</i>
<i>Enterococcus dispar</i>	<i>Yersinia pseudotuberculosis</i>
<i>Enterococcus faecalis</i>	<i>Yersinia rohdei</i>

Table 4 a Reproducibility

c*	Number of total tests	Number of pos. tests	Number of neg. tests	Proportion of positive/negative results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
blank	108	0	108	100	96.6 % – 100 %	96.6 % – 100 %
3x c95	108	108	0	100	96.6 % – 100 %	96.6 % – 100 %
<c100	108	108	0	100	96.6 % – 100 %	96.6 % – 100 %

3x c95 = 95 % predicted positive agreement
 < c100 = < 100 % predicted positive agreement
 (verified with 1,83x10⁵ copies/ml; c100 = 5,48x10⁵ copies/ml)

Table 4 b Repeatability

LOT	No of total tests	No of pos. tests (3x c95)	Proportion of positive/negative results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
1	20	20	100	83.9 % – 100 %	83.2 % – 100 %
2	20	20	100	83.9 % – 100 %	83.2 % – 100 %
3	20	20	100	83.9 % – 100 %	83.2 % – 100 %
Total	60	60	100	94 % – 100 %	94 % – 100 %

3x c95 = 95 % predicted positive agreement

Table 5 – Tested Substances for Interference

No interference detected

For swab samples	For urine samples
Cholesterol; 0.5 % (v/v)	LOPERAMID-ratiopharm acute 2 mg film-coated tablet; 0.2 % (v/v)
Mucus; 3 mg/mL	MAALOXAN 25 mVal Liquid; 0.1 mg/mL
Palmitic acid; 2mg/mL	METRONIDAZOL Aristo 400 mg tablet; 0.5 % (w/v)
Stearic Acid; 2mg/mL	Naproxen-ratiopharm® tablet; 0.5 % (w/v)
Triglyceride; 0.8 % (w/v)	Nonoxynol-9; 0.015 % (w/v)
Whole blood; 3.0 (v/v)	Nystatin STADA®; 1 % (w/v)
Barium Sulfate; 1.3 % (w/v)	Postericrot ointment; 0.5 % (v/v)
GLYCILAX for adults; 1.0 % (w/v)	Claversal 4 g/60 ml clysms; 0.5 % (v/v)
GRÜNWALDER Sennalax film-coated tablet; 0.1 mg/mL	Vancomycin; 1.4 mg/mL
Hemorrhoid ointment with witch hazel; 1 % (v/v)	Vaseline; 1 % (w/v)

Interference was experimentally verified 3xLoD for C. difficile using a spiking approach in eNAT® (COPAN) plus 1000 hPBTEC and 6 % (v/v) Stool Matrix Negative Control (artificial simulated matrix).

Table 6 – Clinical Sensitivity (PPA) [1] and Specificity (NPA) [2] for samples in eNAT® (95 % confidence interval, clinical study)

[1] 96.61 % (88.29 – 99.59 %)*

[2] 95.24 % (86.71 - 99.01 %)**

Clinical study within total 122 valid results, generated at two study sites and gained from 63 negative and 59 positive samples. Clinical performance data of Vivalytic C. difficile in comparison to the reference tests RIDA®GENE Clostridium difficile (R-Biopharm) on Roche LightCycler 480II device Discrepant results were verified using Allplex™ GI-Bacteria(I) Assay (5 samples) and Cepheid Xpert® C. difficile BT (5 samples).

*2 samples were detected as false negative.

** 3 samples were detected as false positive.

Table 7 – Document History

Revision 01	Initial document
Revision 02	Update including new analytical and clinical performance data

For more information see
www.bosch-vivalytic.com