SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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
mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

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Introduction

Bordetella pertussis is a small gram-negative, immobile, encapsulated, aerobic pathogen that produces a variety of toxins and virulence factors, including pertussis toxin (PT), filamentous haemagglutinin (FHA), tracheal cytotoxin, pertactin, heat-labile toxin, and adenylate cyclase haemolysin. *Bordetella* multiplies on the ciliated epithelium of the respiratory mucosa, causing local destruction of the mucosa. Some of the toxins also impair the immune system and cause tissue damage. [1, 2]

B. pertussis is the primary pathogen causing whooping cough. Infections with *B. parapertussis* or *B. holmesii* can also lead to a clinical picture similar to whooping cough, but the course is typically milder and shorter than that of a disease caused by *B. pertussis*.

Whooping cough is highly contagious and is transmitted through respiratory droplets generated by coughing, sneezing, and even talking.

Pertussis can last several weeks to months. The typical initial infection in unvaccinated individuals progresses in three stages:

- Catharral phase (duration 1-2 weeks; interval 5-21 days): This stage is characterized by cold-like symptoms, such as a cold and mild cough, but usually no or only moderate fever.
- Paroxysmal phase (duration 4-6 weeks): In this stage, there are the classic symptoms of bouts of coughing (staccato cough), followed by inspiratory dragging. The coughing attacks are often accompanied by regurgitation of viscous mucus followed by vomiting.
- Convalescent phase (duration 6-10 weeks): The coughing fits gradually subside.

Since pertussis often does not present with the classic symptoms, especially in adults, but also in older children and adolescents, laboratory diagnostic methods are crucial for diagnostics. In patients with cough (regardless of its duration), laboratory diagnostics for pertussis should be performed if there has been contact with a confirmed pertussis case or if classic symptoms such as cough attacks, inspiratory stridor or vomiting after the cough attacks are present; or in the case of a longer persistent cough (>14 days duration), even if these are absent. Patients who have received a pertussis vaccine are not excluded. Even if classical symptoms occur laboratory diagnostics should be performed for confirmation. The laboratory diagnostic method recommended depends on the stage of the disease. During the first 2 to 3 weeks after the onset of coughing, it is strongly advised to use a culture or nucleic acid amplification technique (NAT), usually PCR, to detect these pathogens from deep nasopharyngeal swabs, nasopharyngeal secretions, or material obtained by suction. [1, 2]

Package Contents

15 Vivalytic Bordetella test cartridges for the qualitative detection of *B. pertussis*, *B. parapertussis* and *B. holmesii* specific gene sequences in eNAT® transport medium (eNAT® 608C, COPAN Italia s.p.a.).

The Vivalytic Bordetella test is configured for use with a Vivalytic *one* analyser.

Intended Use

The Vivalytic Bordetella test is an automated, qualitative *in vitro* diagnostic test based on real-time polymerase chain reaction (PCR) for the detection of nucleic acids from *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* from human nasopharyngeal swabs to aid in the diagnosis of pertussis, parapertussis or pertussis like infections of symptomatic individuals.

Results should not be used as the sole basis for diagnosis, treatment, or other patient management decisions. Positive results do not exclude co-infection with other pathogens. The agent(s) detected may not be the definite cause of disease. Negative results do not exclude a *Bordetella* infection or another respiratory tract infection. Results must be clinically correlated with patient history, clinical observations, and epidemiological information. Other diagnostic information is necessary to determine patient infection status. Intended for use with a Vivalytic *one* analyser by healthcare professionals only in laboratory settings such as hospital laboratories and reference laboratories.

Vivalytic Bordetella – Instructions for Use

Safety Information

These Instructions for Use contain test-specific information only. For additional warnings and instructions refer to the Instructions for Use provided with your Vivalytic *one* analyser (chapter device safety information). Only use Vivalytic cartridges and accessories approved for the Vivalytic *one* analyser. Take care to avoid any contamination when handling patient samples and cartridges. When sample was spilled on the cartridge, do not use the cartridge, and dispose it.

For *in vitro* diagnostic use by trained healthcare professional.

WARNING

- Always follow good laboratory practice to ensure the proper performance of this test.
- Make sure to wear appropriate personal protective equipment (PPE).
- Do not use a cartridge if the sealed pouch or the cartridge itself is visibly damaged.
- Do not touch or scratch the detection area of the cartridge.
- Do not reuse a cartridge.
- Do not use expired cartridges. The expiration date can be found on the packaging and the cartridge label.
- Do not wait longer than 15 minutes after opening the cartridge pouch to begin the test. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.
- Do not shake a cartridge that contains a sample.
- Do not turn the cartridge upside down.
- Place the cartridge on a clean and flat surface only.
- Do not use sample types, media and volumes that are not approved for the test.
- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions. Handle potentially infectious patient samples and cartridges according to national laboratory standards and dispose samples and cartridges according to regional and laboratory standards.
- Be compliant with the national safety regulations and practices.

Note: Further information can be found in the safety data sheet (SDS) of the product. Please contact the customer support of your local distributor.

Additional Equipment & Consumables Required but not Provided

- Bosch Vivalytic *one* analyser (reference number F09G300115)
- Pipettor (100–1000 µl)
- Sterile filter pipette tips 100–1000 µl
- Swab collection kits:
Flexible Flocked Swab FLOQSwabs® 553C (COPAN Italia s.p.a.)
- Transport medium:
eNAT® 608C (COPAN Italia s.p.a.)
- Suitable protective clothing

Test Principle

Vivalytic Bordetella is a qualitative real-time PCR based test.

Storage und Usage Conditions

Product is stable until the expiry date if stored at +15 °C to +25 °C. Storage and usage conditions can be taken either from the cartridge, pouch, or box label. Cartridge has to be used at +15 °C to +25 °C, relative humidity < 65 %, within 15 min upon pouch opening. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.

Reagents

All reagents necessary for the sample processing are integrated into the cartridge. The processing includes cell lysis, nucleic acid extraction, DNA amplification and detection.

Reagents are PCR bead, binding buffer, washing buffer, and elution buffer. The PCR bead contains the DNA polymerase, primers, and probes. Binding buffer facilitates binding of nucleic acids during the purification process. Washing buffer is a formulation of different salts and solvents to remove impurities e. g. proteins during the extraction process. Elution Buffer is a low-salt buffer and contains the purified nucleic acids at the end of the extraction process.

Sample Types / Medium

The test is intended to detect nucleic acids of *B. pertussis* (gene target sequence IS481), *B. parapertussis* (gene target sequence pIS1001) and *B. holmesii* (gene target sequences IS481, hIS1001) of nasopharyngeal swab samples (flexible Flocked Swab FLOQSwabs® 553C) in eNAT® transport medium (eNAT® transport medium Ref. 608C, COPAN Italia s.p.a.).

Collect and store samples as indicated in the manufacturer's instructions [3]. In case the sample is not processed immediately after sample collection, nucleic acids will be preserved in eNAT® transport medium for up to 4 weeks at room temperature and at 4°C and up to 6 months at 4 °C and up to 6 months at –20 °C to –80 °C (see Instructions for Use eNAT® transport medium [3]).

Sample Preparation

Use a COPAN FLOQSwabs® flexible swab for collection of nasopharyngeal swab samples. Each swab is transferred into 1 ml eNAT® transport medium. Break the swab and close the tube. Shake the sample tube containing the swab sample in eNAT® medium for homogenization (avoid foaming). Fill 300 µl of homogenized sample in the sample input of the cartridge and close the cartridge lid. Insert the cartridge into the Vivalytic *one* analyser. Do not use viscous samples that are difficult to pipette.

Test Result

After automatic processing of the sample with the Vivalytic *one* analyser the test result is shown on the screen of the Vivalytic *one* analyser. The time to result is about 50 min. For high titer specimens results are available after less than 33 min and the test run can be terminated earlier (see chapter Test Termination).

The sample is classified either as *B. pertussis* positive, *B. parapertussis* positive, *B. holmesii* positive, negative or invalid or as a combination of the aforementioned. In case of a positive detection of one or more targets the test is considered valid even if the Human Control is negative.

Detection of the human cell based whole process control (Human Control) in negative samples shows a successful DNA purification and amplification procedure and excludes an inhibition of the PCR reaction. Interpretation of results is shown in the table below.

Vivalytic Bordetella – Instructions for Use

IS481 ⁵	hIS1001	pIS1001 ⁵	Human Control	Result
+	-	-	+/-	Sample is considered positive for <i>B. pertussis</i> .
+	+	-	+/-	Sample is considered positive for <i>B. holmesii</i> . ¹
-	-	+	+/-	Sample is considered positive for <i>B. parapertussis</i> .
+	-	+	+/-	Sample is considered positive for <i>B. pertussis</i> and <i>B. parapertussis</i> .
+	+	+	+/-	Sample is considered positive for <i>B. parapertussis</i> and <i>B. holmesii</i> . ¹
-	+	+	+/-	Sample is considered positive for <i>B. parapertussis</i> (!). ²
-	+	-	+/-	Sample is considered negative (!). ³
-	-	-	+	Sample is considered negative.
-	-	-	-	Invalid. ⁴

¹In the case of a positive signal for *B. holmesii*, a mixed sample with *B. pertussis* cannot be excluded as the two species share identical sequences within the amplified region of IS481.

²*B. parapertussis* positive (!)

Samples close to LOD might not show a positive signal for all defined targets and are therefore classified as negative for *B. pertussis* and *B. holmesii*. That does not exclude an infection with these pathogens. Retesting is advised.

³Negative (!)

Samples close to LOD might not show a positive signal for all defined targets and are therefore classified as negative for *B. pertussis*, *B. parapertussis* and *B. holmesii*. That does not exclude an infection with these pathogens. The result must not be used. Retesting is advised.

⁴Retesting is recommended.

⁵The IS481 and pIS1001 insertional elements can also be present in *B. bronchiseptica*.

PCR – Curve and C_q Value

Real-time PCR curves (software-modified) are shown and classified as positive or negative by the software. In case of positive curves, the respective C_q value is displayed. Inconclusive results are marked by the software (△). Retesting is advised.

Invalid or Failed Tests

A test is rated as invalid if neither target DNA nor Human Control is detected. Possible reasons for an invalid run might be poor sample quality due to a partial or complete absence of human cellular material in the sample. Pay attention to use the correct sample type, sample collection and storage of the sample and cartridges prior to the test run. If required, repeat the analysis with a new sample.

In case of a failed test, first check for correct operating conditions of the Vivalytic *one* analyser (refer to Vivalytic *one* analyser's Instructions for Use). Restart the Vivalytic *one* analyser. If the problem persists, contact the customer support of your local distributor.

Test Termination (Early Finish)

As soon as a valid, positive result for *B. paraptussis* or *B. holmesii* is shown on the screen, the user has the option to finish the test. Only finish the test if no further targets need to be analyzed. Pending results will not be displayed.

Test Report

In the printed test report, pathogen, results, control and information on user, patient sample and Vivalytic *one* analyser are listed with a signature field.

Quality Control

If required by your local or laboratory standards, quality control testing must be performed. The integrated control (Human Control) covers the processes of sample collection, nucleic acid purification, amplification, and detection.

In case a reference testing method for quality control is needed the following reference materials can be used:

AMPLIRUN® BORDETELLA PERTUSSIS TOTAL CONTROL (NEUTRAL MATRIX) (MC078, Vircell S.L. Spain).

AMPLIRUN® BORDETELLA PARAPTUSSIS TOTAL CONTROL (NEUTRAL MATRIX) (MC077, Vircell S.L. Spain).

AMPLIRUN® BORDETELLA HOLMESII TOTAL CONTROL (NEUTRAL MATRIX) (MC168, Vircell S.L. Spain).

Please follow the manufacturer's instructions. Reconstitute the reference material in Copan eNAT® medium.

In case of unexpected results, repeat the analysis with another sample. If the result of a negative quality control sample e. g. pure eNAT® medium remains positive, the Vivalytic *one* analyser or its environment might be contaminated. Stop using the Vivalytic *one* analyser and clean the device as described in the Instructions for Use of the Vivalytic *one* analyser.

In case of repeated false results for quality control samples please contact the customer support of your local distributor.

Notice to Users in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established. The summary of safety and performance (SSP) will be available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI. If the respective module in Eudamed will be fully functional the SSP will be available at the public website (<https://ec.europa.eu/tools/eudamed>).

Limitations

The results of the Vivalytic Bordetella test must be interpreted by a trained healthcare professional only. The results of the Vivalytic Bordetella test must not be used as the sole parameter for diagnosis.

- A negative result does not exclude pathogens being present in the sample at a level below assay sensitivity or other pathogens being present not covered by this assay.
- There is a risk of false negative or false positive results due to improperly collected, transported, or handled samples.
- In borderline cases atypical PCR characteristics (e. g. flat curve with low or high C_q-value) can occur. In case of atypical characteristics results are not allowed to be used for diagnostic interpretation. Retesting is advised.
- A positive result does not necessarily mean that viable pathogens are present.
- Vivalytic Bordetella is a qualitative real-time PCR test and does not provide a quantitative result.
- A positive IS481 and HIS1001 can either be caused by the presence of *B. holmesii* or by presence of *B. pertussis* and *B. holmesii* in the sample.
- A positive IS481 can be considered as probable *B. pertussis* positive, although *B. bronchiseptica* may also lead to a positive IS481 signal. If needed, additional testing should be performed to differentiate *B. bronchiseptica* and *B. pertussis*.

Vivalytic Bordetella – Instructions for Use

Analytical Performance Evaluation

Analytical Sensitivity (Limit of Detection)

The limit of detection (LoD) of the Vivalytic Bordetella test was determined as the lowest concentration of analyte that can be consistently detected ($\geq 95\%$ of samples tested under routine laboratory conditions using a defined type of sample, [Table 1](#)).

Inclusivity

To evaluate inclusivity, an *in silico* analysis (BLAST alignment) of the genomic sequence of various relevant *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* strains against the sequence of the PCR primers and hydrolysis probe used in the Vivalytic Bordetella test for amplification and detection of the respective pathogens was performed. Inclusivity could be shown for organisms listed in [Table 2](#).

Exclusivity / Analytical Specificity

To exclude cross-reactivity (exclusivity), an *in silico* analysis (BLAST alignment) of the target region of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* against the genomic sequence of various other pathogens representing common respiratory pathogens or closely related species was conducted. For the detection system of *Bordetella pertussis*, sequence matches in the probe and primer area could be detected for the *Bordetella bronchiseptica* strains NCTC10543, MBORD731, NCTC8344, E001, MBORD595, concluding a possible amplification. For the detection system of *Bordetella parapertussis*, sequence matches in the probe and primer area could be detected for the *Bordetella bronchiseptica* strains KM22, BB5, FR3523, FR2017, FR3823, concluding a possible amplification. There was no evidence of an interference for the detection system of *Bordetella holmesii* ([Table 3](#)).

Reproducibility

The reproducibility of the Vivalytic Bordetella test was established using a panel with 3 different concentrations of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii*. At 3 test sites, each mix was tested on the same set of Vivalytic instruments by the same operator with 3 LOTs in 4 replicates on 3 days, respectively. The obtained positivity rates for the different combinations were correlated to the expected positivity rate ([Table 4a](#)).

Repeatability

The repeatability of the Vivalytic Bordetella test was established using a panel with 1 concentration (3x c95) of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii*. At 1 test site, the mix was tested on the same set of Vivalytic instruments by the same operator with 3 LOTs in 20 replicates, respectively, yielding in a total of 60 observations per target pathogen. The obtained positivity rates for the different combinations were correlated to the expected positivity rate ([Table 4b](#)).

Interferences

Interferences were evaluated for endogenous and exogenous substances, that are potentially present in the patient sample. Refer to [Table 5](#) for substances that have the potential to interfere with the test. No interferences were detected.

Clinical Performance Evaluation

Sensitivity and specificity results derived from native nasopharyngeal swab samples. Samples were collected in a clinical setting and compared with results of reference methods.

561 nasopharyngeal swab samples in eNAT were tested at one study site using a screening approach with potentially clinical positive samples. 8 clinical samples were found positive for one of the three pathogens (all *Bordetella parapertussis*) leading to a positive rate within the sample set of 1.9% (8/411). Additionally, 150 negatively pre-characterized samples were spiked with respective reference material (50 sample per pathogen).

In total, 495 results were included in the data set gained from testing with a reference test (Bordetella Speciation Plus Toxin - OSR for BD MAX™) and Vivalytic Bordetella, leading to the performance results shown in [Table 6](#).



Technical Support

If you require any support, technical help or have additional questions, please contact your local distributor or visit the Bosch Vivalytic website at www.bosch-vivalytic.com.

References

- [1] https://www.rki.de/DE/Content/Infekt/EpidBull/Merkblaetter/Ratgeber_Pertussis.html (status 30.11.2023)
- [2] van der Zee A, Schellekens JF, Mooi FR. Laboratory Diagnosis of Pertussis. Clin Microbiol Rev. 2015 Oct; 28(4): 1005-26. doi: 10.1128/CMR.00031-15. PMID: 26354823; PMCID: PMC4575397.
- [3] COPAN Instructions for Use of Copan eNAT® Collection and Preservation System; Version EIFU002R01 Date 2021.06.

Symbols

- | | | | |
|--|---------------------|---|---|
|  | Manufacturer |  | Serial number |
|  | Date of manufacture |  | Temperature limit |
|  | Expiry date |  | Do not use if package is damaged |
|  | Lot number |  | For single use only |
|  | Reference number |  | Consult Instructions for Use |
|  | Contains <n> tests |  | <i>in vitro</i> diagnostic medical device |
-  ₀₁₂₃ CE mark incl. identification number of Notified Body
-  Pipette the indicated sample volume in the sample input of the cartridge as marked by the black triangle.

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WARNING**Hazard component in cartridge**

(binding buffer RB-BB-34a)

Guanidinium chloride; guanidine hydrochloride

Hazard statements

H302	Harmful if swallowed.
H315	Causes skin irritation.
H319	Causes serious eye irritation.

Precautionary statements

P280	Wear protective gloves and eye/face protection.
P301+P312	IF SWALLOWED: Call a POISON CENTER/ doctor if you feel unwell.
P302+P352	IF ON SKIN: Wash with plenty of water and soap.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.
P501	Dispose of contents/container in accordance with local/ regional/ national/ international regulation

Table 1 – Limit of Detection* (LoD)

<i>B. pertussis</i>	200 genome copies/mL
<i>B. parapertussis</i>	200 genome copies/mL
<i>B. holmesii</i>	200 genome copies/mL

*Based on determination by using Vircell material (MC078, MC077, MC168) for a spiking approach in eNAT® (COPAN) supplemented with approximately 10.000 Human Bronchial/Tracheal Epithelial Cells (hPBTEC).

Table 2 – Inclusivity

<i>Bordetella pertussis</i>
<i>Bordetella parapertussis</i>
<i>Bordetella holmesii</i>

Strains validated for inclusivity via BLAST alignment.

Table 3 – Exclusivity

<i>Acinetobacter</i>	<i>Lactobacillus acidophilus</i>
<i>Schaalia odontolytica</i> (<i>Actinomyces odontolyticus</i>)	<i>Lactiplantibacillus plantarum</i>
Mastadenovirus	<i>Legionella</i>
<i>Arcanobacterium haemolyticum</i>	<i>Listeria monocytogenes</i>
<i>Bacillus cereus</i>	Measles morbillivirus
<i>Bacillus subtilis</i>	Metapneumovirus
<i>Bacteroides fragilis</i>	<i>Moraxella catarrhalis</i>
Bocaparvovirus	<i>Morganella morganii</i>
<i>Bordetella avium</i>	Mumps orthorubulavirus
<i>Bordetella bronchiseptica</i>	<i>Mycobacterium avium</i>
<i>Bordetella hinzii</i>	<i>Mycobacterium tuberculosis</i>
<i>Bordetella petrii</i>	Metamycoplasmataceae
<i>Bordetella trematum</i>	<i>Neisseria</i>
<i>Burkholderia</i>	<i>Nocardia asteroides</i>
<i>Candida albicans</i>	<i>Oligella ureolytica</i>
<i>Nakaseomyces glabratus</i> (<i>Candida glabrata</i>)	<i>Pandoraea oxalativorans</i>
<i>Chlamydia pneumoniae</i>	Human respirovirus 1
<i>Chlamydia trachomatis</i>	Human orthorubulavirus 2
<i>Citrobacter freundii</i>	Human respirovirus 3
<i>Clostridioides difficile</i>	Human orthorubulavirus 4
<i>Orthocoronavirinae</i>	<i>Parvimonas micra</i>
<i>Corynebacterium diphtheriae</i>	Parechovirus A
<i>Corynebacterium pseudodiphtheriticum</i>	<i>Peptostreptococcus anaerobius</i>
Cytomegalovirus	<i>Pneumocystis jirovecii</i>
<i>Enterobacter cloacae</i>	<i>Proteus mirabilis</i>
<i>Enterococcus faecalis</i>	<i>Proteus vulgaris</i>
Enterovirus	<i>Pseudomonas aeruginosa</i>
Human gammaherpesvirus 4 (Eppstein-Barr Virus)	<i>Pseudomonas fluorescens</i>
<i>Escherichia coli</i>	<i>Cupriavidus pauculus</i>
<i>Fusobacterium necrophorum</i>	HRSV-A
<i>Fusobacterium nucleatum</i>	HRSV-B
<i>Haemophilus influenzae</i>	<i>Serratia liquefaciens</i>
<i>Haemophilus parainfluenzae</i>	<i>Serratia marcescens</i>
Human alphaherpesvirus 1	<i>Staphylococcus</i>
Human alphaherpesvirus 2	<i>Stenotrophomonas maltophilia</i>
Influenza A	<i>Streptococcus</i>
Influenza B	<i>Ureaplasma urealyticum</i>
<i>Klebsiella</i>	Human alphaherpesvirus 3

Table 4a Reproducibility							
Target	c	No of total tests	No of pos. tests	No of neg. tests	Proportion of positive/negative results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
none	blank	108	0	108	100	96.57 % – 100 %	96.64 % – 100 %
Bordetella pertussis	3x c95	108	108	0	100	96.57 % – 100 %	96.64 % – 100 %
	c100	108	107	1	99.1	94.94 % – 99.84 %	94.95 % – 99.98 %
Bordetella parapertussis	3x c95	108	108	0	100	96.57 % – 100 %	96.64 % – 100 %
	c100	108	108	0	100	96.57 % – 100 %	96.64 % – 100 %
Bordetella holmesii	3x c95	108	104	4	96.3	90.86 % – 98.55 %	90.79 % – 98.98 %
	c100	108	107	1	99.1	94.94 % – 99.84 %	94.95 % – 99.98 %

3x c95 = 95 % predicted positive agreement

c100 = 100 % predicted positive agreement

Table 4 b Repeatability						
Target	LOT	No of total tests	No of pos. tests (3x c95)	Proportion of positive/negative results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
<i>Bordetella pertussis</i>	1	20	20	100.0	83.89 % – 100 %	83.16 % – 100 %
	2	20	20	100.0	83.89 % – 100 %	83.16 % – 100 %
	3	20	19	95.0	76.39 % – 99.11 %	75.13 % – 99.87 %
	Total	60	59	98.3	91.14 % – 99.71 %	91.06 % – 99.96 %
<i>Bordetella para-pertussis</i>	1	20	20	100.0	83.89 % – 100 %	83.16 % – 100 %
	2	20	19	95.0	76.39 % – 99.11 %	75.13 % – 99.87 %
	3	20	20	100.0	83.89 % – 100 %	83.16 % – 100 %
	Total	60	59	98.3	91.14 % – 99.71 %	91.06 % – 99.96 %
<i>Bordetella holmesii</i>	1	20	20	100.0	83.89 % – 100 %	83.16 % – 100 %
	2	20	19	95.0	76.39 % – 99.11 %	75.13 % – 99.87 %
	3	20	20	100.0	83.89 % – 100 %	83.16 % – 100 %
	Total	60	59	98.3	91.14 % – 99.71 %	91.06 % – 99.96 %

3x c95 = 95 % predicted positive agreement

Table 5 – Tested Substances for Interference

EDTA whole blood; 2 % (v/v)

Aspirin 500 mg; 5 mg/mL

Azithromycin HEXAL® 500 mg; 5 mg/mL

Budesonid Easyhaler® 0.2 mg/Dosis; 25 mg/mL

SILOMAT DMP pastille; 0.1 mg/mL

DOLO-DOBENDAN 1.4 mg/10 mg throat lozenge; 2.5 % (w/v)

Em-eukal classic bonbons; 0.26 % (w/v)

IBU ratiopharm 400 mg akut analgesic tablet; 10 mg/mL

Mometason-ratiopharm® hay fever spray; 2.5 mg/mL

Nasivin® Spray; 10 % (v/v)

Otri-Allergie Spray Fluticason; 5 % (v/v)

Paracetamol ratiopharm 500 mg; 10 mg/mL

Rhinomer nasal spray; 15 % (v/v)

Salbuhexal N 200 puff dosing aerosol; 0.1 % (v/v)

Tamiflu® 75 mg; 3.75 mg/mL

TobraZid®; 4 µg/mL

Turixin®; 10 mg/mL

Vaseline; 1 % (w/v)

DOREGRIPPIN tablet; 0.075 mg/mL

Wick cough solvent sirup; 5 % (v/v)

Interference was experimentally verified at 3x c95 for *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* using a spiking approach in eNAT® (COPAN) supplemented with approximately 10.000 Human Bronchial/Tracheal Epithelial Cells (hPBTEC).

Table 6 – Clinical Sensitivity (PPA) and Specificity (NPA) for samples in eNAT® (95 % confidence interval, clinical study)

Pathogen	Sensitivity (PPA) (95%-width of CI)	Specificity (NPA) (95%-width of CI)
<i>Bordetella pertussis</i>	94.0 % (83.45 – 98.75 %)	100.0 % (96.84 – 100.0 %)
<i>Bordetella parapertussis</i>	100.0 % (93.51 – 100.0 %)	95.5 % (89.71 – 98.51 %)
<i>Bordetella holmesii</i>	100.0 % (90.97 – 100.0 %)	97.6 % (93.20 – 99.51 %)
Total	97.9 % (94.03 – 99.57 %)	97.7 % (95.56 – 99.01 %)

Table 7 – Document History

Revision 01	Initial document
Revision 02/03	Pre-release changes

For more information see
www.bosch-vivalytic.com