SZABO SCANDIC

Part of Europa Biosite

Produktinformation

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mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Latest change 2024-08-16 – Rev. 02

BOSCH

For more information see
www.bosch-vivalytic.com

vivalytic

C. auris
PCR Test for detection
of *Candida auris*

Instructions for use

CH REP Decomplex AG
Freiburgstrasse 3
3010 Bern
Switzerland

 Bosch Healthcare Solutions GmbH
Stuttgarter Straße 130
71332 Waiblingen, Germany

REF

F 09G 301 032

F09G301047_02

IVD

CE 0123

Table of contents

English	2
Deutsch	8
Nederlands	14
Français	20
Italiano	26
Norsk	32
Svenska	38
Suomi	44
Dansk	50
Eesti keel	56
Lietuvių	62
Latviešu	68
Polski	74
Čeština	80
Slovenský	86
Magyar	92
Română	98
Hrvatski	104
Slovenščina	110
Български	116
Ελληνικά	122
Español	128
Português	134
Annex	140

Introduction

Candida auris (*C. auris*) is a pathogenic fungus belonging to the *Candida* genus. It was first discovered in Japan in 2009 and has spread globally since then. *C. auris* is a commensal skin inhabitant frequently found in axilla, groin or nasal swab samples.¹ The fungus is notorious for its resistance to many common antifungal drugs, making the treatment of infections challenging. *C. auris* can survive on artificial surfaces and has a high resilience against disinfection strategies, allowing it to cause outbreaks and persist in hospitals and other medical facilities.² The clinical symptoms caused by *C. auris* can vary depending on the site of invasive infection. It can lead to infections of the bloodstream, respiratory tract, urinary tract, and other areas of the body. Symptoms may include fever, chills, fatigue, and respiratory difficulties. The clinical relevance of *C. auris* lies in its ability to cause critical infections in immunocompromised individuals. This particularly affects patients in hospitals and nursing homes.^{3,4}

Package Contents

15 Vivalytic *C. auris* test cartridges for the detection of *Candida auris* (ITS2 gene).

Pathogen List
<i>Candida auris</i> (ITS2 gene)

Intended Use

The Vivalytic *C. auris* test is an automated qualitative *in vitro* diagnostic test based on real-time PCR for the detection of *Candida auris* DNA from human bilateral axilla, groin and nasal combination swabs and nasal swabs to screen individuals at risk for colonization and to aid in the prevention and control of *Candida auris* outbreaks and infections in healthcare settings. Negative results do not exclude a *Candida auris* colonization. Results should not be used as the sole basis for treatment, or other patient management decisions. Results must be clinically correlated with patient history, clinical observations, and epidemiological information. Intended for use with a Vivalytic *one* analyser by healthcare professionals only in a laboratory setting.

Safety Information

These Instructions for Use contain test-specific information only. For additional warnings and instructions refer to the Instructions for Use provided with your Vivalytic *one* analyser (chapter device safety information). Only use Vivalytic cartridges and accessories approved for the Vivalytic *one* analyser. Take care to avoid any contamination when handling patient samples and cartridges. When sample was spilled on the cartridge, do not use the cartridge, and dispose it. For *in vitro* diagnostic use by trained healthcare professional.

Vivalytic *C. auris* – Instructions for use

WARNING

- Always follow good laboratory practice to ensure the proper performance of this test.
- Make sure to wear appropriate personal protective equipment (PPE).
- Do not use a cartridge if the sealed pouch or the cartridge itself is visibly damaged.
- Do not touch or scratch the detection area of the cartridge.
- Do not reuse a cartridge.
- Do not use expired cartridges. The expiration date can be found on the packaging and the cartridge label.
- Do not wait longer than 15 minutes after opening the cartridge pouch to begin the test. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.
- Do not shake a cartridge that contains a sample.
- Do not turn the cartridge upside down.
- Place the cartridge on a clean and flat surface only.
- Do not use sample types, media and volumes that are not approved for the test.
- Biological specimens, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions. Handle potentially infectious patient samples and cartridges according to national laboratory standards and dispose samples and cartridges according to regional and laboratory standards.
- Be compliant with the national safety regulations and practices.

Note: Further information can be found in the material safety data sheet (MSDS) of the product. Please contact the customer support of your local distributor.

Additional Equipment & Consumables Required but not Provided

- Bosch Vivalytic *one* analyser (reference number F09G300115)
- Pipette (100–1000 µl)
- Sterile filter pipette tips 100 –1000 µl
- 2x Copan Liquid Amies Elution Swab (eSwab®) Collection and Transport System (Ref: 480CE, 1 ml liquid Amies with Regular Flocked Swab FLO-QSwabs®)

Test Principle

Vivalytic *C. auris* is a qualitative real-time PCR based test.

Storage und Usage Conditions

Product is stable until the expiry date if stored at +15 °C to +25 °C. Storage and usage conditions can be taken either from the cartridge, pouch, or box label. Cartridge has to be used at +15 °C to +25 °C, relative humidity <65 %, within 15 min upon pouch opening. This maintains hygiene and avoids performance loss due to humidity. Prolonged exposure to humidity has a negative impact on test performance.

Reagents

All reagents necessary for the sample processing are integrated into the cartridge. The processing includes cell lysis, nucleic acid extraction, DNA amplification and detection.

Reagents are PCR bead, binding buffer, washing buffer, and elution buffer. The PCR bead contains the DNA polymerase, primers, and probes. Binding buffer facilitates binding of nucleic acids during the purification process. Washing buffer is a formulation of different salts and solvents to remove impurities, e. g. proteins during the extraction process. Elution Buffer is a low-salt buffer and contains the purified nucleic acids at the end of the extraction process.

Sample Type/Medium

The test is intended for use with a human bilateral axilla, groin and nasal combination swab or nasal swabs in liquid Amies medium (Copan Liquid Amies Elution Swab (eSwab®) Collection and Transport System Ref. 480CE, COPAN Italia s.p.a.).

Vivalytic *C. auris* – Instructions for use

Collect and store samples as indicated in the manufacturer's instructions.

In case the sample is not processed immediately after sample collection, nucleic acids will be preserved in liquid Amies transport medium for up to 5 days at room temperature, 7 days at 4 °C and up to 6 months at –20 °C (see Instructions for Use eSwab®⁵).

Sample Preparation

Use two Regular Flocked Swabs (FLOQSwabs® COPAN Italia s.p.a.) to collect human samples separately in 1 ml liquid Amies. The first swab is used for the bilateral axilla and groin sample and the second swab is used for the nasal sample.

Sample collection procedure description using bilateral axilla and groin swabs

Firmly rub swab tip across left axilla skin surface for at least 3 to 5 times. Repeat procedure with the same swab on the right axilla skin surface.

With the same swab used on the axilla, firmly rub the swab tip across the left groin skin surface for at least 3 to 5 times. Repeat the procedure on the right side. Transfer the swab into 1 ml liquid Amies transport medium, "break the swab" and close the tube.

Sample collection procedure description using nasal swabs

Insert a swab into the nostril parallel to the palate. Leave the swab in place for a few seconds to absorb secretions. Rotate the swab tip against the nasal mucosa 3 to 5 times. Transfer the swab into 1 ml liquid Amies transport medium, "break the swab" and close the tube. Transfer 300 µl with a pipette into the cartridge sample inlet.

Sample collection and procedure description using bilateral axilla, groin and nasal combination swabs

Take the bilateral axilla and groin swab and the nasal swab as mentioned above.

After transferring each swab to 1 ml liquid Amies transport medium gently mix the samples for homogenization (avoid foaming). Take 150 µl with a pipette from either the nasal sample or the bilateral axilla and groin sample and transfer it into the cartridge sample inlet. Again take 150 µl from the matching sample of the remaining sample type and transfer it into the cartridge sample inlet. Gently mix the sample in the cartridge sample inlet by using the pipette (avoid foaming).

Do not use viscous samples that are difficult to pipette.

Test Result

After automatic processing of the sample with the Vivalytic *one* analyser the test result is shown on the screen of the Vivalytic *one* analyser. The time to result is about 60 min. For high titer specimens results are available after less than 35 min and the test run can be terminated earlier (see Chapter Test Termination).

The sample is classified either as *C. auris* positive, *C. auris* negative or invalid. In case of a positive detection of the pathogen, the test is considered valid even if the Human Control is negative.

Detection of the human cell based whole process control (Human Control) in negative samples shows a successful extraction procedure and excludes an inhibition of the PCR reaction. Interpretation of results is shown in the table below.

<i>C. auris</i>	Human Control	Validity	Result
+	+/-	valid	Sample is considered <i>C. auris</i> positive.
-	+	valid	Sample is considered <i>C. auris</i> negative.
-	-	invalid	Not evaluable. ¹

¹Retesting is recommended

Vivalytic *C. auris* – Instructions for use

PCR – Curve and C_q Value

Real-time PCR curves (software-modified) are shown and classified as positive or negative by the software. In case of positive curves, the respective C_q value is displayed. Inconclusive results are marked by the software (Δ). Retesting is advised.

Invalid or Failed Tests

A test is rated as invalid if neither target DNA nor Human Control is detected. Possible reasons for an invalid run might be poor sample quality due to a partial or complete absence of human cellular material in the sample. Results are displayed for an invalid test but are not allowed to be used for diagnostic interpretation. Pay attention to use the correct sample type, sample collection and storage of the sample and cartridges prior to the test run. If required, repeat the analysis with a new sample.

In case of a failed test, first check for correct operating conditions of the Vivalytic *one* analyser (refer to Vivalytic *one* analyser's Instructions for Use). Restart the Vivalytic *one* analyser. If the problem persists, contact the customer support of your local distributor.

Test Termination

As soon as a valid, positive result is shown for the pathogen on the screen, the user has the option to finish the test.

Test Report

In the printed test report, pathogen, results, control and information on user, patient and Vivalytic *one* analyser are listed with a signature field.

Notice to Users in EU

Any serious incident that has occurred in relation to the device, should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Limitations

The result of the Vivalytic *C. auris* test must be interpreted by a trained healthcare professional only. The results of the Vivalytic *C. auris* test must not be used as the sole parameter for diagnosis.

- A negative result does not exclude pathogen being present in the sample at a level below assay sensitivity or other pathogens being present but not covered by this assay.
- There is a risk of false negative or false positive results due to improperly collected, transported, or handled samples.
- In borderline cases atypical PCR characteristics (e. g. flat curve with low or high C_q -value) can occur. In case of atypical characteristics results are not allowed to be used for diagnostic interpretation. Inconclusive results are marked by the software. Retesting is advised.
- Vivalytic *C. auris* is a qualitative real-time PCR test and does not provide a quantitative result.
- A positive result does not necessarily mean that a viable pathogen is present.
- In artificial sample setups an inhibitory effect of aluminium chlorohydrate was observed on PCR as ingredient in deodorants. This was not observed in clinical samples but could still cause invalid results.

Analytical Performance Evaluation

Analytical Sensitivity (Limit of Detection)

The limit of detection of the Vivalytic *C. auris* test was determined as the lowest concentration of analyte that can be consistently detected ($\geq 95\%$ of samples tested under routine laboratory conditions using a defined type of sample) ([Table 1](#)).

Inclusivity

To evaluate inclusivity, an *in silico* analysis (BLAST alignment) of the genomic sequence of various relevant *Candida auris* strains against the sequence of the PCR primers and hydrolysis probe used in the Vivalytic *C. auris* test for amplification and detection of the respective pathogen was performed. Inclusivity could be shown for strains listed in [Table 2](#).

Exclusivity / Analytical Specificity

To exclude cross-reactivity (exclusivity), an *in silico* analysis (BLAST alignment) of the target region of *C. auris* against the genomic sequence of various other pathogens representing common commensal microorganisms or closely related species was conducted. There was no evidence of an interference for the detection system of *C. auris* ([Table 3](#)).

Reproducibility

The reproducibility of the Vivalytic *C. auris* test was established using a panel with 3 different concentrations of *Candida auris*. At 3 test sites, each mix was tested on the same set of Vivalytic instruments by the same operator with 3 LOTs in 4 replicates on 3 days, respectively. The obtained positivity rates for the different combinations were correlated to the expected positivity rate ([Table 4a](#)).

Repeatability

The repeatability of the Vivalytic *C. auris* test was established using a panel with 1 concentration (3x c95) of *Candida auris*. At 1 test site, the mix was tested on the same set of Vivalytic instruments by the same operator with 3 LOTs in 20 replicates, respectively, yielding a total of 60 observations per target pathogen. The obtained positivity rates for the different combinations were correlated to the expected positivity rate ([Table 4b](#)).

Interferences

Interferences were evaluated for endogenous and exogenous substances that are potentially present in the patient sample. Refer to [Table 5](#) for substances that have the potential to interfere with the test.

Clinical Performance evaluation Sensitivity and Specificity

In total, 231 samples with valid results were tested (126 clinical samples, 50 spiked samples and 55 isolates) at two study sites. All samples were prepared according to recommendations of used methods. Sensitivity or Positive Percent Agreement (PPA) was calculated as $100\% \times TP / (TP+FN)$. Specificity or Negative Percent Agreement was calculated as $100\% \times TN / (TN+FP)$. The results of the clinical performance evaluation are shown in [Table 6a](#) and [Table 6b](#).

Technical Support


If you require any support, technical help or have additional questions, please contact your local distributor or visit the Bosch Vivalytic website at www.bosch-vivalytic.com.

References

- ¹ Fernandes L, Ribeiro R, Henriques M, Rodrigues ME. *Candida auris*, a singular emergent pathogenic yeast: Its resistance and," *Eur J Clin Microbiol Infect Dis*, Nr. 12, pp. 1371-1385, 2022;12;1361-1385;
- ² Preda M, Chivu RD, Ditu LM, Popescu O, Manolescu LSC. Pathogenesis, Prophylaxis and Treatment of *Candida auris*. *Biomedicines*. 2024;12(3):561. doi: 10.3390/biomedicines12030561.
- ³ Cristina ML, Spagnolo AM, Sartini M, Carbone A, Oliva M, Schinca E, Boni S, Pontali E. An Overview on *Candida auris* in Healthcare Settings. *J Fungi (Basel)*. 2023 Sep 8;9(9):913. doi: 10.3390/jof9090913
- ⁴ Centers for Disease Control. <https://www.cdc.gov/candida-auris/>
- ⁵ COPAN Liquid Amies Elution Swab (ESwab®) Collection and Preservation System Instructions for use; Version HPC031R08 Date 2021.10

Vivalytic C. auris – Instructions for use

Symbols

	Manufacturer		Serial number
	Date of manufacture		Temperature limit
	Expiry date		Do not use if package is damaged
	Lot number		For single use only
	Reference number		Consult Instructions for Use
	Contains <n> tests		<i>in vitro</i> diagnostic medical device

 ₀₁₂₃ CE mark

 Pipette the indicated sample volume in the sample input of the cartridge as marked by the black triangle.

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WARNING**Hazard components in cartridge**

(Binding buffer LB-05)
Guanidinium thiocyanate

Hazard statements

H302 Harmful if swallowed.
H319 Causes serious eye irritation.

Precautionary statements

P260 Do not breathe dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves and eye protection/face protection/ hearing protection.
P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.
P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER/doctor.

**Hazard statements**

P310 Causes severe skin burns and eye damage.

Performance Data**Table 1 Limit of Detection (LoD)**

Candida auris	933 copies/mL
---------------	---------------

*Determined by using Vircell material (Candida Auris Total Control (Neutral Matrix); MC234) for a spiking approach in Copan Liquid Amies Elution Swab (ESwab®) supplemented with approximately 1000 Human Bronchial/Tracheal Epithelial Cells (hPBTEC).

Table 2 Inclusivity

Candida auris

Strain validated for inclusivity via BLAST alignment

Vivalytic *C. auris* – Annex

Table 3 – Exclusivity

<i>Acinetobacter baumannii</i>	<i>Streptococcus dysgalactiae</i>
<i>Anaerococcus prevotii</i>	<i>Streptococcus mitis</i>
<i>Bacillus cereus</i>	<i>Streptococcus mutans</i>
<i>Bifidobacterium adolescentis</i>	<i>Vibrio cholerae</i>
<i>Campylobacter jejuni</i>	<i>Vibrio parahaemolyticus</i>
<i>Campylobacter rectus</i>	<i>Yersinia enterocolitica</i>
<i>Citrobacter freundii</i>	<i>Entamoeba histolytica</i>
<i>Citrobacter koseri</i>	<i>Giardia lamblia</i>
<i>Clostridioides difficile</i> (<i>Clostridium difficile</i>)	<i>Aspergillus niger</i>
<i>Corynebacterium diphtheriae</i>	<i>Candida albicans</i>
<i>Enterobacter aerogenes</i>	<i>Candida auris</i>
<i>Enterobacter cloacae</i>	<i>Candida dubliniensis</i>
<i>Enterococcus faecalis</i>	<i>Candida duobushaemulonii</i>
<i>Enterococcus faecium</i>	<i>Candida glabrata</i> (<i>Nakaseomyces glabratus</i>)
<i>Escherichia coli</i>	<i>Candida guilliermondii</i>
<i>Fusobacterium necrophorum</i>	(<i>Meyerozyma guilliermondii</i>)
<i>Helicobacter pylori</i>	<i>Candida haemulonii</i>
<i>Klebsiella oxytoca</i>	<i>Candida kefyr</i> (<i>Kluyveromyces marxianus</i>)
<i>Klebsiella pneumoniae</i>	<i>Candida krusei</i> (<i>Pichia kudriavzevii</i>)
<i>Listeria monocytogenes</i>	<i>Candida lipolytica</i>
<i>Moraxella catarrhalis</i>	<i>Candida lusitanae</i> (<i>Clavispora lusitanae</i>)
<i>Morganella morganii</i>	<i>Candida parapsilosis</i>
<i>Mycobacterium avium</i>	<i>Candida sake</i>
<i>Mycobacterium fortuitum</i>	<i>Candida tropicalis</i>
<i>Mycobacteroides abscessus</i>	<i>Fusarium moniliforme</i>
(<i>Mycobacterium abscessus</i>)	<i>Fusarium oxysporum</i>
<i>Mycobacteroides chelonae</i>	<i>Fusarium solani</i>
(<i>Mycobacterium chelonae</i>)	<i>Fusarium verticillioides</i>
<i>Nocardia farcinica</i>	<i>Hortaea werneckii</i>
<i>Peptostreptococcus anaerobius</i>	(<i>Phaeoannellomyces werneckii</i>)
<i>Plesiomonas shigelloides</i>	<i>Kodamaea ohmeri</i>
<i>Propionibacterium acnes</i>	<i>Penicillium citrinum</i>
<i>Proteus mirabilis</i>	<i>Penicillium marneffeii</i>
<i>Providencia alcalifaciens</i>	<i>Saccharomyces cerevisiae</i>
<i>Providencia stuartii</i>	<i>Trichophyton rubrum</i>
<i>Pseudomonas aeruginosa</i>	<i>Trichophyton tonsurans</i>
<i>Salmonella enterica</i>	Enterovirus 69 (Human enterovirus B)
<i>Shigella flexneri</i>	Enterovirus 70 (Human enterovirus D)
<i>Shigella sonnei</i>	Enterovirus 71 (Human enterovirus A)
<i>Staphylococcus aureus</i>	Human Adenovirus subgroup A serotype 31

Table 3 – Exclusivity

<i>Stenotrophomonas maltophilia</i>	Human Adenovirus subgroup F serotype 40
<i>Streptococcus agalactiae</i>	Human Adenovirus subgroup F serotype 41
<i>Streptococcus anginosus</i>	Norovirus genogroup II genotype 4 (GI.4)

Strains validated for exclusivity via BLAST alignment

Table 4 a Reproducibility

Target	c	No of total tests	No of pos. tests	No of neg. tests	Proportion of positive/negative results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
none	blank	108	0	108	100	96.6 % – 100 %	96.6 % – 100 %
<i>C. auris</i>	3x c95	108	107	1	99.1	94.9 % – 100 %	95.0 % – 100 %
	c100	108	108	0	100	96.6 % – 100 %	96.6 % – 100 %

3x c95 = 95 % predicted positive agreement
c100 = 100 % predicted positive agreement

Table 4 b Repeatability

Target	LOT	No of total tests	No of pos. tests (3x c95)	Proportion of positive/negative results (%)	95 % Wilson-Score confidence interval (%)	95 % Pearson-Clopper confidence interval (%)
<i>C. auris</i>	1	20	19	95	76.4 % – 99.1 %	75.1 % – 99.9 %
	2	20	20	100	83.9 % – 100 %	83.2 % – 100 %
	3	20	20	100	83.9 % – 100 %	83.2 % – 100 %
	Total	60	60	98.3	91.1 % – 99.7 %	91.1 % – 100 %

3x c95 = 95 % predicted positive agreement

Vivalytic *C. auris* – Annex

Table 5 – Tested Substances for Interference

No interference detected

Human whole blood; 1 % (v/v)

Mucin; 5 mg/mL

Mupirocin; 0.2 mg/mL

Nasivin® spray; 0.011 mg/mL

Otri-Allergie spray fluticason; 0.005 mg/mL

hysan® Salinspray; 1 % (v/v)

TobraZid®; 4 µg/mL

WICK VapoRub; 1 % (w/v)

Casein; 5 mg/mL

Cornstarch; 1.25 mg/mL

Desitin ointment; 3.5% (w/v)

Eucerin Anti-Transpirant 48 h Roll-on; 0.1 % (v/v)

Vichy DEO Roll-On Mineral 48h; 5 % (v/v)

Vichy DEO Roll-on anti-stain 48h; 5 % (v/v)

Widmer Deo Roll-on unperfumed; 0.1% (v/v)

Olive oil roll-on deodorant green tea; 1% (v/v)

HIDROFUGAL classic sprayer; 0.1 % (v/v)

LINOLA douche and wash; 7 % (v/v)

Glucose; 11 mg/mL

Lactate; 2.2 mg/mL

Interference was experimentally verified at 3xLoD for *Candida auris* using a spiking approach Copan Liquid Amies Elution Swab (ESwab®) supplemented with approximately 1000 Human Bronchial/Tracheal Epithelial Cells (hPBTEC).

Table 6a – Clinical Sensitivity (PPA) [1] and Specificity (NPA) [2] for clinical samples and isolates in comparison to MALDI-TOF (95 % confidence interval, clinical study)

[1] 98.36 % (91.20 – 99.96 %)

[2] 72.22 % (58.36 – 83.54 %)*

*15 samples were found false positive in comparison to MALDI-TOF testing, even though for at least 2 of these samples, initially performed pre-characterization using MALDI-TOF found the samples *Candida auris* positive.

In total, 61 positive and 54 negative samples were included in the dataset against MALDI-TOF as the gold standard method. Sensitivity (PPA) of 98.34% could be reached. Most of the positive results (55) were generated using *Candida auris* isolates. These isolates were collected during several outbreaks in Germany. All available *Candida auris* clades (I, II, III, IV) were included in the testing.

Table 6b – Clinical Sensitivity (PPA) [1] and Specificity (NPA) [2] for clinical samples and isolates in comparison to MALDI-TOF and additional qPCR test (95 % confidence interval, clinical study)

[1] 97.52 % (92.93 – 99.49 %)

[2] 94.79 % (88.26 – 98.29 %)

In total, 121 positive (16 positive clinical samples) and 96 negative samples could be included in the final dataset. Sensitivity (PPA) was determined to be 97.52 % and specificity (NPA) was 94.79 %.

Table 7 – Document History

Revision 01	Initial document
Revision 02	European languages have been added.