



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



ZYMO RESEARCH

DNA
Purification
Made Simple

Quick-DNA™ Fecal/Soil Microbe Midiprep Kit

DNA from fecal, soil, and microbial samples.

Highlights

- Simple, efficient isolation of humic-free, PCR-quality DNA from microbes including Gram positive and Gram negative bacteria, fungi, algae, protozoa, etc. in fecal and soil samples in as little as 25 minutes.
- State-of-the-art, ultra-high density **BashingBeads™** are fracture resistant and chemically inert.
- Omits the use of organic denaturants as well as proteinases.

Catalog Numbers:
D6110



Scan with your smart-phone camera to
view the online protocol/video.



tech@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682

Table of Contents

Product Contents	01
Specifications	02
Product Description	03
Protocol	06
Ordering Information	08
Guarantee	10

Product Contents

Quick-DNA™ Fecal/Soil Microbe Midiprep Kit	D6110 (25 Preps.)	Storage Temperature
ZR BashingBead™ Lysis/Filtration Tubes	25	Room Temp.
BashingBead™ Buffer	150 ml	Room Temp.
Genomic Lysis Buffer ¹	2 x 250 ml	Room Temp.
DNA Pre-Wash Buffer ²	15 ml	Room Temp.
g-DNA Wash Buffer	50 ml	Room Temp.
DNA Elution Buffer	16 ml	Room Temp.
Prep Solution	30 ml	Room Temp.
Zymo-Spin™ V-E Columns w/ Zymo-Midi Filters™	25	Room Temp.
Zymo-Spin™ III-HRC Filters	50	Room Temp.
Collection Tubes	100	Room Temp.
Instruction Manual	1	-

1 For optimal performance, add beta-mercaptoethanol to 0.5% (v/v) *i.e.*, 500 µl per 100 ml.


2 A precipitate may have formed in the DNA Pre-Wash Buffer during shipping. To completely resuspend the buffer, incubate the bottle at 30 – 37 °C for 30 minutes and mix by inversion. DO NOT MICROWAVE.

Specifications

- **Format** – Bead Beating, Spin/Vacuum Filtration, and Spin Column Purification.
- **Sample Sources** – Host, bacterial, fungal, algal, protozoan, viral DNA can be isolated from up to 375 mg of feces or up to 5 g of soil (2.5 g recommended). The amount of soil sample processed will vary depending on the composition of the sample: process more soil material for wet muddy samples and less for dry sandy samples. Additionally, DNA can be isolated directly from pelleted fungi and bacteria.
- **DNA Purity** – High quality, humic/fulvic-free DNA is eluted with **DNA Elution Buffer** making it perfect for PCR. ($A_{260}/A_{280} > 1.8$).
- **DNA Size Limits** – Capable of recovering genomic DNA up to and above 40 kb. In most instances, mitochondrial DNA and viral DNA (if present) will also be recovered.
- **DNA Recovery** – Typically, up to ~125 μg total DNA is eluted into $\geq 150 \mu\text{l}$ DNA Elution Buffer per sample.
- **Equipment** – Centrifuge, vacuum source and manifold, microcentrifuge, cell disrupter or pulverizer w/ 50 ml tube adapter


Product Description

The **Quick-DNA™ Fecal/Soil Microbe Midiprep Kit** is designed for the simple, rapid isolation of inhibitor-free, PCR-quality DNA from a variety of fecal (including humans, birds, rats, mice, cattle, etc.) and soil (including clay, sandy, silty, peaty, chalky, and loamy soils) samples. The procedure is easy and can be completed in as little as 25 minutes: fecal samples (≤ 375 mg each) or soil samples (≤ 5 g) are added directly to a **ZR BashingBead™ Lysis/Filtration Tube**, where microbes are rapidly and efficiently lysed by bead beating without the use of organic denaturants or proteinases. The DNA is then isolated and purified using our Zymo-Spin™ Technology, which is subsequently filtered to remove humic acids/polyphenols that inhibit PCR. The entire procedure can be performed in as little as 25 minutes, and there is no need for organic denaturants or proteinases. A schematic of the **Quick-DNA™ Fecal/Soil Microbe Midiprep Kit** procedure is shown below.




* Vacuum is the preferred method for DNA binding to the Zymo-Spin™ V-E column.

Fecal DNA Isolation





Comparison of DNA yields from rat feces using the **Quick-DNA™ Fecal/Soil Microbe Kit** and kits from suppliers A and B. Equivalent amounts of feces were processed using each kit and then equal volumes of eluted DNA were analyzed in a 0.8% (w/v) agarose/ethidium bromide gel. Samples were processed in triplicate.




PCR of DNAs from rat and human fecal samples isolated with the **Quick-DNA™ Fecal/Soil Microbe Kit**. Panels A and B show the results of PCR with DNA isolated from rat and human fecal samples, respectively, using primers specific for prokaryotic 16S rRNA. Panel C shows the results of PCR of DNA isolated from human feces with and without the addition of *E. coli* containing pUC19 plasmid DNA (indicated at the top of the image) using primers specific for the pUC19 sequence. In each case, amplicons were analyzed in a 1.5% (w/v) agarose / ethidium bromide gel using a UV imager. Numbers above each lane of the gel images are the volumetric equivalent (in μl) of eluted DNA (100 μl) used for PCR. Arrows mark the relative migration of amplicons in the gels, and M is a 100 bp DNA ladder (NEB).

Soil Microbe DNA Isolation



The **Quick-DNA™ Fecal/Soil Microbe Kit** can be used to isolate high quality DNA from a variety of soil types which yields robust products following PCR. **Panel A:** Physical characteristics of sampled soils (1-5) (Ref. 1). **Panel B:** Microbial DNA was isolated from soil samples (1-5) using the **Quick-DNA™ Fecal/Soil Microbe Kit**. Approximately 10% of the eluted DNA was then separated in a 0.8% (w/v) agarose/ethidium bromide gel. **Panels C and D** show the results of PCR of microbial DNA isolated from the samples with primers specific for prokaryotic 16S rRNA (**C**) or eukaryotic rRNA (**D**). In the figures, the 1 kb size marker (NEB) is as indicated and the arrows show the prokaryotic 16S rRNA and eukaryotic rRNA PCR products.



DNA isolated from *Saccharomyces cerevisiae* (strain TMY18) using the **Quick-DNA™ Fecal/Soil Microbe Kit** is high-quality and structurally intact. Equivalent amounts of yeast were processed using the **Quick-DNA™ Fecal/Soil Microbe Kit** or the kits from suppliers Q and E. Equal volumes of eluted DNA were then analyzed in a 0.8% (w/v) agarose/ethidium bromide gel. The size marker is a 1 kb ladder (NEB).

References:

1. Soil and Plant Laboratory, Inc. P.O. Box 11744, Santa Ana, California 92711

Protocol

For optimal performance, add beta-mercaptoethanol (user supplied) to the **Genomic Lysis Buffer** to a final dilution of 0.5% (v/v) *i.e.*, 2.5 ml per 500 ml.

1. Add 2.5 grams (5 g max.)¹ of soil sample or up to 375 mg of fecal samples to the bead/filter chamber of a **ZR BashingBead™ Lysis/Filtration Tube**. Add 6 ml **BashingBead™ Buffer** to the sample, cap tube tightly, and process.

Note: To prevent the **BashingBead™ Buffer** from leaking into the bottom of the 50 ml tube, place the **ZR BashingBead™ Lysis/Filtration Tube** on its side prior to processing. Alternatively, add 250-500 mg (wet weight) fungal and/or bacterial cells that have been resuspended in 6 ml of **BashingBead™ Buffer** to a **ZR BashingBead™ Lysis/Filtration Tube**.

2. Secure in a bead beater fitted with a 50 ml tube holder assembly (see page 9) to process samples. Optimization of processing time/speed will be necessary for complete sample lysis.

Note: Processing times may be as little as one minute when using high-speed bead beaters (*e.g.*, 2010 GenoGrinder®). See manufacturer's literature for specific operating information.

3. Centrifuge the **ZR BashingBead™ Lysis/Filtration Tube** in a centrifuge at $\geq 3,000 \times g$ (5,000 $\times g$ max.) for 5 minutes.
4. Remove bead/filter chamber from the top of the **ZR BashingBead™ Lysis/Filtration Tube** and transfer supernatant² from the bottom of the tube to a clean 50 ml tube (not provided).
5. Add 18 ml of **Genomic Lysis Buffer** to the supernatant. Mix well.
6. Filter the entire mixture from Step 5 using a **Zymo-Spin™ V-E Column/Zymo-Midi Filter™** assembly mounted on a vacuum manifold³ (see diagram on page 3) with a vacuum source set at ≥ 600 mm Hg.

¹ Although 2.5 g is recommended for most applications, the amount of sample will vary depending on its composition: process more material for wet muddy samples and less for dry sandy samples.

² Be careful to avoid the pelleted material at the bottom of the tube when transferring the supernatant.

³ Alternatively, the **Zymo-Spin™ V-E Column/Zymo-Midi Filter™** assembly can be placed in a 50 ml tube and centrifuged at 2,000 $\times g$ max. for 5 minutes. Filtration of the entire mixture will require several spins. Empty the flow through from the tube after each spin. **CAUTION:** Make sure the connection between the column and filter is secure (finger tight) prior to centrifugation!

7. Disconnect the Zymo-Spin™ V-E Column/Zymo-Midi Filter™ assembly and transfer the **Zymo-Spin™ V-E Column** to a **Collection Tube**. Spin the column at 10,000 x *g* for 1 minute in a microcentrifuge⁴, then add 300 µl **DNA Pre-Wash Buffer** to the column and spin at 10,000 x *g* for 1 minute. Discard the flow through.
 8. Add 400 µl **g-DNA Wash Buffer** to the column and centrifuge at 10,000 x *g* for 1 minute. Discard flow through and repeat wash step.
 9. Transfer the **Zymo-Spin™ V-E Column** to a 1.5 ml microcentrifuge tube and add 150 µl **DNA Elution Buffer** directly to the column matrix⁵. Wait for 1 minute and then centrifuge at 10,000 x *g* for 1 minute to elute the DNA⁶.
- Note:** If fungal or bacterial cultures were sampled, the DNA is now suitable for PCR as well as other downstream applications.
10. Place a **Zymo-Spin™ III-HRC Filter** in a clean Collection Tube and add 600 µl **Prep Solution**. Centrifuge at 8,000 x *g* for 3 minutes.
 11. Transfer the eluted DNA to a prepared Zymo-Spin™ III-HRC Filter in a clean 1.5 ml microcentrifuge tube and centrifuge at exactly 16,000 x *g* for 3 minutes. The filtered DNA is now suitable for PCR and other downstream applications.

⁴ Leave the rotor cover off the microcentrifuge if clearance with the column tops is a problem.

⁵ DNA yields can be increased by performing a second elution and pooling the eluates.

⁶ In some cases a brown-colored pellet may form at the bottom of the tube after centrifugation. Avoid this pellet when collecting the eluted DNA.

Ordering Information

Product Description	Catalog No.	Size
Quick-DNA™ Fecal/Soil Microbe Microprep Kit	D6012	50 Preps.
Quick-DNA™ Fecal/Soil Microbe Miniprep Kit	D6010	50 Preps.
Quick-DNA™ Fecal/Soil Microbe Midiprep Kit	D6110	25 Preps.
Quick-DNA™ Fecal/Soil 96 Kit	D6011	2 x 96 Preps.

Individual Kit Components	Catalog No.	Amount
ZR BashingBead™ Lysis/Filtration Tubes (50 ml) w/ 0.5 mm Beads	S6010	25 Pack
Genomic Lysis Buffer	D3004-1-100	100 ml
BashingBead™ Buffer	D6001-3-150	150 ml
DNA Pre-Wash Buffer	D3004-5-15	15 ml
g-DNA Wash Buffer	D3004-2-50	50 ml
DNA Elution Buffer	D3004-4-16	16 ml
Prep Solution	D6035-1-30	30 ml
Zymo-Spin™ V-E Columns w/ Zymo-Midi Filters™	C1021-25	25 Pack
OneStep™ PCR Inhibitor Removal Kit	D6030	50 Preps.
Collection Tubes	C1001-50 C1001-500 C1001-1000	50 Pack 500 Pack 1,000 Pack

Lysis Instruments	Catalog No.	Amount
FastPrep™ 24	S6005	1 Unit
BigPrep™ Adapter (2 x 50ml Tubes) for FastPrep-24	S6005-4	1 Unit
2010 Geno/Grinder	S6006	1 Unit
50 ml Tube Holder/Cryo Block Assembly for Geno/Grinder	S6006-3	2 Blocks



**100% satisfaction guarantee on all Zymo Research products,
or your money back.**

Zymo Research is committed to simplifying your research with quality products and services. If you are dissatisfied with this product for any reason, please call 1(888) 882-9682.

Integrity of kit components is guaranteed for up to one year from date of purchase.
Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

TM Trademarks of Zymo Research Corporation
Other trademarks: Vortex Genie[®] (Scientific Industries, Inc.), Precellys[®] (Bertin Instruments), Fastprep-24TM (MP Biomedicals), Bead Ruptor[®] (Omni International).



ZYMO RESEARCH

The **BEAUTY** of **SCIENCE** is to Make Things **SIMPLE**®



tech@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682