



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product Information

Nerve Terminal Staining Kit II (A)

Catalog Number: 70031

Kit Contents

Component	Size
70024 AM1-43	1 mg
70029-1 ADVASEP-7	100 mg

Storage and Handling

Store desiccated at 4°C or below. Protect AM1-43 from light, especially in solution. Components are stable for at least 12 months from date of receipt when stored as recommended. AM1-43 and ADVASEP-7 are soluble in water. To prepare 10 mM AM1-43, dissolve 1 mg in 179 μ L dH₂O. To prepare 100 mM ADVASEP-7, dissolve 100 mg in 462 μ L dH₂O. Stock solutions can be stored at 4°C or -20°C for six months or longer.

Molecular information

AM1-43

MW: 559.5

Formula: C₂₉H₄₉Cl₃N₄

See Fig. 1

ADVASEP-7

Average MW: 2163

Formula: C₄₂H_{70-n}O₃₅(C₄H₈SO₃Na)_n•hydrate where n = Average Degree of Substitution 6.5

Spectral Properties

AM1-43

Abs/Em: 510/625 nm (in MeOH); 480/598 nm (in membranes) (Fig. 2)




Figure 1. Structure of AM1-43.




Figure 2. Absorption and emission spectra of SynaptoGreen C4 (also known as FM1-43) in liposomes. The spectra for AM1-43, AM1-44, AM2-10 and other SynaptoGreen dyes are similar.

Product Description

Nerve terminal probes are a series of fluorescent cationic styryl dyes developed to follow synaptic activities at neuromuscular junctions or synapses. These dyes typically have a lipophilic tail (two carbon chains) at one end and a highly hydrophilic, cationically charged head group at the other end. These nerve terminal probes were originally called FM® dyes, and are available from Biotium under the trademark names of SynaptoGreen™ and SynaptoRed™. SynaptoGreen probes are dyes with a single double bond (n = 1) while SynaptoRed probes are dyes with three double bonds (n = 3). A nerve terminal dye is named as either SynaptoGreen or SynaptoRed followed by a carbon number designating the length of the lipophilic tail.

Cationic styryl dyes are believed to function by staining synaptic vesicles in an activity-dependent fashion. In the presence of cells or tissue preparations, the dyes partition between the aqueous phase, where the dyes are virtually non-fluorescent, and the outer leaflet of the cell surface membranes, where the dyes insert the lipophilic end into the membranes and become intensely fluorescent. During endocytosis following nerve stimulation, the dyes become trapped inside the vesicles. Thus, after washing off the dyes on the cell surface, the fluorescent signal is proportional to the number of newly formed vesicles. On the other hand, during exocytosis, the dyes are released from the vesicles along with neurotransmitters, causing a decrease in fluorescent signal. As a result, the change in fluorescent intensity reflects the amount of endocytosis/exocytosis or synaptic activity. The rate of fluorescence increase during endocytosis, the "on-rate", and the rate of fluorescence decrease during exocytosis, the "off-rate", vary from dye to dye. In general, dyes with longer lipophilic tails and more double bonds have a higher affinity toward membrane and thus a higher on-rate and lower off-rate.

Biotium developed AM1-43 (1) as a fixable version of SynaptoGreen C4 (FM®1-43). AM1-43 has the same absorption/emission spectra and the same lipophilic tail as SynaptoGreen C4 does except that the former has an additional amine group that makes the dye aldehyde-fixable.

When using nerve terminal dyes, one frequent problem researchers encounter is background fluorescence due to nonspecific membrane staining. Although most of the background fluorescence can be removed by repeated washing, the problem is still significant with dyes that have a longer tail or more double bonds, particularly when the dyes are used in tissue preparations. Washing with sulfonated beta-cyclodextrin ADVASEP-7 helps to remove nerve terminal dyes from the cell surface and reduces background (2).

Biotium offers additional nerve terminal staining kits with other pairings of nerve terminal dyes and background reducing agents (see Related Products).

References

1. Renger, JJ, et al. Neuron 29, 469(2001).
2. Kay, AR, et al. Neuron 24, 809 (1999).

Continued next page

Assay Protocol

The following is an example of a protocol for nerve terminal staining of cultured neurons on coverslips. Nerve terminal dyes also can be used to label endocytic vesicles in non-neuronal cell types. Staining can be performed at 4°C for selective labeling of the plasma membrane; at room temperature or 37°C, endocytosis of the dye generally occurs within 10 minutes. Buffers other than Tyrode solution may be used. The addition of the sodium channel blocker tetrodotoxin (TTX) is optional, its purpose is to block action potentials and prevent synaptic vesicle release after staining. Optimal protocols for specific applications may need to be determined by the user; see reference 2 for examples of protocols for staining brain slices and other tissue samples.

1. Dilute AM1-43 stock solution to a final concentration of 4 µM in 50 mM Tyrode solution (for example, 1 µL 10 mM dye per 2.5 mL solution). Place the coverslip with your cells in this solution for 1 minute at room temperature. Use enough solution to completely submerge the cells.
2. Transfer the coverslip to Tyrode + 0.5 µM tetrodotoxin (TTX, catalog no. 00061) solution for 1 minute at room temperature.
3. Wash for 5 minutes at room temperature with Tyrode + 0.5 µM TTX + ADVASEP-7. The typical concentration of ADVASEP-7 working solution is 1 mM (for example, 10 µL of 100 mM ADVASEP-7 per mL solution).
4. Rinse the coverslip twice with Tyrode + 0.5 µM TTX.
5. Transfer the coverslip to fixation solution (4% formaldehyde, 4% sucrose, 1 µM TTX in PBS) for 20 minutes at room temperature.
6. Transfer the coverslip directly to pre-chilled 0.01% Triton X-100 in PBS for 12 minutes at 4°C.
7. Wash 3 x 1 minute with cold PBS.
8. Stain with primary antibody in 10% serum/PBS for 3 hours at 4°C. The concentration of the antibody should be double of what you would use for regular immunofluorescence staining.
9. Wash 3 x 1 minute with cold PBS.
10. Stain the preparation with secondary antibody at the concentration you would normally use for immunofluorescence staining in 2% serum/PBS for 40 minutes at 4°C.
11. Wash as in step 8.
12. Mount the coverslip in PBS and image.

Related Products

Catalog number	Product
70042	SynaptoGreen™ C1
70044	SynaptoGreen™ C2 (equivalent to FM®2-10)
70023	SynaptoGreen™ C3
70020	SynaptoGreen™ C4 (equivalent to FM®1-43)
70046	SynaptoGreen™ C5 (equivalent to FM®4-84)
70048	SynaptoGreen™ C18 (equivalent to FM®1-84)
70040	SynaptoRed™ C1
70021	SynaptoRed™ C2 (equivalent to FM®4-64)
70028	SynaptoRed™ C2M (equivalent to FM®5-95)
70024	AM1-43
70038	AM1-44
70036	AM2-10
70051	AM3-25
70025	AM4-64
70039	AM4-65
70050	AM4-66
70029	ADVASEP-7
70037	SCAS
80101	Sulforhodamine 101
70030	Nerve Terminal Staining Kit I 5 x 1 mg SynaptoGreen™ C4 and 250 mg ADVASEP-7
70031-1	Nerve Terminal Staining Kit II (B) 1 mg AM1-43 and 100 mg SCAS
70032	Nerve Terminal Staining Kit III 5 x 1 mg SynaptoGreen™ C4 and 100 mg Sulforhodamine 101
70034	Nerve Terminal Staining Kit V 5 x 1 mg SynaptoRed™ C2 and 250 mg ADVASEP-7
00060	Tetrodotoxin, citrate-free
00061	Tetrodotoxin, with citrate
00010	α-Bungarotoxin
00019	β-Bungarotoxin

Please visit our website at www.biotium.com for information on our life science research products, including fluorescent CF™ dye bungarotoxins, antibodies, and other conjugates, calcium and other ion indicator dyes, apoptosis detection reagents, and other fluorescent probes and kits for cell biology research.

FM is a registered trademark of Thermo Fisher Scientific.

Materials from Biotium are sold for research use only, and are not intended for food, drug, household, or cosmetic use.