



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Anti-Fluoroquinolone [F9] Bulk Size Ab04000-30.11-BT

This is an scFv fragment with a His tag.

This is a reformatted human scFv antibody, based on the original human scFv format, created for improved compatibility with existing reagents, assays and techniques.

Isotype and Format: scFv fragment (His), ScFv

Clone Number: F9

Alternative Name(s) of Target: Quinoline monocarboxylic acid; antibiotic

UniProt Accession Number of Target Protein:

Published Application(s): ELISA

Published Species Reactivity: Species independent

Immunogen: The original antibody was isolated from a human scFv phage display library by panning against biotinylated fluoroquinolone antibiotic.

Specificity: This antibody is a general fluoroquinolone binder that likely binds to the antigen at the common pyridone-carboxylic acid pharmacophore.

Application Notes: This antibody was capable of binding the biotin analogue of fluoroquinolone in an ELISA. The binding affinity of the original scFv antibody to norfloxacin-biotin analogue was reported to be $K_d = 1.58 \mu\text{M}$. This antibody was used in the development of an ScFv biosensor for qualitative detection of antibiotics in livestock drinking water (PMID: 21142055),

Antibody First Published in: Young Cha et al. Pharmacophore-based strategy for the development of general and specific scFv biosensors for abused antibiotics. *Bioconjug Chem.* 2011 Jan 19;22(1):88-94.

[PMID:21142055](#)

Note on publication: Describes the development of a fluorescent biosensor system using antibodies for detection of fluoroquinolone antibiotics.

Product Form

Size: 500 µg Purified antibody in bulk size.

Purification: Purified by Immobilized Metal Affinity Chromatography

Supplied In: PBS only.

Storage Recommendation: Store at 4°C for up to 3 months. Note, this antibody is provided without added preservatives, it is therefore recommended this antibody be handled under sterile conditions. For longer

storage, aliquot and store at -20°C.

Concentration: 1 mg/ml.

Important note - This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals.