

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

## Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

#### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in







# Anti-Envelope glycoprotein H; Envelope glycoprotein L [Fab-RC] Standard Size Ab04452-23.0

Isotype and Format: Rabbit IgG, Kappa

Clone Number: Fab-RC

Alternative Name(s) of Target: gH; gL; ORF60; ORF37

UniProt Accession Number of Target Protein: Q775J3; Q9J3N1 Published Application(s): crystallography, IP, neutralization, SPR, IHC

Published Species Reactivity: Varicella-zoster virus

**Immunogen:** The original antibody was raised from the white blood cells of a human individual immunized

with zoster vaccine.

**Specificity:** This antibody is specific for the combination of gH and gL and does not detect the proteins in isolation.

**Application Notes:** The original antibody was used for surface plasmon resonance on VZV gHgL. This antibody was used for a neutralization assay. Here is showed it was able to inhibit gB/gHgL-mediated membrane fusion. The structure of this antibody was determined using x-ray crystallography (Xing et al., 2015; PMID:25918416). This antibody was used for immunohistochemistry on fixed and unfixed VZV-infected and control HFL cells. This antibody was also used for a plaque reduction assay on VZV-infected HFL cells. Further this antibody was used for immunoprecipitation on VZV-infected and uninfected HFL cells (Birlea et al., 2013; PMID:23077312).

**Antibody First Published in:** Birlea et al. Human Anti-Varicella-Zoster Virus (VZV) Recombinant Monoclonal Antibody Produced after Zostavax Immunization Recognizes the gH/gL Complex and Neutralizes VZV Infection J Virol. 2013 Jan; 87(1): 415–421. PMID:23077312

**Note on publication:** One week after immunizing a 59-year-old VZV-seropositive man with Zostavax, the original paper sorted his circulating blood plasma blasts and amplified expressed immunoglobulin variable domain sequences by single-cell PCR.

#### **Product Form**

**Size:** 100 µg Purified antibody.

**Purification:** Protein A affinity purified **Supplied In:** PBS with 0.02% Proclin 300.

Storage Recommendation: Store at 4°C for up to 3 months. For longer storage, aliquot and store at -

