



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### EGFR/CEN 7 CISH Probe

**Catalog Number:** CG0034

**Regulatory Status:** For research use only (RUO)

**Product Description:** EGFR/CEN 7 CISH Probe is designed for the qualitative detection of human EGFR gene amplifications as well as the detection of chromosome 7 alpha satellites in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH).

**Interpretation of Result:**

**Probe Position:**

**Regulatory Status:** For research use only (RUO)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** The product is ready-to-use. No reconstitution, mixing, or dilution is required. Bring probe to room temperature (18-25°C) and mix briefly before use.

**Supplied Product:** Reagent Provided:

This Probe is composed of:

1. Digoxigenin-labeled polynucleotides, which target sequences mapping in 7p11.2\* (chr7:55,082,262-55,278,647) harboring the EGFR gene region.
2. Dinitrophenyl-labeled polynucleotides, which target sequences mapping in 7p11.1-q11.1 specific for the alpha satellite centromeric region D7Z1 of chromosome 7.
3. Formamide based hybridization buffer.

**Storage Instruction:** Store at 2-8°C in an upright position. Return to storage conditions immediately after use.