



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### DNA (+) Control CISH Probe

**Catalog Number:** CO0008

**Regulatory Status:** For research use only (RUO)

**Product Description:** DNA (+) Control CISH Probe is designed to be used for the detection of human Alu repetitive sequences in formalin-fixed, paraffin-embedded tissue or cells by chromogenic *in situ* hybridization (CISH).

**Regulatory Status:** For research use only (RUO)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** Bring probe to hybridization temperature before use.

**Supplied Product:** Reagent Provided:

Digoxigenin-labeled oligonucleotides with GC contents of 40-70% without known consensus to any naturally occurring sequences.

**Storage Instruction:** Store at 2-8°C in an upright position. Return to storage conditions immediately after use.