



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### CytoChipNano Combo Pack

**Catalog Number:** DM0002

**Regulatory Status:** For research use only (RUO)

**Product Description:** CytoChipNano Combo Pack contains CytoChipNano, Plugs and SCx™ Spiral Chamber for separation of CRC on CytoQuest™ CR.

**Instrument Requirement:** [CytoQuest™ CR](#)

**Regulatory Status:** CE-IVD  
(EDMA Code: 21.09.10.01.00)

**Applications:** IF-CTC  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Supplied Product:** This pack contains 20 combos. Each combo includes:  
CytoChipNano  
SCx™ Spiral Chamber  
Plugs

\*1 Clean Chip is included for both Prime and Clean purposes.

\*1 Dummy Chip is included for Tension Alignment purpose if needed.

**Storage Instruction:** Store CytoChipNano at 2-8°C.  
Leave CytoChipNano in their sealed pouches until use.  
Store Clean Chip, Dummy Chip, SCx™ Spiral Chamber and Plugs at room temperature.