



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### CEN5p (Cy5) FISH Probe

**Catalog Number:** FC0040

**Regulatory Status:** For research use only (RUO)

**Product Description:** Labeled FISH probes as chromosomal markers using Fluorescent In Situ Hybridization Technique. ([Technology](#))

**Applications:** FISH-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Recommend Usage:** The probe is provided in 5x concentrated format, to allow mixing of up to 5 chromosome FISH probes in a single hybridization assay. When used alone, it should be diluted to 1x with FISH Hybridization Buffer (Cat # [U0028](#) or [U0029](#)) before use.

**Supplied Product:** FISH Hybridization Buffer (80 uL)

**Storage Instruction:** Store at 4°C in the dark.

#### References:

1. SUN2 Modulates HIV-1 Infection and Latency through Association with Lamin A/C To Maintain the Repressive Chromatin. Sun WW, Jiao S, Sun L, Zhou Z, Jin X, Wang JH. MBio. 2018 May 1;9(3). pii: e02408-17.