



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

mutaFISH™ CEP1wt CEP7wt CEP8wt CEP17wt DNA Probes

Catalog Number: FP0009

Regulatory Status: For research use only (RUO)

Product Description: mutaFISH™ CEP1wt CEP7wt CEP8wt CEP17wt DNA Probes is designed to identify human CEP1, CEP7, CEP8, CEP17 amplification on dsDNA in cells using padlock probe and *in situ* rolling-circle amplification technology.

Applications: mutaFISH-Ce
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Supplied Product: Content:

1. mutaFISH™ CEP1wt DNA Probe
2. mutaFISH™ CEP7wt DNA Probe
3. mutaFISH™ CEP8wt DNA Probe
4. mutaFISH™ CEP17wt DNA Probe
5. Detection Probe-FITC
6. Detection Probe-6-HEX
7. Detection Probe-Texas Red X
8. Detection Probe-Aqua 431

Storage Instruction: Store at -20°C.
Aliquot to avoid repeated freezing and thawing.