



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

ERBB2/CEN 17 CISH Probe Kit

Catalog Number: KA6901

Regulation Status: For research use only (RUO)

Product Description: ERBB2/CEN 17 CISH Probe Kit is designed for the qualitative detection of amplifications involving the human ERBB2 gene as well as the detection of chromosome 17 alpha satellites in formalin-fixed, paraffin-embedded specimens by chromogenic *in situ* hybridization (CISH).

Interpretation of Result:

Probe Position:

Regulatory Status: For research use only (RUO)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Liquid

Recommend Usage: 20x Wash Buffer TBS (WB5) is to be prepared according to the instructions in Assay procedure. All other kit reagents are ready-to-use. No reconstitution, mixing, or dilution is required. Bring probe to room temperature (18-25°C) and mix briefly before use.

Supplied Product: Reagent Provided:

This Probe is composed of:

1. Digoxigenin-labeled polynucleotides, which target sequences mapping in 17q12 (chr17:37,725,661-37,973,541) harboring the ERBB2 gene region.
2. Dinitrophenyl-labeled polynucleotides, which target sequences mapping in 17p11.1-q11.1 specific for the alpha satellite centromeric region D17Z1 of chromosome 17.
3. Formamide based hybridization buffer.

Storage Instruction: Store at 2-8°C in an upright position. Return to storage conditions immediately after use.