



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### AP-Red Solution Set

**Catalog Number:** KA6902

**Regulation Status:** For research use only (RUO)

**Product Description:** AP-Red Solution Set is intended to be used as a substrate for an AP-conjugated antibody in chromogenic *in situ* hybridization (CISH) applications.

**Regulatory Status:** For research use only (RUO)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Recommend Usage:** Refer to the instructions for use of the respective CISH probe and implementation kit.

**Storage Instruction:** Store at 2-8°C in an upright position. Return to storage conditions immediately after use.