



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

AAV9 Titration ELISA Kit

Catalog Number: KA6921

Regulation Status: For research use only (RUO)

Product Description: AAV9 Titration ELISA Kit is a sandwich ELISA to determine the titer of AAV9 capsids in the test sample.

Calibration Range: 2.03×10^7 - 1.30×10^9 capsids/mL

Limit of Detection: 1.02×10^7 capsids/mL

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Instruction: Store kit components at 2-8°C for 12 months from the production date. After reconstitution, store AAV9 standard at 2-8°C for 2 weeks; for long-term storage, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.