



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### AAV5 Titration ELISA Kit

**Catalog Number:** KA6922

**Regulation Status:** For research use only (RUO)

**Product Description:** AAV5 Titration ELISA Kit is a sandwich ELISA to determine the titer of AAV5 capsids in the test sample.

**Calibration Range:**  $3.13 \times 10^7$  -  $2 \times 10^9$  capsids/mL

**Limit of Detection:**  $1.56 \times 10^7$  capsids/mL

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Instruction:** Store kit components at 2-8°C for 12 months from the production date. After reconstitution, store AAV5 standard at 2-8°C for 2 weeks; for long-term storage, store at -20°C or lower.

Aliquot to avoid repeated freezing and thawing.