



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### FITC Labeling Kit

**Catalog Number:** KA6932

**Regulation Status:** For research use only (RUO)

**Product Description:** FITC Labeling Kit is designed for production scale, provide a convenient and efficient method for labeling large volumes of antibodies with iFluor dyes, Alexa Fluor, and various other labels.

**Suitable Sample:** Monoclonal antibody, Polyclonal antibody, Proteins

**Excitation (Max):** 491

**Emission (Max):** 516

**Regulation Status:** For research use only (RUO)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Instruction:** Store FITC (Component A) at -20°C, and avoid from light. Store other components at 4°C, and Reaction Buffer (Component B) avoid from light.