



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### Inorganic Pyrophosphatase (Yeast)

**Catalog Number:** P9372

**Regulatory Status:** For research use only (RUO)

**Product Description:** Inorganic Pyrophosphatase (Yeast) is derived from a recombinant *E.Coli* strain carrying *Saccharomyces cerevisiae* ppa gene. The enzyme catalyzes the hydrolysis of inorganic pyrophosphate to form orthophosphate.  $P_2O_7^{-4} + H_2O + PPase \rightarrow 2HPO_4^{-2}$  Inorganic Pyrophosphatase (yeast) could hydrolyze the inorganic pyrophosphate generated by nucleic acid amplification, in vitro transcription and other reactions, removing the inhibition of the inorganic pyrophosphate generated on the reaction system. The removal of pyrophosphate can shift the reaction equilibrium to the end of product formation, which is beneficial to increase the yield of synthetic product. RNA and DNA synthesis are examples of reactions that can be pulled far in the synthesis direction by the action of inorganic pyrophosphatase.

**Theoretical MW (kDa):** 71

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Concentration:** 0.1 U/?L

**Recommend Usage:** In Vitro Transcription  
The optimal working dilution should be determined by the end user.

**Storage Buffer:** 20 mM Tris-HCl (pH 8.0), 100 mM KCl, 0.1 mM EDTA, 1 mM dithiothreitol, 50% glycerol

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.