



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### SETD7 recombinant monoclonal antibody, clone RAB-C220

**Catalog Number:** RAB03843

**Regulatory Status:** For research use only (RUO)

**Product Description:** Rabbit recombinant monoclonal antibody raised against recombinant antibody (rAb) phage display recognizing SETD7 protein under non-denaturing conditions.

**Clone Name:** RAB-C220

**Immunogen:** Original antibody is raised against recombinant protein corresponding to recombinant antibody (rAb) phage display recognizing SETD7 protein under non-denaturing conditions.

**Antibody Species:** Rabbit

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Conjugation:** Unconjugated

**Concentration:** batch dependent

**Isotype:** IgG

**Recommend Usage:** ChIP

ELISA

Flow cytometry

Immunofluorescence

Immunoprecipitation-MS

spiked-Immunoprecipitation

Western Blot

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS with 0.02% Proclin 300

**Storage Instruction:** Store at 4°C for 3 months. For long term storage store at -20°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 80854

**Gene Symbol:** SETD7

**Gene Alias:** FLJ21193, KIAA1717, KMT7, SET7, SET7/9, SET9