



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### **Gallbladder Tumor (adenocarcinoma, poorly differentiated) FFPE tissue slide**

**Catalog Number:** S1009

**Regulation Status:** For research use only (RUO)

**Product Description:** Gallbladder Tumor  
(adenocarcinoma, poorly differentiated) FFPE tissue  
slide contains tissue section from a single human donor.

**Tissue:** Gallbladder

**Host:** Human

**Applications:** IHC-P, Profiling  
(See our web site product page for detailed applications  
information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product  
page for detailed protocols

**Form:** Tissue sections are 5  $\mu$ m in thickness, and  
mounted on positively charged glass slides.

**Recommend Usage:** Please bake slides at 60°C for 30  
minutes before use.

**Storage Instruction:** Store at 4°C.