



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Datasheet

### Frontal Lobe (Alzheimer's disease) FFPE tissue slide

**Catalog Number:** S1207

**Regulation Status:** For research use only (RUO)

**Product Description:** Frontal Lobe (Alzheimer's disease) FFPE tissue slide contains tissue section from a single human donor.

**Tissue:** Brain: Frontal lobe

**Host:** Human

**Applications:** IHC-P, Profiling  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Tissue sections are 5 um in thickness, and mounted on positively charged glass slides.

**Recommend Usage:** Please bake slides at 60°C for 30 minutes before use.

**Storage Instruction:** Store at 4°C.