



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Temporal Lobe (progressive supranuclear palsy) FFPE tissue slide

Catalog Number: S1257

Regulation Status: For research use only (RUO)

Product Description: Temporal Lobe (progressive supranuclear palsy) FFPE tissue slide contains tissue section from a single human donor.

Tissue: Brain: Temporal lobe

Host: Human

Applications: IHC-P, Profiling

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Form: Tissue sections are 5 μ m in thickness, and mounted on positively charged glass slides.

Recommend Usage: Please bake slides at 60°C for 30 minutes before use.

Storage Instruction: Store at 4°C.