



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Datasheet

Protein Staining Solution

Catalog Number: TM0008

Regulation Status: For research use only (RUO)

Product Description: Protein Staining Solution stains proteins in polyacrylamide gel after SDS gel electrophoresis (detection limit \geq 8 ng of protein) with *CoomassieR-G250. All processes including washing and destaining processes can be performed in 1.5 hr. After staining of gel (20 min), proteins can be seen before destaining. Overnight staining and overnight destaining is allowable.

Applications: Electrophoresis
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Instruction: Store at 4°C.